



**CIDA DIGITAL THERMOMETER**

Accurately records body temperature in 60 seconds. Liquid crystal display readout displays temperature within 1/10th of a degree. Audible beep signals peak temperature. A replaceable battery provides 200 hours of use. Each package comes complete with thermometer, 5 probe covers, battery, storage case, instructions, and one year limited warranty. Can be used orally, rectally or under the arm. Latex free.

52-0410 Digital Thermometer w/beeper, °F

**ADC® ADTEMP™ I DIGITAL THERMOMETER**

An electronic thermometer designed and priced to be a viable replacement for mercury filled glass thermometers, which are gradually being phased out of healthcare. This new model features a striking translucent housing with integral case. It is dual-scale (C/F) switchable, water resistant and offers a last reading memory. It does not have a beeper, but it will lock in the reading once measurement is complete. This makes an ideal item as a true mercury thermometer replacement. Features a storage case, long life battery and auto-off feature.

412-00 Single Patient Use Digital Thermometer, Bulk, 20/pkg  
 412K Single Patient Use Digital Thermometer SPU Kit

**ADC® ADTEMP™ II DIGITAL THERMOMETER**

Features a water resistant design, last reading memory and is dual scale (C/F) switchable. Replaces the 410 and 6001 series thermometers.

413B Oral Digital Thermometer  
 413B-00 Oral Digital Thermometer, 25/pkg  
 413BK Oral Digital Thermometer SPU Kit  
 413R Rectal Digital Thermometer  
 413R-00 Rectal Digital Thermometer, 25/pkg  
 413RK Rectal Digital Thermometer SPU Kit

**ADC® ADTEMP™ III DIGITAL THERMOMETER**

An electronic thermometer with a flexible probe designed and priced to be a viable replacement for mercury filled glass thermometers, which are gradually being phased out of healthcare. The flexible probe provides greater patient safety and is ideal for pediatric applications. This new model features a striking translucent housing with integral case. It is dual-scale (C/F) switchable, water resistant and offers a last reading memory. Beeps when measurement is complete. Features a storage case, long life battery and auto-off feature.

414 Digital Thermometer Kit, Flex Tip, Retail Boxed  
 414K Digital Thermometer Kit, Flex Tip, Boxed w/5 sheaths


**ADC® ADTEMP™ IV DIGITAL THERMOMETER**

415 Digital Thermometer  
 415K Digital Thermometer SPU Kit

**ADC® ADTEMP™ V DIGITAL THERMOMETER**

Represents state-of-the-art in micro digital fever thermometers. This new unit features a stylish European design, extra large display and flexible tip. It is water resistant, dual scale (C/F) switchable and has a last reading memory. The response time is 10 seconds rectal and 18 seconds oral. Beeps when measurement is complete. Features a storage case, long life battery and auto-off feature.

418 Super Fast Read Digital Thermometer, Flex Tip, Boxed w/5 sheaths



**ADC® ADTEMP™ THERMOMETER SHEATHS**

Disposable design eliminates cross contamination. Cost effective. Easy to apply. Suitable for oral or rectal use. Digital model designed for use with ADC® and most comparably sized digital thermometers. Glass model designed for use on all glass thermometers. Made in the U.S.A.

416-100 ADTEMP™ Digital Sheaths, 100/bx  
 416-1000 ADTEMP™ Digital Sheaths, 1000/bx  
 416-50 ADTEMP™ Digital Sheaths, 50/bx

**BD THERMOMETERS**

522862 Basic Digital Thermometer, 12/bx

**Accessories**

522870 Probe Covers For Oral Thermometer, 50/sp, 20 sp/cs  
 522871 Probe Covers For Rectal Thermometer, 50/sp, 20 sp/cs



**OMRON COMPACT DIGITAL THERMOMETER**

No fear of broken glass or mercury poisoning. Fast, accurate reading in 60 seconds. Premium long life battery (up to 400 hours of continuous use). Waterproof for easy cleaning.

MC-102 Digital Thermometer (Spanish/Portuguese) Centigrade  
 MC-110B Digital Rectal Thermometer, 20/bx



MC-130MO



MC-130MR

**OMRON DIGITAL FEVER THERMOMETER**

Oral or rectal available. Quick accurate temperature in 60 seconds. Beeps at peak temperature. Memory recalls last temperature. Meets ASTM (E1112) standard. Waterproof for easy cleaning.

MC-130MO Digital Oral Thermometer  
 MC-130MR Digital Rectal Thermometer



**OMRON DIGITAL FLEXIBLE THERMOMETER**

Flexible tip for added comfort and safety. Beeps at peak temperature. Memory recalls last temperature. 20 second reading. Meets ASTM (E1112) standard.

MC-206 Flexible Digital Thermometer

**OMRON DIGITAL RIGID THERMOMETER**

3 times faster than traditional digital thermometers. Features include jumbo display for easy reading, memory and long-life 3-volt battery. For use oral, underarm and rectal locations. Safe, no fear of glass breakage or mercury. Includes 5 probe covers.

MC-106N Rigid Digital Thermometer

## OMRON DIGITAL THERMOMETER & PROBE COVERS

Individually wrapped disposable probe covers fit securely over the sensor portion on any Omron or Marshall thermometer.

|             |                                 |
|-------------|---------------------------------|
| 0007PROB    | Probe Covers, 50/bx             |
| 0009PROB    | Probe Covers, 100/bx            |
| 600RPROBE   | Probe Covers For MC-600, 25/bx  |
| 601MPROB500 | Probe Covers For MC-601, 500/bx |

## ADC® ADTEMP™ 422 VETERINARY DIGITAL THERMOMETER

Auto-off function conserves battery life. Range 90°F ± 2°F or 32°C - 43.9°C ± 1°C depending upon scale selection. Replaceable 1.55v (LR41) type battery provides up to 1,500 measurements. Specially designed housing with lanyard attachment for enhanced removal. Integral carry case.

422 Adtemp™ 422 Veterinary Digital Thermometer

## ADC® ADTEMP™ HYPOTHERMIA THERMOMETER

The first digital extended range thermometer. Auto off function conserves battery life. Replaceable 1.55v type battery provides up to 1,500 measurements. Clear two-piece plastic storage case.

419 Digital Hypothermia Thermometer



## OMRON BASAL DIGITAL THERMOMETER

Identify ovulation timing. Includes temperature tracking chart. Beeps at peak temperature. Waterproof for easy cleaning. Measures in 60 seconds and is three times faster than basal glass thermometers. Meets ASTM (E1112) standard.

MC-301N Basal Digital Thermometer



## TIDI SANITHERM THERMOMETER SHEATHS

|       |  |
|-------|--|
| 20521 | Mercury Oral Thermometer Sheaths, 100/bx, 50 bx/cs   |
| 20524 | Mercury Oral Thermometer Sheaths, 500/bx, 10 bx/cs   |
| 20525 | Mercury Oral Thermometer Sheaths, 5000/cs  |
| 20538 | Mercury Oral Thermometer Sheaths, Bulk, 11000/cs   |
| 20631 | White Digital Oral Thermometer Sheaths, Bulk, 11000/cs                                     |
| 20633 | Digital Oral Thermometer Sheaths, 100/bx, 50 bx/cs   |
| 20634 | Digital Oral Thermometer Sheaths, 500/bx, 10 bx/cs   |
| 20636 | Digital Oral Thermometer Sheaths, 5000/cs  |
| 20657 | Digital Oral Thermometer Sheaths, Packed in Sheets of 10, Indexed Every 50, Bulk, 10000/cs |
| 20710 | Prelubricated Mercury Rectal Thermometer Sheaths, 100/bx, 50 bx/cs                         |
| 20713 | Prelubricated Mercury Rectal Thermometer Sheaths, 500/bx, 10 bx/cs                         |
| 20714 | Prelubricated Mercury Rectal Thermometer Sheaths, 5000/cs                                  |
| 20775 | Digital Rectal Thermometer Sheaths, 100/bx, 50 bx/cs                                       |
| 20777 | Digital Rectal Thermometer Sheaths, 500/bx, 10 bx/cs                                       |



8884221000



8884309800

## KENDALL FILAC F-1500 ELECTRONIC THERMOMETER

Large, clear LCD temperature display. Automatic display. Auto Off feature. 30-35 second temperature reading. Audible signal. In-house two point calibration. Two lithium batteries. Lightweight.

|            |                          |
|------------|--------------------------|
| 8884309800 | Thermometer w/oral probe |
| 8884103650 | Oral Probe Only          |
| 8884103660 | Rectal Probe Only        |

### Accessories

|            |   |
|------------|---|
| 8884221000 | Probe Covers, 20/pkg, 25 pkg/bx, 10 bx/cs |
|------------|---|



## KENDALL FILAC FASTEMP ELECTRONIC THERMOMETER

4-10 second reading. Isolation chambers. Global iconography. Oversized backlit LCD display. Fahrenheit/Celsius button. Last temperature recall. 30 second pulse timer. Temperature range of 85°F-110°F/30C-43C. Electronic anti-theft system. ASTM and CEN compliance.

|        |   |
|--------|---|
| 202000 | Complete Oral/Ax System Including Electric Thermometer, 4 Lang Ops Man, Blue Oral/Ax Iso Chamber, Blue Probe w/4' cord, Wall Mountable Base, Mounting Hardware                        |
| 202018 | Fastemp Monitor, Including 4 Language Operating Manuals   |
| 202026 | Fastemp Blue Oral/Axillary Probe w/blue 4' cord   |
| 202027 | Fastemp Blue Oral/Axillary Probe w/blue 9' cord   |
| 202028 | Fastemp Blue Oral/Axillary Isolation Chamber  |
| 202034 | Fastemp Red Rectal Probe w/red 4' cord  |
| 202036 | Fastemp Red Rectal Isolation Chamber  |
| 202037 | Fastemp Red Rectal Probe w/red 9' cord  |
| 202038 | Fastemp Red Rectal Isolation Chamber, Red Rectal Probe w/red 4' cord, Wall Mountable BS, Hardware   |
| 202099 | Fastemp Calibration Plug  |
| 2020PC | Primary Care Kit Includes: FastTemp Thermometer, Blue Oral/Axillary Isolation Chamber, Blue Oral/Axillary Probe w/4' cord, Wall Mountable Base, Mounting Hardware, (500) Probe Covers |




## KENDALL FASTEMP PROBE COVERS

Disposable, single-use. Smooth, rigid plastic. Dispenser package fits in the isolation chamber. Contains no latex or natural rubber.

|        |   |
|--------|---|
| 202020 | Probe Covers for FASTEMP Thermometer, 20/pk, 24 pk/bx, 10 bx/cs |
|--------|---|



01690-200



01690-400

**WELCH ALLYN SURETEMP® PLUS  
ELECTRONIC THERMOMETER**

- 01690-200 Model 690 Electronic Thermometer, Oral Probe w/oral probe well, 2-Year Limited Warranty
- 01690-201 Model 690 Electronic Thermometer, Rectal Probe w/rectal probe well, 2-Year Limited Warranty
- 01690-300 Model 690 Electronic Thermometer, One Per Room, 9 ft Oral Probe w/oral probe well, 2-Year Limited Warranty (must be locked to the wall)
- 01690-301 Model 690 Electronic Thermometer, One Per Room, 9 ft Rectal Probe
- 01690-400 Model 690 Electronic Thermometer, Wall Mount, Oral Probe w/oral probe well, 2-Year Limited Warranty
- 01690-700 Model 690 Electronic Thermometer, Wall Mount, 9 ft Oral Probe w/oral probe well, rolling stand, spare probe & well holde, 2-Year Limited Warranty
- 01692-200 Model 692 Electronic Thermometer, Wall Mount, Security System w/ID location field, Oral Probe w/oral probe well, 3-Year Limited Warranty
- 01692-201 Model 692 Electronic Thermometer, Wall Mount, Security System w/ID location field, Rectal Probe w/rectal probe well, 3-Year Limited Warranty
- 01692-300 Model 692 Electronic Thermometer, One Per Room, Wall Mount, 9 ft Oral Probe w/oral probe well, 3-Year Limited Warranty (must be locked to the wall)
- 01692-700 Model 692 Electronic Thermometer, Wall Mount, 9 ft Oral Probe w/oral probe well, rolling stand, spare probe & well holder, 3-Year Limited Warranty



**WELCH ALLYN SURETEMP® THERMOMETER**

Oral response time: approx. 4 seconds. Axillary response time: approx. 10 seconds. Rectal response time: approx. 15 seconds. Continuous temperature reading. Measurement range: 84.0°F to 108.0°F. Unit dimensions: 2.25" x 3" x 7". Weight: 11.5 oz. Three AA alkaline batteries included. Three-year limited warranty on instrument. One-year limited warranty on probe.

- 01692-301 Model 692 SureTemp Plus Thermometer, One Per Room, Wall Mount, 9 ft Rectal Probe w/rectal probe well (must be locked to the wall)

**Contact your  
representative if you  
need more information.**



02678-000



05031-101

**WELCH ALLYN SURETEMP® THERMOMETER ACCESSORIES**

- 02678-000 4 ft Oral/Axillary Probe
- 02678-100 9 ft Oral/Axillary Probe
- 02679-000 4 ft Rectal Probe
- 02679-100 9 ft Rectal Probe
- 05031-101 Disposable Probe Covers, 1000/cs
- 06133-000 Calibration Key For M600
- 06137-000 Calibration Key For 767T, M678, M679
- 08043-200 Spare Probe Holder (does not include probe)
- 08070-000 Stand For SureTemp Units
- 08071-000 Infared Stand For SureTemp Units
- 08072-000 Thermometer Stand For Model M678
- 21004-0000 Additional Wall Mount
- 21004-1000 SureTemp Wall Mount w/security magnet
- 21004-2000 Wall Holder Assembly
- 21012-0000 Neck & Wrist Straps
- 21013-0000 Wall Holder For 250 Probe Covers

**ADC® ADTEMP™ TYMPANIC THERMOMETER**

Measures core temperature instantly from the ear reliably and accurately. Patented technology scans 512 times per second. Probe cover free design is less costly and more convenient to use. Large 6mm display is easy to read. Ultra compact ergonomic design fits comfortably in operator's hand and patient's ear. Fever alarm at 100.4°F. Complete with lithium battery, zippered nylon storage/travel case, and operating instructions.

- 421 Tympanic Thermometer



8884810055



8884830004

**KENDALL GENIUS MODEL 3000A TYMPANIC THERMOMETER**

Takes an accurate, non-invasive temperature. Peak Select System records and analyzes measurements. Scan button instantly activates thermometer. Release button ejects disposable probe. LCD menu prompt window. Incorporates 15-second pulse timer.

- 8884830004 Genius Model 3000A Thermometer
- Accessories**
- 8884810055 Cassette of 21 Probe Covers, 5 cassettes/pkg, 10 pkg/bx, 2 bx/cs
- 8884813711 Side Utility Basket for Mobile Thermometer Stand, 1/cs
- 8884814065 Wall Mount
- 8884830301 Mobile Thermometer Stand, 1/cs

## OMRON MC-514 EAR THERMOMETER

Instant, one second reading (takes 250 scans per second). Convenient storage case keeps thermometer and lens covers together. Engineered button placement ensures, safe comfortable use. Sleek design fits comfortably in the ear and in the palm of any size hand. Meets accuracy standards of ASTM (E 1965-98). Memory feature recalls last temperature reading. Measuring mode reads in Fahrenheit and Celsius. Lens cover detection for hygienic use. Five-year warranty.

MC-514 Ear Thermometer



## OMRON INSTANT EAR THERMOMETER LENS COVER

MC-509NLC Thermometer Lens Covers, 40/bx

## OMRON LENS COVER FOR GENTLE TEMP® INSTANT EAR THERMOMETER

MC-505LC Lens Covers For Ear Thermometer, 40/bx



## WELCH ALLYN PRO4000 EAR THERMOMETER

- 04000-200 Braun ThermoScan Pro4000 Ear Thermometer w/3 yr limited warranty
- 04000-300 Braun ThermoScan Pro4000 Ear Thermometer w/security base station and 3 yr limited warranty
- 04000-400 Braun ThermoScan Pro4000 Ear Thermometer w/roll stand and 3 yr limited warranty

### Accessories

- 01802-110 Braun Black Body Calibration Device
- 24001-0000 Rolling Thermometer Stand and Rechargeable Base Station
- 24001-1000 Recharging Base Station
- 53009-000 Replacement Alkaline Battery
- 53020-0000 Braun ThermoScan Pro4000 Rechargeable Battery Pack
- 80199-0000 Braun ThermoScan Pro4000 Replacement Power Cord for Base Station
- 85378-0000 Rolling Stand Mounting Bracket (Mounting Bracket for WA Universal Stand)


## WELCH ALLYN EAR THERMOMETER ACCESSORIES

- 05075-005 Probe Covers, 5000/cs
- 05075-800 Probe Covers, 800/cs
- 08075-000 Pro3000 Security Rolling Stand w/security base
- 23001-0000 Wall Holder & Probe
- 23001-1000 Wall Holder
- 80195-0000 Security Tether

## WELCH ALLYN 9020 TYMPANIC THERMOMETER ACCESSORIES

While the 9020 Tympanic Thermometer is no longer available, certain accessories for it are.

05051-001 Probe Covers, 1000/bx




## 3M™ TEMPA-DOT™ AX THERMOMETERS

Single-use and disposable with an adhesive strip to hold itself in place in the underarm area when measuring axillary body temperature. Individually wrapped.

5001 Fahrenheit Thermometer, 80/bx, 20 bx/cs

5002 Celsius Thermometer, 80/bx, 20 bx/cs




## 3M™ TEMPA-DOT™ PLUS THERMOMETERS

Single-use and disposable pre-assembled in a sheath for use in measuring rectal body temperature.

5192 Fahrenheit Thermometer, 100/bx, 20 bx/cs

5502 Celsius Thermometer, 100/bx, 20 bx/cs



## 3M™ TEMPA-DOT™ SINGLE-USE THERMOMETERS

Single-use and disposable thermometer that can measure oral or axillary body temperature.

- 5122 Fahrenheit Thermometer, Sterile, Individually Wrapped, 100/bx, 20 bx/cs
- 5125 Fahrenheit Thermometer, Sterile, Individually Wrapped, 100/bx, 5 bx/cs
- 5126 Fahrenheit Thermometer, Sterile, Individually Wrapped, 24/evn, 84 env/cs
- 5147 Fahrenheit Thermometer, Non-Sterile, Unwrapped, 250/bx, 24 bx/cs
- 5532 Celsius Thermometer, Sterile, Individually Wrapped, 100/bx, 20 bx/cs
- 5557 Celsius Thermometer, Non-Sterile, Unwrapped, 250/bx, 24 bx/cs





**CIDA ANEROID SPHYGMOMANOMETER**

Black enamel 300mmHg no-pin stop manometer. Range markings to facilitate selection of correct cuff size. Artery label and gauge holder. Latex inflation bladder and bulb. Chrome plated deflation valve. Complete with leatherette zippered carrying case and operating instructions. Assembled, inspected and packaged in the U.S.A. Ten year calibration warranty.

- 37-0400 Adult Aneroid Sphygmomanometer, Blue
- 37-0405 Adult Aneroid, Gray



**CIDA POCKET ANEROID SPHYGMOMANOMETER**

This is our standard sphygmomanometer. Features are the same as the premium pocket series, except for a standard bulb/valve and bladder combination. Latex. (Latex free available with LF suffix.)

- 37-0407 Adult Aneroid, Black
- 37-0408 Child Aneroid, Black
- 37-0409 Infant Aneroid, Black
- 37-0410 Newborn Aneroid, Black
- 37-0411 Thigh Aneroid, Black
- 37-0412 Large Adult Aneroid, Black
- 37-0413 Adult Aneroid, Gray
- 37-0417 Adult Aneroid, Green
- 37-0418 Adult Aneroid, Red
- 37-0419 Adult Aneroid, Orange
- 37-0420 Adult Aneroid, Teal
- 37-0421 Adult Aneroid, Royal Blue
- 37-0422 Adult Aneroid, Purple
- 37-0424 Adult Aneroid, Magenta
- 37-0607 Adult Aneroid, Black, Latex Free
- 37-0608 Child Aneroid, Black, Latex Free
- 37-0612 Large Adult Aneroid, Black, Latex Free
- 37-0648 Adult Aneroid, Dark Green
- 37-0652 Adult Aneroid, Burgundy



**CIDA PREMIUM POCKET ANEROID SPHYGMOMANOMETER**

This is our premier sphygmomanometer. Features lifetime warranted gauge with red needle and top quality inflation system, including a nylon Velcro cuff and highest quality bulb and valve. Complete with leatherette case. Latex. (Latex free available with LF suffix.)

- 37-0248 Adult Aneroid, Dark Green
- 37-0252 Adult Aneroid, Burgundy
- 37-0450 Adult Aneroid, Teal
- 37-0451 Adult Aneroid, Royal Blue
- 37-0452 Adult Aneroid, Purple
- 37-0454 Adult Aneroid, Magenta
- 37-0461 Adult Aneroid, Black
- 37-0463 Adult Aneroid, Gray
- 37-0464 Adult Aneroid, Green
- 37-0465 Adult Aneroid, Red
- 37-0466 Adult Aneroid, Orange
- 37-0470 Infant Aneroid, Black
- 37-0473 Newborn Aneroid, Black
- 37-0474 Child Aneroid, Black
- 37-0483 Thigh Aneroid, Black
- 37-0484 Large Adult Aneroid, Black
- 37-0661 Adult Aneroid, Black, Latex Free
- 37-0674 Child Aneroid, Black, Latex Free
- 37-0684 Large Adult Aneroid, Black, Latex Free



**CIDA SPHYG/SPRAGUE KIT**

Black enamel 300mmHg no-pin stop manometer. Color matched sprague-type stethoscope. Calibrated nylon cuff. Latex bladder and bulb. Chrome plated valve. Color matched nylon carrying case that is extra large with one large zippered compartment, outside Velcro® pocket, end strap and handy key ring. Case measures 10" x 6" x 3". Assembled, inspected and packaged in the U.S.A. Two year inflation system warranty. 20 year calibration warranty.

- 37-0700 Silver Sphyg/Sprague Kit
- 37-0734 Teal Sphyg/Sprague Kit
- 37-0736 Purple Sphyg/Sprague Kit
- 37-0740 Dark Green Sphyg/Sprague Kit
- 37-0750 Magenta Sphyg/Sprague Kit
- 37-0752 Burgundy Sphyg/Sprague Kit



## ADC® E-SPHYG DIGITAL ANEROID SPHYGMOMANOMETER

Virtually indestructible design can withstand repeated drops of 5 feet or more. Extra large backlit LCD with 7/8" digits, is easy to read in virtually any light. Splash resistant housing for use in harsh environments. Swivel cuff clip allows the manometer to be re-positioned for easy viewing. Auto off function conserves battery strength. Dual low battery warnings. Recalibration reminder after 10,000 cycles (recalibration is free for the life of the instrument). Systolic assist alerts practitioner to the onset of Korotkoff sounds - makes a great training aid. Deflation rate display promotes proper technique in accordance with AHA guidelines. Pulse rate display at completion of measurement. Manometer works with virtually every cuff on the market. Requires familiar and trusted manual auscultation to obtain blood pressure readings (use of stethoscope). 2 "AAA" batteries (included). May be ordered as a complete sphygmomanometer or gauge only for use with existing inflation systems. Lifetime calibration warranty.

|          |   |
|----------|---|
| 7002     | Digital Aneroid, Adult, Black                   |
| 7002C    | Digital Aneroid, Child, Black                   |
| 7002CLF  | Digital Aneroid, Child, Latex Free, Black       |
| 7002DB   | Digital Aneroid, Adult, Navy                    |
| 7002DBLF | Digital Aneroid, Adult, Latex Free, Navy        |
| 7002I    | Digital Aneroid, Infant, Black                  |
| 7002LF   | Digital Aneroid, Adult, Latex Free, Black       |
| 7002N    | Digital Aneroid, Newborn, Black                 |
| 7002OR   | Digital Aneroid, Adult, Orange                  |
| 7002ORLF | Digital Aneroid, Adult, Latex Free, Orange      |
| 7002RB   | Digital Aneroid, Adult, Royal Blue              |
| 7002RBLF | Digital Aneroid, Adult, Latex Free, Royal Blue  |
| 7002V    | Digital Aneroid, Adult, Purple                  |
| 7002VLF  | Digital Aneroid, Adult, Latex Free, Purple      |
| 7002X    | Digital Aneroid, Large Adult, Black             |
| 7002XLF  | Digital Aneroid, Large Adult, Latex Free, Black |
| 8002     | Digital Aneroid, Gauge Only                     |



## ADC® GENERAL PRACTICE MULTICUFF KIT

Adult, child, and large adult ADCUFF™ calibrated nylon cuffs with one tube latex inflation bladders and palm style manometer in a compact storage tray. Rugged gauge features precision crafted German mechanism housed in a patented double shell ABS impact polymer casing. Ambidextrous trigger release air valve. Quick release luer slip connectors. Four section smoke gray plexiglass storage tray features a small 14" x 3" footprint to save on counter space or may be wall mounted with included hardware. Assembled, inspected, and packed in the U.S.A. from foreign components. Two year inflation system warranty. Lifetime calibration warranty.

|        |                             |
|--------|-----------------------------|
| 705GPK | Multicuff Kit, Black, Latex |
|--------|-----------------------------|



## ADC® SYSTEM 5™ BLOOD PRESSURE SYSTEM

A multipocketed, heavy duty 1000D nylon attaché style case that keeps each component securely in place yet instantly accessible. Dual zipper closure permits rapid access to the contents, protects against the elements. Measures a compact 10" x 9" x 4" when closed, weighs just 1 1/2 lbs. Compact enough to fit in most emergency bags or trauma boxes. Reinforced carry handle doubles as a hang tab for wall mounting in the open position. Includes infant, child, adult, large adult and thigh ADCUFF™ nylon cuffs featuring the Size Guide™ marking system with a one tube latex inflation bladder. Cuffs equipped with quick release luer slip connector to permit rapid and secure selection of the appropriate cuff. Both systems utilize the revolutionary DIAGNOSTIX™ 804 palm style manometer with precision crafted German mechanism housed in a patented double shell ABS impact polymer casing, Ambidextrous trigger release air valve. Case and cuffs available in color coordinated navy or orange. Assembled, inspected and packaged in the U.S.A. Two year inflation system warranty. Lifetime calibration warranty.

|       |  |
|-------|--|
| 740   | 5 Cuff Blood Pressure System, Navy               |
| 740LF | 5 Cuff Blood Pressure System, Navy, Latex Free   |
| 746   | 5 Cuff Blood Pressure System, Orange             |
| 746LF | 5 Cuff Blood Pressure System, Orange, Latex Free |



## ADC® DIAGNOSTIX™ 703 SERIES SPHYGMOMANOMETER

Revolutionary palm style manometer (#804) with integral latex free bulb and trigger style release valve for one hand convenience. Precision crafted 300mmHg no-pin stop German mechanism tested to 200,000 cycles. Mechanism housed in a patented double shell ABS impact polymer casing. Ambidextrous infinite control trigger release air valve. Quick release luer slip connector. Six clinical cuff sizes available. ADCUFF™ nylon cuff with Size Guide™ marking system ensures proper cuff positioning and use of correct cuff size. Dip molded inflation bladder withstands heavy use. Complete with leatherette zippered carrying case with nylon zipper. Assembled, inspected, and packaged in the U.S.A. from foreign components. Two year inflation system warranty. Lifetime calibration warranty. Aneroids completely latex free have a LF in their number.

|     |                             |
|-----|-----------------------------|
| 703 | Adult Aneroid, Black, Latex |
|-----|-----------------------------|



760GARCQ

**ADC® PROSPHYG™ 760 SERIES ANEROID**

- 760BOOCQ Child, Aneroid, Betty Boop
- 760GARCQ Child, Aneroid, Garfield
- 760PEACQ Child, Aneroid, Peanuts
- 760BOOQ Adult, Aneroid, Betty Boop
- 760GARQ Adult, Aneroid, Garfield
- 760PEAQ Child, Aneroid, Peanuts
- 760BOOXQ Large Adult, Aneroid, Betty Boop
- 760GARXQ Large Adult, Aneroid, Garfield
- 760PEAXQ Large Adult, Aneroid, Peanuts
- 768-641BOO Adult, Aneroid, Betty Boop
- 768-641GAR Adult, Aneroid, Garfield
- 768-641PEA Adult, Aneroid, Peanuts

**ADC® PRO'S COMBO II ANEROID**

- 768-641BOO Adult, Aneroid, Betty Boop
- 768-641GAR Adult, Aneroid, Garfield
- 768-641PEA Adult, Aneroid, Peanuts



**OMRON ECONOMY ANEROIDS**

Calibrated nylon or cotton cuff with artery label and manometer tab. Bold, black gradations on white faceplate. Gray or blue enamel, no stop pin, 300 mm Hg gauge. Latex inflation system. Standard air release valve. Zippered carrying case. 10 year calibration warranty > 3 mm Hg.

- 11-200 Adult Aneroid w/nylon cuff, Blue, Blue Enamel
- 11-200C Child Aneroid w/nylon cuff, Blue, Blue Enamel
- 11-200X Large Adult Aneroid w/nylon cuff, Blue, Blue Enamel



**OMRON MARSHALL® SPHYGMOMANOMETER WITH COTTON CUFF**

Long-life cotton cuff with artery label and manometer tab. Bold, black gradations on white faceplate. Gray enamel, no stop pin, 300 mm Hg gauge. Latex inflation system. Standard air release valve. Zippered carrying case. Lifetime calibration warranty > 3 mm Hg.

- 0115M Adult Sphygmomanometer, Gray
- 0115MC Child Sphygmomanometer, Gray
- 0115MO Large Adult Sphygmomanometer, Gray



**OMRON MARSHALL® SPHYGMOMANOMETER WITH NYLON CUFF**

Long-life, calibrated nylon cuff, with artery label and manometer tab. Bold, black gradations on white faceplate. Gray enamel, no stop pin, 300 mm Hg gauge. Latex inflation system. Deluxe air release valve. Zippered, leatherette carrying case. Lifetime calibration warranty > 3 mm Hg.

- 0108M Adult Sphygmomanometer, Black
- 0108MC Child Sphygmomanometer, Black
- 0108MI Infant Sphygmomanometer, Black
- 0108ML Thigh Sphygmomanometer, Black
- 0108MNL Adult Sphygmomanometer, Black, Latex Free
- 0108MO Large Adult Sphygmomanometer, Black



**OMRON SELF-TAKING BLOOD PRESSURE KIT**

Long-lasting cotton D-ring cuff. One-hand operation. Easy-to-read gauge. Professional stethoscope attached to cuff. Includes vinyl carrying case, illustrated instructions, and blood pressure chart. One year warranty.

- 0104 Adult Blood Pressure Kit, Black
- 0104MAJ Large Adult Blood Pressure Kit, Blue
- 0116 Adult Blood Pressure Kit, Two-Party, Black
- HEM-18 Self Taking Manual Blood Pressure
- HEM-18S Self Taking Manual Blood Pressure Kit w/vinyl bladder





5098-02



5098-18



5098-20

## WELCH ALLYN CLASSIC HAND ANEROID

Unsurpassed accuracy and durability in the traditional classic design. For more than 30 years, the Tycos Classic Hand Aneroid has delivered performance and reliability in worldwide clinical use. Precision, jewelled bearings resist wear for consistently accurate readings. Every dial face is individually laser-engraved for exact calibration. Crisp white graphics on large charcoal-gray scale for better readability. Unique phosphor-bronze aneroid bellows ensure consistent pressure measurement. Recessed dial brings needle and dial markings close together and eliminates parallax. Single-tube leuc connector allows for easy interchange of cuff sizes. Lifetime calibration and manufacturing warranty. All latex free.

- 5098-02 Aneroid w/Adult Cuff, Latex Free
- 5098-05 Aneroid w/Child Cuff, Latex Free
- 5098-06 Aneroid w/Infant Cuff, Latex Free
- 5098-11 Aneroid w/Adult Cuff, w/o Case, Latex Free
- 5098-18 Pediatric Blood Pressure Kit w/case, Includes Newborn, Infant & Child Size Cuffs, Latex Free
- 5098-20 Family Practice Kit w/case, Includes Child, Adult & Large Adult Cuffs, Latex Free
- 5098-23 Family Practice Kit w/case, Includes Child Print Cuff, Adult & Large Adult, Latex Free

### Accessories

- 5098-03 Gauge & Bulb Only, Latex Free



## WELCH ALLYN CLASSIC POCKET ANEROID

Provides protection to the manometer during the rigors of clinical use. Precision, jewelled bearings resist wear for consistently accurate readings. Gauge attached to cuff protects instrument and prevents cuff and gauge from being easily separated. Laser-engraved dial for precise calibration. Unique phosphor-bronze bellows ensure consistent pressure measurement. Crisp, white graphics on large charcoal-gray scale for better readability. Lifetime warranty on calibration and manufacturing quality. All latex free.

- 5090-01 Aneroid w/Adult Hook Cuff, Latex Free
- 5090-02 Aneroid w/Adult Cuff, Latex Free
- 5090-06 Aneroid w/Child Cuff, Latex Free
- 5090-14 Aneroid w/Adult Cuff, w/o case, Latex Free
- 5090-22 Aneroid w/adult conductive cuff, Latex Free
- 5090-37 Aneroid w/Thigh Cuff w/o case, Latex Free
- 5090-41 Aneroid w/Large Adult Cuff, Latex Free

### Accessories

- 5090-03 Gauge Only, Latex Free



## WELCH ALLYN TYCOS® DS66 HAND ANEROID

More accurate blood pressure readings & lower lifetime ownership cost due to award-winning design that is gear-free & shock-resistant. Welch Allyn's line of aneroid sphygmomanometers featuring gear-free, shock-resistant DuraShock technology offer practical benefits compared to traditional sphygmomanometers w/gears. This is especially true in busy medical environments where instruments may be handled roughly or dropped. With DuraShock technology, Welch Allyn sphygmomanometers will remain in calibration longer, providing more accurate readings than traditional sphygmomanometers w/gears. They won't need to be replaced as often as traditional gauges w/gears, providing lower lifetime ownership cost. Only Welch Allyn aneroids feature gear-free, shock resistant DuraShock technology to deliver unprecedented durability. The only gauge capable of passing AAMI Shock Resistant standard ensures your sphygmomanometer will remain accurately calibrated longer. Lower lifetime cost of ownership since DuraShock gauges don't need replacing as often due to drops & rough handling. Only Welch Allyn offers a jeweled movement w/an industrial sapphire at the key friction point for long product life. Only Welch Allyn calibrates by computer & laser engraves a unique, custom-matched dial face for unprecedented accuracy over the instrument life. Lighter weight than traditional sphygmomanometers so it won't weigh down your lab coat. 1-piece cuff design is more durable than traditional cuff & bladder systems & it eliminates the need to remove or replace bladders. Latex-free cuff is easy to clean, disinfect, or sterilize for safety. 15-year calibration warranty for overall value.

- 5098-27 Aneroid w/Adult Cuff in Zippered Case, Latex Free Cuff
- 5098-28 Aneroid w/Large Adult Cuff in Zippered Case, Latex Free Cuff
- 5098-29 Aneroid w/Child Cuff in Zippered Case, Latex Free Cuff
- 5098-30 Multi-Cuff Kit, Includes Adult, Large Adult & Child-Print Inflation Systems in Zippered Case, Latex Free Cuff



5086-05



5098-70



5098-73

## WELCH ALLYN TYCOS® TR-2 HAND ANEROID

Lightweight and compact, the Tycos TR-2 hand aneroid provides portability in an affordably priced, hand-held design. The trigger-release air valve is positioned for optional left or right-handed use and allows for easy fingertip control of cuff deflation. Large easy-to-read dial face with bright graphics. Balanced palm design allows for comfortable use over repeated measurements. Meets AAMI accuracy standard of  $\pm 3$ mmHg. Comes in a zippered case. Five year warranty. All latex free.

- 5098-70 Aneroid w/Adult Cuff, Latex Free
- 5098-71 Aneroid w/Large Adult Cuff, Latex Free
- 5098-72 Aneroid w/Child Print Cuff, Latex Free
- 5098-73 Aneroid Multi-Cuff Kit, Includes Adult, Large Adult & Child Print Cuff, Latex Free

### Accessories

- 5082-63 Child Print Cuff, Latex Free
- 5086-05 Bulb Assembly, Latex Free
- 5098-79 Crystal, Latex Free






**WELCH ALLYN TYCOS® TR-2 HOME ANEROID**

Designed for simple and accurate personal blood pressure monitoring at home or work. Built-in stethoscope chestpiece allows for clear detection of blood pressure sounds with minimal training. Easy-to-use trigger deflation valve provides smooth release of air pressure. D-ring cuff for convenient self-application. Large & bright dial face for easy viewing of numbers. Instructional audio cassette and blood pressure recording forms. Meets AAMI accuracy standard of  $\pm 3$ mmHg. Five year warranty. All latex free.

- 5098-75 Aneroid w/Adult Cuff, Latex Free
- 5098-76 Aneroid w/Large Adult Cuff, Latex Free

**Accessories**

- 5082-67 Adult D-Ring Cuff & Bag, Latex Free
- 5082-68 Large Adult D-Ring Cuff & Bag, Latex Free
- 5098-74 Gauge & Bulb Only, Latex Free



750CW



752CM



752M

**ADC® DIAGNOSTIX™ 750 SERIES ANEROID**

Large, easy to read 6" dial with white scale on contrasting black face for high visibility. Luminescent dial for visibility in low light conditions. Lifetime calibration warranty. Contemporary Euro-style black casing complements any decor. Rear mounting bracket doubles as an inflation system storage compartment. ADCUFF™ inflation system with black nylon cuff, Size Guide™ marking system, premium latex bladder and bulb and 8' (when extended) coiled tubing with luer connectors. Includes all necessary hardware and instructions. Wall bracket swivels through 110° arc for enhanced visibility. Includes an ABS cuff storage compartment in main unit.

- 750CW Child Wall-Mounted Aneroid, Black, Latex
- 750CWLF Child Wall-Mounted Aneroid, Black, Latex Free
- 750IW Infant Wall-Mounted Aneroid, Black, Latex
- 750IWL Infant Wall-Mounted Aneroid, Black, Latex Free
- 750NW Newborn Wall-Mounted Aneroid, Black, Latex
- 750TW Thigh Wall-Mounted Aneroid, Black, Latex
- 750TWLF Thigh Wall-Mounted Aneroid, Black, Latex Free
- 750W Adult Wall-Mounted Aneroid, Black, Latex
- 750WLF Adult Wall-Mounted Aneroid, Black, Latex Free
- 750XW Large Adult Wall-Mounted Aneroid, Black, Latex
- 750XWLF Large Adult Wall-Mounted Aneroid, Black, Latex Free
- 752CM Child Mobile Aneroid, Black, Latex
- 752CMLF Child Mobile Aneroid, Black, Latex Free
- 752IM Infant Mobile Aneroid, Black, Latex
- 752IMLF Infant Mobile Aneroid, Black, Latex Free
- 752M Adult Mobile Aneroid, Black, Latex
- 752MLF Adult Mobile Aneroid, Black, Latex Free
- 752NM Newborn Mobile Aneroid, Black, Latex
- 752TM Thigh Mobile Aneroid, Black, Latex
- 752TMLF Thigh Mobile Aneroid, Black, Latex Free
- 752XM Large Adult Mobile Aneroid, Black, Latex
- 752XMLF Large Adult Mobile Aneroid, Black, Latex Free

**Can't find what  
you are looking for?  
Give us a call for fast,  
personal service.**



7670-01



7670-03

## WELCH ALLYN 767 SERIES WALL & MOBILE ANEROIDS

The Tycos 767 Wall and Mobile Aneroids offer all the dependability and accuracy of our classic aneroids, along with an updated look. Precision, jeweled-movement manometer that is certified to be accurate within one percent of scale ( $\pm 3\text{mmHg}$ ). Laser-engraved dial face for precision calibration. Built-in basket for securely storing the inflation system. High-impact housing and bezel to protect the movement and crystal. Large dial face (6 1/4") and bright orange pointer for easy reading. 40° side-to-side swivel for clear viewing at all angles. The Mobile aneroid has adjustable height stand and locking collar, 4' coiled tubing and 360° manometer swivel for ease of reading. Lifetime calibration and manufacturing warranty.

- 7670-01 Wall Aneroid w/Adult Cuff, Latex Free
- 7670-02 Wall Aneroid Manometer Only w/8 ft tubing, Latex Free
- 7670-03 Mobile Aneroid w/Adult Cuff, Latex Free
- 7670-04 Mobile Aneroid w/weighted base & adult cuff, Latex Free

### Accessories

- 7670-06 Mobile Stand & Basket Assembly Only (base weight not included), Latex Free
- 7670-07 Inflation System Basket, White, Latex Free
- 7670-08 Crystal, Latex Free



5091-13



5091-38



5091-41



5091-47

## WELCH ALLYN TYCOS® 509 SERIES WALL & MOBILE ANEROIDS

Tycos 509 Series aneroids provide all the accuracy and reliability you expect from Tycos - at an affordable price. Lifetime calibration and manufacturing warranty. Precision, jeweled bearings resist wear for consistently accurate readings. Unique phosphor-bronze aneroid bellows ensure consistent pressure measurement. Visual check of accuracy on the dial. Strong die-cast zinc internal housing protects mechanism. The Wall style has an 8' coiled tubing with 90° swivel for flexible reading, 10° forward tilt to reduce glare. The Mobile style has 360° swivel with adjustable height stand, and 4' coiled tubing.

- 5091-23 Burnt Orange Wall Aneroid Manometer Only w/8 ft tubing, Latex Free
- 5091-38 Charcoal Wall Aneroid w/Adult Cuff, Latex Free
- 5091-39 Charcoal Wall Aneroid Manometer Only w/8 ft tubing, Latex Free
- 5091-41 Charcoal Mobile Aneroid w/Adult Cuff, Latex Free

### Accessories

- 5091-13 Inflation System Basket, Small, Stainless Steel, Latex Free
- 5091-47 Inflation System Basket, Large, Black, Latex Free
- 5091-48 Mobile Aneroid Stand & Basket Assembly Only, Latex Free
- 7670-05 Mobile Stand Base Weight Only, Latex Free

## ADC® ADVANTAGE™ 6012 SEMI-AUTOMATIC BP MONITOR

Compact and easy-to-use. Measures systolic, diastolic pressure and pulse. Manual inflation. Self adjusting latex free D-bar cuff. Automatic deflation valve. Auto off.

- 6012 Advantage™ 6012 BP Monitor

## ADC® ADVANTAGE™ 6013 AUTOMATIC BP MONITOR

Smart Logic™ Technology determines ideal inflation level. Measures systolic, diastolic pressure, and pulse. Automatic inflation and deflation control. Self adjusting latex free D-bar cuff. Recalls last 14 readings with date and time stamp. Auto off. Optional AC adapter (sold separately).

- 6013 Advantage™ 6013 BP Monitor, 9" to 13" Limb Range
- 6013SA Advantage™ 6013 BP Monitor, 7.5" to 9.4" Limb Range
- 6013X Advantage™ 6013 BP Monitor, 13" to 17" Limb Range

### Accessories

- 850-6013 Cuff & Bladder, Adult For 6013 BP Monitor, Navy, Latex Free
- 850-6013SA Cuff and Bladder, Adult for 6013SA BP Monitor, Navy, Latex Free
- 850-6013X Cuff & Bladder, Large Adult For 6013 BP Monitor, Navy, Latex Free

## ADC® ADVANTAGE™ 6014 ADVANCED BP MONITOR

Smart Logic™ Technology determines ideal inflation level. Measures systolic, diastolic pressure, and pulse. Automatic inflation and deflation control. Average mode automatically averages three readings for more meaningful results. Self adjusting latex free D-bar cuff. Advanced two zone memory with date and time stamp recalls last 30 readings for two different people. Auto off. Optional thermal printer (sold separately). Optional AC adapter (sold separately).

- 6014 Advantage™ 6014 BP Monitor

### Accessories

- 6013ZAC AC Adapter For 6013, 6014, 6014P
- 6014P Printer For 6014 & 6016 Monitors
- 6014PAPER Paper For 6014P Printer, 5/pkg

## ADC® ADVANTAGE™ 6016 DIGITAL WRIST BP MONITOR

Smart Logic™ Technology determines ideal inflation level. Measures systolic, diastolic pressure, and pulse. New 8 bit microprocessor and advanced algorithm for more accurate readings. Average mode automatically averages three readings for more meaningful results. Latex free wrist cuff fits 5" to 7" wrists. Advanced two zone memory with date and time stamp recalls last 30 readings for two different people. Auto off. Compact storage case. Optional thermal printer (sold separately).

- 6016 Advantage™ 6016 BP Monitor

## OMRON AUTOMATIC BLOOD PRESSURE MONITOR WITH COMFIT™ CUFF

Comfit™ cuff is pre-formed for quick and proper fit (9" to 17" arms). Clinically proven accuracy. Technology offers personalized inflation for maximum comfort. 90 memory recall with date and time stamp. Plastic case for convenient storage included. Blood pressure and pulse rate displayed on a large digital panel. Large, clear digital display. AC adapter and 4 AA batteries included. Five-year warranty.

- HEM-780 Automatic Blood Pressure Monitor w/Comfit™ Cuff



## OMRON DIGITAL BP MONITOR WITH EZ WRAP CUFF

In the clinical setting, the EZ Wrap Cuff (9"-17") allows the unit to be used on a wider range of arm sizes, often negating the need for separate, larger sized cuffs. For home use, the EZ Wrap Cuff makes applying the cuff much easier for an individual without assistance. The IntelliSense technology allows for comfortable, personalized inflation and accurate results. Two person, 21 memory recall with date and time stamped for each person. AC adapter and 4 "AA" batteries included. Five-year warranty.

- HEM-773AC Digital BP Monitor w/EZ Wrap Cuff



**OMRON DIGITAL WRIST BP MONITOR WITH APS**

APS (Advanced Positioning Sensor) technology means that unless the wrist is held in the proper position, slightly above the heart, it will not provide a reading, eliminating mixed results due to user technique. 90 memory recall. IntelliSense technology offers personalized inflation for maximum comfort, graphing capabilities with date and time stamp. Five-year warranty.

HEM-637 Digital Wrist BP Monitor w/APS



**OMRON INTELLISENSE™ AUTOMATIC PLUS BLOOD PRESSURE MONITOR**

IntelleSense™ technology offers automated level setting and electronic deflation control. Completely automatic, operates with the push of a button. Large, clear digital display. 14 memory storage with AC adapter and illustrated instruction manual included. Quickly and accurately measures blood pressure and pulse. Tailored arm cuff is durable and fits for arms 9"-13" in circumference (large cuff, 13" to 17" available). Five year warranty.

HEM-711AC Automatic Inflation Monitor w/fuzzy logic & AC adapter



HEM-422CRLC



HEM-432CN

**OMRON REGENCY SERIES MANUAL INFLATION BLOOD PRESSURE MONITOR**

Large digital panel with 14 reading memory. Easy squeeze inflation bulb quickly fills cuff with air. Complete with zippered carrying case. Operates on 4 AA batteries (not included).

HEM-422CRLC Manual Inflation Monitor w/large adult cuff  
HEM-432CN Manual Inflation Monitor



**ADC® DIAGNOSTIX™ 922 SERIES SPHYGMOMANOMETER**

Rugged diecast aluminum housing with textured black baked enamel finish complements any decor. Fold down cover with spring lock mechanism. Graduated unbreakable plastic cartridge tube features 4.00mm inside bore for unrestricted, friction-free mercury flow. High visibility red graduations. Filter system keeps mercury clean for reliable performance. Locking reservoir secures mercury during storage, transport, or maintenance. ADCUFF™ black nylon cuff with ADC's proprietary Size Guide™ marking system ensures proper cuff positioning and use of correct cuff size. Dip molded latex inflation bladder withstands heavy use. Latex inflation bulb fitted with filter screen end valve reduces dust build-up. ADFLOW™ chrome plated brass air release valve with filter screen protection permits precise deflation control and years of reliable service. Comprehensive owner's manual. Assembled, inspected, and packaged in the U.S.A. from foreign components. Two year inflation system warranty. Lifetime warranty on cartridge tube.

- 922 Adult Sphygmomanometer, Black, Latex
- 922C Child Sphygmomanometer, Black, Latex
- 922CLF Child Sphygmomanometer, Black, Latex Free
- 922I Infant Sphygmomanometer, Black, Latex
- 922ILF Infant Sphygmomanometer, Black, Latex Free
- 922LF Adult Sphygmomanometer, Black, Latex Free
- 922N Newborn Sphygmomanometer, Black, Latex
- 922T Thigh Sphygmomanometer, Black, Latex
- 922X Large Adult Sphygmomanometer, Black, Latex
- 922XLF Large Adult Sphygmomanometer, Black, Latex Free
- 922TLF Thigh Sphygmomanometer, Black, Latex Free

**Accessories**

- 922-00 Manometer Only For 922, Latex Free
- 922-126 Plastic Cartridge Tube For 922, Latex Free

**ADC® DIAGNOSTIX™ 932 SERIES SPHYGMOMANOMETER**

Rugged housing constructed from a combination of virtually indestructible ABS and diecast aluminum. Striking black finish complements any decor. Contemporary upright design with weighted base and integral wire inflation system basket. Graduated unbreakable plastic cartridge tube features a massive 5.00mm inside bore for unrestricted, friction-free mercury flow. White oversized numerals on contrasting black scale combined with high visibility red graduations for easy reading. EZ-Tube™ system allows safe, simplified maintenance of the cartridge tube or top cap filters. Additional safety provided by the locking reservoir which secures mercury for transport or maintenance. ADCUFF™ black nylon cuff with ADC's proprietary Size Guide™ marking system ensures proper cuff positioning and use of correct cuff size. Dip molded latex inflation bladder withstands heavy use. Latex inflation bulb fitted with filter screen end valve reduces dust build-up. ADFLOW™ chrome plated brass air release valve with filter screen protection permits precise deflation control. Comprehensive owner's manual. Assembled, inspected, and packaged in the U.S.A. from foreign components. Two year inflation system warranty. Lifetime warranty on cartridge tube.

- 932 Adult Sphygmomanometer, Black, Latex
- 932C Child Sphygmomanometer, Black, Latex
- 932CLF Child Sphygmomanometer, Black, Latex Free
- 932I Infant Sphygmomanometer, Black, Latex
- 932ILF Infant Sphygmomanometer, Black, Latex Free
- 932LF Adult Sphygmomanometer, Black, Latex Free
- 932N Newborn Sphygmomanometer, Black, Latex
- 932T Thigh Sphygmomanometer, Black, Latex
- 932TLF Thigh Sphygmomanometer, Black, Latex Free
- 932X Large Adult Sphygmomanometer, Black, Latex
- 932XLF Large Adult Sphygmomanometer, Black, Latex Free





**ADC® DIAGNOSTIX™ 952B SERIES SPHYGMOMANOMETER**

Rugged diecast aluminum housing. Striking black finish complements any decor. Graduated unbreakable plastic cartridge tube (recessed for safety) features 5.00mm inside bore for unrestricted, friction-free mercury flow. Wall bracket swivels through 110° arc for enhanced visibility and viewer comfort. EZ-Tube™ system allows safe, simplified maintenance of the cartridge tube or top cap filters. Locking reservoir secures mercury during transport or maintenance. ADCUFF™ nylon cuff with ADC®'s proprietary Size Guide™ marking system ensures proper cuff positioning and use of correct cuff size. Complete with vinyl coated wire cuff basket, 8 ft length (when extended) of coiled tubing, all mounting hardware and comprehensive owner's manual. Assembled, inspected and packaged in the U.S.A. Two year inflation system warranty. Lifetime warranty on cartridge tube.

- 952B Adult Sphygmomanometer, Black, Latex
- 952BC Child Sphygmomanometer, Black, Latex
- 952BCLF Child Sphygmomanometer, Black, Latex Free
- 952BI Infant Sphygmomanometer, Black, Latex
- 952BILF Infant Sphygmomanometer, Black, Latex Free
- 952BLF Adult Sphygmomanometer, Black, Latex Free
- 952BN Newborn Sphygmomanometer, Black, Latex
- 952BT Thigh Sphygmomanometer, Black, Latex
- 952BTLF Thigh Sphygmomanometer, Black, Latex Free
- 952BX Large Adult Sphygmomanometer, Black, Latex
- 952BXLF Large Adult Sphygmomanometer, Black, Latex Free

**Accessories**

- 952-025 Cuff Storage Basket For 952
- 952-129 Plastic Cartridge Tube For 952
- 952B-00 Manometer Only

**ADC® ADVANCED DIGITAL BP MONITOR WITH PC LINK**

- 6017 Advanced Digital BP w/PC link (software not included), Navy
- 6017PC PC Link Software & USB Cable

**Call for  
current pricing!**



H-003D



HEM-907-STAND

**OMRON DIGITAL BLOOD PRESSURE PARTS & ACCESSORIES**

- 004IC Interface Cable RS232C For Model HEM-7471C
- 0080AGRAY D-Ring Contour Cuff For HEM-412C, HEM-712C & HEM-703CP, Gray
- 0090TRP Replacement Roll of Thermal Paper For HEM-705CP, 5 rl/bx
- 0094AGRAY D-Ring Curved Contour Cuff, Gray
- 0097IR Ink Roller For HEM-703CP & 97
- 0097RP Replacement Roll Paper For HEM-703CP & 97, 5/bx
- 0146GRAY Gray Bulb & Valve For HEM-412C & HEM-425C
- H-003D Large Adult Cuff, Gray
- H-003DS Small Size Cuff For IntelliSense Monitors
- HEM-907-BLDRL Bladder Large Cuff For HEM-907 & 907XL
- HEM-907-BLDRM Bladder Medium Cuff For HEM-907 & 907XL
- HEM-907-BLDRS Bladders Small Cuff For HEM-907 & 907XL
- HEM-907-CL19 Cuff & Bladder Set, Large 32-42cm For HEM-907 & 907XL
- HEM-907-CR19 Cuff & Bladder Set, Medium 22-32cm For HEM-907 & 907XL
- HEM-907-CS19 Cuff & Bladder Set, Small 17-22cm For HEM-907 & 907XL
- HEM-907-CUFFL Cuff For Large Bladder 32-42cm For HEM-907 & 907XL
- HEM-907-CUFFM Cuff For Medium Bladder 22-32cm For HEM-907 & 907XL
- HEM-907-CUFFS Cuff For Small Bladder 17-22cm For HEM-907 & 907XL
- HEM-907-PBAT Battery Pack For HEM-907 & 907XL BP Unit
- HEM-907-PKIT Pole Mount For HEM-907 & 907XL BP Unit
- HEM-907-STAND Stand For HEM-907 & 907XL BP Unit
- HEM-907-WKIT Wall Mount Kit For HEM-907 & 907XL BP Unit
- HEM-AC-J HEM-907 & 907XL Adpater
- HEM-ADPT1 AC Adapter For Omron & Marshall Auto Inflate Monitors
- HEM-TUBE-100 Air Tube 1.0m WL For HEM-907 & 907XL
- HEM-TUBE-130 Air Tube 1.3m WL For HEM-907 & 907XL

**ADC® ADCUFF™ 1-TUBE BLADDERS**

Dip molded latex or latex free. 55cm tubing length. Size and model number printed on bladder. For use with most palm style aneroids. Individually boxed. Two year warranty.

- 820 One-Tube Bladder, Latex, Adult
- 820C One-Tube Bladder, Latex, Child
- 820CLF One-Tube Bladder, Latex Free, Child
- 820I One-Tube Bladder, Latex, Infant
- 820ILF One-Tube Bladder, Latex Free, Infant
- 820LF One-Tube Bladder, Latex Free, Adult
- 820N One-Tube Bladder, Latex, Newborn
- 820T One-Tube Bladder, Latex, Thigh
- 820TLF One-Tube Bladder, Latex Free, Thigh
- 820X One-Tube Bladder, Latex, Large Adult
- 820XLF One-Tube Bladder, Latex Free, Large Adult



**ADC® ADCUFF™ 2-TUBE BLADDERS**

Dip molded latex or latex free. 55cm tubing length. Size and model number printed on bladder. For use with most aneroid and mercurial blood pressure instruments. Individually boxed. Two year warranty.

|         |  |
|---------|--|
| 810     | Two-Tube Bladder, Latex, Adult                                   |
| 810C    | Two-Tube Bladder, Latex, Child                                   |
| 810CH   | Two-Tube Bladder, Latex, Hospital Pack, Child, 10/pkg            |
| 810CHLF | Two-Tube Bladder, Latex Free, Hospital Pack, Child, 10/pkg       |
| 810CLF  | Two-Tube Bladder, Latex Free, Child                              |
| 810E    | Two-Tube Bladder, Latex, Adult, Standard                         |
| 810H    | Two-Tube Bladder, Latex, Hospital Pack, Adult, 10/pkg            |
| 810HLF  | Two-Tube Bladder, Latex Free, Hospital Pack, Adult, 10/pkg       |
| 810I    | Two-Tube Bladder, Latex, Infant                                  |
| 810IH   | Two-Tube Bladder, Latex, Hospital Pack, Infant, 5/pkg            |
| 810IHLF | Two-Tube Bladder, Latex Free, Hospital Pack, Infant, 5/pkg       |
| 810ILF  | Two-Tube Bladder, Latex Free, Infant                             |
| 810LF   | Two-Tube Bladder, Latex Free, Adult                              |
| 810N    | Two-Tube Bladder, Latex, Newborn                                 |
| 810NH   | Two-Tube Bladder, Latex, Hospital Pack, Newborn, 5/pkg           |
| 810T    | Two-Tube Bladder, Latex, Thigh                                   |
| 810TH   | Two-Tube Bladder, Latex, Hospital Pack, Thigh, 5/pkg             |
| 810THLF | Two-Tube Bladder, Latex Free, Thigh, 5/pkg                       |
| 810TLF  | Two-Tube Bladder, Latex Free, Thigh                              |
| 810X    | Two-Tube Bladder, Latex, Large Adult                             |
| 810XH   | Two-Tube Bladder, Latex, Hospital Pack, Large Adult, 10/pkg      |
| 810XHLF | Two-Tube Bladder, Latex Free, Hospital Pack, Large Adult, 10/pkg |
| 810XLF  | Two-Tube Bladder, Latex Free, Large Adult                        |

**CIDA INFLATION BULBS & VALVES**

Dip molded latex or latex free. For use with most aneroid and mercurial blood pressure instruments. End valve included. Individually boxed. Two year warranty.

|         |                              |
|---------|------------------------------|
| 38-0105 | Standard Bulb & Valve, Latex |
|---------|------------------------------|

**CIDA 2-TUBE BLADDERS**

Dip molded latex or latex free. 55cm tubing length. Size and model number printed on bladder. For use with most aneroid and mercurial blood pressure instruments. Individually boxed. Two year warranty.

|         |                                |
|---------|--------------------------------|
| 38-0630 | Two-Tube Bladder, Latex, Adult |
|---------|--------------------------------|

**CIDA CALIBRATED NYLON CUFF**

Featuring marking system. For use with most aneroid or mercurial blood pressure instruments. Two year warranty.

|         |                         |
|---------|-------------------------|
| 38-0212 | Adult Nylon Cuff, Black |
|---------|-------------------------|

**CIDA CUFFS & BLADDERS**


For use with most aneroid or mercurial blood pressure instruments. Individually boxed.

|         |   |
|---------|---|
| 38-0121 | Adult Cuff & Bladder, 2-Tube, Black             |
| 38-0160 | Adult Cuff & Bladder, 1-Tube, Black             |
| 38-0621 | Adult Cuff & Bladder, 2-Tube, Black, Latex Free |
| 38-0660 | Adult Cuff & Bladder, 1-Tube, Black, Latex Free |

**CIDA INFLATION SYSTEM**

Includes nylon cuff, deluxe 2-tube bladder, bulb and valve. For use with most aneroid or mercurial blood pressure instruments. Individually boxed.

|         |   |
|---------|---|
| 38-0400 | Adult Inflation System, Black, Latex      |
| 38-0600 | Adult Inflation System, Black, Latex Free |



**ADC® ADCUFF™ & BLADDERS, 1-TUBE**

Durable nylon construction. Exclusive Size Guide™ marking system. Latex or latex free inflation bladder with one 55cm tubes. For use with most palm style aneroid instruments. Individually boxed. Two-year warranty.

|          |  |
|----------|--|
| 850      | ADCUFF™ & Bladder, 1-Tube, Latex, Black, Adult             |
| 850BD    | ADCUFF™ & Bladder, 1-Tube, Latex, Burgundy, Adult          |
| 850BDLF  | ADCUFF™ & Bladder, 1-Tube, Latex Free, Burgundy, Adult     |
| 850C     | ADCUFF™ & Bladder, 1-Tube, Latex, Black, Child             |
| 850CLF   | ADCUFF™ & Bladder, 1-Tube, Latex Free, Black, Child        |
| 850DB    | ADCUFF™ & Bladder, 1-Tube, Latex, Navy, Adult              |
| 850DBC   | ADCUFF™ & Bladder, 1-Tube, Latex, Navy, Child              |
| 850DBCLF | ADCUFF™ & Bladder, 1-Tube, Latex Free, Navy, Child         |
| 850DBI   | ADCUFF™ & Bladder, 1-Tube, Latex, Navy, Infant             |
| 850DBILF | ADCUFF™ & Bladder, 1-Tube, Latex Free, Navy, Infant        |
| 850DBLF  | ADCUFF™ & Bladder, 1-Tube, Latex Free, Navy, Adult         |
| 850DBN   | ADCUFF™ & Bladder, 1-Tube, Latex, Navy, Newborn            |
| 850DBT   | ADCUFF™ & Bladder, 1-Tube, Latex, Navy, Thigh              |
| 850DBTLF | ADCUFF™ & Bladder, 1-Tube, Latex, Navy, Thigh, Latex-Free  |
| 850DBX   | ADCUFF™ & Bladder, 1-Tube, Latex, Navy, Large Adult        |
| 850DBXLF | ADCUFF™ & Bladder, 1-Tube, Latex Free, Navy, Large Adult   |
| 850DG    | ADCUFF™ & Bladder, 1-Tube, Latex, Dark Green, Adult        |
| 850DGLF  | ADCUFF™ & Bladder, 1-Tube, Latex Free, Dark Green, Adult   |
| 850I     | ADCUFF™ & Bladder, 1-Tube, Latex, Black, Infant            |
| 850ILF   | ADCUFF™ & Bladder, 1-Tube, Latex Free, Black, Infant       |
| 850LF    | ADCUFF™ & Bladder, 1-Tube, Latex Free, Black, Adult        |
| 850M     | ADCUFF™ & Bladder, 1-Tube, Latex, Magenta, Adult           |
| 850MLF   | ADCUFF™ & Bladder, 1-Tube, Latex Free, Magenta, Adult      |
| 850N     | ADCUFF™ & Bladder, 1-Tube, Latex, Black, Newborn           |
| 850RB    | ADCUFF™ & Bladder, 1-Tube, Latex, Royal Blue, Adult        |
| 850RBLF  | ADCUFF™ & Bladder, 1-Tube, Latex Free, Royal Blue, Adult   |
| 850T     | ADCUFF™ & Bladder, 1-Tube, Latex, Black, Thigh             |
| 850TL    | ADCUFF™ & Bladder, 1-Tube, Latex, Teal, Adult              |
| 850TLF   | ADCUFF™ & Bladder, 1-Tube, Latex Free, Black, Thigh        |
| 850TLLF  | ADCUFF™ & Bladder, 1-Tube, Latex Free, Teal, Adult         |
| 850V     | ADCUFF™ & Bladder, 1-Tube, Latex, Purple, Adult            |
| 850VLF   | ADCUFF™ & Bladder, 1-Tube, Latex Free, Purple, Adult       |
| 850X     | ADCUFF™ & Bladder, 1-Tube, Latex, Purple, Large Adult      |
| 850XLF   | ADCUFF™ & Bladder, 1-Tube, Latex Free, Purple, Large Adult |
| 851      | ADCUFF™ & Bladder, 1-Tube, Latex, Gray, Adult              |
| 851LF    | ADCUFF™ & Bladder, 1-Tube, Latex Free, Gray, Adult         |
| 854      | ADCUFF™ & Bladder, 1-Tube, Latex, Green, Adult             |
| 854LF    | ADCUFF™ & Bladder, 1-Tube, Latex Free, Green, Adult        |
| 855      | ADCUFF™ & Bladder, 1-Tube, Latex, Red, Adult               |
| 855LF    | ADCUFF™ & Bladder, 1-Tube, Latex Free, Red, Adult          |
| 856      | ADCUFF™ & Bladder, 1-Tube, Latex, Orange, Adult            |
| 856C     | ADCUFF™ & Bladder, 1-Tube, Latex, Orange, Child            |
| 856CLF   | ADCUFF™ & Bladder, 1-Tube, Latex Free, Orange, Child       |
| 856I     | ADCUFF™ & Bladder, 1-Tube, Latex, Orange, Infant           |
| 856ILF   | ADCUFF™ & Bladder, 1-Tube, Latex Free, Orange, Infant      |
| 856LF    | ADCUFF™ & Bladder, 1-Tube, Latex Free, Orange, Adult       |
| 856T     | ADCUFF™ & Bladder, 1-Tube, Latex, Orange, Thigh            |
| 856TLF   | ADCUFF™ & Bladder, 1-Tube, Latex Free, Orange, Thigh       |
| 856X     | ADCUFF™ & Bladder, 1-Tube, Latex, Orange, Large Adult      |
| 856XLF   | ADCUFF™ & Bladder, 1-Tube, Latex Free, Orange, Large Adult |



**ADC® ADCUFF™ & BLADDERS, 2-TUBE**

Durable nylon construction. Exclusive Size Guide™ marking system. Latex or latex free inflation bladder with two 55cm tubes. For use with most aneroid or mercurial blood pressure instruments. Individually boxed. Two-year warranty.

- 840 ADCUFF™ & Bladder, 2-Tube, Latex, Black, Adult
- 840BD ADCUFF™ & Bladder, 2-Tube, Latex, Burgundy, Adult
- 840BDLF ADCUFF™ & Bladder, 2-Tube, Latex Free, Burgundy, Adult
- 840C ADCUFF™ & Bladder, 2-Tube, Latex, Black, Child
- 840CH ADCUFF™ & Bladder, 2-Tube, Latex, Black, Child, 10/pkg
- 840CHLF ADCUFF™ & Bladder, 2-Tube, Latex Free, Black, Child, 10/pkg
- 840CLF ADCUFF™ & Bladder, 2-Tube, Latex Free, Black, Child
- 840DB ADCUFF™ & Bladder, 2-Tube, Latex, Navy, Adult
- 840DBC ADCUFF™ & Bladder, 2-Tube, Latex, Navy, Child
- 840DBCLF ADCUFF™ & Bladder, 2-Tube, Latex Free, Navy, Child
- 840DBI ADCUFF™ & Bladder, 2-Tube, Latex, Navy, Infant
- 840DBILF ADCUFF™ & Bladder, 2-Tube, Latex Free, Navy, Infant
- 840DBLF ADCUFF™ & Bladder, 2-Tube, Latex Free, Navy, Adult
- 840DBN ADCUFF™ & Bladder, 2-Tube, Latex, Navy, Newborn
- 840DBT ADCUFF™ & Bladder, 2-Tube, Latex, Navy, Thigh
- 840DBTLF ADCUFF™ & Bladder, 2-Tube, Latex Free, Navy, Thigh
- 840DBX ADCUFF™ & Bladder, 2-Tube, Latex, Navy, Large Adult
- 840DBXLF ADCUFF™ & Bladder, 2-Tube, Latex Free, Navy, Large Adult
- 840DG ADCUFF™ & Bladder, 2-Tube, Latex, Dark Green, Adult
- 840DGLF ADCUFF™ & Bladder, 2-Tube, Latex Free, Dark Green, Adult
- 840H ADCUFF™ & Bladder, 2-Tube, Latex, Black, Adult, 10/pkg
- 840HLF ADCUFF™ & Bladder, 2-Tube, Latex Free, Black, Adult, 10/pkg
- 840I ADCUFF™ & Bladder, 2-Tube, Latex, Black, Infant
- 840IH ADCUFF™ & Bladder, 2-Tube, Latex, Black, Infant, 5/pkg
- 840IHLF ADCUFF™ & Bladder, 2-Tube, Latex Free, Black, Infant, 5/pkg
- 840ILF ADCUFF™ & Bladder, 2-Tube, Latex Free, Black, Infant
- 840LF ADCUFF™ & Bladder, 2-Tube, Latex Free, Black, Adult
- 840M ADCUFF™ & Bladder, 2-Tube, Latex, Magenta, Adult
- 840MLF ADCUFF™ & Bladder, 2-Tube, Latex Free, Magenta, Adult
- 840N ADCUFF™ & Bladder, 2-Tube, Latex, Black, Newborn
- 840NH ADCUFF™ & Bladder, 2-Tube, Latex, Black, Newborn, 5/pkg
- 840RB ADCUFF™ & Bladder, 2-Tube, Latex, Royal Blue, Adult
- 840RBLF ADCUFF™ & Bladder, 2-Tube, Latex Free, Royal Blue, Adult
- 840T ADCUFF™ & Bladder, 2-Tube, Latex, Black, Thigh
- 840TH ADCUFF™ & Bladder, 2-Tube, Latex, Turquoise, Adult, 10/pkg
- 840THLF ADCUFF™ & Bladder, 2-Tube, Latex Free, Black, Thigh, 5/pkg
- 840TL ADCUFF™ & Bladder, 2-Tube, Latex, Teal, Adult
- 840TLF ADCUFF™ & Bladder, 2-Tube, Latex Free, Black, Thigh
- 840TLLF ADCUFF™ & Bladder, 2-Tube, Latex Free, Teal, Adult
- 840V ADCUFF™ & Bladder, 2-Tube, Latex, Purple, Adult
- 840VLF ADCUFF™ & Bladder, 2-Tube, Latex Free, Purple, Adult
- 840X ADCUFF™ & Bladder, 2-Tube, Latex, Black, Large Adult
- 840XH ADCUFF™ & Bladder, 2-Tube, Latex, Black, Large Adult, 10/pkg
- 840XHLF ADCUFF™ & Bladder, 2-Tube, Latex Free, Black, Large Adult, 10/pkg
- 840XLF ADCUFF™ & Bladder, 2-Tube, Latex Free, Black, Large Adult
- 841 ADCUFF™ & Bladder, 2-Tube, Latex, Gray, Adult
- 841LF ADCUFF™ & Bladder, 2-Tube, Latex Free, Gray, Adult
- 844 ADCUFF™ & Bladder, 2-Tube, Latex, Green, Adult
- 844LF ADCUFF™ & Bladder, 2-Tube, Latex Free, Green, Adult
- 845 ADCUFF™ & Bladder, 2-Tube, Latex, Red, Adult
- 845LF ADCUFF™ & Bladder, 2-Tube, Latex Free, Red, Adult
- 846 ADCUFF™ & Bladder, 2-Tube, Latex, Orange, Adult
- 846C ADCUFF™ & Bladder, 2-Tube, Latex, Orange, Child

**ADC® ADCUFF™ & BLADDERS, 2-TUBE (Continued)**


- 846CLF ADCUFF™ & Bladder, 2-Tube, Latex Free, Orange, Child
- 846I ADCUFF™ & Bladder, 2-Tube, Latex, Orange, Infant
- 846ILF ADCUFF™ & Bladder, 2-Tube, Latex Free, Orange, Infant
- 846LF ADCUFF™ & Bladder, 2-Tube, Latex Free, Orange, Adult
- 846T ADCUFF™ & Bladder, 2-Tube, Latex, Orange, Thigh
- 846TLF ADCUFF™ & Bladder, 2-Tube, Latex Free, Orange, Thigh
- 846X ADCUFF™ & Bladder, 2-Tube, Latex, Orange, Large Adult
- 846XLF ADCUFF™ & Bladder, 2-Tube, Latex Free, Orange, Large Adult



875PLF-10



8971-10



8972-10

**ADC® ADCUFF™ BLOOD PRESSURE CUFF ACCESSORIES**

- 611-20-10 Y Adaptor, No Fittings For Use w/two tube cuffs, 10/pkg
- 875PLF-10 Bulb & Valve, S.P.U., 10/pkg
- 891 Luer Connector Set, Metal, Includes 891F & 891M
- 891F Luer Connector, Female, Aluminum
- 891F-10 Female Luer Slip Connector w/barbed end for use w/ADC & Baum (metal), 10/pkg
- 891M Luer Connector, Male Metal
- 891M-10 Male Luer Slip Connector w/barbed end for use w/ADC & Baum (metal), 10/pkg
- 8971-10 Screw Adaptor w/barbed end for use w/Dynamap, 10/pkg
- 8972-10 Male luer slip connector w/barbed end (plastic), 10/pkg
- 8973-10 Female Locking Luer Connector w/barbed end (plastic), 10/pkg
- 8974-10 Bayonet Connector For Use w/Hewlett Packard & Siemens Monitors, 10/pkg
- 8975-10 Neonatal Cuff Connector, Male Luer Slip Fitting For 3/32" ID Tubing, 10/pkg
- 920-011 Double Male Connector
- 920-011-10 Double Male Connector For 5/32" ID Tubing (plastic), 10/pkg
- 920-039 Filter/Diaphragm/Gasket Set



860



860BD

**ADC® ADCUFF™ INFLATION SYSTEMS**

Includes ADCUFF™ nylon cuff featuring the Size Guide™ marking system. Latex or latex free inflation bladder with two 55cm tubes. ADFLOW™ bulb and valve with filter screen protection. For use with most aneroid or mercurial blood pressure instruments. Individually boxed. Two-year warranty.

- 860 Inflation System, Latex, Black, Adult
- 860BD Inflation System, Latex, Burgundy, Adult
- 860BDF Inflation System, Latex Free, Burgundy, Adult
- 860C Inflation System, Latex, Black, Child
- 860CLF Inflation System, Latex Free, Black, Child
- 860DB Inflation System, Latex, Navy, Adult
- 860DBC Inflation System, Latex, Navy, Child
- 860DBCLF Inflation System, Latex Free, Navy, Child
- 860DBI Inflation System, Latex, Navy, Infant
- 860DBILF Inflation System, Latex Free, Navy, Infant
- 860DBLF Inflation System, Latex Free, Navy, Adult
- 860DBN Inflation System, Latex, Navy, Newborn
- 860DBT Inflation System, Latex, Navy, Thigh
- 860DBTLF Inflation System, Latex Free, Navy, Thigh
- 860DBX Inflation System, Latex, Navy, Large Adult
- 860DBXLF Inflation System, Latex Free, Navy, Large Adult
- 860DG Inflation System, Latex, Dark Green, Adult
- 860DGLF Inflation System, Latex Free, Dark Green, Adult
- 860I Inflation System, Latex, Black, Infant
- 860ILF Inflation System, Latex Free, Black, Infant
- 860LF Inflation System, Latex Free, Black, Adult
- 860M Inflation System, Latex, Magenta, Adult
- 860MLF Inflation System, Latex Free, Magenta, Adult
- 860N Inflation System, Latex, Black, Newborn
- 860RB Inflation System, Latex, Royal Blue, Adult
- 860RBLF Inflation System, Latex Free, Royal Blue, Adult
- 860T Inflation System, Latex, Black, Thigh
- 860TL Inflation System, Latex, Teal, Adult
- 860TLF Inflation System, Latex Free, Black, Thigh
- 860TLLF Inflation System, Latex Free, Teal, Adult
- 860V Inflation System, Latex, Purple, Adult
- 860VLF Inflation System, Latex Free, Purple, Adult
- 860X Inflation System, Latex, Black, Large Adult
- 860XLF Inflation System, Latex Free, Black, Large Adult
- 861 Inflation System, Latex, Gray, Adult
- 861LF Inflation System, Latex Free, Gray, Adult
- 864 Inflation System, Latex, Green, Adult
- 864LF Inflation System, Latex Free, Green, Adult
- 865 Inflation System, Latex, Red, Adult
- 865LF Inflation System, Latex Free, Red, Adult
- 866 Inflation System, Latex, Orange, Adult
- 866C Inflation System, Latex, Orange, Child
- 866CLF Inflation System, Latex Free, Orange, Child
- 866I Inflation System, Latex, Orange, Infant
- 866ILF Inflation System, Latex Free, Orange, Infant
- 866LF Inflation System, Latex Free, Orange, Adult
- 866T Inflation System, Latex, Orange, Thigh
- 866TLF Inflation System, Latex Free, Orange, Thigh
- 866X Inflation System, Latex, Orange, Large Adult
- 866XLF Inflation System, Latex Free, Orange, Large Adult

**ADC® ADCUFF™ NYLON CUFF**

Featuring the Size Guide™ marking system. For use with most aneroid or mercurial blood pressure instruments. Two year warranty.

- 830 Nylon Cuff, Black, Adult
- 830BD Nylon Cuff, Burgundy, Adult
- 830C Nylon Cuff, Burgundy, Child
- 830CH Nylon Cuff, Black, Child, 10/pkg
- 830DB Nylon Cuff, Navy, Adult
- 830DBC Nylon Cuff, Navy, Child
- 830DBI Nylon Cuff, Navy, Infant
- 830DBN Nylon Cuff, Navy, Newborn
- 830DBT Nylon Cuff, Navy, Thigh
- 830DBX Nylon Cuff, Navy, Large Adult
- 830DG Nylon Cuff, Dark Green, Adult
- 830H Nylon Cuff, Black, Adult, 10/pkg
- 830I Nylon Cuff, Black, Infant
- 830IH Nylon Cuff, Black, Infant, 5/pkg
- 830M Nylon Cuff, Magenta, Adult
- 830N Nylon Cuff, Black, Newborn
- 830NH Nylon Cuff, Black, Newborn, 5/pkg
- 830RB Nylon Cuff, Royal Blue, Adult
- 830T Nylon Cuff, Black, Thigh
- 830TH Nylon Cuff, Black, Thigh, 5/pkg
- 830TL Nylon Cuff, Teal, Adult
- 830V Nylon Cuff, Purple, Adult
- 830X Nylon Cuff, Black, Large Adult
- 830XH Nylon Cuff, Black, Large Adult, 10/pkg
- 831 Nylon Cuff, Gray, Adult
- 834 Nylon Cuff, Green, Adult
- 835 Nylon Cuff, Red, Adult
- 836 Nylon Cuff, Orange, Adult
- 836C Nylon Cuff, Orange, Child
- 836I Nylon Cuff, Orange, Infant
- 836T Nylon Cuff, Orange, Thigh
- 836X Nylon Cuff, Orange, Large Adult

**ADC® INFLATION BULBS & VALVES**

Dip molded latex or latex free. For use with most aneroid or mercurial blood pressure instruments. End valve included. Individually boxed. Two-year warranty.

- 870 ADFLOW™ Bulb w/filter screen end valve, Latex
- 870-1 Filter Screen End Valve Only
- 870H ADFLOW™ Bulb w/filter screen end valve, Latex, 10/pkg
- 870HLF ADFLOW™ Bulb w/filter screen end valve, Latex Free, 10/pkg
- 870L Large ADFLOW™ Bulb w/filter screen end valve, Latex
- 870LF ADFLOW™ Bulb w/filter screen end valve, Latex Free
- 870LH Large ADFLOW™ Bulb w/filter screen end valve, Latex, 10/pkg
- 872 ADFLOW™ Bulb & Valve (870 & 871), Latex
- 872H ADFLOW™ Bulb & Valve (870 & 871), Latex, 10/pkg
- 872HLF ADFLOW™ Bulb & Valve (870 & 871), Latex Free, 10/pkg
- 872L ADFLOW™ Bulb & Valve (870 & 871), Large, Latex
- 872LF ADFLOW™ Bulb & Valve (870 & 871), Latex Free
- 872LH ADFLOW™ Bulb & Valve (870 & 871), Large, Latex, 10/pkg
- 873 Standard Bulb w/regular end valve, Latex
- 873-1 Standard End Valve
- 873LF Standard Bulb w/regular end valve, Latex Free
- 875 Standard Bulb & Valve, Latex
- 875LF Standard Bulb & Valve, Latex Free



## OMRON CUFF AND BLADDER REPLACEMENTS

|          |   |
|----------|---|
| 21-071   | Calibrated Nylon Cuff, Blue, Adult                      |
| 21-063   | Cotton Cuff, Adult, Gray                                |
| 21-063C  | Cotton Cuff, Child, Gray                                |
| 21-063X  | Cotton Cuff, Large Adult, Gray                          |
| 21-071C  | Calibrated Nylon Cuff, Blue, Child                      |
| 21-071I  | Calibrated Nylon Cuff, Blue, Infant                     |
| 21-071L  | Calibrated Nylon Cuff, Blue, Thigh                      |
| 21-071X  | Calibrated Nylon Cuff, Blue, Large Adult                |
| 21-073   | Calibrated Nylon Cuff, Gray, Adult                      |
| 21-073C  | Calibrated Nylon Cuff, Gray, Child                      |
| 21-073I  | Calibrated Nylon Cuff, Gray, Infant                     |
| 21-073L  | Calibrated Nylon Cuff, Gray, Thigh                      |
| 21-073X  | Calibrated Nylon Cuff, Gray, Large Adult                |
| 21-081   | Double Tube Bladder, Adult                              |
| 21-081NL | Double Tube Bladder, Adult, Latex Free (Neoprene)       |
| 21-082   | Double Tube Bladder, Infant                             |
| 21-082NL | Double Tube Bladder, Infant, Latex Free (Neoprene)      |
| 21-083   | Double Tube Bladder, Child                              |
| 21-083NL | Double Tube Bladder, Child, Latex Free (Neoprene)       |
| 21-084   | Double Tube Bladder, Large Adult                        |
| 21-084NL | Double Tube Bladder, Large Adult, Latex Free (Neoprene) |
| 21-085   | Double Tube Bladder, Thigh                              |
| 21-085NL | Double Tube Bladder, Thigh, Latex Free (Neoprene)       |
| 21-091   | Single Tube Bladder, Adult                              |
| 21-092   | Single Tube Bladder, Infant                             |
| 21-093   | Single Tube Bladder, Child                              |
| 21-094   | Single Tube Bladder, Large Adult                        |
| 21-095   | Single Tube Bladder, Thigh                              |

## OMRON CUFF AND BLADDER SETS

|           |  |
|-----------|--|
| 21-031    | Calibrated Blue Nylon Cuff w/double-tube bladder, Adult  |
| 21-031C   | Calibrated Blue Nylon Cuff w/double-tube bladder, Child  |
| 21-031CNL | Calibrated Blue Nylon Cuff w/double-tube bladder, Child, Latex Free (Neoprene)                   |
| 21-031I   | Calibrated Blue Nylon Cuff w/double-tube bladder, Infant   |
| 21-031INL | Calibrated Blue Nylon Cuff w/double-tube bladder, Infant, Latex Free (Neoprene)                  |
| 21-031L   | Calibrated Blue Nylon Cuff w/double-tube bladder, Thigh  |
| 21-031LNL | Calibrated Blue Nylon Cuff w/double-tube bladder, Thigh, Latex Free (Neoprene)                   |
| 21-031NB  | Calibrated Blue Nylon Cuff w/double-tube bladder, Newborn  |
| 21-031NL  | Calibrated Blue Nylon Cuff w/double-tube bladder, Adult, Latex Free (Neoprene)                   |
| 21-031X   | Calibrated Blue Nylon Cuff w/double-tube bladder, Large Adult                                    |
| 21-031XNL | Calibrated Blue Nylon Cuff w/double-tube bladder, Large Adult, Latex Free (Neoprene)             |
| 21-033    | Calibrated Gray Nylon Cuff w/double-tube bladder, Adult  |
| 21-033C   | Calibrated Gray Nylon Cuff w/double-tube bladder, Child  |
| 21-033I   | Calibrated Gray Nylon Cuff w/double-tube bladder, Infant   |
| 21-033X   | Calibrated Gray Nylon Cuff w/double-tube bladder, Large Adult                                    |
| 21-037    | Calibrated Black Nylon Cuff w/double-tube bladder, Adult   |
| 21-037C   | Calibrated Black Nylon Cuff w/double-tube bladder, Child   |
| 21-037X   | Calibrated Black Nylon Cuff w/double-tube bladder, Large Adult                                   |
| 21-053    | Calibrated Gray Nylon Cuff w/single-tube bladder, Adult  |
| 21-053C   | Calibrated Gray Nylon Cuff w/single-tube bladder, Child  |
| 21-053I   | Calibrated Gray Nylon Cuff w/single-tube bladder, Infant   |
| 21-053L   | Calibrated Gray Nylon Cuff w/single-tube bladder, Thigh  |
| 21-053X   | Calibrated Gray Nylon Cuff w/single-tube bladder, Large Adult                                    |
| 21-201    | Calibrated Blue Nylon Cuff w/double-tube bladder, inflation bulb, air release valve, Adult       |
| 21-201C   | Calibrated Blue Nylon Cuff w/double-tube bladder, inflation bulb, air release valve, Child       |
| 21-201I   | Calibrated Blue Nylon Cuff w/double-tube bladder, inflation bulb, air release valve, Infant      |
| 21-201L   | Calibrated Blue Nylon Cuff w/double-tube bladder, inflation bulb, air release valve, Thigh       |
| 21-201X   | Calibrated Blue Nylon Cuff w/double-tube bladder, inflation bulb, air release valve, Large Adult |

## OMRON INFLATION BULBS

|          |   |
|----------|---|
| 21-122R  | Inflation Bulb For 107M Palm Aneroid & R-1077 |
| 21-125   | Inflation Bulb w/end valve only               |
| 21-125NL | Inflation Bulb, Latex Free w/end valve only   |
| 21-140   | Inflation Bulb w/standard air release valve   |
| 21-141   | Inflation Bulb w/deluxe air release valve     |

## OMRON MANOMETER PARTS & ACCESSORIES

|        |  |
|--------|--|
| 0179   | Calibration Tool For Standard Aneroids |
| 0180   | Gauge Guard, Blue                      |
| 148    | Standard Aneroid Gauge, Gray           |
| 156    | Marshall® Professional Gauge For 108M  |
| 21-001 | Exactus II Gauge, Blue                 |
| 21-150 | Crystal For Standard Aneroids          |

## OMRON MERCURIAL SPHYGMOMANOMETER PARTS

|          |   |
|----------|---|
| 21-321   | 8 ft Latex Coiled Tubing                  |
| 21-321NL | Coiled Tubing 8 ft, Latex Free (Neoprene) |
| 21-347   | BMS Wall Basket                           |

## TRIMLINE BLADDERLESS SOFT CUFFS, LATEX FREE

Economical, soft fabric blood pressure cuffs are comfortable for patients, even during lengthy procedures. Each size is printed in a different color print for easy identification. Cuffs feature the exclusive Pre-Gaged® BALANCED™ design with an exterior sizing guide and a "Use Larger Cuff" area to help prevent miscuffing, and a longer bladder length to facilitate right or left arm use. Cuffs and tubing are completely latex free. Can be used with all manual and NIBP monitors; connectors/adapters available. Conforms to American Heart Association and AAMI recommendations. Made in the USA.


|        |   |
|--------|---|
| 8602XL | Two Tubes Cuff w/bulb, valve & clamp, Large Adult, Extra Long, Red Print, 5/bx          |
| 8900   | Two Tubes cuff w/no connectors (Add "S" to # For Single Tube), Thigh, Brown Print, 5/bx |



5080-01



5082-01H



5082-43



5086-01



5087-01



5090-65

## WELCH ALLYN ANEROID ACCESSORIES & PARTS

|          |                                 |
|----------|---------------------------------|
| 5080-01  | Adult, Non-Conductive Hook Cuff |
| 5082-01  | Adult Non-Conductive Cuff       |
| 5082-01H | Adult Cuff, Latex Free, 12/bx   |
| 5082-02  | Child Non-Conductive Cuff       |
| 5082-03  | Infant Cuff w/1-tube bag        |
| 5082-06  | Adult Non-Conductive Cuff, Blue |
| 5082-07  | Newborn Cuff w/1-tube bag       |
| 5082-08  | Infant Cuff w/2-tube bag        |
| 5082-11  | Thigh Non-Conductive Cuff       |



**WELCH ALLYN ANEROID ACCESSORIES & PARTS**

**(Continued)**

- 5082-16 Large Adult Non-Conductive Cuff
- 5082-16H Large Adult Cuff, Latex Free, 6/bx
- 5082-21 Child Inflation System w/2 tube bag, Latex Free
- 5082-21H Inflation System, Child, Latex Free, 6/bx
- 5082-22 Adult Inflation System w/2 tube bag, Latex Free
- 5082-22H Inflation System, Adult, Latex Free, 6/bx
- 5082-23 Large Adult Inflation System w/2 tube bag, Latex Free
- 5082-23H Inflation System, Large Adult, Latex Free, 6/bx
- 5082-24 Child Cuff w/2-tube bag
- 5082-25 Adult Cuff w/2-tube bag
- 5082-257 Integrated Gauge Adapter Only
- 5082-26 Large Adult Cuff w/2-tube bag
- 5082-42 Child Cuff w/1-tube bag
- 5082-43 Adult Cuff w/1-tube bag
- 5082-44 Large Adult Cuff w/1-tube bag
- 5082-45 Thigh Inflation System
- 5085-07 Sphyg Carrying Case, Polyester
- 5085-08 Two-Sided Sphyg Carrying Case, Polyester
- 5086-01 Latex Inflation Bulb
- 5086-01H Latex Inflation Bulb, Latex Free, 12/bx
- 5086-03H Large Inflation Bulb, Latex Free, 6/bx
- 5087-01 Air Release Valve
- 5087-01H Air Release Valve, Latex Free, 12/bx
- 5087-03 Valve Thumb Screw For Hand Model (pre-1979)
- 5087-04 Valve Assembly For Hand Model (pre-1979)
- 5087-08 Valve Assembly For Hand Model (after Jan. 1, 1979)
- 5087-09 Valve Thumb Screw For Hand Model (after Jan. 1, 1979)
- 5088-01 Air Release Valve & Bulb
- 5088-01H Bulb & Valve Combination, Latex Free, 12/bx
- 5089-01 Adult Inflation Bags, Non-Conductive, 2 Tubes
- 5089-01H Inflation Bag, Adult, 2 Tubes, Latex Free, 12/bx
- 5089-02 Child Inflation Bags, Non-Conductive, 2 Tubes
- 5089-02H Inflation Bag, Child, 2 Tubes, Latex Free, 12/bx
- 5089-03 Adult Inflation Bags, Non-Conductive, 1 Tube
- 5089-04 Child Inflation Bags, Non-Conductive, 1 Tube
- 5089-09 Metal Y Adaptor
- 5089-11 Metal Twist Fit Connector Set
- 5089-12 Tubing 27" Length
- 5089-13H Coil Tubing 8 ft Length, Latex Free, 12/bx
- 5089-18 Thigh Inflation Bags, Non-Conductive, 1 Tube
- 5089-19 Thigh Inflation Bags, Non-Conductive, 2 Tubes
- 5089-19H Inflation Bag, Thigh, 2 Tubes, Latex Free, 6/bx
- 5089-21 Large Adult Inflation Bags, Non-Conductive, 1 Tube
- 5089-22 Large Adult Inflation Bags, Non-Conductive, 2 Tubes
- 5089-22H Inflation Bag, Large Adult, 2 Tubes, Latex Free, 6/bx
- 5089-25 Infant Inflation Bags, Non-Conductive, 1 Tube
- 5089-26 Infant Inflation Bags, Non-Conductive, 2 Tubes
- 5090-29 Bulb End Valve
- 5090-30 Screen Assembly
- 5090-44 Front Ring Bezel For Aneroid
- 5090-45 Shatterproof Crystal For Aneroid
- 5090-64 Clip Assembly For Pocket Aneroid
- 5090-65 Sphyg Repair Kit
- 5091-16 Crystal For Wall Gauge, Square
- 5091-33 Cuff Clip
- 5097-15 Filter & Gasket Set For Mercurial Gauge
- 5097-16 Tube Cap, Wall Mercurial
- 5097-25 Incline Tube, Mobile Mercurial
- 5097-32 Vertical Tube, Wall Mercurial, Marked 4.2 mm
- 5097-33 Incline Tube, Mobile Mercurial, Marked 4.2 mm
- 5098-07 Luer Lock (set)
- 5098-17 Housing Screws For Hand Gauge (Packet of 12)
- 7052-25 Printer Paper For 7052-23 Only, 5 rolls/pack
- DS44A Model DS44 Gauge w/integrated gauge adapter only
- DS45A Model DS45 Gauge w/integrated gauge adapter only



**WELCH ALLYN TYCOS® DS66 HAND ANEROID ACCESSORIES**

- 5082-18 Child Cuff, Print, Latex Free
- 5082-56 Adult Cuff & Black Inflation Bag, Latex Free
- 5082-58 Large Adult Cuff & Black Inflation Bag, Latex Free
- 5082-59 Adult Cuff, Black, Latex Free
- 5082-61 Large Adult Cuff, Black, Latex Free
- 5082-62 Child Print Cuff & Black Inflation Bag, Latex Free
- 5086-04 Bulb Assembly, Latex Free
- 5098-31 Gauge, Latex Free
- 5098-40 Crystal, Latex Free

**BCI MINI-TORR PLUS® NON-INVASIVE BLOOD PRESSURE MONITOR**

Provides a base to build on for your monitoring needs. Offers cost-effective, non-invasive blood pressure monitoring. Options which can be added include pulse oximetry, a built-in printer, oral or rectal temperature, and/or a lightweight roll stand or IV pole mount. The ability to combine the options needed in your clinical environment at an affordable price, makes this monitor an important addition to your equipment needs.

- 6004-000 Mini-Torr Plus® NIBP Monitor
- 6004-001 Mini-Torr Plus® NIBP/Oximeter
- 6004-004 Mini-Torr Plus® NIBP/Electronic Thermometry
- 6004-005 Mini-Torr Plus® NIBP/Printer
- 6004-006 Mini-Torr Plus® NIBP/Oximeter/Printer
- 6004-009 Mini-Torr Plus® NIBP/Electronic Thermometry/Printer
- 6004-011 Mini-Torr Plus® NIBP/SpO2/Electronic Thermometry
- 6004-013 Mini-Torr Plus® NIBP/SpO2/Thermometry/Printer



143265



143268



90121-1

**BURDICK AMBULATORY BLOOD PRESSURE SYSTEMS**

Provides unmatched accuracy for identifying and managing your hypertensive patients.

- 143265 Ultralite ABP Monitor w/software
- 143266 Ultralite ABP Monitor w/report generator
- 143268 ABP Monitor w/report generator
- 90121-1 ABP Report Management System Base Station Interface (for Windows), Features Remote Communication Capability For Data Retrieval & Monitor Programming via Modem
- 90239A ABP Report Generator/Printer
- 143267 ABP Monitor w/software



## BURDICK AMBULATORY BLOOD PRESSURE SYSTEM ACCESSORIES

- 000-0027-02 Patient Diary For ABP
- 006-0210-02 Thermal Paper for 90239A ABP Report Generator, 57mm x 32, 6 rl/bx
- 012-0034-03 90239-80 Interface Cable (90239-80 to PC Serial Port)
- 012-0096-01 Modem Interface Cable (ABP Monitor to Modem)
- 012-0097-02 ABP Interface Cable (ABP Monitor to PC)
- 012-0105-02 Base Station Interface Cable (PC to Modem)
- 015-0070-00 Cuff Support Harness For All Cuffs
- 016-0077-01 Large Adult Cuff, 32-42cm
- 016-0080-00 Waist Belt for ABP
- 016-0109-01 XL Adult Cuff, 38-50cm (includes cuff support harness)
- 016-0262-00 Strap For Use w/vinyl pouch for model 90207 ABP
- 016-0340-00 Cloth Pouch for Model 90217 ABP
- 016-0410-00 Cloth Pouch for Model 90207 ABP (SL logo)
- 040-0549-00 Data Management Option Kit, Permits Storage of up to 40 Reports on the Report Generation Data Card
- 146-0008-00 Lithium Battery, 3-Volt, For Use w/90207 ABP
- 672-0167-00 Report Generation Data Card



90207-1Q



90217-1Q

## BURDICK AMBULATORY BLOOD PRESSURE SYSTEMS WITH QUICK CONNECT

Our ABP monitors are among the smallest and lightest weight on the market today, with advanced reporting software that allows rapid analysis of patient information. Extensive testing by independent clinicians has validated ABPM accuracy time and time again, making them widely recognized and the "gold standard" for ABP. Both the 90217 Ultralite and 90207 ABP monitors feature a compact, lightweight design that optimizes patient comfort—a key factor in successful ABP monitoring. Comfort helps ensure patient compliance and, consequently, may decrease the number of artifactual readings. Collects measurements for a period of 24 hours or up to 8 days. Lightweight recorder operates quietly for patient comfort. Oscillometry measures systolic, diastolic and mean blood pressure. No chest electrodes or microphones reduce operator costs. Uses a real-time clock to facilitate diary notations.

- 90207-1Q ABP Monitor w/quick Connect, English
- 90217-1Q ABP Ultralite Monitor w/quick Connect, English

## BURDICK CUFFS WITH QUICK CONNECT

- 015-0067-01Q Small Adult 17-26cm
- 015-0068-02Q Standard Adult 24-32cm
- 015-0118-01Q Child 12-20cm



## WELCH ALLYN AMBULATORY BLOOD PRESSURE MONITOR ACCESSORIES

- 5100-30 Microphone/Pneumatic Tubing Assembly
- 5100-31 Large Cuff & Pneumatic System Assembly
- 5100-32 Medium Cuff & Pneumatic System Assembly
- 5100-34 Large Cuff & Bladder
- 5100-35 Medium Cuff & Bladder
- 5100-36 Small Cuff & Bladder
- 5100-37 Large Cuff
- 5100-38 Medium Cuff
- 5100-41E Patient Hook-Up Kit, 25/bx
- 5100-42E Patient Diary, 50/bx
- 5100-43 Microphone Anchoring Pads, 50/bx
- 5100-44 Cuff Anchoring Pads, 50/bx
- 5100-45 Cuff Anchoring Straps, 5/bx
- 5100-46 Carrying Pouch
- 5100-47 Waist Belt
- 5100-48 Shoulder Strap
- 5100-49 Calibration T-Connector
- 5100-51 Small Bladder
- 5100-52 Medium Bladder
- 5100-53 Large Bladder
- 5100-61 Q-Trak PC Software (includes PC cable)

## WELCH ALLYN ATLAS MONITOR

- 622NP-E1 Monitor with Printer



4200B-E1



4200B-E1-M

## WELCH ALLYN SPOT VITAL SIGNS

Standard non-invasive blood pressure (NIBP) as well as pulse rate and MAP. Large LCD display with backlight. Extremely user-friendly and easy to use. Uses the oscillometric method to determine systolic and diastolic pressure, MAP and pulse rate. Lightweight and portable: weights just 4 pounds. Training video included.

- 4200B-E1 Device w/BP, Domestic/English
- 4200B-E1-M Device w/BP & Mobile Stand, Domestic/English
- 420TB-E1 Device w/BP & Temperature, Domestic/English
- 420TB-E1-M Device w/BP & Temperature & Mobile Stand, Domestic/English
- 42N0B-E1 Device w/BP & Pulse Oximeter, Domestic/English
- 42N0B-E1-M Device w/BP & Pulse Oximeter, Mobile Stand, Domestic/English
- 42NTB-E1 Device w/BP, Pulse Oximeter & Temperature, Domestic/English
- 42NTB-E1-M Device w/BP, Pulse Oximeter, Temperature & Mobile Stand, Domestic/English

**WELCH ALLYN SPOT VITAL SIGNS ACCESSORIES**

- 4200-100 Carrying Case
- 4200-145E Technical Service Manual
- 4200-60 Mobile Stand
- 4200-62 Wall Mount
- 4200-64 IV Pole Mount
- 4200-70 Anti-Theft Kit
- 4200-84 Lead Acid Battery, 6V
- 5200-01 Adult Cuff Assembly, Latex Free (Cuff, Bladder & Connector)
- 5200-02 Large Adult Cuff Assembly, Latex Free (Cuff, Bladder & Connector)
- 5200-03 Child Print Cuff Assembly, Latex Free (Cuff, Bladder & Connector)
- 5200-10 Thigh Cuff Assembly, Latex Free (Cuff, Bladder & Connector)
- D-YSPD Pedicheck Pediatric Spot Check




**GE MEDICAL DINAMAP® PROCARE PATIENT MONITORS**

More than 3 million times every single day critical patient care decisions are made based on blood pressure determinations using DINAMAP® technology. With the ProCare Series, you can take BP measurements that are fast, accurate and comfortable. Get reliable temperature readings in seconds and get oxygen saturation readings with your choice of leading SpO2 technologies.

- DPC120X-US DINAMAP ProCare 120 Monitor, NIBP, Heart Rate, & Printer; Typical Accessories Included At No Additional Charge; Two (2) Year Parts & Labor Warranty
- DPC220X-US DINAMAP ProCare 220 Monitor, NIBP, Heart Rate, IVAC® TurboTemp®, & Printer; Typical Accessories Included At No Additional Charge; Two (2) Year Parts & Labor Warranty
- DPC320N-US DINAMAP ProCare 320 Monitor w/Nellcor SpO2, NIBP, Heart Rate, Nellcor® Oximax® SpO2 & Printer; Typical Accessories Included At No Additional Charge; Two (2) Year Parts & Labor Warranty
- DPC420N-US DINAMAP ProCare 420 Monitor w/Nellcor SpO2, NIBP, Heart Rate, IVAC® TurboTemp®, Nellcor® Oximax® SpO2 & Printer; Typical Accessories Included At No Additional Charge; Two (2) Year Parts & Labor Warranty

**Accessories**


- 633178CR 12 Volt Lead Acid Battery For Use w/DPC420N-US



**GE MEDICAL DINAMAP TEMPERATURE PROBES & PROBE COVERS**

Reusable temperature probes. Accurate temperature reading. Probe covers reduce potential for cross-contamination between patients. Latex free.


- 088012 Oral Temperature Probe, Blue
- 088015 Electronic Probe Cover, 20/bx, 250 bx/cs
- 088016 Rectal Temperature Probe, Red
- 2008774-001 TurboTemp Oral Probe, Blue
- 2008775-001 TurboTemp Rectal Probe, Red



1026BOO




760GARCC



768-641PEA

**ADC® CARTOON CHARACTER PRODUCTS**

- 1026BOO Nurse Bag, Betty Boop
- 1026GAR Nurse Bag, Garfield
- 1026HK Nurse Bag, Hello Kitty
- 1026PEA Nurse Bag, Peanuts
- 760BOOCQ Child, Aneroid, Betty Boop
- 760GARCC Child, Aneroid, Garfield
- 760HKCC Child, Aneroid, Hello Kitty
- 760PEACQ Child, Aneroid, Peanuts
- 760BOOC Adult, Aneroid, Betty Boop
- 760GARQ Adult, Aneroid, Garfield
- 760HKQ Adult, Aneroid, Hello Kitty
- 760PEAQ Child, Aneroid, Peanuts
- 760BOOXQ Large Adult, Aneroid, Betty Boop
- 760GARXQ Large Adult, Aneroid, Garfield
- 760PEAXQ Large Adult, Aneroid, Peanuts
- 768-641BOO Adult, Aneroid, Betty Boop
- 768-641GAR Adult, Aneroid, Garfield
- 768-641HK Adult, Aneroid, Hello Kitty
- 768-641PEA Adult, Aneroid, Peanuts
- 888BOO Nylon Zipper Case, Betty Boop
- 888GAR Nylon Zipper Case, Garfield
- 888HK Nylon Zipper Case, Hello Kitty
- 888PEA Nylon Zipper Case, Peanuts



**ADC® ADSCOPE™ 665 DISPOSABLE STETHOSCOPE**

Nurse style diaphragm only single patient use scope. Acrylesin construction throughout is super-lightweight. Non-toxic. ETO sterilizable. May be safely incinerated. Weighs about 1.5 oz for comfort, yet is remarkably durable. Ultra sensitive diaphragm for outstanding acoustic performance. Individually packaged in a sealed polybag. Available in yellow and royal blue. Latex free.

- 665RB-10 Single Patient Use Scope, Royal Blue, 10/ctn
- 665RB-100 Single Patient Use Scope, Royal Blue, 100/cs
- 665Y-10 Single Patient Use Scope, Yellow, 10/ctn
- 665Y-100 Single Patient Use Scope, Yellow, 100/cs

**Accessories**

- 698BK Stethoscope ID Tag
- 697BK Premium Stethoscope ID Tag

**ADC PROSCOPE™ 664Y DISPOSABLE STETHOSCOPE**

Diaphragm only disposable ABS chestpiece with non-chill rim for patient comfort. 22" PVC Yellow Y-tubing. Chrome-plated brass binaurals. Weighs 2.5 oz, measures 32".

- 664Y Proscope™ 664Y, Disposable Scope



**WELCH ALLYN/TYCOS UNISCOPE DISPOSABLE STETHOSCOPE**

- 17460 Adult Stethoscope, Green
- 17461 Adult Stethoscope, Yellow
- 17462 Adult Stethoscope, Dk. Blue
- 17462P Pediatric Stethoscope, Lt. Blue
- 17463 Adult Stethoscope, Orange
- 17464 Adult Stethoscope, Lt. Blue
- 17465 Adult Stethoscope, Green, Anesthesia



**OMRON DISPOSABLE NURSE STETHOSCOPE**

Lightweight plastic scopes are ideal for single-patient use in contagious areas. Sound transmission is comparable to a standard nurse scope.

- 0300BLUE Blue Nurse Stethoscope
- 0300RED Red Nurse Stethoscope
- 0300YELLOW Yellow Nurse Stethoscope
- 13-176H Yellow Nurse Stethoscope w/aluminum binaurals, 50/cs
- HSS-490 Lightweight Nurse Stethoscope



**CIDA SINGLE HEAD NURSE STETHOSCOPES**

Diaphragm only aluminum chestpiece. Adjustable chrome-plated binaurals. Flexible 22" PVC Y tubing. White plastic eartips. Overall length 31". Latex free.

- 37-1000 Gray Stethoscope
- 37-1010 Navy Stethoscope
- 37-1014 Light Blue Stethoscope
- 37-1017 Royal Blue Stethoscope
- 37-1020 Gold Stethoscope
- 37-1034 Teal Stethoscope
- 37-1040 Green Stethoscope
- 37-1047 Neon Green Stethoscope
- 37-1048 Dark Green Stethoscope
- 37-1050 Red Stethoscope
- 37-1052 Burgundy Stethoscope
- 37-1055 Pink Stethoscope
- 37-1056 Purple Stethoscope
- 37-1057 Neon Pink Stethoscope
- 37-1060 White Stethoscope
- 37-1070 Black Stethoscope
- 37-1095 Neon Orange Stethoscope
- 37-1100 Orchid Stethoscope
- 37-1110 Lavender Stethoscope
- 37-1145 Yellow Stethoscope
- 37-1147 Neon Yellow Stethoscope
- 37-1150 Magenta Stethoscope



**OMRON LIGHTWEIGHT SINGLE HEAD NURSE STETHOSCOPES**

Match a favorite uniform with one of eleven brightly colored nurse stethoscopes. Silver chestpiece. These colorful stethoscopes are useful for color coding departments or identifying contagious areas.

- 411BLK Nurse Stethoscope w/black tubing
- 411BLUE Nurse Stethoscope w/blue tubing
- 411RED Nurse Stethoscope w/red tubing



**WELCH ALLYN SPECTRUM™ SINGLE-HEAD STETHOSCOPES**

Excellent for taking blood pressures and general auscultation.

- 5079-230 38" Adult Black Stethoscope, Anesthesiology w/taper lock

**MILTEX FORD STETHOSCOPE**

Shallow bell and plastic ear tips. Chrome.

- 1-4 Ford Stethoscope



**A Strong Partner  
in Healthcare  
Distribution**






**OMRON BOWLES STETHOSCOPE**

Most economical professional stethoscope features a chrome-plated binaural and black PVC tubing. Standard single stem.

0407 Black Bowles Stethoscope



**3M™ LITTMANN® MASTER CLASSIC II STETHOSCOPE**

The excellent acoustic sensitivity of a patented tunable diaphragm is the hallmark of the low-profile Littmann Master Classic II. Award-winning design combines high performance with portable convenience in this single-sided stethoscope. Features patented Littmann Snap Tight soft-sealing eartips, patient friendly, nonchill rim and diaphragm. Three-year warranty, plus lifetime service and repair policy.

- 2141 Master Classic II Stethoscope, 27" Black Plated Chestpiece & Eartubes
- 2142G Master Classic II Stethoscope, 27" Gold-Plated Chestpiece & Eartubes
- 2143 Master Classic II Stethoscope, 27" Purple Tubing
- 2144L Master Classic II Stethoscope, 27" Black Tubing
- 2146 Master Classic II Stethoscope, 27" Burgundy Tubing
- 2147 Master Classic II Stethoscope, 27" Navy Blue Tubing
- 2630 Master Classic II Stethoscope, 27" Caribbean Blue Tubing
- 2632 Master Classic II Stethoscope, 27" Hunter Green Tubing
- 2633 Master Classic II Stethoscope, 27", Ceil Blue Tubing, ea



**3M™ LITTMANN® SELECT STETHOSCOPES**

A convenient tunable diaphragm in a single-sided coated alloy chestpiece providing reliable and clear acoustics.

- 2290 Select Stethoscope, 28" Black Tubing
- 2291 Select Stethoscope, 28" Caribbean Blue
- 2293 Select Stethoscope, 28" Burgundy Tubing
- 2294 Select Stethoscope, 28" Purple Tubing
- 2296 Select Stethoscope, 28" Raspberry Tubing
- 2298 Select Stethoscope, 28" Royal Blue Tubing
- 2302 Select Stethoscope, 28", Seafoam Green Tubing, ea
- 2303 Select Stethoscope, 28", Lilac Tubing, ea



**CIDA DUAL-HEAD STETHOSCOPES**

Combination diaphragm/bell chestpiece for wide frequency response. Bell and diaphragm equipped with non-chill ring for patient comfort. Lightweight, durable construction. Adjustable metal binaurals. Flexible 22" PVC tubing. Plastic eartips. Overall length 31 1/2". Latex free.

- 37-2000 Gray Stethoscope
- 37-2010 Navy Stethoscope
- 37-2014 Light Blue Stethoscope
- 37-2017 Royal Blue Stethoscope
- 37-2020 Gold Stethoscope
- 37-2034 Teal Stethoscope
- 37-2040 Green Stethoscope
- 37-2047 Neon Yellow Stethoscope
- 37-2048 Dark Green Stethoscope
- 37-2050 Red Stethoscope
- 37-2052 Burgundy Stethoscope
- 37-2055 Pink Stethoscope
- 37-2056 Purple Stethoscope
- 37-2057 Neon Pink Stethoscope
- 37-2060 White Stethoscope
- 37-2070 Black Stethoscope
- 37-2095 Neon Orange Stethoscope
- 37-2100 Orchid Stethoscope
- 37-2110 Lavender Stethoscope
- 37-2145 Yellow Stethoscope
- 37-2150 Magenta Stethoscope



**ADC® ADSCOPE™ 603 ADULT STETHOSCOPE**

Combination adult size chestpiece precisely machined to exacting tolerances from surgical stainless steel for outstanding performance and rugged durability. Extra large bell (1 3/8" diameter) for unsurpassed low frequency response. Ultra sensitive diaphragm (1 7/8" diameter) for greater amplification and crisper high frequency response. Color coordinate non-chill bell and diaphragm retaining rim for patient comfort. Stainless steel binaurals with double leaf internal spring fixed at 15° angle for greater comfort. Reinforcing yoke molded into flexible 22" PVC "Y" tubing. ADSOFT™ deluxe PVC eartips. Spare diaphragm and scope ID tag included. Weights 6.5 oz. Overall length 31 1/2". Lifetime quality performance guarantee. Latex free.

- 603BD Burgundy Stethoscope
- 603BK Black Stethoscope
- 603DG Dark Green Stethoscope
- 603FG Frosted Glacier Stethoscope
- 603FL Frosted Lilac Stethoscope
- 603FS Frosted Seafoam Stethoscope
- 603G Gray Stethoscope
- 603LV Lavender Stethoscope
- 603N Navy Stethoscope
- 603R Red Stethoscope
- 603RB Royal Blue Stethoscope
- 603TL Teal Stethoscope

## ADC® ADSCOPE™ 603 ADULT STETHOSCOPE ACCESSORIES

|          |  |
|----------|--|
| 603-01BD | Burgundy Chestpiece                          |
| 603-01BK | Black Chestpiece                             |
| 603-01DG | Dark Green Chestpiece                        |
| 603-01FG | Frosted Glacier Chestpiece                   |
| 603-01FL | Frosted Lilac Chestpiece                     |
| 603-01FS | Frosted Seafoam Chestpiece                   |
| 603-01G  | Gray Chestpiece                              |
| 603-01LV | Lavender Chestpiece                          |
| 603-01N  | Navy Chestpiece                              |
| 603-01R  | Red Chestpiece                               |
| 603-01RB | Royal Blue Chestpiece                        |
| 603-01TL | Teal Chestpiece                              |
| 603-04BD | Burgundy Non-Chill Ring For 603 & 609        |
| 603-04BK | Black Non-Chill Ring For 603 & 609           |
| 603-04DG | Dk. Green Non-Chill Ring For 603 & 609       |
| 603-04FG | Frosted Glacier Non-Chill Ring For 603 & 609 |
| 603-04FL | Frosted Lilac Non-Chill Ring For 603 & 609   |
| 603-04FS | Frosted Seafoam Non-Chill Ring For 603 & 609 |
| 603-04G  | Gray Non-Chill Ring For 603 & 609            |
| 603-04LV | Lavender Non-Chill Ring For 603 & 609        |
| 603-04N  | Navy Non-Chill Ring For 603 & 609            |
| 603-04R  | Red Non-Chill Ring For 603 & 609             |
| 603-04RB | Royal Blue Non-Chill Ring For 603 & 609      |
| 603-04TL | Teal Non-Chill Ring For 603 & 609            |
| 603-05BD | Burgundy Binaural/Tubing                     |
| 603-05BK | Black Binaural/Tubing                        |
| 603-05DG | Dk. Green Binaural/Tubing                    |
| 603-05FG | Frosted Glacier Binaural/Tubing              |
| 603-05FL | Frosted Lilac Binaural/Tubing                |
| 603-05FS | Frosted Seafoam Binaural/Tubing              |
| 603-05G  | Gray Binaural/Tubing                         |
| 603-05LV | Lavender Binaural/Tubing                     |
| 603-05N  | Navy Binaural/Tubing                         |
| 603-05R  | Red Binaural/Tubing                          |
| 603-05RB | Royal Blue Binaural/Tubing                   |
| 603-05TL | Teal Binaural/Tubing                         |



## ADC® ADSCOPE™ 609 ADSCOPE-LITE™ STETHOSCOPE

Combination adult chestpiece precisely machined to exacting tolerances from lightweight aluminum is remarkable durable, yet extraordinarily comfortable. Extra large 1 3/8" diameter bell. High sensitivity 1 3/4" diaphragm. Non-chill bell and diaphragm retaining rim for patient comfort. Double leaf internal spring binaurals are fixed at 15° angle for greater comfort. Custom molded 22" PVC tubing. ADSOFT™ PVC eartips and spare diaphragm included. Weights 4.5 oz. Overall length 31 1/2". Five year quality performance guarantee. Latex free.

|       |                             |
|-------|-----------------------------|
| 609BD | Burgundy Stethoscope        |
| 609BK | Black Stethoscope           |
| 609DG | Dk. Green Stethoscope       |
| 609FG | Frosted Glacier Stethoscope |
| 609FL | Frosted Lilac Stethoscope   |
| 609FS | Frosted Seafoam Stethoscope |
| 609G  | Gray Stethoscope            |
| 609LV | Lavender Stethoscope        |
| 609N  | Navy Stethoscope            |
| 609R  | Red Stethoscope             |
| 609RB | Royal Blue Stethoscope      |
| 609TL | Teal Stethoscope            |

## ADC® ADSCOPE™ 609 ADSCOPE-LITE™ STETHOSCOPE ACCESSORIES

|          |                                 |
|----------|---------------------------------|
| 609-01BD | Burgundy Chestpiece             |
| 609-01BK | Black Chestpiece                |
| 609-01DG | Dark Green Chestpiece           |
| 609-01FG | Frosted Glacier Chestpiece      |
| 609-01FL | Frosted Lilac Chestpiece        |
| 609-01FS | Frosted Seafoam Chestpiece      |
| 609-01G  | Gray Chestpiece                 |
| 609-01LV | Lavender Chestpiece             |
| 609-01N  | Navy Chestpiece                 |
| 609-01R  | Red Chestpiece                  |
| 609-01RB | Royal Blue Chestpiece           |
| 609-01TL | Teal Chestpiece                 |
| 609-03   | Diaphragm For All               |
| 609-05BD | Burgundy Binaural/Tubing        |
| 609-05BK | Black Binaural/Tubing           |
| 609-05DG | Dark Green Binaural/Tubing      |
| 609-05FG | Frosted Glacier Binaural/Tubing |
| 609-05FL | Frosted Lilac Binaural/Tubing   |
| 609-05FS | Frosted Seafoam Binaural/Tubing |
| 609-05G  | Gray Binaural/Tubing            |
| 609-05LV | Lavender Binaural/Tubing        |
| 609-05N  | Navy Binaural/Tubing            |
| 609-05R  | Red Binaural/Tubing             |
| 609-05RB | Royal Blue Binaural/Tubing      |
| 609-05TL | Teal Binaural/Tubing            |

## ADC® ADSCOPE™ 645 STETHOSCOPE

Extraordinary design for outstanding acoustic performance and rugged durability. Chestpiece drum, retaining rings, tubing clip and binaurals generously plated with 5 microns of 18 karat gold, then polished to a brilliant shine. Complete with an accessory pouch containing two extra pairs of eartips, spare diaphragms and the three bell chestpieces. Weighs 8 oz. Overall length 30". Inspected and packaged in the U.S.A. One year warranty on finish, lifetime on performance. Latex free.

645 ADCSCOPE™ 645 18 Karat Gold Plated Stethoscope

### Accessories

|         |                            |
|---------|----------------------------|
| 645-01  | Chestpiece                 |
| 645-01D | Chestpiece Drum            |
| 645-02  | Retaining Ring, Adult      |
| 645-02S | Retaining Ring, Pediatric  |
| 645-05  | Binaurals, Adult           |
| 645-05S | Binaural Spring, Pediatric |
| 645-14A | Diaphragm Assembly         |
| 645-15A | Diaphragm Assembly         |
| 645-16  | Tubing Clip Assembly       |



## OMRON LIGHTWEIGHT DUAL HEAD STETHOSCOPES

Lightweight stethoscopes feature adjustable binaurals for a comfortable fit, one-piece molded PVC tubing for good sound transmission, and built-in rubber non-chill rings. Includes extra pair of eartips and diaphragm.

|          |   |
|----------|---|
| 412BLK   | Dual Head Stethoscope, Silver Chestpiece w/black tubing |
| 412BLUE  | Dual Head Stethoscope, Blue Chestpiece w/blue tubing    |
| 412GRAY  | Dual Head Stethoscope, Silver Chestpiece w/gray tubing  |
| 412GREEN | Dual Head Stethoscope, Green Chestpiece w/green tubing  |
| 412RED   | Dual Head Stethoscope, Red Chestpiece w/red tubing      |



**ADC® ADSCOPE™ 615 PLATINUM PROFESSIONAL EDITION**

Combines the acoustic response of a traditional bell and diaphragm into a convenient one-sided chestpiece. New gasketed diaphragm enables the scope to achieve multi-frequency response simply by altering pressure. Light pressure produces bell-type low frequency response. Firm pressure will produce diaphragm-type high frequency performance. New ovoid shaped chestpiece improves ergonomics. Chestpiece made from lightweight zinc alloy plated with an attractive satin finish. Non-chill diaphragm retaining rim for patient comfort. Large bore chrome-plated brass binaurals in satin finish. Double leaf internal spring fixed at 15° angle for greater comfort. Reinforcing yoke molded into flexible 22" PVC "Y" tubing. Adsoft™ deluxe PVC eartips in small and large size for the ultimate in wearing comfort and acoustic seal. Spare diaphragm included. Scope ID Tag included. Weighs 5.4 oz. Overall length 30.5". Lifetime Warranty.

- 615BD Multifrequency Scope, Burgundy
- 615BK Multifrequency Scope, Black
- 615DG Multifrequency Scope, Dark Green
- 615N Multifrequency Scope, Navy
- 615RB Multifrequency Scope, Royal Blue
- 615ST Multifrequency Scope, Ninja
- 615V Multifrequency Scope, Purple



**WELCH ALLYN ELITE® STETHOSCOPE**

Stainless steel chestpiece and binaurals designed for optimal weight and durability. Dual-channel PVC tubing provides separate sound pathways to each ear. Triple-leaf binaural spring encased in polyurethane permits size adjustment while preventing breakage. Ten year warranty.

- 5079-122 22" Stethoscope, Black
- 5079-125 25" Stethoscope, Black
- 5079-125P 22" Stethoscope, Pediatric, Black
- 5079-270 25" Stethoscope, Burgundy
- 5079-271 25" Stethoscope, Navy
- 5079-284 25" Stethoscope, Forest Green



5079-135



5079-149

**WELCH ALLYN SPECTRUM™ DOUBLE-HEAD STETHOSCOPES**

The full-range Spectrum stethoscope is for professionals who demand the best in auscultation and comfort. Superior acoustics mean accurate auscultation of heart, lung and blood pressure sounds. Binaurals adjust to angle of ear. Stainless steel chestpiece. Portable and durable. Five year warranty.

- 5079-135 28" Adult Black Stethoscope
- 5079-137 28" Adult Navy Stethoscope
- 5079-139 28" Adult Burgundy Stethoscope
- 5079-145 28" Pediatric Black Stethoscope
- 5079-147 28" Pediatric Navy Stethoscope
- 5079-149 28" Pediatric Burgundy Stethoscope
- 5079-285 28" Adult Forest Green Stethoscope
- 5079-287 28" Pediatric Forest Green Stethoscope
- 5079-289 28" Adult Blue Stethoscope
- 5079-291 28" Pediatric Blue Stethoscope



**3M™ LITTMANN® CLASSIC II S.E. STETHOSCOPE**

Featuring a two-sided chestpiece, this classic stethoscope offers high acoustic sensitivity for superior performance. A patented tunable diaphragm on one side augments the traditional bell function of the opposite side. Features patented Littmann Snap Tight soft-sealing eartips plus patient friendly nonchill rim and diaphragm. Three-year warranty plus lifetime service and repair policy.

- 2201 S.E. Stethoscope, 28" Black Tubing
- 2203 S.E. Stethoscope, 28" Gray Tubing
- 2205 S.E. Stethoscope, 28" Navy Blue Tubing
- 2206 S.E. Stethoscope, 28" Caribbean Blue Tubing
- 2208 S.E. Stethoscope, 28" Hunter Green Tubing
- 2209 S.E. Stethoscope, 28" Purple Tubing
- 2210 S.E. Stethoscope, 28" Raspberry
- 2211 S.E. Stethoscope, 28" Burgundy Tubing
- 2812 S.E. Stethoscope, 28", Orange Tubing
- 2814 S.E. Stethoscope, 28", Seafoam Green Tubing
- 2815 S.E. Stethoscope, 28", Pink Tubing

## 3M™ LITTMANN® LIGHTWEIGHT II S.E. STETHOSCOPES

A general-purpose stethoscope that offers a traditional bell and diaphragm with an aluminum chestpiece to provide acoustic sensitivity.

- 2450 Lightweight Stethoscope, 28" Black Tubing
- 2451 Lightweight Stethoscope, 28" Burgundy Tubing
- 2452 Lightweight Stethoscope, 28" Caribbean Blue Tubing
- 2453 Lightweight Stethoscope, 28" Lilac Tubing
- 2454 Lightweight Stethoscope, 28" Ceil Blue Tubing
- 2455 Lightweight Stethoscope, 28" Sea Foam Green Tubing

## 3M™ LITTMANN® VETERINARIAN STETHOSCOPES

- 1392 Vet Stethoscope, Master Classic II, Single-Sided Chestpiece, 32" Navy Tubing



## DRG PURETONE™ MENTOR II TRADITIONAL STETHOSCOPES

- PT14B27 Mentor II Traditional Stethoscope, 27" Black Tube
- PT14N27 Mentor II Traditional Stethoscope, 27" Navy Tube



PT3B27



PT4B27

## DRG PURETONE™ TRADITIONAL STETHOSCOPES

The original stethoscope with the patented ENR™ (External Noise Reduction) tube technology is truly an innovation of design and function. Significantly reduces the amount of external noise while optimizing internal acoustics. 40% less weight as a result of titanium chestpiece construction; also available in stainless steel. Both with conventional diaphragm functionality. GELseal™ eartips and Anti-Slip™ rings are included. Three year warranty.

- PT3B27 Traditional Stethoscope, 27" Black Tube
- PT3N27 Traditional Stethoscope, 27" Navy Tube
- PT4B27 Traditional Ti-Lite Stethoscope, 27" Black Tube
- PT4N27 Traditional Ti-Lite Stethoscope, 27" Navy Tube



## CIDA PEDIATRIC STETHOSCOPES

Pediatric sized, combination diaphragm/bell chestpiece for wide frequency response. 1 1/4" diaphragm/bell equipped with non-chill ring for patient comfort. Lightweight aluminum construction. Adjustable chrome-plated binaurals. Flexible 22" PVC Y tubing. White plastic eartips. Overall length 31 1/2". Latex free system.

- 37-3000 Gray Stethoscope
- 37-3010 Navy Stethoscope
- 37-3014 Light Blue Stethoscope
- 37-3017 Royal Blue Stethoscope
- 37-3020 Gold Stethoscope
- 37-3034 Teal Stethoscope
- 37-3040 Green Stethoscope
- 37-3048 Dark Green Stethoscope
- 37-3050 Red Stethoscope
- 37-3052 Burgundy Stethoscope
- 37-3055 Pink Stethoscope
- 37-3056 Purple Stethoscope
- 37-3057 Neon Pink Stethoscope
- 37-3060 White Stethoscope
- 37-3070 Black Stethoscope
- 37-3095 Neon Orange Stethoscope
- 37-3100 Orchid Stethoscope
- 37-3110 Lavender Stethoscope
- 37-3145 Yellow Stethoscope
- 37-3147 Neon Yellow Stethoscope
- 37-3150 Magenta Stethoscope



## ADC® PROSCOPE™ PEDIATRIC STETHOSCOPES

Pediatric sized, combination diaphragm/bell chestpiece for wide frequency response. 1 1/4" diaphragm 3/4" bell equipped with non-chill ring for patient comfort. Non-chill rim. Lightweight aluminum construction. Adjustable chrome-plated binaurals. Flexible 22" PVC Y tubing. White plastic eartips. Overall length 30". Latex free.

- 675BD Burgundy Stethoscope
- 675BK Black Stethoscope
- 675DG Dark Green Stethoscope
- 675G Gray Stethoscope
- 675GD Gold Stethoscope
- 675GR Green Stethoscope
- 675LB Light Blue Stethoscope
- 675LV Lavender Stethoscope
- 675M Magenta Stethoscope
- 675N Navy Stethoscope
- 675NO Neon Orange Stethoscope
- 675NP Neon Pink Stethoscope
- 675NY Neon Yellow Stethoscope
- 675OD Orchid Stethoscope
- 675P Pink Stethoscope
- 675R Red Stethoscope
- 675RB Royal Blue Stethoscope
- 675TL Teal Stethoscope
- 675V Purple Stethoscope
- 675W White Stethoscope
- 675Y Yellow Stethoscope



## ADC® PROSCOPE™ PEDIATRIC STETHOSCOPES ACCESSORIES

### Accessories

|        |                                  |
|--------|----------------------------------|
| 675-01 | Chestpiece For 675               |
| 675-02 | Diaphragm Retaining Ring For 675 |
| 675-03 | Diaphragm For 675                |
| 675-04 | Non-Chill Ring For 675           |



## 3M™ LITTMANN® CLASSIC II PEDIATRIC & INFANT STETHOSCOPES

Two-sided models with traditional bell/diaphragm chestpieces are designed, sized and acoustically precise for children and infants. Feature patented Littmann Snap Tight soft-sealing eartips plus patient-friendly nonchill rim. Three-year warranty plus a lifetime service and repair policy.

|       |  |
|-------|--|
| 2113  | Pediatric Stethoscope, 28" Black Tubing          |
| 2113R | Pediatric Stethoscope, 28" Red Tubing            |
| 2114  | Infant Stethoscope, 28" Black Tubing             |
| 2114R | Infant Stethoscope, 28" Red Tubing               |
| 2119  | Pediatric Stethoscope, 28" Caribbean Blue Tubing |
| 2122  | Pediatric Stethoscope, 28" Raspberry Tubing      |
| 2124  | Infant Stethoscope, 28" Caribbean Blue Tubing    |
| 2125  | Infant Stethoscope, 28" Raspberry Tubing         |



## ADC® ADSCOPE™ 604 PEDIATRIC STETHOSCOPE

Combination pediatric chestpiece precisely machined to exacting tolerances from surgical stainless steel for outstanding performance and rugged durability. Extra large bell (1 3/8" diameter) for unsurpassed low frequency response. Ultra sensitive diaphragm (1 3/4" diameter) for greater amplification and crisper high frequency transmission. Color coordinated non-chill bell and diaphragm retaining rim for patient comfort. Reinforcing yoke molded into flexible 22" PVC Y tubing. ADSOFT™ deluxe PVC eartips. Spare diaphragm included. Weights 5.5 oz. Overall length 31". Lifetime quality performance guarantee. Latex free.

|       |                                |
|-------|--------------------------------|
| 604BK | Black Pediatric Stethoscope    |
| 604G  | Gray Pediatric Stethoscope     |
| 604LB | Lt. Blue Pediatric Stethoscope |
| 604P  | Pink Pediatric Stethoscope     |

## ADC® ADSCOPE™ 604 PEDIATRIC STETHOSCOPE ACCESSORIES

|          |                           |
|----------|---------------------------|
| 604-01BK | Black Chestpiece          |
| 604-01G  | Gray Chestpiece           |
| 604-01LB | Lt. Blue Chestpiece       |
| 604-01P  | Pink Chestpiece           |
| 604-02BD | Burgundy Diaphragm Rim    |
| 604-02BK | Black Diaphragm Rim       |
| 604-02DG | Dark Green Diaphragm Rim  |
| 604-02G  | Gray Diaphragm Rim        |
| 604-02LB | Lt. Blue Diaphragm Rim    |
| 604-02N  | Navy Diaphragm Rim        |
| 604-02P  | Pink Diaphragm Rim        |
| 604-03   | Diaphragm For All         |
| 604-04BD | Burgundy Non-Chill Ring   |
| 604-04BK | Black Non-Chill Ring      |
| 604-04DG | Dark Green Non-Chill Ring |
| 604-04G  | Gray Non-Chill Ring       |
| 604-04LB | Lt. Blue Non-Chill Ring   |
| 604-04N  | Navy Non-Chill Ring       |
| 604-04P  | Pink Non-Chill Ring       |



PT5B27



PT6B27



PT9B27

## DRG PURETONE™ PEDIATRIC & INFANT STETHOSCOPES

Designed specifically for smaller patients. Features the superior acoustic response of the breakthrough ENR™ tube technology found on the Classic. Patented ENR™ (External Noise Reduction) tube technology significantly reduces the amount of external noise while optimizing internal acoustics. 40% less weight as a result of titanium chestpiece construction; also available in stainless steel. Both with conventional diaphragm functionality sized appropriately for use on smaller patients. GELseal™ eartips and Anti-Slip™ rings add comfort and style. Three year warranty.

|        |   |
|--------|---|
| PT5B27 | Pediatric Stethoscope, 27" Black Tube         |
| PT5N27 | Pediatric Stethoscope, 27" Navy Tube          |
| PT6B27 | Infant Stethoscope, 27" Black Tube            |
| PT6N27 | Infant Stethoscope, 27" Navy Tube             |
| PT7B27 | Pediatric Ti-Lite Stethoscope, 27" Black Tube |
| PT7N27 | Pediatric Ti-Lite Stethoscope, 27" Navy Tube  |
| PT9B27 | Infant Ti-Lite Stethoscope, 27" Black Tube    |
| PT9N27 | Infant Ti-Lite Stethoscope, 27" Navy Tube     |



## CIDA INFANT STETHOSCOPES

Infant sized, combination diaphragm/bell chestpiece for wide frequency response. 25mm diaphragm/bell equipped with non-chill ring for patient comfort. Chrome-plated brass chestpiece. Adjustable chrome-plated binaurals. Flexible 22" PVC Y tubing. White plastic eartips. Overall length 31". Latex free.

- 37-4000 Gray Stethoscope
- 37-4010 Navy Stethoscope
- 37-4014 Light Blue Stethoscope
- 37-4017 Royal Blue Stethoscope
- 37-4020 Gold Stethoscope
- 37-4034 Teal Stethoscope
- 37-4040 Green Stethoscope
- 37-4048 Dark Green Stethoscope
- 37-4050 Red Stethoscope
- 37-4052 Burgundy Stethoscope
- 37-4055 Pink Stethoscope
- 37-4056 Purple Stethoscope
- 37-4057 Neon Pink Stethoscope
- 37-4060 White Stethoscope
- 37-4070 Black Stethoscope
- 37-4095 Neon Orange Stethoscope
- 37-4100 Orchid Stethoscope
- 37-4110 Lavender Stethoscope
- 37-4145 Yellow Stethoscope
- 37-4147 Neon Yellow Stethoscope
- 37-4150 Magenta Stethoscope



## ADC® ADSCOPE™ 605 INFANT STETHOSCOPE

Combination infant chestpiece precisely machined to exacting tolerances from surgical stainless steel for outstanding performance and rugged durability. 7/8" diameter bell provides unsurpassed low frequency response. Ultra sensitive diaphragm (1 1/8" diameter) for greater amplification and crisper high frequency transmission. Non-chill bell and diaphragm retaining rim for patient comfort. Reinforcing yoke molded into flexible 22" PVC Y tubing. ADSOFT™ deluxe PVC eartips. Spare diaphragm and scope ID tag included. Weights 5 oz. Overall length 31". Lifetime quality performance guarantee. Latex free.

- 605BK Black Infant Stethoscope
- 605G Gray Infant Stethoscope
- 605LB Lt. Blue Infant Stethoscope
- 605P Pink Infant Stethoscope

### Accessories

- 605-01BK Black Chestpiece
- 605-01G Gray Chestpiece
- 605-02BK Black Diaphragm Rim
- 605-02G Gray Diaphragm Rim
- 605-03 Diaphragm For All
- 605-04BK Black Non-Chill Ring
- 605-04G Gray Non-Chill Ring



## ADC® PLATINUM EDITION ADSCOPE-LITE™

Unique design that combines the acoustic response of a traditional bell and diaphragm into a convenient one-sided chestpiece. Bell-type low frequency response is achieved by simply holding the chestpiece in place with light pressure. Firm pressure will produce diaphragm-type high frequency performance. A rubber gasket recessed into the chestpiece alters the diaphragm permitting the scope to achieve multi-frequency response. Chestpiece made from lightweight zinc alloy plated with an attractive satin finish. Non-chill diaphragm retaining rim for patient comfort. Large bore chrome plated brass binaurals with matching satin finish. Double leaf internal spring fixed at 15° angle for greater comfort. Reinforcing yoke molded into flexible 22" PVC Y tubing. ADSOFT™ deluxe PVC eartips in small and large size for the ultimate in wearing comfort and acoustic seal. Spare diaphragm included. Scope ID tag included. Weights 6 oz. Overall length 31". Lifetime warranty.

- 614BK Black, Child Adscope-Lite™
- 614G Gray, Child Adscope-Lite™
- 614LB Light Blue, Child Adscope-Lite™
- 614P Pink, Child Adscope-Lite™

## MILTEX PINARD FETAL HEART STETHOSCOPE

Aluminum.

- 30-2200 Fetal Heart Stethoscope

## OMRON FETAL STETHOSCOPE

Designed for superior auscultation of the fetal heart. Stethoscope is equipped with chrome-plated brass binaurals and stem. 22" tubing can be easily shortened to facilitate use of the headrest.

- 13-570 Black Fetal Stethoscope

## OMRON NEONATAL STETHOSCOPE


Chrome-plated brass chestpiece with screw-on retaining ring and slip-on non-chill ring. Quality sound transmission. Chrome-plated binaurals and 22" Y-tubing. Overall length, 30".

- 13-310 Gray Neonatal Stethoscope
- 13-311B Light Blue Neonatal Stethoscope
- 13-315P Pink Neonatal Stethoscope

## MILTEX DE LEE OBSTETRICAL STETHOSCOPES

Chrome. Complete with binaurals.

- 30-2205 Stethoscope w/padded head support, Rigid Metal Arm
- 30-2210 Stethoscope w/padded head support, Flexible Arm




### CIDA SPRAGUE RAPPAPORT STETHOSCOPES

Threaded chestpiece drum of chrome plated zinc permits the use of five interchangeable chestpiece fittings. 22" tubing length; overall length 30". Latex free, double tube configuration for improved acoustic separation. Proprietary valve mechanism eliminates acoustic leakage. Adjustable chrome plated brass binaurals. Complete with an accessory pouch containing the bell chestpieces, chestpiece key, 2 spare diaphragms, one pair of mushroom eartips and one pair of silicon eartips. Weighs 7.5 oz. Inspected and packaged in the U.S.A. Five year warranty.

|         |                         |
|---------|-------------------------|
| 37-0513 | Navy Stethoscope        |
| 37-0514 | Royal Blue Stethoscope  |
| 37-0515 | Black Stethoscope       |
| 37-0517 | Gray Stethoscope        |
| 37-0523 | Green Stethoscope       |
| 37-0525 | Light Blue Stethoscope  |
| 37-0527 | Lavender Stethoscope    |
| 37-0529 | Magenta Stethoscope     |
| 37-0531 | Neon Green Stethoscope  |
| 37-0533 | Neon Orange Stethoscope |
| 37-0534 | Teal Stethoscope        |
| 37-0535 | Neon Pink Stethoscope   |
| 37-0537 | Neon Yellow Stethoscope |
| 37-0538 | Neon Blue Stethoscope   |
| 37-0542 | Dark Green Stethoscope  |
| 37-0543 | Pink Stethoscope        |
| 37-0545 | Red Stethoscope         |
| 37-0547 | White Stethoscope       |
| 37-0552 | Burgundy Stethoscope    |
| 37-0556 | Purple Stethoscope      |
| 37-0558 | Boysenberry Stethoscope |

Everything you  
need for your  
state-of-the-art  
facility.



### ADC® ADSCOPE™ SPRAGUESCOPES

Threaded chestpiece drum of chrome plated zinc permits the use of five interchangeable chestpiece fittings. 22" tubing length; overall length 30" (except #640 has 16" tubing length with overall length 24"). Latex free, double tube configuration for improved acoustic separation. Proprietary valve mechanism eliminates acoustic leakage. Adjustable chrome plated brass binaurals. Complete with an accessory pouch containing the bell chestpieces, chestpiece key, 2 spare diaphragms, one pair of mushroom eartips and one pair of silicon eartips. Weighs 6 oz. Inspected and packaged in the U.S.A. Five year warranty. Latex free.

|        |                              |
|--------|------------------------------|
| 640BK  | 16", Black Spraguescope      |
| 641BB  | Boysenberry Spraguescope     |
| 641BD  | Burgundy Spraguescope        |
| 641BK  | Black Spraguescope           |
| 641CA  | Caribbean Blue Spraguescope  |
| 641FG  | Frosted Glacier Spraguescope |
| 641FL  | Frosted Lilac Spraguescope   |
| 641FS  | Frosted Seafoam Spraguescope |
| 641G   | Gray Spraguescope            |
| 641GR  | Green Spraguescope           |
| 641LB  | Light Blue Spraguescope      |
| 641LV  | Lavender Spraguescope        |
| 641M   | Magenta Spraguescope         |
| 641N   | Navy Spraguescope            |
| 641NB  | Neon Blue Spraguescope       |
| 641NGR | Neon Green Spraguescope      |
| 641NO  | Neon Orange Spraguescope     |
| 641NP  | Neon Pink Spraguescope       |
| 641NY  | Neon Yellow Spraguescope     |
| 641P   | Pink Spraguescope            |
| 641R   | Red Spraguescope             |
| 641RB  | Royal Blue Spraguescope      |
| 641TL  | Teal Spraguescope            |
| 641V   | Purple Spraguescope          |

### ADC® SPRAGUE RAPPAPORT STETHOSCOPE ACCESSORIES

|           |   |
|-----------|---|
| 640-01    | Chestpiece For 640 & 641                    |
| 640-01D   | Chestpiece Drum For 640 & 641               |
| 640-02    | Retaining Ring For 640 & 641, Adult         |
| 640-02S   | Retaining Ring For 640 & 641, Pediatric     |
| 640-03    | Diaphragm For All Spragues, Adult           |
| 640-03S   | Diaphragm For All Spragues, Pediatric       |
| 640-05    | Binaurals For 640 & 641                     |
| 640-05S   | Binaural Spring For 640 & 641               |
| 640-05TBK | Tubing For 640, 16" Black                   |
| 640-06    | Tubing Caps For 640 & 641, Black            |
| 640-09    | Eartips, Black, For 640 & 641               |
| 640-10    | Eartips, Clear                              |
| 640-11    | Bell, 1 1/4", For 640, 641 & 645            |
| 640-12    | Bell, 1", For 640, 641 & 645                |
| 640-13    | Bell, 3/4", For 640, 641 & 645              |
| 640-14    | Aluminum Disc For 640 & 641, Adult          |
| 640-14A   | Diaphragm Assembly For 640 & 641            |
| 640-15    | Aluminum Disc For 640 & 641, Pediatric      |
| 640-15A   | Diaphragm Assembly For 640 & 641, Pediatric |
| 640-16    | Tubing Clip For 640 & 641                   |
| 640-17    | Accessory Kit For 640 & 641                 |
| 640-18    | Accessory Pouch For 640 & 641               |

## ADC® SPRAGUE RAPPAPORT STETHOSCOPE ACCESSORIES (Continued)

|            |  |
|------------|--|
| 640-19     | Chestpiece Key For All Spragues                |
| 641-05TBB  | Boysenberry 22" Tubing For 641 Stethoscope     |
| 641-05TBD  | Burgundy 22" Tubing For 641 Stethoscope        |
| 641-05TBK  | Black 22" Tubing For 641 Stethoscope           |
| 641-05TCA  | Caribbean Blue 22" Tubing For 641 Stethoscope  |
| 641-05TDG  | Dark Green 22" Tubing For 641 Stethoscope      |
| 641-05TFG  | Frosted Glacier 22" Tubing For 641 Stethoscope |
| 641-05TFL  | Frosted Lilac 22" Tubing For 641 Stethoscope   |
| 641-05TFS  | Frosted Seafoam 22" Tubing For 641 Stethoscope |
| 641-05TG   | Gray 22" Tubing For 641 Stethoscope            |
| 641-05TGR  | Green 22" Tubing For 641 Stethoscope           |
| 641-05TLB  | Light Blue 22" Tubing For 641 Stethoscope      |
| 641-05TLV  | Lavender 22" Tubing For 641 Stethoscope        |
| 641-05TM   | Magenta 22" Tubing For 641 Stethoscope         |
| 641-05TN   | Navy 22" Tubing For 641 Stethoscope            |
| 641-05TNB  | Neon Blue 22" Tubing For 641 Stethoscope       |
| 641-05TNGR | Neon Green 22" Tubing For 641 Stethoscope      |
| 641-05TNO  | Neon Orange 22" Tubing For 641 Stethoscope     |
| 641-05TNP  | Neon Pink 22" Tubing For 641 Stethoscope       |
| 641-05TNY  | Neon Yellow 22" Tubing For 641 Stethoscope     |
| 641-05TP   | Pink 22" Tubing For 641 Stethoscope            |
| 641-05TR   | Red 22" Tubing For 641 Stethoscope             |
| 641-05TRB  | Royal Blue 22" Tubing For 641 Stethoscope      |
| 641-05TTL  | Teal 22" Tubing For 641 Stethoscope            |
| 641-05TV   | Purple 22" Tubing For 641 Stethoscope          |

## ADC® NURSE COMBO

Combines the Pocket Pal II™ Kit complete with three instruments and choice of 8 Velcro® tag colors with an ADSCOPE™ 641 Sprague Rappaport scope with 22" tubing color coordinated to the Velcro® tags and a nine piece accessory kit.

|            |                          |
|------------|--------------------------|
| 116-641BDQ | Nurse Combo, Burgundy    |
| 116-641BKQ | Nurse Combo, Black       |
| 116-641DGQ | Nurse Combo, Dark Green  |
| 116-641NQ  | Nurse Combo, Navy        |
| 116-641NYQ | Nurse Combo, Neon Yellow |
| 116-641RBQ | Nurse Combo, Royal Blue  |
| 116-641TLQ | Nurse Combo, Teal        |
| 116-641VQ  | Nurse Combo, Purple      |



## OMRON SPRAGUE RAPPAPORT-TYPE STETHOSCOPES

5-in-1 multipurpose Quintscope™ has approximately 22" long tubing (except #416) and silver chestpiece. Accessories include three different size open bells, two different size replacement diaphragms, and two pairs of eartips in a vinyl storage case.

|             |                        |
|-------------|------------------------|
| 416-22-BLK  | Black Stethoscope      |
| 416-22-DB   | Dark Blue Stethoscope  |
| 416-22-DKGN | Dark Green Stethoscope |



## 3M™ LITTMANN® CARDIOLOGY III STETHOSCOPE

Outstanding acoustic performance and exceptional versatility characterize the Littmann Cardiology III stethoscope. The innovative design offers a patented tunable diaphragm on each side of the chestpiece. The large side can be used for adult patients, while the small side is especially useful for pediatric or thin patients. Features patented Littmann Snap Tight soft-sealing eartips, patient-friendly, nonchill rim and diaphragms. Five-year warranty on chestpiece and binaural, lifetime service and repair policy.

|      |   |
|------|---|
| 3127 | Cardiology Stethoscope, 22" Black Tubing        |
| 3128 | Cardiology Stethoscope, 27" Black Tubing        |
| 3129 | Cardiology Stethoscope, 27" Burgundy Tubing     |
| 3130 | Cardiology Stethoscope, 27" Navy Blue Tubing    |
| 3134 | Cardiology Stethoscope, 27" Hunter Green Tubing |
| 3135 | Cardiology Stethoscope, 27", Plum Tubing        |

## 3M™ LITTMANN® CARDIOLOGY S.T.C. STETHOSCOPES

Outstanding acoustics are combined with a unique and easy to grasp chestpiece that allows easy movement between auscultation sites. Single-sided chestpiece provides both bell and diaphragm frequency responses. Dual lumen design helps to eliminate noise artifact. Comfortably angled, anatomically correct headset. Patient-friendly, nonchill rim. Six-year warranty.

|      |  |
|------|--|
| 4470 | 22" Black Cardiology S.T.C. Stethoscope        |
| 4471 | 27" Black Cardiology S.T.C. Stethoscope        |
| 4472 | 27" Burgundy Cardiology S.T.C. Stethoscope     |
| 4473 | 27" Navy Blue Cardiology S.T.C. Stethoscope    |
| 4474 | 27" Hunter Green Cardiology S.T.C. Stethoscope |
| 4475 | 27" Plum Cardiology S.T.C. Stethoscope         |



## 3M™ LITTMANN® MASTER CARDIOLOGY STETHOSCOPE

A patented tunable diaphragm and superb craftsmanship offer the ultimate performance in a superior-single sided stethoscope meets the exacting standards of health care professionals who demand the highest acoustic response. Features patented Littmann Snap Tight soft-sealing eartips, patient-friendly, nonchill rim and diaphragm and special procedures adapter. Seven-year warranty, plus lifetime service and repair policy.

|      |  |
|------|--|
| 2159 | Cardiology Stethoscope, 22" Black Tubing                       |
| 2160 | Cardiology Stethoscope, 27" Black Tubing                       |
| 2161 | Cardiology Stethoscope, 27" Black Plated Chestpiece & Eartubes |
| 2163 | Cardiology Stethoscope, 27" Burgundy Tubing                    |
| 2164 | Cardiology Stethoscope, 27" Navy Blue Tubing                   |
| 2165 | Cardiology Stethoscope, 27" Hunter Green Tubing                |
| 2167 | Cardiology Stethoscope, 27" Plum Tubing                        |





### ADC® ADSCOPE™ 600 CARDIOLOGY STETHOSCOPE

Unique design that combines the acoustic response of a traditional bell and diaphragm into a convenient one-sided chestpiece. Bell-type low frequency response is achieved by simply holding the chestpiece in place with light pressure. Firm pressure will produce diaphragm-type frequency performance. A rubber gasket recessed into the chestpiece alters the diaphragm permitting the scope to achieve multi-frequency response. Custom molded extra thick bilumen 19" PVC tubing stays soft and flexible. Extra thick tubing walls minimize extraneous noise. Includes one pair of plastic eartips and two pairs of ADSOFT™ PVC eartips (small and large). Overall length 28". Weight 8.5 oz. Scope ID tag included. Lifetime quality performance guarantee. Latex free.

|       |   |
|-------|---|
| 600BD | ADSCOPE™ 600 Cardiology Stethoscope, Burgundy           |
| 600BK | ADSCOPE™ 600 Cardiology Stethoscope, Black              |
| 600DG | ADSCOPE™ 600 Cardiology Stethoscope, Dark Green         |
| 600GP | ADSCOPE™ 600 Cardiology Stethoscope, Gold Plated, Black |
| 600N  | ADSCOPE™ 600 Cardiology Stethoscope, Navy               |

### ADC® ADSCOPE™ 600 CARDIOLOGY STETHOSCOPE ACCESSORIES

|          |                                |
|----------|--------------------------------|
| 600-01BD | Chestpiece For 600, Burgundy   |
| 600-01BK | Chestpiece For 600, Black      |
| 600-01DG | Chestpiece For 600, Dark Green |
| 600-01GP | Chestpiece For 600GP, Black    |
| 600-01N  | Chestpiece For 600, Navy       |
| 600-03   | Diaphragm For 600              |
| 600-05GP | Binaural/Tubing For 600GP      |
| 600-06BK | Eartips, Black                 |

### ADC® ADSCOPE™ 601 CARDIOLOGY STETHOSCOPE

Threaded chestpiece drum allows rapid conversion from a pediatric diaphragm to an adult bell. Chestpiece precisely machined to exacting tolerance from surgical stainless steel for outstanding performance and rugged durability. Extra deep bell for unsurpassed low frequency response. Ultra-sensitive diaphragms for greater amplification and crisper high frequency transmission. Non-chill bell and diaphragm retaining rims for patient comfort. Large internal bore stainless steel binaurals are fixed at a 15° angle for greater comfort. Custom molded bilumen 19" PVC tubing stays soft and flexible. Latex free. Extra thick tubing walls minimize extraneous noise. Double leaf internal spring construction reduces bulk and prevents kinking when folded. Complete with one pair of PVC eartips and ADSOFT™ deluxe PVC eartips in small and large sizes for the ultimate in wearing comfort and acoustic seal. Complete with an accessory case containing the adult bell, 2 spare diaphragms, and two sets of eartips. Scope ID tag included. Weighs 9 oz (adult bell) and 10 oz (pediatric diaphragm). Overall length 28". Lifetime warranty.

|       |                        |
|-------|------------------------|
| 601BD | Burgundy Stethoscope   |
| 601BK | Black Stethoscope      |
| 601DG | Dark Green Stethoscope |
| 601N  | Navy Stethoscope       |

### ADC® ADSCOPE™ 601 CARDIOLOGY STETHOSCOPE ACCESSORIES

|          |   |
|----------|---|
| 601-01BD | Chestpiece For 601, Burgundy              |
| 601-01BK | Chestpiece For 601, Black                 |
| 601-01DG | Chestpiece For 601, Dark Green            |
| 601-01N  | Chestpiece For 601, Navy                  |
| 601-05BD | Binaural/Tubing For 600 & 601, Burgundy   |
| 601-05BK | Binaural/Tubing For 600 & 601, Black      |
| 601-05DG | Binaural/Tubing For 600 & 601, Dark Green |
| 601-05N  | Binaural/Tubing For 600 & 601, Navy       |
| 601-11BD | Bell For 601, Burgundy                    |
| 601-11BK | Bell For 601, Black                       |
| 601-11DG | Bell For 601, Dark Green                  |
| 601-11N  | Bell For 601, Navy                        |
| 601-15BD | Diaphragm Assembly For 601, Burgundy      |
| 601-15BK | Diaphragm Assembly For 601, Black         |
| 601-15DG | Diaphragm Assembly For 601, Dark Green    |
| 601-15N  | Diaphragm Assembly For 601, Navy          |
| 601-18   | Accessory Case For 600 & 601              |



### ADC® ADSCOPE™ 602 CARDIOLOGY STETHOSCOPE

Stainless steel construction for rugged durability. Combination chestpiece precisely machined to exacting tolerances for outstanding acoustic performance. Extra deep bell for unsurpassed low frequency response. Non-chill bell and diaphragm ring for patient comfort. Custom molded 16" bilumen PVC tubing improves acoustic separation. Internal spring construction reduces bulk. Binaurals fixed at 15° angle for greater comfort. Includes one pair of plastic eartips and two pair of PVC eartips with threaded metal inserts for the ultimate in wearing comfort and acoustic seal. Weighs 9 oz. Overall length 25". Latex free. Lifetime warranty.

|       |                        |
|-------|------------------------|
| 602BD | Burgundy Stethoscope   |
| 602BK | Black Stethoscope      |
| 602DG | Dark Green Stethoscope |

## ADC® ADSCOPE™ 602 CARDIOLOGY STETHOSCOPE ACCESSORIES

- 600-17 Accessory Case For 602, 603, 604, 605, 612, 655
- 602-01BD Chestpiece For 602, Burgundy
- 602-01BK Chestpiece For 602, Black
- 602-01DG Chestpiece For 602, Dark Green
- 602-02BD Burgundy Diaphragm Rim For 600, 601, 602, 603, 609, 612
- 602-02BK Black Diaphragm Rim For 600, 601, 602, 603, 609, 612
- 602-02DG Dark Green Diaphragm Rim For 600, 601, 602, 603, 609, 612
- 602-02FG Frosted Glacier Diaphragm Rim For 600, 601, 602, 603, 609, 612
- 602-02FL Frosted Lilac Diaphragm Rim For 600, 601, 602, 603, 609, 612
- 602-02FS Frosted Seafoam Diaphragm Rim For 600, 601, 602, 603, 609, 612
- 602-02G Gray Diaphragm Rim For 602, 603, 609, 612
- 602-02LV Lavender Diaphragm Rim For 603, 609, 612
- 602-02N Navy Diaphragm Rim For 600, 601, 602, 603, 609, 612
- 602-02R Red Diaphragm Rim For 603, 609, 612
- 602-02RB Royal Blue Diaphragm Rim For 602, 603, 609, 612
- 602-02TL Teal Diaphragm Rim For 603, 609, 612
- 602-03 Diaphragm For 600GP, 601, 602, 603, 612
- 602-04BD Burgundy Non-Chill Rings For 601, 602
- 602-04BK Black Non-Chill Rings For 601, 602
- 602-04DG Dark Green Non-Chill Rings For 601, 602
- 602-04N Navy Non-Chill Rings For 601, 602
- 602-05BD Binaural/Tubing For 602, Burgundy
- 602-05BK Binaural/Tubing For 602, Black
- 602-05DG Binaural/Tubing For 602, Dark Green
- 602-07BK ADSOFT Eartips, Black, Large
- 602-07G ADSOFT Eartips, Gray, Large
- 602-08BK ADSOFT Eartips, Black, Small
- 602-08G ADSOFT Eartips, Gray, Small
- 602-09 ADSOFT Eartips, White
- 602-02V Purple Diaphragm Rim For 600, 601, 602, 603, 609, 612



PT2B27



PT8B27

## DRG PURETONE™ CARDIOLOGY STETHOSCOPES

Provides the very best acoustics available through a dual lumen tube and a double-sided chestpiece. Patented dual lumen ENR™ (External Noise Reduction) tube technology significantly reduces the amount of external noise while optimizing internal acoustics. 40% less weight as a result of titanium chestpiece construction; also available in stainless steel. Both with conventional adult- and pediatric-sized diaphragm functionality and the added convenience of a traditional bell. GELseal™ eartips and Anti-Slip™ rings add to the comfort and style. Five year warranty.

- PT2B27 Cardiology Stethoscope, 27" Black Tube
- PT2N27 Cardiology Stethoscope, 27" Navy Tube
- PT8B27 Cardiology Ti-Lite Stethoscope, 27" Black Tube
- PT8N27 Cardiology Ti-Lite Stethoscope, 27" Navy Tube



## DRG PURETONE™ MENTOR III CARDIOLOGY STETHOSCOPES

- PT11B27 Mentor III Cardiology Stethoscope, 27" Black Tube
- PT11N27 Mentor III Cardiology Stethoscope, 27" Navy Tube



## DRG PURETONE™ STAFF CARDIOLOGY STETHOSCOPES

- PT10B27 Staff Cardiology Stethoscope, 27" Black Tube
- PT10N27 Staff Cardiology Stethoscope, 27" Navy Tube
- PT12B27 Staff Cardiology Ti-Lite Stethoscope, 27" Black Tube
- PT12N27 Staff Cardiology Ti-Lite Stethoscope, 27" Navy Tube



5079-321



5079-325

## WELCH ALLYN HARVEY DELUXE DOUBLE AND TRIPLE HEAD STETHOSCOPES

The Harvey™ DLX Stethoscope is our most prestigious and durable model with dual-lumen PVC tubing, brass chestpieces and binaurals, and triple-leaf binaural springs. Optimized for listening to higher frequency and critically important sounds like murmurs, clicks, and ejection sounds. Provides clear transmission of the first and second heart sounds. Ideally sized bell bridges intercostal space between the ribs. Tightly constrained diaphragm provides superior response at large frequencies. Additional corrugated diaphragm for detecting mid-range and low frequencies.

- 5079-321 Harvey Deluxe Triple Head Stethoscope, Black
- 5079-322 Harvey Deluxe Triple Head Stethoscope, Burgundy
- 5079-323 Harvey Deluxe Triple Head Stethoscope, Navy
- 5079-324 Harvey Deluxe Triple Head Stethoscope, Forest Green
- 5079-325 Harvey Deluxe Double Head Stethoscope, Black
- 5079-326 Harvey Deluxe Double Head Stethoscope, Burgundy
- 5079-327 Harvey Deluxe Double Head Stethoscope, Navy
- 5079-328 Harvey Deluxe Double Head Stethoscope, Forest Green



## WELCH ALLYN ORIGINAL HARVEY DOUBLE & TRIPLE HEAD STETHOSCOPE ACCESSORIES

- 5079-05 Corrugated Diaphragm Only
- 5079-09 Flat Pediatric Chestpiece
- 5079-102 Tubing, 2 Tubed Joined For Use w/two-piece spindle
- 5079-103 Name Clip For Use w/5079-102 & 105
- 5079-104 Flat Diaphragm Only
- 5079-105 Tubing, 2 Tubes Joined For Use w/single-piece spindle
- 5079-11 Pediatric Diaphragm Only
- 5079-20 Accessory Kit, Includes 1 Pediatric Chestpiece, 1 Flat Diaphragm, 2 Corrugated Diaphragms, 1 Pair Soft Rubber Eartips, 1 Pair Small Delrin Eartips
- 5079-315 Accessory Kit
- 5079-318 Bell Chestpiece Assembly
- 5079-32 Eartips, Rubber
- 5079-33 Eartips, Small, White, Pair
- 5079-60 Binaurals
- 5079-61 Corrugated Adult Chestpiece
- 5079-62 Flat Adult Chestpiece
- 5079-63 Tubing, 2 Separate Tubes For Use w/two-piece spindle
- 5079-65 Bell Chest Piece



5079-27



5079-28

### WELCH ALLYN ORIGINAL HARVEY DOUBLE & TRIPLE HEAD STETHOSCOPES

Original Harvey Stethoscopes come with solid brass chestpieces and dual-lumen latex tubes. Features three easily rotated chestpieces (bell, corrugated, diaphragm, and flat diaphragm). Clearly detects the full range of important sounds from low-frequency gallops to the faint high-frequency blow of aortic insufficiency. Sound detection capability of the three chestpieces allows for the most complete auscultatory examination.

- 5079-27 Triple Head Stethoscope
- 5079-28 Double Head Stethoscope



### ADC® ADSCOPE™ 655 AMPLIFYING STETHOSCOPE

Diaphragm only type chestpiece fabricated from lightweight acrylic provides up to twice the amplification of conventional acoustic scopes. No electronics to malfunction, no batteries to wear out. Stainless steel binural with double leaf internal spring fixed at 15° angle for greater comfort. Reinforcing yoke molded into flexible 22" PVC Y tubing. ADSOFT™ deluxe PVC eartips in small and large size. Scope ID tag included. Weighs 5 oz. Overall length 31". Lifetime quality performance guarantee. Latex free.

- 655BD Burgundy Stethoscope
- 655BK Black Stethoscope
- 655DG Dark Green Stethoscope
- 655N Navy Stethoscope
- 655R Red Stethoscope
- 655RB Royal Blue Stethoscope

#### Accessories

- 655-01 Chestpiece For 655
- 655-02 Retaining Ring For 655
- 655-03K Diaphragm Kit For 655



### ADC® ADSCOPE™ ADVANCED ELECTRONIC SCOPE

Combines traditional scope appearance with state-of-the-art electronics for the ultimate in acoustic performance and ease of use. Amplifies up to 18x over traditional acoustic scopes. 8 level volume control. Three frequency modes. Bell (20-230Hz) for cardiac sounds. Diaphragm (100-800Hz) for lung sounds. Extended Range (20-800Hz). Electronic transmission of sounds to separate speakers near earpiece for improved performance. LCD with iconic display of heart and lung for more intuitive use, and volume level and battery strength. Impact muting switch protects clinician's ears from high intensity impact noises. Battery warning signals when replacement needed. Binaurals fixed at 15° angle for snug aural fit. Specially constructed Adsoft™ PVC eartips for the ultimate in wearing comfort and acoustic seal. Lightweight and compact - weighs under 7.5 oz., measures just 36". Auto shut off extends battery life. Latex Free. Uses 2 "AA" batteries. 1 year warranty.

- 657BK Advanced Electronic Scope, Black



### WELCH ALLYN MEDITRON™ ELECTRONIC STETHOSCOPE SYSTEM

Welch Allyn has teamed up with Meditron AS of Norway to provide the first computer power and Auscultation combination. A unique, patented listening system that dramatically enhances sound clarity. Ideal educational tool. Records sounds and images and speeds consultant review and future analysis (based on Microsoft™ wav files). Should be used by any physical examination where additional amplification is required, including doctor's offices, EMS sites, veterinarian's office, tele-medicine site, hospitals, and medical schools. Reimbursement is the same as regular physical examination (multiple codes). 200 hours or four years of use from the battery. Excellent for use with difficult patients (weak sounds, obese patients, children). May be cleaned with a lightly moistened cotton wipe using chlorhexidine or alcohol. Does not contain any latex.

- 5079-400 Electronic Stethoscope

#### Accessories

- 5079-401 Accessory Package
- 5079-414 Headset For Additional Users
- 5079-430 Analyzer Software Accessory



### DRG ECHO™ AMPLIFIED STETHOSCOPES

Highlighted by the benefits of both a conventional and digital stethoscope in one lightweight, portable design. Digital amplification assists in the identification of even the most subtle cardiopulmonary events, while its conventional bell and diaphragm modes provide clinicians with the security against unexpected battery failure during critical clinical interactions. 400% increase in amplification when in digital mode over typical stethoscopes, a heightened frequency response range, and batteries warranted to last up to 100 hours of use. Switches easily between three modes of operation: digital, diaphragm, and bell. Conventional acoustics accompany the bell and diaphragm modes. Patented dual lumen ENR™ (External Noise Reduction) tube technology significantly reduces the amount of external noise while optimizing internal acoustics. GELseal™ eartips and Anti-Slip™ rings are included. One year warranty.

- EC3B27 Amplified Stethoscope, 27" Black Tube




### 3M™ LITTMANN® TEACHING STETHOSCOPES

A stethoscope with two headsets which allow instructor and student to perform simultaneous auscultation. Patented tunable diaphragms.

- 2138 Classic II S.E. Teaching Stethoscope, Bell & Diaphragm w/2 headsets, 40" 1 Tube Black, 1 Tube Burgundy
- 2139 Master Classic II Teaching Stethoscope, Single-Sided Chestpiece w/2 headsets, 40" 1 Tube Black, 1 Tube Burgundy





## ADC® PROSCOPE™ TEACHING STETHOSCOPE

Choice of diaphragm only (660 Series) or combination diaphragm/bell (670 Series) PROSCOPE™ chestpiece. Two sets of adjustable chrome plated brass binaurals. Two sets of flexible 22" and one straight 10" length of black PVC tubing. Two spare pair of black mushroom eartips included. Overall length 39". Inspected and packaged in the U.S.A. One year quality performance guarantee. Latex free.

- 661 PROSCOPE™ 661 Teaching Scope, Nurse Style
- 671 PROSCOPE™ 671 Teaching Scope, Dualhead

### Accessories

- 661-02 Diaphragm Retaining Ring For 660, 661, 670, 671
- 661-03 Diaphragm For 660, 661, 670, 671



## OMRON DUAL-HEAD TRAINING SCOPE

Two-channel stethoscope provides accurate sound transmission to both student and instructor. Dual-head chest piece converts from diaphragm to an open bell head with a single turn of the stem.

- 0405 Dual-Head Training Scope

## TRIMLINE TEACHING STETHOSCOPE

- SK20358 Teaching Stethoscope w/Two Sets of Binaurals



## DRG ECHO+PLUS™ DIGITAL RECORDING STETHOSCOPES

Versatility and digitally-enhanced amplification plus the diagnostic and teaching options of digital record and replay. Allows clinicians to instantly capture critical, unexpected heart sounds for future analysis or instant consultation with colleagues. Assists instructors in teaching students the fundamentals of auscultation. 400% increase in amplification when in digital mode over typical stethoscopes, a heightened frequency response range, and batteries warranted to last up to 100 hours of use. Switches easily between three modes of operation: digital, diaphragm, and bell. Allows fifteen seconds of recording time and two replay speeds. Conventional acoustics accompany the bell and diaphragm modes. Patented dual lumen ENR™ (External Noise Reduction) tube technology significantly reduces the amount of external noise while optimizing internal acoustics. GELseal™ eartips and Anti-Slip™ rings are included. One year warranty.

- EC2B27 Digital Recording Stethoscope, 27" Black Tube

## DRG ECHO+PLUS™ TEACHING MODULE

Effective, yet affordable instrumentation to assist in teaching the fundamentals of auscultation as well as the identification of more complicated cardiac pathologies. The instructor can record patient sounds directly at the bedside using the DataLink and replay them through the teaching module to students who are using their own stethoscope in a simulated setting. Ten individual pushbutton sound storage locations. Nine separate listening stations and universal volume control. Self-contained design is effective, portable, and a valuable asset. Universal 110V power requirements. One year warranty.

- TM100 Teaching Module



## DRG ECHO+PLUS™ WITH DATALINK™ STETHOSCOPES

The push-button record, playback, and digital amplification features allows the ability to quickly download recorded patient data. By plugging in directly to a standard PC, the DataLink feature enable .wav file outputs for patient tracking and archiving, case study slide presentations and unlimited educational applications. A simple device to more effectively communicate with your students and remain current on your patient's progress. 400% increase in amplification when in digital mode over typical stethoscopes, a heightened frequency response range, and batteries warranted to last up to 100 hours of use. Switches easily between three modes of operation: digital, diaphragm, and bell. Allows fifteen seconds of recording time and two replay speeds. Conventional acoustics accompany the bell and diaphragm modes. Patented dual lumen ENR™ (External Noise Reduction) tube technology significantly reduces the amount of external noise while optimizing internal acoustics. GELseal™ eartips and Anti-Slip™ rings are included. One year warranty.

- EC1B27 DataLink Recording Stethoscope, 27" Black Tube



## 3M™ DOUBLE STICK DISCS

Pre-cut, hypoallergenic discs that are ideal for securing 3M™ Littmann® stethoscope chestpieces. The discs have an 11/32" aperture and a 1 1/4" diameter.

- 2181 Double Stick Discs w/11/32" aperture, 102/pk, 6 pk/cs



Easy to apply. Easy to live with.  
3M™ Scotchcast™ Wet or Dry Padding.



3M Medical

*from the 3M Health Care family*

© 2005 3M All Rights Reserved



From now on, you only need one cast padding for either wet or dry use.

3M™ Scotchcast™ Wet or Dry Padding is a hypoallergenic, synthetic cast padding that allows patients to get their cast wet.

Plus, the padding can be applied using the same techniques used for synthetic or cotton cast padding. And, there are no special techniques required for removal of the cast.

For more information, visit us at [www.3m.com/healthcare](http://www.3m.com/healthcare).

**3M**

## Could you be missing more than you think?

Why did cardiologists miss hearing a Grade 2 aortic regurgitation murmur 5 times more often using a conventional, cardiology-type stethoscope than they did with the new 3M™ Littmann® Electronic Stethoscope Model 3000?\*

Because the Littmann Model 3000 with proprietary Ambient Noise Reduction technology reduces distracting noise by an average of 75% (-12dB) without eliminating critical body sounds. So difficult body sounds may become easier to hear. Get the complete story at [www.littmann.com](http://www.littmann.com).

**Don't miss the sounds you need to hear.**

\*100 cardiologists listening to a pre-recorded heart sound with and without a Grade I or Grade II aortic regurgitation murmur, presented in random order in the presence of 70-75 dB ambient noise, using a 3M Littmann Electronic Stethoscope Model 3000 and an acoustic (non-electronic) cardiology-type stethoscope.

3M Medical

*from the 3M Health Care family*

© 2005 3M All Rights Reserved



**3M**



## Your one-stop shop for office setups

Look to Welch Allyn for all of these high quality products:

**AEDs**  
**Defibrillators**  
**Otoscopes**  
**Thermometers**

**Audiometers**  
**ECGs**  
**Lighting Products**  
**Vag Specs**

**Blood Pressure**  
**Laryngoscopes**  
**Stethoscopes**  
**Vital Signs Monitors**

**Colposcopes**  
**Ophthalmoscopes**  
**Tympanometric Devices**  
**Headlights**

---

Welch Allyn's easy-to-use, intuitive products were developed with busy clinicians in mind and are designed to work together. For example, our Integrated Diagnostic System combines all the standard physical exam products—ophthalmoscope, otoscope, blood pressure gauge, thermometer and ear tips—in one easy-to-order, easy-to-install package. This flexible product can be configured to fit any doctor's need, and is a great way to start every office setup.

To learn more about Welch Allyn products, visit our website at [www.welchallyn.com](http://www.welchallyn.com)

© 2005 Welch Allyn MC 3197

**WelchAllyn®**

Advancing Frontline Care™

Alko Enterprises




**3M™ LITTMANN® STETHOSCOPE BINAURALS**

Includes eartips, eartubes and tubing.

|             |  |
|-------------|--|
| 12237612523 | Binaurals For Select & Lightweight, 28" Light Blue   |
| 1223768527  | Binaurals For Classic II S.E., Classic II & Classic II Pediatric & Infant, 28" Black                           |
| 12237685289 | Binaurals For Classic II S.E., Classic II & Classic II Pediatric & Infant, 28" Gray                            |
| 12237685305 | Binaural Assembly for Littmann® Classic II Pediatric, Navy Blue Tubing, 28"                                    |
| 12237685321 | Binaurals For Classic II S.E., Classic II & Classic II Pediatric & Infant, 28" Red                             |
| 12237685339 | Binaurals For Classic II S.E., Classic II & Classic II Pediatric & Infant, 28" Burgundy                        |
| 7880555776  | Binaurals For Classic II S.E., Classic II & Classic II Pediatric & Infant, 28" Pink                            |
| 78806297442 | Binaurals For Master Classic II, 27" Navy  |
| 78806297459 | Binaurals For Master Classic II, 27" Burgundy  |
| 78806299273 | Binaurals For Master Classic II, 27" Black   |
| 78806311375 | Binaurals For Master Classic II, 27" Black, Black Plated Eartubes  |
| 78806311425 | Binaurals For Master Classic II, 27" Black, 23 Karat Gold-Plated Eartubes                                      |
| 78807802844 | Binaurals For Teaching Stethoscope, 40" Black w/2 headsets For Classic II S.E.                                 |
| 78807802851 | Binaurals For Teaching Stethoscope, 40" Black w/2headsets For Master Classic II                                |
| 78807812066 | Binaurals For Select & Lightweight, 28" Black  |
| 78807812074 | Binaurals For Select & Lightweight, 28" Pink   |
| 78807833831 | Binaural Assembly For Littmann Master Classic II Veterinarian Stethoscope, 32" Navy Blue                       |
| 78807845165 | Binaurals For Classic II S.E., Classic II & Classic II Pediatric & Infant, 28" Royal Blue                      |
| 78807845181 | Binaurals For Select & Lightweight, 28" Burgundy   |
| 78807845207 | Binaurals For Select & Lightweight, 28" Royal Blue   |
| 78807858390 | Binaurals For Select & Lightweight, 28" Purple   |
| 78807858465 | Binaurals For Master Classic II, 27" Purple  |
| 78807865023 | Binaurals For Classic II S.E., Classic II & Classic II Pediatric & Infant, 28" Hunter Green                    |
| 78807865031 | Binaurals For Classic II S.E., Classic II & Classic II Pediatric & Infant, 28" Purple                          |
| 78807865825 | Binaurals For Master Cardiology, Cardiology III & Cardiology II S.E., 22" Black, Short                         |
| 78807865833 | Binaurals For Master Cardiology, Cardiology III & Cardiology II S.E., 27" Black                                |
| 78807865858 | Binaurals For Master Cardiology, Cardiology III & Cardiology II S.E., 27" Navy                                 |
| 78807865866 | Binaurals For Master Cardiology, Cardiology III & Cardiology II S.E., 27" Hunter Green                         |
| 78807865874 | Binaurals For Master Cardiology, Cardiology III & Cardiology II S.E., 27" Hunter Green                         |
| 78807866211 | Binaurals For Master Cardiology, Cardiology III & Cardiology II S.E., 27" Black, Black Plated Eartubes         |
| 78807866229 | Binaurals For Master Cardiology, Cardiology III & Cardiology II S.E., 27" Black, 23 Karat Gold-Plated Eartubes |
| 78807866359 | Binaurals For Classic II S.E., Classic II & Classic II Pediatric & Infant, 28" Raspberry                       |
| 78807866367 | Binaurals For Classic II S.E., Classic II & Classic II Pediatric & Infant, 28" Caribbean Blue                  |
| 78807866508 | Binaurals For Master Classic II, 27" Caribbean Blue  |
| 78807866516 | Binaurals For Master Classic II, 27" Hunter Green  |
| 78807866581 | Binaural Assembly For Littmann Select & Lightweight Stethoscope, 28" Raspberry                                 |
| 7880786659  | Binaurals For Select & Lightweight, 28" Caribbean Blue   |
| 7880789464  | Binaurals For Cardiology S.T.C., 27" Burgundy  |
| 7880789465  | Binaurals For Cardiology S.T.C., 27" Navy Blue   |
| 7880789466  | Binaurals For Cardiology S.T.C., 27" Black   |
| 7880789467  | Binaurals For Cardiology S.T.C., 27" Hunter Green  |
| 7880789513  | Binaurals For Cardiology S.T.C., 22" Black   |
| 78808303578 | Binaural Assembly for Littmann® Master Cardiology and Cardiology III, Plum Tubing, 27"                         |
| 78808303586 | Binaural Assembly for Littmann® Master Classic II, Ceil Blue Tubing, 27"                                       |
| 78808303594 | Binaural Assembly for Littman® Classic II S.E., Ceil Blue Tubing, 28"  |
| 78808303602 | Binaural Assembly for Littmann® Classic II S.E., Seafoam Green Tubing, 28"                                     |
| 78808303610 | Binaural Assembly for Littmann® Classic II Infant, Lilac Tubing, 28"   |
| 78808303628 | Binaural Assembly for Littmann® Classic II S.E., Orange Tubing, 28"  |
| 78808303636 | Binaural Assembly for Littmann® Select, Ceil Blue Tubing, 28"  |
| 78808303644 | Binaural Assembly for Littmann® Select, Seafoam Green Tubing, 28"  |
| 78808303651 | Binaural Assembly for Littmann® Select, Lilac Tubing, 28"  |
| 78808303909 | Binaural Assembly for Littman® Cardiology S.T.C., Plum Tubing, 27"   |

**3M™ LITTMANN® STETHOSCOPE DIAPHRAGMS**

|       |   |
|-------|---|
| 36559 | Stethoscope Diaphragms For the Cardiology II, Cardiology, Classic II, Classic, Lightweight, Blood Pressure & Differential Stethoscope, Printed w/Littmann Logo, 10/bg |
| 36560 | Stethoscope Diaphragms For Classic II Pediatric Stethoscope, Printed w/Littmann Logo, 10/bg   |
| 36561 | Stethoscope Diaphragms For Classic II Infant Stethoscope, Printed w/Littmann Logo, 10/bg  |


**3M™ LITTMANN® STETHOSCOPE EARTIPS**

|       |   |
|-------|---|
| 37808 | Snap-Tight Stethoscope Eartips, Firm, Gray, Large, 10 pr/bg           |
| 37809 | Snap-Tight Stethoscope Eartips, Soft-Sealing, Black, Small, 10 pr/ctn |
| 37810 | Snap-Tight Stethoscope Eartips, Soft-Sealing, Gray, Small, 10 pr/ctn  |
| 37811 | Snap-Tight Stethoscope Eartips, Soft-Sealing, Black, Large, 10 pr/ctn |
| 37812 | Snap-Tight Stethoscope Eartips, Soft-Sealing, Gray, Large, 10 pr/ctn  |
| 38742 | Snap-Tight Stethoscope Eartips, Firm, Gray, Small, 10 pr/bg           |
| 51161 | Soft-Sealing Eartips, Light Brown, Large, 10 pr/ctn                   |

**3M™ LITTMANN® STETHOSCOPE IDENTIFICATION TAGS**

Fits all sizes.

|      |  |
|------|--|
| 2170 | Stethoscope Identification Tags, Black, 5/cs |
| 2171 | Stethoscope Identification Tags, Gray, 5/cs  |

**3M™ LITTMANN® STETHOSCOPE NONCHILL BELL SLEEVES**

For use with traditional two-sided chestpieces on 3M™ Littmann™ stethoscopes.

|       |  |
|-------|--|
| 36544 | Stethoscope Nonchill Bell Sleeves For Classic II S.E., Classic II & Lightweight, Black, 10/bg          |
| 36545 | Stethoscope Nonchill Bell Sleeves For Classic II S.E., Classic II & Lightweight, Gray, 10/bg           |
| 36546 | Stethoscope Nonchill Bell Sleeves For Classic II Pediatric, Black, 5/bg                                |
| 36547 | Stethoscope Nonchill Bell Sleeves For Classic II Infant, Black, 5/bg                                   |
| 36568 | Stethoscope Nonchill Bell Sleeves For Classic II Pediatric, Gray, 5/bg                                 |
| 36569 | Stethoscope Nonchill Bell Sleeves For Classic II Infant, Gray, 5/bg                                    |
| 36570 | Stethoscope Nonchill Bell Sleeves For Cardiology III, Cardiology II S.E. & Cardiology II, Black, 10/bg |
| 36571 | Stethoscope Nonchill Bell Sleeves For Cardiology III, Cardiology II S.E. & Cardiology II, Gray, 10/bg  |
| 51162 | Stethoscope Nonchill Bell Sleeve For Littmann Lightweight II S.E. Stethoscope, Black, 10/bg            |
| 51163 | Stethoscope Nonchill Bell Sleeve For Littmann Lightweight II S.E. Stethoscope, Light Brown, 10/bg      |

## 3M™ LITTMANN® STETHOSCOPE REPAIR KIT CASE

44525 Stethoscope Repair Kit Case

## 3M™ LITTMANN® STETHOSCOPE SNAP-ON RIMS

- 36548 Stethoscope Snap-On Rims For Cardiology III (large side), Cardiology II S.E., Classic II S.E., Select & Lightweight, Black, 10/bg
- 36549 Stethoscope Snap-On Rims For Cardiology III (large side), Cardiology II S.E., Classic II S.E., Select & Lightweight, Gray, 10/bg
- 36550 Stethoscope Snap-On Rims For Cardiology III (small side) & Classic II Pediatric, Black, 5/bg
- 36551 Stethoscope Snap-On Rims For Cardiology III (small side) & Classic II Pediatric, Gray, 5/bg
- 36552 Stethoscope Snap-On Rims For Classic II Infant, Black, 5/bg
- 36553 Stethoscope Snap-On Rims For Classic II Infant, Gray, 5/bg



## 3M™ LITTMANN® STETHOSCOPE SPECIAL PROCEDURE ADAPTORS

- 36563 Stethoscope Adaptors For Cardiology II S.E. & Cardiology II, Black, 3/bg
- 36564 Stethoscope Adaptors For Cardiology II S.E. & Cardiology II, Gray, 3/bg
- 36565 Stethoscope Adaptors For Master Cardiology, Black, 3/bg
- 36566 Stethoscope Adaptors For Master Cardiology, Gray, 3/bg



## 3M™ LITTMANN® TUNABLE DIAPHRAGM & RIM ASSEMBLY

- 36554 Diaphragm & Rim Assembly For Master Cardiology, Black Rim, 5/bg
- 36555 Diaphragm & Rim Assembly For Master Cardiology, Gray Rim, 5/bg
- 36556 Diaphragm & Rim Assembly For Cardiology III (large side), Cardiology II S.E., Master Classic II, Classic II S.E. & Select, Black Rim, 5/bg
- 36557 Diaphragm & Rim Assembly For Cardiology III (large side), Cardiology II S.E., Master Classic II, Classic II S.E. & Select, Gray Rim, 5/bg
- 36572 Diaphragm & Rim Assembly For Cardiology III (small side), Black Rim, 5/bg
- 36573 Diaphragm & Rim Assembly For Cardiology III (small side), Gray Rim, 5/bg
- 51159 Diaphragm & Rim Assembly For Littmann Lightweight II S.E. Stethoscope, Black Rim, 5/bg
- 51160 Diaphragm & Rim Assembly For Littmann Lightweight II S.E. Stethoscope, Light Brown Rim, 5/bg
- 51185 Diaphragm & Rim Assembly For Littmann Cardiology S.T.C., Black Rim, 5/bg



GET100-1



SS50-4CL-4C



SS50-4PK-4C-1

## DRG STETHOSCOPE DIAPHRAGMS FOR ECHO™, ECHO+PLUS™ & PURETONE™

Durable, chemically inert and acoustically sensitive, and provide a level of protection unmatched in the industry. Add AgION™ Antimicrobial Compound to help prevent the growth of life-threatening microbes. Latex free, C-Flex, medical-grade synthetic elastomer. Chemically inert, non-allergenic, and resistant to ripping or hardening over time. Warm against patient's skin, unlike traditional anti-chill rings. Disposable diaphragms are single use (packaged clean); antimicrobial diaphragms with AgION can be used for up to one month before replacement. Standard diaphragms available in white; AgION and disposable diaphragm colors are: red, green, orange, royal blue and maroon.

- BR100-C-BK Bell Ring, Cardiology, Black
- D100-A Stethoscope Hard Diaphragm, Adult
- D100-I Stethoscope Hard Diaphragm, Infant
- D100-P Stethoscope Hard Diaphragm, Pediatric
- DR100-BK Diaphragm Ring, Black
- GET100-1 GELSeal Eartips, Add Color Code to Part #: BK-Black, BL-Blue, GR-Green, RD-Red, or YL-Yellow
- HSAD150-A Diaphragm Anti-Microbial, Adult, Specify Color: Green, Orange, Blue, Red, 50/bx
- HSAD150-I Diaphragm Anti-Microbial, Infant, Specify Color: Green, Orange, Blue, Red, 50/bx
- HSAD150-P Diaphragm Anti-Microbial, Pediatric, Specify Color: Green, Orange, Blue, Red, 50/bx
- HUM100-S Hummer Relex Hammer & Tuning Fork
- ID100-BK ID Tag
- SS50-4CL-4C SafeSeal, 4 Column 4 Color, Adult, 200/bx
- SS50-4CL-4C-I SafeSeal, 4 Color, Infant, 200/bx
- SS50-4CI-4C-P SafeSeal, 4 Color, Pediatric, 200/bx
- SS50-4PK-4C-I SafeSeal, 4 Color (Blue, Green, Orange, Red), Infant, 200/bx
- SS50-4PK-4C-P SafeSeal, 4 Color (Blue, Green, Orange, Red), Pediatric, 200/bx
- SS50-4PK-PAT SafeSeal, Patriotic (2 Red, 2 Blue), Adult, 200/bx
- BR100-CL-BK Bell Ring, Traditional, Black
- SS50-4PK-4C SafeSeal, 4 Color (Blue, Green, Orange, Red), Adult, 200/bx

## OMRON STETHOSCOPE REPLACEMENT PARTS

- 0418 Binaurals For Sprague Rappaport-Type
- 0426 Chestpiece For Bowles 407
- 0431 Plastic Eartips For Lightweight Binaurals, Black
- 0431CLEAR Soft Clear Eartips For Sprague Rappaport-Types
- 0437 Plastic Eartips For Lightweight Binaurals, White
- 0463 Retaining Rings For 412 Series
- 0485 Chestpiece Sprague Rappaport-Type
- 0487 Mushroom Eartips For Lightweight Binaurals, Black
- 423 Black PVC "Y" Tubing
- 423BLUE Blue PVC "Y" Tubing
- 423GRAY Gray PVC "Y" Tubing
- 423GREN Green PVC "Y" Tubing
- 423RED Red PVC "Y" Tubing





**WELCH ALLYN ELITE® STETHOSCOPE ACCESSORIES**

- 5079-108 Nametag
- 5079-115 Pediatric Diaphragm Assembly, Black
- 5079-116 Adult Diaphragm Assembly, Black
- 5079-117 Pediatric Bell, Non-Chill Rim, Black
- 5079-118 Adult Bell, Non-Chill Rim, Black
- 5079-119 22" Tubing, Black
- 5079-120 25" Tubing, Black
- 5079-121 Binaural/Spring Assembly, Black
- 5079-123 Adult Diaphragm Disc Only
- 5079-124 Pediatric Diaphragm Disc Only
- 5079-126 Adult Diaphragm, Non-Chill Rim, Black
- 5079-127 Pediatric Diaphragm, Non-Chill Rim, Black
- 5079-129 Black Tubing Assembly
- 5079-184 Adult Diaphragm, Non-Chill Rim, Gray
- 5079-185 Pediatric Diaphragm, Non-Chill Rim, Gray
- 5079-232 Accessory Kit, Black, Includes a Pediatric Diaphragm, Bell Ring & 3 Additional Set of Eartips
- 5079-233 Comfort Sealing Eartips, Gray
- 5079-272 25" Tubing, Burgundy
- 5079-273 25" Tubing, Navy
- 5079-274 Binaural/Spring Assembly, Burgundy
- 5079-275 Binaural/Spring Assembly, Navy
- 5079-276 Pediatric Bell, Non-Chill Rim, Gray
- 5079-277 Adult Bell, Non-Chill Rim, Gray
- 5079-278 Pediatric Diaphragm Assembly, Gray
- 5079-279 Adult Diaphragm Assembly, Gray
- 5079-280 Elite Accessory Kit, Gray, Includes Pediatric Diaphragm, Bell Ring & 3 Additional Sets of Eartips
- 5079-305 25" Tubing, Forest Green
- 5079-307 Binaural/Spring Assembly, Forest Green
- 5079-314 Comfort Sealing Eartips, Small
- 5079-66 Hard Eartips, Large

**WELCH ALLYN SPECTRUM™ ACCESSORIES & PARTS**

- 5079-170 Comfort Sealing Eartips, Large
- 5079-172 Hard Eartips, Medium
- 5079-173 Hard Eartips, Large
- 5079-180 Adult Bell, Non-Chill Rim, Black
- 5079-181 Adult Bell, Non-Chill Rim, Gray
- 5079-182 Pediatric Bell, Non-Chill Rim, Black
- 5079-183 Pediatric Bell, Non-Chill Rim, Gray
- 5079-195 Binaural/Spring Assembly & Tubing, 28", Black
- 5079-196 Binaural/Spring Assembly & Tubing, 28", Gray
- 5079-197 Binaural/Spring Assembly & Tubing, 28", Navy
- 5079-199 Binaural/Spring Assembly & Tubing, 28", Burgundy
- 5079-210 Double-Head Chestpiece Assembly, Adult, Gray
- 5079-211 Double-Head Chestpiece Assembly, Adult, Black
- 5079-214 Double-Head Chestpiece Assembly, Pediatric, Black
- 5079-215 Double-Head Chestpiece Assembly, Pediatric, Gray
- 5079-313 Comfort Sealing Eartips, Small



**SUMMIT DOPPLERS**

Interchangeable probes. Superior sound quality. Extended battery life.

- L150ASD2 No Display Doppler w/2MHz OB probe & recorder
- L150ASD3 No Display Doppler w/3MHz OB probe & recorder
- L150ASD4 No Display Doppler w/4MHz vascular probe & recorder
- L150ASD5 No Display Doppler w/5MHz vascular probe & recorder
- L150ASD8 No Display Doppler w/8MHz vascular probe & recorder
- L150ASDW No Display Doppler w/OB waterproof probe & recorder
- L150RSD2 No Display Doppler w/2MHz OB probe & recharger
- L150RSD3 No Display Doppler w/3MHz OB probe & recharger
- L150RSD4 No Display Doppler w/4MHz vascular probe & recharger
- L150RSD5 No Display Doppler w/5MHz vascular probe & recharger
- L150RSD8 No Display Doppler w/8MHz vascular probe & recharger
- L150RSDW No Display Doppler w/OB waterproof probe & recharger
- L150SD2 No Display Doppler w/2MHz OB probe
- L150SD3 No Display Doppler w/3MHz OB probe
- L150SD4 No Display Doppler w/4MHz vascular probe
- L150SD5 No Display Doppler w/5MHz vascular probe
- L150SD8 No Display Doppler w/8MHz vascular probe
- L150SDW No Display Doppler w/OB waterproof probe
- L250ARSD2 Display Doppler w/2MHz OB probe, recharger & recorder
- L250ARSD3 Display Doppler w/3MHz OB probe, recharger & recorder
- L250ARSD4 Display Doppler w/4MHz vascular probe, recharger & recorder
- L250ARSD5 Display Doppler w/5MHz vascular probe, recharger & recorder
- L250ARSD8 Display Doppler w/8MHz vascular probe, recharger & recorder
- L250ARSDW Display Doppler w/OB waterproof probe, recharger & recorder
- L250RSD2 Display Doppler w/2MHz OB probe & recharger
- L250RSD3 Display Doppler w/3MHz OB probe & recharger
- L250RSD4 Display Doppler w/4MHz vascular probe & recharger
- L250RSD5 Display Doppler w/5MHz vascular probe & recharger
- L250RSD8 Display Doppler w/8MHz vascular probe & recharger
- L250RSDW Display Doppler w/OB waterproof probe & recharger
- L250SD2 Display Doppler w/2MHz OB probe
- L250SD3 Display Doppler w/3MHz OB probe
- L250SD4 Display Doppler w/4MHz vascular probe
- L250SD5 Display Doppler w/5MHz vascular probe
- L250SD8 Display Doppler w/8MHz vascular probe
- L250SDW Display Doppler w/OB waterproof probe
- L350RSD2 Table-Top Doppler w/2MHz OB probe & recharger
- L350RSD3 Table-Top Doppler w/3MHz OB probe & recharger
- L350RSD4 Table-Top Doppler w/4MHz Vascular probe & recharger
- L350RSD5 Table-Top Doppler w/5MHz Vascular probe & recharger
- L350RSD8 Table-Top Doppler w/8MHz Vascular probe & recharger
- L350RSDW Table-Top Doppler w/2MHz OB waterproof probe & recharger
- SD2 2MHz OB Probe Only
- SD3 3MHz OB Probe Only
- SD4 4MHz Vascular Probe Only
- SD5 5MHz Vascular Probe Only
- SD8 8MHz Vascular Probe Only
- SDW 2MHz OB Waterproof Probe Only
- L250ABI Reimbursable Vascular Doppler w/printer & 8 MHz Bi-Directional Probe
- Accessories
- K150 ABI Kit Includes Aneroid, 2 Cuffs, 50 Forms, Chart & Video

## SUMMIT DOPPLER ACCESSORIES

|      |   |
|------|---|
| B150 | Rechargeable Battery, 3/pk  |
| B155 | Alkaline Battery, 3/pk  |
| C150 | Transeptic Cleansing Solution   |
| G150 | Aquasonic Clean Ultrasound Gel, 60 gm   |
| G155 | Aquasonic Clean Ultrasound Gel, 250 gm  |
| H150 | Headphones  |
| K160 | ABI 2 Cuff Package w/trigger release aneroid, Includes: 1-10cm and 1-12cm Cuff  |
| K170 | ABI 4 Cuff Package w/trigger release aneroid, Includes: 2-10cm and 2-12cm Cuffs |
| K180 | Printer Label Paper, 5rolls/pk  |
| K185 | Printer Paper, 5 rolls/pk   |
| K200 | LifeDop Hand-Held Stand w/tray, 5 Castor  |
| K220 | Pinnacle Hand-Held Stand w/tray, 5 Castor                                       |
| K230 | LifeDop Hand-Held Stand w/out tray, 5 Castor                                    |
| K240 | Pinnacle Hand-Held Stand w/out tray, 5 Castor                                   |
| K250 | LifeDop Universal Clamp   |
| K260 | LifeDop Carrying Case   |



P200



P500



P800

## NICOLET IMEX POCKET-DOP II™ DOPPLER SYSTEM

Hand-held ultrasonic Doppler features interchangeable probes and built-in speaker for vascular and obstetrical applications. Supplied in soft sided case with 1 or 2 probes, recharger, walkman style headset, gel, batteries and CD ROM of pulse sounds. Includes a five-year parts and one-year labor warranty.

|      |                                  |
|------|----------------------------------|
| P200 | Unit w/2 MHz OB Probe            |
| P230 | Unit w/2 & 3 MHz Probes          |
| P250 | Unit w/2 & 5 MHz Probes          |
| P280 | Unit w/2 & 8 MHz Probes          |
| P2WP | Unit w/2 MHz Waterproof OB Probe |
| P300 | Unit w/3 MHz OB Probe            |
| P350 | Unit w/3 & 5 MHz Probes          |
| P380 | Unit w/3 & 8 MHz Probes          |
| P500 | Unit w/5 MHz Vascular Probe      |
| P580 | Unit w/5 & 8 MHz Probes          |
| P800 | Unit w/8 MHz Vascular Probe      |
| PAD5 | Unit w/5 MHz Probe & ABI Kit     |
| PAD8 | Unit w/8 MHz Probe & ABI Kit     |



## NICOLET IMEX STETHODOP™ DOPPLER SYSTEM

Mounts to most bell-style stethoscopes using a custom non-chill ring. Easy one-hand operation. Momentary on/off actuator operates only when depressed to save battery life. Sliding volume control to adjust to individual user comfort level. Super sensitive Doppler. Custom pager style carrying case. Includes one-year labor and parts warranty.

|      |   |
|------|---|
| SD50 | Imex StethoDop w/5 MHz Probe                  |
| SD52 | Imex StethoDop w/Cardiology Style Stethoscope |
| SS10 | Imex Cardiology Stethoscope                   |
| SS20 | Imex Classic Stethoscope                      |



ED20



EN2W



EN2WR

## NICOLET IMEX ELITE™

A fully configurable ultrasound Doppler used to determine fetal viability and to assist in monitoring peripheral arterial and venous blood flow. Options include interchangeable OB, vascular and waterproof probes, display or non-display models and alkaline or rechargeable batteries.

|       |   |
|-------|---|
| ED20  | Elite 200 Doppler, Digital Display w/alkaline batteries & 2 MHz OB Probe              |
| ED20R | Elite 200R Doppler, Digital Display w/rechargeable batteries & 2 MHz OB Probe         |
| ED23  | Elite 200 Doppler, Digital Display w/alkaline batteries & 2 & 3 MHz OB Probes         |
| ED23R | Elite 200R Doppler, Digital Display w/rechargeable batteries & 2 & 3 MHz OB Probes    |
| ED25  | Elite 200 Doppler, Digital Display w/alkaline batteries & 2 & 5 MHz Probes            |
| ED25R | Elite 200R Doppler, Digital Display w/rechargeable batteries & 2 & 5 MHz Probes       |
| ED28  | Elite 200 Doppler, Digital Display w/alkaline batteries & 2 & 8 MHz Probes            |
| ED28R | Elite 200R Doppler, Digital Display w/rechargeable batteries & 2 & 8 MHz Probes       |
| ED2W  | Elite 200 Doppler, Digital Display w/alkaline batteries & 2 MHz Waterproof Probe      |
| ED2WR | Elite 200R Doppler, Digital Display w/rechargeable batteries & 2 MHz Waterproof Probe |
| ED30  | Elite 200 Doppler, Digital Display w/alkaline batteries & 3 MHz OB Probe              |
| ED30R | Elite 200R Doppler, Digital Display w/rechargeable batteries & 3 MHz OB Probe         |
| ED35  | Elite 200 Doppler, Digital Display w/alkaline batteries & 3 & 5 MHz Probes            |
| ED35R | Elite 200R Doppler, Digital Display w/rechargeable batteries & 3 & 5 MHz Probes       |
| ED38  | Elite 200 Doppler, Digital Display w/alkaline batteries & 3 & 8 MHz Probes            |
| ED38R | Elite 200R Doppler, Digital Display w/rechargeable batteries & 3 & 8 MHz Probes       |

**NICOLET IMEX ELITE™ (Continued)**

|       |  |
|-------|--|
| ED50  | Elite 200 Doppler, Digital Display w/alkaline batteries & 5 MHz Vascular Probe               |
| ED50R | Elite 200R Doppler, Digital Display w/rechargeable batteries & 5 MHz Vascular Probe          |
| ED58  | Elite 200 Doppler, Digital Display w/alkaline batteries & 5 & 8 MHz Probes                   |
| ED58R | Elite 200R Doppler, Digital Display w/rechargeable batteries & 5 & 8 MHz Vascular Probes     |
| ED80  | Elite 200 Doppler, Digital Display w/alkaline batteries & 8 MHz Vascular Probe               |
| ED80R | Elite 200R Doppler, Digital Display w/rechargeable batteries & 8 MHz Probe                   |
| ED8A  | Elite 200 Doppler, Digital Display w/alkaline batteries, 8 MHz Probe & ABI Kit               |
| EN20  | Elite 100 Doppler, Non-Digital Display w/alkaline batteries & 2 MHz OB Probe                 |
| EN20R | Elite 100R Doppler, Non-Digital Display w/rechargeable batteries & 2 MHz OB Probe            |
| EN23  | Elite 100 Doppler, Non-Digital Display w/alkaline batteries & 2 & 3 MHz OB Probes            |
| EN23R | Elite 100R Doppler, Non-Digital Display w/rechargeable batteries & 2 & 3 MHz OB Probes       |
| EN25  | Elite 100 Doppler, Non-Digital Display w/alkaline batteries & 2 & 5 MHz Probes               |
| EN25R | Elite 100R Doppler, Non-Digital Display w/rechargeable batteries & 2 & 5 MHz Probes          |
| EN28  | Elite 100 Doppler, Non-Digital Display w/alkaline batteries & 2 & 8 MHz Probes               |
| EN28R | Elite 100R Doppler, Non-Digital Display w/rechargeable batteries & 2 & 8 MHz Probes          |
| EN2W  | Elite 100 Doppler, Non-Digital Display w/alkaline batteries & 2 MHz Waterproof Probe         |
| EN2WR | Elite 100R Doppler, Non-Digital Display w/rechargeable batteries & 2 MHz Waterproof Probe    |
| EN30  | Elite 100 Doppler, Non-Digital Display w/alkaline batteries & 3 MHz OB Probe                 |
| EN30R | Elite 100R Doppler, Non-Digital Display w/rechargeable batteries & 3 MHz OB Probe            |
| EN35  | Elite 100 Doppler, Non-Digital Display w/alkaline batteries & 3 & 5 MHz Probes               |
| EN35R | Elite 100R Doppler, Non-Digital Display w/rechargeable batteries & 3 & 5 MHz Probes          |
| EN38  | Elite 100 Doppler, Non-Digital Display w/alkaline batteries & 3 & 8 MHz Probes               |
| EN38R | Elite 100R Doppler, Non-Digital Display w/rechargeable batteries & 3 & 8 MHz Probes          |
| EN50  | Elite 100 Doppler, Non-Digital Display w/alkaline batteries & 5 MHz Vascular Probe           |
| EN50R | Elite 100R Doppler, Non-Digital Display w/rechargeable batteries & 5 MHz Vascular Probe      |
| EN58  | Elite 100 Doppler, Non-Digital Display w/alkaline batteries & 5 & 8 MHz Vascular Probes      |
| EN58R | Elite 100R Doppler, Non-Digital Display w/rechargeable batteries & 5 & 8 MHz Vascular Probes |
| EN5A  | Elite 100 Doppler, Non-Digital Display w/alkaline batteries, 5 MHz Probe & ABI Kit           |
| EN80  | Elite 100 Doppler, Non-Digital Display w/alkaline batteries & 8 MHz Vascular Probe           |
| EN80R | Elite 100R Doppler, Non-Digital Display w/rechargeable batteries & 8 MHz Vascular Probe      |
| EN8A  | Elite 100 Doppler, Non-Digital Display w/alkaline batteries, 8 MHz Probe & ABI Kit           |
| EN8AR | Elite 100R Doppler, Non-Digital Display w/rechargeable batteries, 8 MHz Probe & ABI Kit      |

**NICOLET IMEX POCKET-DOP 3™ DOPPLER SYSTEM**

A valuable aid in detecting pregnancy and monitoring the fetal heartbeat. Features interchangeable probes: the 3 MHz probe is for early fetal heartbeat detection and the 2 MHz probe is for use in later pregnancy or with larger patients. Also features auto-correlation which is a pattern recognition technique that quickly acquires and accurately tracks the complex signal of the fetal heart rate. Comes with a convenient recharging stand to be placed on a counter or mounted to a wall. Covered by a 5-year parts and 1-year labor warranty.

|      |                               |
|------|-------------------------------|
| J200 | Doppler w/2 MHz OB Probe      |
| J230 | Doppler w/2 & 3 MHz OB Probes |
| J300 | Doppler w/3 MHz OB Probe      |



CT20+



CT30+



CT50+



CT80+


**NICOLET IMEXDOP CT+™ DOPPLER SYSTEM**

Counter top/wall mount Doppler designed for hospitals, clinics and heavy volume offices. Features interchangeable probes, convenient recharging stand and a fetal heart rate display in easy-to-read digital numbers. Supplied with recharging stand, gel, battery pack, monaural earphone and mounting kit. Includes a five-year parts, one-year labor warranty.

|       |                                  |
|-------|----------------------------------|
| CT20+ | Unit w/2 MHz OB Probe            |
| CT23+ | Unit w/2 & 3 MHz Probes          |
| CT25+ | Unit w/2 & 5 MHz Probes          |
| CT28+ | Unit w/2 & 8 MHz Probes          |
| CT30+ | Unit w/3 MHz OB Probe            |
| CT35+ | Unit w/3 & 5 MHz Probes          |
| CT38+ | Unit w/3 & 8 MHz Probes          |
| CT50+ | Unit w/5 MHz Vascular Probe      |
| CT58+ | Unit w/5 & 8 MHz Probes          |
| CT80+ | Unit w/8 MHz Vascular Probe      |
| CTWP  | Unit w/2 MHz Waterproof OB Probe |



VLS48



VLS80

**NICOLET IMEX VERSALAB BI-DIRECTIONAL DOPPLER SYSTEM**

Continuous wave doppler transducers optimized for detecting blood flow throughout peripheral parts of the body. Optional color display feature allows quantification, documentation and qualification of blood flow using a color frequency distribution spectral analysis for accurate determination of direction and velocity of blood flow. Fast, built-in printer for quality documentation of all information obtained. A variety of industry standard calculations and indices calculated to complement the diagnostic capabilities of the VersaLab.

|       |   |
|-------|---|
| VLL40 | Unit w/4 MHz LE Bi-Directional Probe      |
| VLL48 | Unit w/4 & 8 MHz LE Bi-Directional Probes |
| VLL80 | Unit w/8 MHz LE Bi-Directional Probe      |
| VLS40 | Unit w/2 MHz LE Bi-Directional Probe      |
| VLS48 | Unit w/4 & 8 MHz LE Bi-Directional Probes |
| VLS80 | Unit w/8 MHz LE Bi-Directional Probe      |



## NICOLET IMEX VERSALAB ACCESSORIES

|       |                                       |
|-------|---------------------------------------|
| A220  | Personal Stereo Headset               |
| A350  | Printer Pack                          |
| A355  | VersaLab Printer Paper, 5 r/l/cs      |
| A430  | Soft-Sided Carrying Case              |
| BD40  | 4 MHz Bi-Directional Probe            |
| BD80  | 8 MHz Bi-Directional Probe            |
| C633  | Battery Pack                          |
| C644  | Power Supply 110-250V                 |
| C644I | Power Supply 220V                     |
| C645  | Power Cord, 110V                      |
| C645I | Power Cord 220V                       |
| DV100 | VersaLab APM/APM2 Event Marker        |
| DV400 | VersaLab Toco Probe Covers, 5/pk      |
| ST30  | Tilt Stand                            |
| TV100 | VersaLab APM/APM2 Toco Probe          |
| TV180 | VersaLab APM2, 1.8 MHz Doppler Probe  |
| TV200 | VersaLab APM/APM2, 2MHz Doppler Probe |
| V200  | VersaLab IV Pole Mount                |



MVAPM



MVAPT

## NICOLET IMEX VERSALAB ANTEPARTUM MONITOR™

Easy-to-use antepartum monitor designed for the office or the antepartum testing area of the hospital. Pulsed, directional 1.5 MHz ultrasound locks into the fetal heart rate; sensitive toco transducer tracks maternal contractions. Trace prints out on Z-fold paper for easy scanning and storage. Supplied with chart paper, belt, plastic membrane for toco, gel and patient event marker. Includes a one-year parts and labor warranty.

|       |  |
|-------|--|
| MVAPM | Single APM Unit w/toco & 2 MHz Doppler Transducer    |
| MVAPT | Twin APM Unit w/toco 1.8 & 2 MHz Doppler Transducers |

## NICOLET IMEX ANTEPARTUM MONITOR ACCESSORIES

|      |                      |
|------|----------------------|
| D200 | Z-Fold Chart Paper   |
| D300 | Belts, Pre-Cut, 2/pk |



FD2PUSA/EZ8



FD2PUSA/OP2HS

## HUNTLEIGH FETAL DOPPLEX® II POCKET DOPPLER

Advanced pocket Doppler offers three Fetal Heart Rate (FHR) modes, LCD display, auto-correlation, active noise reduction, digital interface to Dopplex® Printa™ to produce FHR tracings. Supplied with one transducer of choice.

|                |  |
|----------------|--|
| FD2PUSA/EZ8    | Fetal Dopplex® II w/EZ8 Widebeam Probe |
| FD2PUSA/OP2HS  | Fetal Dopplex® II w/2 MHz Probe        |
| FD2PUSA/OP3HS  | Fetal Dopplex® II w/3 MHz Probe        |
| FD2PUSA/VP10HS | Fetal Dopplex® II w/10 MHz Probe       |
| FD2PUSA/VP4HS  | Fetal Dopplex® II w/4 MHz Probe        |
| FD2PUSA/VP5HS  | Fetal Dopplex® II w/5 MHz Probe        |
| FD2PUSA/VP8HS  | Fetal Dopplex® II w/8 MHz Probe        |



SD2PUSA/EZ8



SD2PUSA/OP2HS

## HUNTLEIGH SUPER DOPPLEX® II DOPPLER

Sophisticated pocket Doppler with bi-directional indication of blood flow velocity, LC display, true separated audio signals, clear visual indication of forward and reverse flow, 4 level gain control, active noise reduction, connects to EKG with custom cable. Supplied with one transducer of choice.

|                |                                      |
|----------------|--------------------------------------|
| SD2PUSA/EZ8    | Super Dopplex 2 w/EZ8 Widebeam Probe |
| SD2PUSA/OP2HS  | Super Dopplex 2 w/2 MHz Probe        |
| SD2PUSA/OP3HS  | Super Dopplex 2 w/3 MHz Probe        |
| SD2PUSA/VP10HS | Super Dopplex 2 w/10 MHz Probe       |
| SD2PUSA/VP4HS  | Super Dopplex 2 w/4 MHz Probe        |
| SD2PUSA/VP5HS  | Super Dopplex 2 w/5 MHz Probe        |
| SD2PUSA/VP8HS  | Super Dopplex 2 w/8 MHz Probe        |



## HUNTLEIGH FETAL DOPPLEX® I PLUS DOPPLER

Fetal dopplex with LCD display of fetal heart rate, fixed 2 MHz transducer and built-in loudspeaker. Waterproof.

|            |                       |
|------------|-----------------------|
| FD1+-P-USA | Fetal Dopplex® I Plus |
|------------|-----------------------|



## HUNTLEIGH AUDIO DOPPLEX® DOPPLER (WITHOUT DISPLAY)

Dedicated fetal Doppler with fixed 2 MHz transducer and integral loudspeaker. Waterproof. Includes carry case, headphones, alkaline battery (9V) and gel.

|            |                            |
|------------|----------------------------|
| D920-P-USA | Audio Dopplex® w/o display |
|------------|----------------------------|






**HUNTLEIGH TRANSDUCERS**

For FD2, MD2, SD2, and D900.

- EZ8-HS EZ8 Widebeam Probe
- OP2-HS 2 MHz Probe For Fetal Heart Detection, High Sensitivity
- OP3-HS 3 MHz Probe For Fetal Heart Detection, High Sensitivity
- VP10-HS 10 MHz Probe For Very Superficial Vascular Applications, High Sensitivity
- VP4-HS 4 MHz Probe For Very Deep Vessels, High Sensitivity
- VP5-HS 5 MHz Probe Moderately Deep Vessels, High Sensitivity
- VP8-HS 8 MHz Transducer For Peripheral/Superficial Vessels, High Sensitivity



**HUNTLEIGH BABY DOPPLEX® DOPPLER 3000 MARK 2**

The advanced electronic design enables accurate data to be produced using autocorrelation signal processing techniques, and displayed and recorded on the thermal paper trace. The pulsed Doppler, 7 Crystal Array Ultrasound Transducer provides excellent fetal heart pick-up and wrap-around bumpers provide additional impact protection. The standard monitor also includes: ultrasound channel, contractions measurement (Toco), automatic fetal movement detection with adjustable sensitivity, patient event marker, clinical event marker and annotation facility, clear LED indication of FHR and UA, high quality audio FHR reproduction, signal quality indicator (dual for twins), integral thermal printer.

BD3000MK2 Baby Doppler® 3000 Mark 2

**HUNTLEIGH TWIN BABY DOPPLEX® DOPPLER 3002**

Having all of the features of the 3000 Mark 2 version, the 3002 has a unique and extremely useful feature. Ultrasound FHR traces are produced on "Double Width" paper allowing two completely separate traces to be viewed. This avoids possible confusion which can occur, if conventional twins monitors which superimpose one trace on top of another, are used. A carry case is available to transport the products from clinic to clinic.

BD3002 Baby Doppler® 3002

**HUNTLEIGH BABY DOPPLEX® DOPPLER 3000 MARK 2 ACCESSORIES**

- ACC-15 Printer Paper, Plain, No Grid, 10/bx
- ACC-30 Carry Case/Shoulder Bag
- ACC-31 Wall Bracket
- ACC-32 Desk Top Tray
- ACC-53 Trolley Cart



**HUNTLEIGH FLEXI DOPPLEX® DOPPLER**

Dedicated fetal Doppler with fixed 3 MHz transducer and integral loudspeaker. Waterproof. Includes carry case, headphones, alkaline battery (9V) and gel.


D930-P-USA Flexi Doppler®



**HUNTLEIGH MINI DOPPLEX® DOPPLER (WITHOUT DISPLAY)**

Non-directional pocket Doppler with integral loudspeaker. Supplied with one transducer of choice, carry case, headphones, alkaline battery (9V), gel and operating instructions.

- D900USA/EZ8 Mini Doppler w/EZ8 Widebeam Probe
- D900USA/OP2HS Mini Doppler w/2 MHz Probe
- D900USA/OP3HS Mini Doppler w/3 MHz Probe
- D900USA/VP10HS Mini Doppler w/10 MHz Probe
- D900USA/VP4HS Mini Doppler w/4 MHz Probe
- D900USA/VP5HS Mini Doppler w/5 MHz Probe
- D900USA/VP8HS Mini Doppler w/8 MHz Probe



**HUNTLEIGH MULTI DOPPLEX® II DOPPLER**

Sophisticated pocket Doppler with bi-directional indication of blood flow velocity, 4 level gain control, LCD display, outputs and true separated audio signals, active noise reduction, digital interface to Dopplex® Printa/EKG/or Dopplex® Reporter software package for printout of waveforms (RS232). Supplied with one transducer of choice.

- MD2PUSA/EZ8 Multi Doppler II Pocket w/EZ8 Widebeam Probe
- MD2PUSA/OP2HS Multi Doppler II Pocket w/2 MHz Probe
- MD2PUSA/OP3HS Multi Doppler II Pocket w/3 MHz Probe
- MD2PUSA/VP10HS Multi Doppler II Pocket w/10 MHz Probe
- MD2PUSA/VP4HS Multi Doppler II Pocket w/4 MHz Probe
- MD2PUSA/VP5HS Multi Doppler II Pocket w/5 MHz Probe
- MD2PUSA/VP8HS Multi Doppler II Pocket w/8 MHz Probe



**HUNTLEIGH DOPPLER ACCESSORIES**

- ACC-14 Thermal Dopplex® Printa™ Paper, Plain, No Grid
- ACC-67 Ultrasonic Gel, 12/bx
- DP100 Dopplex® Stand, Locking Device To Help Prevent The Doppler From Being Mishandled, Dropped or Stolen



**PARKER LABS AQUASONIC® 100 ULTRASOUND TRANSMISSION GEL**

The most widely used gel for diagnostic and therapeutic medical ultrasound. Recommended for all procedures where a viscous gel is required. Completely aqueous, will not stain clothing or damage transducers. Hypoallergenic, bacteriostatic, non-sensitizing and non-irritating. No formaldehyde. Not a spermicide.

- 01-02 Ultrasound Transmission Gel, Non-Sterile, 60g Doppler® Size Tube, 12/bx, 12 dz/cs
- 01-08 Ultrasound Transmission Gel, Non-Sterile, 0.25 Liter Bottle, 6 dz/cs
- 01-34 Ultrasound Transmission Gel, Non-Sterile, 1 Liter dispenser bottle, 6/bx, 12/cs
- 01-50 Ultrasound Transmission Gel, Non-Sterile, 5 Liter Sonicpac® w/dispenser bottle, 4/cs



**PARKER LABS AQUASONIC® CLEAR ULTRASOUND GEL**

Based on the same "can't be copied" formula as Aquasonic 100 Ultrasound Transmission Gel, without fragrance or color. Acoustically correct for a broad range of frequencies. Hypoallergenic, non-irritating and bacteriostatic. Water soluble. No formaldehyde. Not a spermicide.

- 03-02 Ultrasound Gel, 60g Doppler® Size Tube, 12/bx, 12 dz/cs
- 03-08 Ultrasound Gel, .25 Liter Bottle, 6 dz/cs
- 03-34 Ultrasound Gel, 1 Liter w/dispenser cap, 6/bx, 2 bx/cs
- 03-50 Ultrasound Gel, 5 Liter SONICPAC® w/disposable bottle, 4/cs
- 03-54 Ultrasound Gel, 5 Liter ECONOPAC w/2 dispensers, 4/cs



**PARKER LABS STERILE AQUASONIC® 100 ULTRASOUND TRANSMISSION GEL**

Recommended for diagnostic and therapeutic ultrasound applications when sterility is indicated. Overwrapped, sterilized foil pouches. Acoustically correct for the broad range of frequencies used in medical ultrasound. Non-irritating, hypoallergenic. Completely aqueous, will not stain clothing or damage equipment. Not a spermicide.

- 01-01 Ultrasound Transmission Gel, Sterile, 20g Overwrapped Sterilized Foil Pouches, 48/bx, 6 bx/cs



**PARKER LABS THERMASONIC® GEL WARMER**

Warms three bottles of ultrasound gel to body temperature. Durable wall mountable construction with warranty. Durable case. Low electrical demand. UL and CSA listed.

- 82-04 Three-Bottle Gel Warmer (plus 3 dispensers of ultrasound gel), 4/cs

**PARKER LABS TRANSEPTIC® CLEANSING SOLUTION**

Multi-purpose spray for professional use in cleansing ultrasound transducer/probe surfaces between procedures. Will not cause transducer membrane swelling. Does not contain phenol or glutaraldehyde.

- 09-25 Cleansing Solution, 250ml (8.5 fl oz) Clear Spray Bottle, 12/bx, 4 bx/cs



**CHESTER LIQUASONIC® ULTRASOUND GEL**

LiquaSonic Ultrasound Transmission Gel has a thicker, no-slip formula that remains viscous even when applied to salty skin. Our salt-free, formaldehyde free formula eliminates the air between the transducer and skin, offering exceptional coupling efficiency for more effective treatment/diagnosis. LiquaSonic Gel has been formulated to be acoustically correct for the broad range of frequencies used in medical ultrasound applications. Our formula is water soluble and non-greasy preventing damage to expensive transducers and clothing.


- 001205 8 oz bottle w/mustard cap, Ultrasound Gel, 12/cs
- 001222 5 Liter Collapsible Cubitainer w/empty bottle, Ultrasound Gel, 4/cs
- 001297 128 oz Multi-Purpose, Medium Viscosity Ultrasound Gel, 4/cs
- 001298 8 oz Bottle w/mustard cap, Multi-Purpose, Medium Viscosity Ultrasound Gel, 12/cs
- 001305 8 oz Bottle w/mustard cap, Clear Ultrasound Gel, 12/cs
- 001333 128 oz Clear Ultrasound Gel, 4/cs



**CHESTER LIQUASONIC® ULTRASOUND LOTION**

Similar in use to the ultrasound gel, the lotion is often preferred by medical professionals. The aqueous formulation is pleasing to the patient and is acoustically correct for the large spectrum of frequencies used in medical ultrasound applications. The lotion is hypoallergenic, bacteriostatic, salt free, formaldehyde free and will not harm clothing or equipment.

- 001267 128 oz Ultrasonic Lotion, 4/cs
- 001268 8 oz Bottle w/mustard cap, Ultrasonic Lotion, 12/cs



**CIDA PHYSICIAN'S PENLIGHT**

Penlight has high-intensity prefocused bulb powered by two AAA batteries. Aluminum contact strip assures good contact for brighter light. Extruded case has a unique V-clip. 1-year shelf life. Useful life is 3- to 6-months with intermittent use.

- 6666C Dr. Pack Disposable Penlight, 6/pk



**AARON DIAGNOSTIC LIGHT EXAM KITS**

Disposable. Two AAA batteries. Cobalt blue with Wratten #47A cobalt filter. For use with fluorescein dye. One-year shelf life.


- 6303 Penlights w/cobalt filter, 3/pk



**AARON FLUORO-DOT™ PENLIGHT**

Battery-operated penlight designed for use with fluorescein dye. Excellent for diagnosing abrasions and fitting contact lenses. Constant on/off switch. Contains two AA alkaline batteries. Penlight has a 4-year shelf life and a 4-hour burn time.

- 6803 Fluoro Slit™ Penlight
- 6802 Fluoro Dot™ Penlight w/wratten cobalt blue filter



222K, 7156

**AARON PENLIGHTS**

- 1125 Disposable Penlight
- 1135 Reusable Penlight

**Accessories**

- 222 Replacement Bulb For 1135
- 7156 Cobalt Filter Only



**AARON PHYSICIAN'S PENLIGHT**

Penlight has high-intensity prefocused bulb powered by two AAA batteries. Aluminum contact strip assures good contact for brighter light. Extruded case has a unique V-clip. 1-year shelf life. Useful life is 3- to 6-months with intermittent use.

- 6666 Dr. Pack Disposable Penlight, 6/pk
- 666R Doctors' Replaceable Battery Penlight, 6/pk
- 66RN RN Pack Disposable Penlight w/pupil gauge, 6/pk



**AARON SPECIALTY LIGHTS**

Individually packaged. Two AA batteries. All lights have a 4-year shelf life.

- 6800 Tech Light™ w/concentric beam
- 6801 White Slit Beam For Cataract Screening
- 6804 1mm Light (for muscle tests)



**AARON WOODS LIGHT**

- UV59 Hand-Held Ultra-Violet Light

**Accessories**

- UV59B Replaceable Black Lamp For UV59



### ADC® ADLITE™ DISPOSABLE PENLIGHT

Sealed concentric beam activated by depressing pocket clip. Available with or without pupil gauge. 5" long, 1/2" diameter.

- 351 Disposable Penlight, 6/pk
- 351P Disposable Penlight w/pupil gauge, 6/pk



### ADC® ADLITE II™ PENLIGHT

Reusable penlight with attractive tapered barrel. Activated by depressing pocket clip. Two AAA batteries included. Measures 5" long, 1/2" diameter.

- 354 Reuseable Penlight



### ADC® ADLITE™ PLUS PENLIGHT

Sealed concentric beam, disposable diagnostic penlight. Activated by depressing pocket clip. Measures 5" long, 1/2" diameter. Sold individually or in display cylinder of 22.

- 356BK ADLITE™ Disposable Penlight, Black
- 356OR ADLITE™ Disposable Penlight, Orange
- 356R ADLITE™ Disposable Penlight, Red
- 356TL ADLITE™ Disposable Penlight, Teal
- 356V ADLITE™ Disposable Penlight, Purple
- 356W ADLITE™ Disposable Penlight, White
- 356WP ADLITE™ Disposable Penlight, White w/pupil gauge
- 358AST Disposable Penlight, Assorted Colors, 22/cylinder
- 358TVB Disposable Penlight, Black, Purple & Teal, 22/cylinder
- 358W Disposable Penlight, White, 22/cylinder
- 358WP Disposable Penlight, White w/pupil gauge, 22/cylinder



### ADC® METALITE™ PENLIGHT

Reusable with brushed aluminum housing. Activated by depressing pocket clip. Two AAA batteries included. Measures 6" long, 1/2" diameter.

- 352 Penlight

### ADC® POCKET PAL III™ ORGANIZER AND KIT

An indispensable uniform accessory that organizes while it protects. The organizer features washable 420D nylon pocket sewn with nylon thread. Five slide-in pockets and one covered pocket. Key ring. Two matched closures secure instruments and double as ID tags. Sewn loop permits belt wearing. Measures 6 1/4" x 4 3/4" x 1/4". Weighs 1 oz. The kit includes: Pocket Pal™ Organizer (216), 5 1/2" Lister Bandage scissors, 5" straight Mosquito forceps, 3 color pen and ADLITE™ disposable penlight. All components are latex free.

- 117BK Pocket Pal III™ Kit, Black
- 117TL Pocket Pal III™ Kit, Teal
- 117V Pocket Pal III™ Kit, Purple

### ADC® TOGGELITE™ PENLIGHT

Constructed of chrome plated brass for rugged durability. Activated by a locking toggle switch. Two AA batteries included. Measures 6" long, 3/4" diameter.

- 353 Penlight

### OMRON HIGH INTENSITY DIAGNOSTIC PENLIGHTS

Focused beam lens-end lamps. Each light has a clip to prevent accidental discharge of the batteries when carried in a pocket. Lights weigh less than one ounce.

- 24-100 Reusable Pen Light, 6/pk
- 24-150 Reusable Pen Light, 100/bx



### WELCH ALLYN PROFESSIONAL PENLITE

Halogen light for true tissue color and consistent, long lasting illumination. Brass with black matte finish for durability.

- 76600 PenLite





21504



24031



25020

**WELCH ALLYN 2.5V/3.5V HALOGEN DIAGNOSTIC OTOSCOPE**

Halogen illumination with fiber optic light transmission provides brilliant, cool, distal illumination without visual obstruction or specular reflections. Wide angle viewing lens slides left or right for instrumentation under magnification.

- 24031 2.5V Halogen Otoscope w/4 polypropylene specula (sizes 2, 3, 4 & 5 mm)
- 24032 2.5V Halogen Otoscope w/insufflator bulb, tube w/tip & 4 polypropylene specula (sizes 2, 3, 4 & 5 mm)
- 24035 2.5V Halogen Otoscope w/o specula
- 25020 3.5V Halogen Otoscope Head w/4 polypropylene specula (sizes 2, 3, 4, 5 mm)
- 25021 Same as No. 25020, Including Insufflator Bulb & Tube

**Accessories**

- 200055-502 Lens Holder Assembly
- 21504 Insufflator Bulb & Tube w/tip



21600



21700

**WELCH ALLYN 2.5V/3.5V HALOGEN OPERATING OTOSCOPE**

Open system for procedures with fully rotatable 5/16" (7.9 mm) diameter magnifying lens and rotatable speculum. Halogen illumination with fiber optic light transmission.

- 21600 2.5V Otoscope w/5 polypropylene specula (sizes 2, 3, 4, 5 & 9mm)
- 21601 2.5V Otoscope w/o specula
- 21700 3.5V Otoscope w/5 polypropylene specula (sizes 2, 3, 4, 5, 9 mm)
- 21701 Same as No. 21700, w/o specula

**Accessories**

- 209026-501 Lens and Attachable Arm




20100



20200



20201



21501

**WELCH ALLYN 2.5V/3.5V HALOGEN PNEUMATIC OTOSCOPE**

Closed system, with insufflator fitting for pneumatic otoscopy. Large magnifying lens can be easily rotated to facilitate use of swabs, curettes, etc. Halogen illumination with fiber optic light transmission.

- 20100 2.5V Otoscope w/5 polypropylene specula (sizes 2, 3, 4, 5 & 9 mm)
- 20101 2.5V Otoscope w/o specula
- 20200 3.5V Otoscope w/5 polypropylene specula (sizes 2, 3, 4, 5, 9 mm)
- 20201 Same as No. 20200, w/o specula

**Accessories**

- 21501 Insufflator Bulb & Tube w/o tip

**WELCH ALLYN 3.5V HALOGEN PNEUMATIC OTOSCOPE/THROAT ILLUMINATOR**

Dual-purpose instrument allows practitioner to examine the ears, then detach otoscope head and have an instant, high-intensity throat illuminator. Wide angle viewing lens slides left or right. Halogen illumination with fiber optic light transmission.

- 25257 Illuminator Section Only



**WELCH ALLYN MACROVIEW™ OTOSCOPE**

30% higher magnification over standard otoscopes. Nearly a full view of the tympanic membrane. Ability to focus to fine tune the view. Greater tip retention and easier tip removal.

- 23810 3.5V Otoscope w/specula (head only)
- 23814 3.5V Otoscope (23810) w/insufflation bulb
- 23820 3.5V Otoscope w/throat illuminator
- 23824 3.5V Otoscope (23820) w/insufflation bulb

**WELCH ALLYN MACROVIEW™ OTOSCOPE ACCESSORIES**

- 05215-M Soft Case For 3.5V Set
- 05258-M Hard Case For PanOptic Ophthalmoscope 3.5V Set
- 05259-M Hard Case For 3.5V Set
- 05815-M Soft Case For PanOptic Ophthalmoscope 3.5V Set
- 06500-U 3.5V Halogen Lamp For Otoscope
- 23804 Insufflation Bulb For Otoscope
- 23857 3.5V Throat Illuminator Section Only For Otoscope
- 26038 Bivalve Nasal Speculum For MacroView™ Complete
- 26538 Nasal Illuminator For MacroView™ Complete



22820



22821

**WELCH ALLYN 2.5V POCKETSCOPE™ OTOSCOPE/THROAT ILLUMINATOR**

Halogen light for true tissue color and consistent, long lasting illumination. Fiber optics for cool light with no reflection, no obstruction. Wide angle viewing lens. Sealed system for pneumatic otoscopy. Built-in throat illuminator/penlight for added versatility.

- 21111 PocketScope Otoscope/Throat Illuminator w/o handle
- 22800 PocketScope Otoscope/Throat Illuminator w/handle & rechargeable battery
- 22811 PocketScope Otoscope/Throat Illuminator w/handle, rechargeable battery & soft case
- 22820 PocketScope Otoscope/Throat Illuminator w/AA alkaline battery handle
- 22821 PocketScope Otoscope/Throat Illuminator w/AA alkaline battery handle & soft case

**Accessories**

- 211005 Replacement Lens

**ADC® STANDARD & VETERINARY OTOSCOPE**

German 2.7V lamp for bright illumination. Bayonet locking otoscope head. Removable 4X magnifying lens. Rheostatic on/off switch. Autoclavable, polypropylene ear specula in sizes 2.5mm, 3.5mm and 4.5mm (one of each). Veterinary model includes 4mm, 5mm, and 7mm specula sizes. Knurled brass handle with a beautiful satin finish. Ergonomically designed for comfort. Compact instrument measures just under 7 1/2" and weighs approximately 7 oz. Powered by 2 C cell batteries (not included). Comes with zippered leatherette storage pouch that holds the otoscope and included speculas. Inspected and packaged in the U.S.A. from foreign components. One year quality performance guarantee.

- 5211 Standard Otoscope
- 5211V Veterinary Otoscope

**ADC® STANDARD & VETERINARY OTOSCOPE ACCESSORIES**

- 5211-1 Standard Otoscope Specula, 2.5mm
- 5211-2 Standard Otoscope Specula, 3.5mm
- 5211-3 Standard Otoscope Specula, 4.5mm
- 5211-4 Standard Otoscope Replacement Lamp, 2.7V
- 5211-7 Veterinary Otoscope Specula, 4.0mm
- 5211-8 Veterinary Otoscope Specula, 5.0mm
- 5211-9 Veterinary Otoscope Specula, 7.0mm
- 5211-AD Specula Adapter For 5211

**MILTEX BRUENING MAGNIFIER OTOSCOPE**

Chrome. 3.5x magnification with 3 specula, small, medium and large sizes. Complete with inflation bulb and tubing in a fitted case.

- 19-70 Otoscope Set
- 19-71 Set of 3 Specula: Small, Medium & Large Sizes For Otoscope
- 19-71-1 Small Speculum
- 19-71-2 Medium Speculum
- 19-71-3 Large Speculum

**MILTEX SIEGEL MAGNIFIER OTOSCOPE**

Chrome. 3.15x magnification with 3 specula: small, medium and large sizes. Complete with inflation bulb and tubing.

- 19-75 Rubber Bulb w/tubing For Otoscope



22002



24305

**WELCH ALLYN POLYPROPYLENE REUSABLE SPECULA**

Material is lightweight and durable. Can be autoclaved, boiled or cleaned with standard germicides.

- 22002 2 mm Speculum, For Use With Pneumatic, Operating & Consulting Otoscopes, Dark Green
- 22003 3 mm Speculum, For Use With Pneumatic, Operating & Consulting Otoscopes, Dark Green
- 22004 4 mm Speculum, For Use With Pneumatic, Operating & Consulting Otoscopes, Dark Green
- 22005 5 mm Speculum, For Use With Pneumatic, Operating & Consulting Otoscopes, Dark Green
- 22009 9 mm Speculum (nasal), For Use With Pneumatic, Operating & Consulting Otoscopes, Dark Green
- 22100 Set of 5 Specula (2, 3, 4, 5, 9mm), For Use With Pneumatic, Operating & Consulting Otoscopes, Dark Green
- 24302-U 2 mm Speculum, Universal Poly Reusable
- 24303-U 3 mm Speculum, Universal Poly Reusable
- 24304-U 4 mm Speculum, Universal Poly Reusable
- 24305 5 mm Speculum, For Use With Diagnostic Otoscopes, Black
- 24305-U 5 mm Speculum, Universal Poly Reusable
- 24400-U Set of 4 Speculum, Universal Poly Reusable



**WELCH ALLYN KLEENSPEC® DISPOSABLE OTOSCOPE SPECULA**

Make of nontoxic plastic; to be used once then discarded. No adapter required. Black.

- 52133 3 mm Specula, For Use With Pneumatic, Operating & Consulting Otoscopes, 500/bg, 10 bg/cs
- 52134 4 mm Specula, For Use With Pneumatic, Operating & Consulting Otoscopes, 500/bg, 10 bg/cs
- 52135 5 mm Specula, For Use With Pneumatic, Operating & Consulting Otoscopes, 500/bg, 10 bg/cs
- 52432-U 2.75mm Universal KleenSpec Disposable Otoscope Specula, 850/bg, 10 bg/cs
- 52434-U 4.25mm Universal KleenSpec Disposable Otoscope Specula, 850/bg, 10 bg/cs



**WELCH ALLYN SOFSPEC™ OTOSCOPE SPECULA**

Polypropylene base with soft, flexible tip provides improved seal for pneumatic otoscopy. Can be autoclaved, boiled or cleaned with standard germicides.

- 22023 3 mm Speculum, For Use With Pneumatic, Operating and Consulting Otoscopes
- 22025 5 mm Speculum, For Use With Pneumatic, Operating and Consulting Otoscopes
- 22027 7 mm Speculum, For Use With Pneumatic, Operating and Consulting Otoscopes
- 22120 Set of Above 3 Specula, For Use With Pneumatic, Operating and Consulting Otoscopes
- 24323 3 mm Speculum, For Use With Diagnostic Otoscopes
- 24325 5 mm Speculum, For Use With Diagnostic Otoscopes
- 24327 7 mm Speculum, For Use With Diagnostic Otoscopes
- 24420 Set of Above 3 Specula, For Use With Diagnostic Otoscopes

**SPECLINC OTREC® DISPOSABLE OTOSCOPE SPECULA**

Complete speculum is disposable, providing a clean tip for each use. No covers or adapters are needed. Black color for less reflection and greater clarity.

- 7225 2.5 mm Otoscope Specula, 500/pk, 10 pk/cs (5000/cs)
- 7250 Pediatric, 2.5 mm Otoscope Specula, (cross reference #52432), 1000/pk, 10 pk/cs (10,000/cs)
- 7300 Youth, 3 mm Otoscope Specula, 1000/pk, 10 pk/cs (10,000/cs)
- 7333 3 mm Otoscope Specula (cross reference #52133) , 500/pk, 10 pk/cs (5000/cs)
- 7400 Med. Adult, 4 mm Otoscope Specula (cross reference #52434), 1000/pk, 10 pk/cs (10,000/cs)
- 7444 4 mm Otoscope Specula (cross reference #51234), 500/pk, 10 pk/cs (5000/cs)
- 7500 Lg. Adult, 5 mm Otoscope Specula, 1000/pk, 10 pk/cs (10,000/cs)
- 7555 5 mm Otoscope Specula (cross reference #52135), 500/pk, 10 pk/cs (5000/cs)
- 7666 6 mm Otoscope Specula, 500/pk, 10 pk/cs (5000/cs)
- 7999 9 mm Otoscope Specula, 500/pk, 10 pk/cs (5000/cs)

**SPECLINC OTREC DISPOSABLE OTOSCOPE SPECULA WALL DISPENSERS**

- 7000 Wall Dispenser for Diagnostic: 7250, 7300, 7400, 7500 5/pk (includes Wall Bracket), 10pk/cs
- 7777 Wall Dispenser for Operating and Consulting: 7225, 7333, 7555, 7666 and 7999, 5/pk (includes Wall Bracket), 10pk/cs



**WELCH ALLYN 2.5V/3.5V HALOGEN OPHTHALMOSCOPE**

Features high-intensity halogen illumination and mirror optical system. 28 lenses with illuminated lens dial. Apertures (microspot, small spot, large spot, fixation, red-free, and slit) for all diagnostic applications.

- 11470 2.5V Halogen Ophthalmoscope, Head Only
- 11710 3.5V Halogen Ophthalmoscope, Head Only



**WELCH ALLYN HALOGEN COAXIAL AUTOSTEP OPHTHALMOSCOPE**

Provides 68 viewing lenses with continuous single diopter steps or automatic progression from -30 to +38. Wide range of diopter steps facilitates precise focusing. Lens values fully displayed, no mathematical calculations required. Advanced halogen coaxial optics for high intensity illumination.

- 11730 3.5-V Halogen Coaxial Autostep Ophthalmoscope



**WELCH ALLYN HALOGEN COAXIAL OPHTHALMOSCOPE**

23 lenses. Illuminated lens dial. Sealed optics provide better dust protection. Coaxial optical system eliminates shadows and facilitates exam through virtually any size pupil or vitreous disorder. Special contoured mirror and matching lamp filament allow axis of illumination to be practically coincidental with axis of physician's vision. Apertures include microspot, small spot, large spot, fixation, cobalt blue, and slit.

- 11720 3.5V Halogen Coaxial Ophthalmoscope, Head Only



**WELCH ALLYN PANOPTIC™ OPHTHALMOSCOPE**

Patented new Axial PointSource™ Optics. Illumination system converges the light to a point at the cornea, allowing easy entry into small, undilated pupils. The illumination pathway then diverges to the retina, illuminating a very wide area of the fundus. Viewing system enables the operator to view the illuminated area on virtually the same axis, creating the widest field-of-view attainable. Provides a five-times larger view of the fundus than ever before achieved with a standard ophthalmoscope in an undilated eye. Increases magnification by 26% over a standard ophthalmoscope. Patented Glare Extinguishment System prevents interference from unwanted glare and reflections. Works with all existing Welch Allyn 3.5v power sources (battery handle or wall transformer) saving additional investment. The patient eyecup establishes and maintains proper viewing distance, provides stabilization for the view, is a pivoting point for leverage in panning around the retina, occludes ambient light, and helps with orientation. The dynamic focusing wheel adjusts the focus in a continuous, smooth action for more precise control and the optimum view with a focusing range of -20 to +20 diopters. The aperture dial has micro, small, and large spot sizes, a slit aperture, and red-free filter, and cobalt blue filter with add-on corneal magnifying lens (Model 11820 only).

- 11810 PanOptic Ophthalmoscope w/o Blue Filter
- 11820 PanOptic Ophthalmoscope w/Blue Filter & Lens

**Accessories**

- 03800-U PanOptic Lamp
- 05258-MBX PanOptic Diagnostic Case
- 05815-MBX PanOptic Diagnostic Soft Case
- 11870 Patient Eye Cup, 5/set
- 11875 Corneal Lens





## WELCH ALLYN 2.5V HALOGEN STREAK RETINOSCOPIES

Continuously rotating external focusing sleeve is easy to grip and manipulate regardless of hand size. Magnetic fixation cards. Rubber brow rest prevents scratching of eyeglasses. Halogen light for bright, consistent, long-lasting illumination.

17200 2.5V Halogen Streak Retinoscope



18320

## WELCH ALLYN 3.5V HALOGEN RETINOSCOPE SETS

- 18310 18200 Retinoscope, 11710 Ophthalmoscope, 71000 Handle, 05820 Hard Case
- 18311 18200 Retinoscope, 11710 Ophthalmoscope, 71000 Handle, 05180 Soft Case
- 18320 18200 Retinoscope, 11720 Ophthalmoscope, 71000 Handle, 05820 Hard Case
- 18321 18200 Retinoscope, 11720 Ophthalmoscope, 71000 Handle, 05180 Soft Case
- 18330 18200 Retinoscope, 11730 Ophthalmoscope, 71000 Handle, 05820 Hard Case
- 18331 18200 Retinoscope, 11730 Ophthalmoscope, 71000 Handle, 05180 Soft Case
- 18340 18200 Retinoscope, 71000 Handle, 05820 Hard Case
- 18341 18200 Retinoscope, 71000 Handle, 05180 Soft Case
- 18410 18300 Retinoscope, 11710 Ophthalmoscope, 71000 Handle, 05820 Hard Case
- 18411 18300 Retinoscope, 11710 Ophthalmoscope, 71000 Handle, 05180 Soft Case
- 18420 18300 Retinoscope, 11720 Ophthalmoscope, 71000 Handle, 05820 Hard Case
- 18421 18300 Retinoscope, 11720 Ophthalmoscope, 71000 Handle, 05180 Soft Case
- 18430 18300 Retinoscope, 11730 Ophthalmoscope, 71000 Handle, 05820 Hard Case
- 18440 18300 Retinoscope, 71000 Handle, 05820 Hard Case
- 18441 18300 Retinoscope, 71000 Handle, 05180 Soft Case



18200



18250



18300

## WELCH ALLYN 3.5V RETINOSCOPIES

External focusing sleeve. Brighter halogen light. Crossed linear polarizing filter. Magnetic, age appropriate targets. One-hand operation for streak focus and 360° streak rotation. Rubber brow rest. Interchangeable streak retinoscope can be converted to a spot retinoscope by simple changing the lamp.

- 18200 3.5V Streak Retinoscope
- 18300 3.5V Spot Retinoscope

### Accessories

- 18250 Complete Set of Fixation Cards



## WELCH ALLYN 2.5V POCKETSCOPE™ OPHTHALMOSCOPE

Halogen light for true tissue color and consistent, long lasting illumination. Five aperture selections. 48 lenses for better resolution. Polarized filter to reduce corneal reflection. Compact size fits in pocket.

- 12800 PocketScope Ophthalmoscope w/handle & rechargeable battery
- 12811 PocketScope Ophthalmoscope w/handle, rechargeable battery & soft case
- 12820 PocketScope Ophthalmoscope w/AA alkaline battery handle
- 12821 PocketScope Ophthalmoscope w/AA alkaline battery handle & soft case
- 13000 PocketScope Ophthalmoscope w/o handle

### Accessories

- 05280-1 Compact Zipper Case
- 05928-U Soft Case



## WELCH ALLYN STRABISMOSCOPE™ UNIDIRECTIONAL OCCLUDER

Test for both heterophorias (latent Strabismus) and heterotropias (manifest Strabismus) at the same time. One-way mirror allows direct observation of occluded eye. View movement of both covered and uncovered eyes simultaneously. Allows for confident diagnosis of Dissociated Vertical Deviation. Rugged ABS housing withstands drops and falls. Lightweight and easy to clean. Fits all Welch Allyn 3.5V power sources. Handles must be ordered separately.

- 12400 StrabismoScope Head Only
- 12451 3.5-V Halogen StrabismoScope Set Includes StrabismoScope, Rechargeable Handle & Soft Carrying Case

## ADC® COMPLETE 2.5V INSTRUMENT SET

Five different instruments in a fitted case. Features: bayonet locking mounts for each instrument head. C cell powered, satin finish chrome-plated battery handle with knurled finish. Rheostatic on/off switch. Complete in zippered case with fitted foam insert that is scuff and scratch resistant and cleans easily

- 5215 Complete Diagnostic Set

## ADC® COMPLETE 2.5V INSTRUMENT SET ACCESSORIES

- 5215-M3 #3 Laryngeal Mirror
- 5215-M4 #4 Laryngeal Mirror
- 5215-NS Nasal Speculum
- 5215-TD Metal Tongue Depressor
- 5215-TH Illuminated Tongue Blade Holder
- 5215-TO Tongue Depressor Attachment For Disposable
- 5215-ZC Black Zipper Case
- 5220 Otoscope Head Only
- 5240 Ophthalmoscope Head Only
- 5260 Battery Handle
- 5312 Dermoscope



**ADC® POCKET OTOSCOPE/OPHTHALMOSCOPE SET**

Ophthalmoscope features five aperture selections for general and specialist use; 17 lenses from -20 to +20 diopters; illuminated dial and built-in lens dust cover. Otoscope features fiberoptics for cool, obstruction free light transmission; wide angle, removable lens and disposable polypropylene ear specula. Both units feature a 2.5V halogen light source, interchangeable metal handles and a lifetime quality performance guarantee on optics. Includes fitted case. Latex free system.

- 5110 Pocket Otoscope/Ophthalmoscope Set
- 5111 Pocket Otoscope Set
- 5112 Pocket Ophthalmoscope Set

**ADC® POCKET OTOSCOPE/OPHTHALMOSCOPE SET ACCESSORIES**

- 5120 Otoscope Head Only For 5111
- 5121 Insufflator For Otoscope
- 5122 Insufflator Bulb For Otoscope
- 5140 Ophthalmoscope Head Only For 5112
- 5160 Battery Handle For Pocket Instruments
- 925110-00 Hard Case Only For 5110
- 925111-00 Hard Case Only For 5111
- 925112-00 Hard Case Only For 5112
- 5182 2.75mm Disposable Specula, 850/pkg
- 5185 4.25mm Disposable Specula, 850/pkg

**ADC® STANDARD 2.5V DIAGNOSTIC SET**

Ophthalmoscope head with lens wheel features 25 corrective lenses from -20 to +20 diopters. Otoscope head with removable 4x magnifying lens and insufflator fitting. Bayonet head locking mounts for secure connection to the battery handle. C cell powered battery handle with rheostatic on/off switch. Auto-clavable, polypropylene ear specula in 2.5mm, 3.5mm, and 4.5mm sizes. Complete in zippered case with fitted foam insert that is scuff and scratch resistant and cleans easily. Batteries not included.

- 5210 Standard Diagnostic Set



72600



92821

**WELCH ALLYN 2.5V HALOGEN POCKETSCOPE DIAGNOSTIC SETS**

Complete PocketScope Set with ophthalmoscope, otoscope/throat illuminator, two rechargeable battery handles and desk charger (110-130V, 50/60 Hz) in hard case.

- 92820 PocketScope Set Includes Ophthalmoscope, Otoscope/Throat Illuminator, 2 "AA" Alkaline Battery Handles, Hard Case (w/o charger)
- 92821 Same as No. 92820 in Soft Case

**Accessories**

- 72600 2.5V Rechargeable Battery
- 72801 PocketScope Set Handle w/o batteries
- 72800 PocketScope Set Handle w/rechargeable battery

**WELCH ALLYN 3.5V COAXIAL MACROVIEW OTOSCOPE/OPHTHALMOSCOPE SETS**

- 97200-M 23820 Otoscope, 11710 Ophthalmoscope, 71000 Handle and Hard Case
- 97200-MC 11720 Ophthalmoscope, 23820 Otoscope w/hard case and Convert. Handle



05215-MBX



97120



97220

**WELCH ALLYN 3.5V HALOGEN DIAGNOSTIC SETS**

Portable 3.5-V halogen diagnostic sets feature ophthalmoscopes and otoscopes. Optional nasal illuminator. With 71000 handle.

- 97120 Standard Set, 11710 Ophthalmoscope, 20200 Pneumatic Otoscope w/hard case
- 97121 Standard Set, 11710 Ophthalmoscope, 20200 Otoscope w/soft case
- 97170 Standard Set, 11710 Ophthalmoscope, 21700 Otoscope w/hard case
- 97171 Standard Set, 11710 Ophthalmoscope, 21700 Otoscope w/soft case
- 97220 Coaxial Set, 11720 Ophthalmoscope, 20200 Otoscope w/hard case
- 97221 Coaxial Set, 11720 Ophthalmoscope, 20200 Otoscope w/soft case
- 97270 Coaxial Set, 11720 Ophthalmoscope, 21700 Otoscope w/hard case
- 97271 Coaxial Set, 11720 Ophthalmoscope, 21700 Otoscope w/soft case
- 97320 AutoStep Set, 11730 Ophthalmoscope, 20200 Otoscope w/hard case
- 97321 AutoStep Set, 11730 Ophthalmoscope, 20200 Otoscope w/soft case
- 97370 AutoStep Set, 11730 Ophthalmoscope, 21700 Otoscope w/hard case
- 97371 AutoStep Set, 11730 Ophthalmoscope, 21700 Otoscope w/soft case

**Accessories**

- 05215-MBX Soft Case



**WELCH ALLYN 3.5V INSTRUMENT SETS**

- 93100 11710 Ophthalmoscope, 20000 Otoscope, 28100 Tongue Blade Holder, 27000 Larynx Illuminator, 43300 Transilluminator, 71000 Handle, 05510 Case
- 93120 11710 Ophthalmoscope, 20200 Otoscope, 28100 Tongue Blade Holder, 27000 Larynx Illuminator, 43300 Transilluminator, 71000 Handle, 05510 Case
- 93150 11710 Ophthalmoscope, 25020 Otoscope, 28100 Tongue Blade Holder, 27000 Larynx Illuminator, 43300 Transilluminator, 71000 Handle, 05520 Case
- 93170 11710 Ophthalmoscope, 21700 Otoscope, 28100 Tongue Blade Holder, 27000 Larynx Illuminator, 43300 Transilluminator, 71000 Handle, 05510 Case
- 93200 11720 Ophthalmoscope, 20000 Otoscope, 28100 Tongue Blade Holder, 27000 Larynx Illuminator, 43300 Transilluminator, 71000 Handle, 05510 Case
- 93220 11720 Ophthalmoscope, 20200 Otoscope, 28100 Tongue Blade Holder, 27000 Larynx Illuminator, 43300 Transilluminator, 71000 Handle, 05510 Case
- 93250 11720 Ophthalmoscope, 25020 Otoscope, 28100 Tongue Blade Holder, 27000 Larynx Illuminator, 43300 Transilluminator, 71000 Handle, 05520 Case
- 93270 11720 Ophthalmoscope, 21700 Otoscope, 28100 Tongue Blade Holder, 27000 Larynx Illuminator, 43300 Transilluminator, 71000 Handle, 05510 Case

**Accessories**

05510-U Hard Case



**WELCH ALLYN 719 SERIES POWER HANDLE DIAGNOSTIC SETS**

Our coaxial ophthalmoscope uses a halogen HPX™ lamp that provides 30% more light output for true tissue color and consistent, long-lasting illumination. It has 28 lenses (except the Prestige Set has 68 lenses) from -25 to +40 diopters for better resolution. The coaxial optical system produces a shadow-free spot and easier entry into undilated pupils. The crossed-linear polarizing filter reduces corneal reflection. The red-free filter may be used with any aperture. It has 18 possible aperture combinations for greater versatility. The otoscope uses a halogen HPX™ lamp and fiber optics for cool light with no reflections, no obstruction. It has a wide-angle viewing lens for instrumentation under magnification. The sturdy Cycology® resin handle is half the weight of current 3.5v handles and provides twice the on-time.

- 97310-SM Prestige Plus Set, Otoscope w/Throat Illuminator & Nasal Illuminator, Coaxial-Plus Ophthalmoscope, Lithium Ion Handle, Specula & Hard Case



**WELCH ALLYN COMPACSET™ PORTABLE DIAGNOSTIC SET**

Integrates both a full-featured ophthalmoscope and fiber-optic otoscope into one durable, easy-to-use instrument. Small enough to fit in any pocket, weighs just 8.5 ounces. Fully protected by a durable ABS plastic housing. Halogen illumination provides brighter, whiter, longer-lasting light. Instrument heads swivel to allow access from any position.

- 92000 2.5V CompactSet w/AA batteries
- 92100 2.5V CompactSet w/rechargeable batteries

**Accessories**


- 72610 2.5V Rechargeable Battery



**WELCH ALLYN 2.5V/3.5V DESK CHARGERS WITH RECHARGEABLE HANDLES**

Compact units provide constant, automatic recharging. Can not be overcharged. Heads must be ordered separately.

- 71500 2.5V Rechargeable Handle
- 71670 Rechargeable 3.5V Handle Only



**WELCH ALLYN 719 SERIES POWER HANDLE**

Half the weight of current 3.5v handles. Average 110 minutes of on-time between charging. Smallest diameter of any 3.5v power handle. Made from Cycology® resin from GE plastics. Low-battery indicator. Lithium ion battery technology.

- 71900 719 Series Power Handle





76520



76710



76720



76730



79520

**WELCH ALLYN 767 WALL TRANSFORMERS**

Puts everything needed for fast, accurate diagnoses at arm's length. Convenient, modular package. Handle-based rheostats for quick light-intensity adjustments. Comfortable, contoured handles are easy to hold and maneuver. Modular coiled cords are guaranteed for 10 years. High-impact ABS housing is durable and easy to clean. Improved cradle design featuring OptiSense™, a reliable optical sensor that automatically turns handles on when they are removed from the cradles. New Halogen HPX™ Lamp provides 30% more light output. Optional digital clock and pulse timer. KleenSpec® Plus Dispenser with a large, 200-specula capacity and integral storage area. Optional third handle module makes expansion a snap.

|           |  |
|-----------|--|
| 76510     | 76710 Wall Transformer & 52401 KleenSpec Plus Dispenser  |
| 76520     | 76720 Wall Transformer & 52401 KleenSpec Plus Dispenser  |
| 76710     | 3.5V Wall Transformer w/o clock  |
| 76710-L   | 76710 Wall Transformer & 76752 SureTemp Thermometer w/back light   |
| 76710-LA  | 76710 Wall Transformer, 76752 SureTemp Thermometer w/back light & 767 Wall Aneroid                                 |
| 76710-LAK | 76710 Wall Transformer, 76752 SureTemp Thermometer w/back light, 767 Wall Aneroid & 52401 KleenSpec Plus Dispenser |
| 76710-LAP | 76710 Wall Transformer, 76752 SureTemp Thermometer w/back light, 767 Wall Aneroid & 52101 KleenSpec Plus Dispenser |
| 76710-LK  | 76710 Wall Transformer, 76752 SureTemp Thermometer w/back light & 52401 KleenSpec Plus Dispenser                   |
| 76710-LP  | 76710 Wall Transformer, 76752 SureTemp Thermometer w/back light & 52101 KleenSpec Plus Dispenser                   |
| 76710-T   | 76710 Wall Transformer & 76751 SureTemp Thermometer  |
| 76710-TA  | 76710 Wall Transformer, 76751 SureTemp Thermometer & 767 Wall Aneroid  |
| 76710-TAK | 76710 Wall Transformer, 76751 SureTemp Thermometer, 767 Wall Aneroid & 52401 KleenSpec Plus Dispenser              |
| 76710-TAP | 76710 Wall Transformer, 76751 SureTemp Thermometer, 767 Wall Aneroid & 52101 KleenSpec Plus Dispenser              |

**WELCH ALLYN 767 WALL TRANSFORMERS (Continued)**

|           |  |
|-----------|--|
| 76710-TK  | 76710 Wall Transformer, 76751 SureTemp Thermometer & 52401 KleenSpec Plus Dispenser                                |
| 76710-TP  | 76710 Wall Transformer, 76751 SureTemp Thermometer & 52101 KleenSpec Plus Dispenser                                |
| 76720     | 3.5V Wall Transformer w/clock  |
| 76720-L   | 76720 Wall Transformer & 76752 SureTemp Thermometer w/back light   |
| 76720-LA  | 76720 Wall Transformer, 76752 SureTemp Thermometer w/back light & 767 Wall Aneroid                                 |
| 76720-LAK | 76720 Wall Transformer, 76752 SureTemp Thermometer w/back light, 767 Wall Aneroid & 52401 KleenSpec Plus Dispenser |
| 76720-LAP | 76720 Wall Transformer, 76752 SureTemp Thermometer w/back light, 767 Wall Aneroid & 52101 KleenSpec Plus Dispenser |
| 76720-LK  | 76720 Wall Transformer, 76752 SureTemp Thermometer w/back light & 52401 KleenSpec Plus Dispenser                   |
| 76720-LP  | 76720 Wall Transformer, 76752 SureTemp Thermometer w/back light & 52101 KleenSpec Plus Dispenser                   |
| 76720-T   | 76720 Wall Transformer & 76751 SureTemp Thermometer  |
| 76720-TA  | 76720 Wall Transformer, 76751 SureTemp Thermometer & 767 Wall Aneroid  |
| 76720-TAK | 76720 Wall Transformer, 76751 SureTemp Thermometer, 767 Wall Aneroid & 52401 KleenSpec Plus Dispenser              |
| 76720-TAP | 76720 Wall Transformer, 76751 SureTemp Thermometer, 767 Wall Aneroid & 52101 KleenSpec Plus Dispenser              |
| 76720-TK  | 76720 Wall Transformer, 76751 SureTemp Thermometer & 52401 KleenSpec Plus Dispenser                                |
| 76720-TP  | 76720 Wall Transformer, 76751 SureTemp Thermometer & 52101 KleenSpec Plus Dispenser                                |
| 76730     | Third Handle Module  |
| 76910     | 76710 Wall Transformer & 767 Wall Aneroid  |
| 76920     | 76720 Wall Transformer & 767 Wall Aneroid  |
| 79510     | 76710 Wall Transformer, 52401 KleenSpec Plus Dispenser & 767 Wall Aneroid  |
| 79520     | 76720 Wall Transformer, 52401 KleenSpec Plus Dispenser & 767 Wall Aneroid  |



70000



70500



70700



71000-A



71050




71000-C

**WELCH ALLYN HANDLE FOR 3.5-V HALOGEN INSTRUMENTS**

Rechargeable, high-performance handles provide constant 3.5V output. Use rechargeable nickel-cadmium batteries. Heads must be ordered separately.

|         |   |
|---------|---|
| 70000   | Standard Battery Handle "D" Size                          |
| 70500   | Standard Battery Handle "C" Size                          |
| 70700   | 2.5-V Handle  |
| 71000-A | 3.5V Handle   |
| 71050   | Extra Bottom Section Only For 71000, 110-130V, 50/60 Hz   |
| 71050-C | Extra Bottom Section Only For 71000-C, 110-130V, 50/60 Hz |
| 71000-C | 3.5-V Handle, Can Convert to C Battery Use in Emergencies |



**WELCH ALLYN REPLACEMENT RECHARGEABLE BATTERIES**

Nickel-cadmium rechargeable batteries for Welch Allyn handles. Color coded by voltage.

- 72000 2.5V Rechargeable Batteries For Nos. 70700, 70710, 70715, 70720, 70750, 70751, 70754, 70755, 71500 & Discontinued 60700 Handles, Color-Coded Red
- 72100 2.5V Rechargeable Batteries For Discontinued Nos. 617, 715 or 715 Handles, Color-Coded Green
- 72200 3.5V Rechargeable Batteries For Nos. 71000, 71010, 71015, 71020, 71050, 71051, 71054, 71055 & 71670 Handles, Color-Coded Black

**WELCH ALLYN UNIVERSAL CHARGER**

Accepts all Welch Allyn 2.5V, 3.5V, and PocketScope™ (requires adapter) well rechargeable handles. Has a green LED indicating power on. Has a yellow LED under each well indicating the handles are seated and charging. Sports an updated design to fit with our old and new products. Meets all new Agency standards.

- 71630 Universal Charger w/o instrument heads, Two 71670 NiCad Handles
- 71640 Universal Charger w/o instrument heads, Two 71900 L-ion Handles
- 71641-M Universal Charger, (2) 71670 NiCad Handles, 11720 Ophthalmoscope & 23810 MacroView Otoscope
- 71641-MS Universal Charger, (2) 71900 L-ion Handles, 11720 Ophthalmoscope & 23810 MacroView Otoscope
- 71811-MP Universal Charger, (2) 71670 NiCad Handles, 11810 Ophthalmoscope & 23810 MacroView Otoscope
- 71811-MPS Universal Charger, (2) 71900 L-ion Handles, 11810 Ophthalmoscope & 23810 MacroView Otoscope


**Accessories**

- 71140 Universal Charger Only
- 71249 Pocket Scope Adapter Sleeve Only
- 71310 Specula Tray
- 71420 Wall Mount Kit

**ADC® DIAGNOSTIC REPLACEMENT LAMPS**

Manufactured to provide extra long life and reliable illumination. Unconditionally guaranteed to be free from defects and give trouble free service. Rated for 20 hours continuous use. Packaged individually in ziplock poly bags.

- 4500-1 Lamp, WA#001100
- 4501-1 Lamp, WA#00400
- 4502-1 Lamp, WA#00200
- 4504-1 Lamp, WA#00800
- 4507-1 Lamp, WA#01300
- 4513-1 Lamp, WA#03000
- 4514-1 Lamp, WA#03100
- 4516-1 Lamp, WA#03400
- 4518-1 Lamp, WA#03900
- 4522-1 Lamp, WA#04900
- 4523-1 Lamp, WA#06000
- 4530-1 Lamp, 2.5V Halogen For 5120
- 4531-1 Lamp, 2.5V Halogen For 5121
- 4532-1 Lamp, 2.5V Halogen For 4065F to 4067F




**WELCH ALLYN REPLACEMENT LAMPS**

- 00200-U Vacuum Replacement Lamp For 41000
- 00900-U 2.5V Vacuum Lamp for OBS17600
- 01300-U 2.5V Vacuum Replacement Lamp
- 01400-U 2.5V Vacuum Replacement Lamp For 77800
- 02000-U 7V Vacuum Replacement Lamp For 45003
- 02500-U Replacement Lamp
- 02600-U Replacement Lamp
- 02900-U 2.5V Halogen Replacement Lamp For 290S
- 03000-U Halogen Replacement Lamp
- 03100-U Halogen 3.5V Replacement Lamp
- 03300-U 2.5V Halogen Replacement Lamp
- 03400-U 2.5V Halogen Replacement Lamp
- 03600-U Vacuum Replacement Lamp For 60500
- 03700-U 3.5V Halogen Replacement Lamp For 18100
- 03900-U Halogen Replacement Lamp For 12810
- 04050-U 150-W Halogen Replacement Lamp
- 04100-U Halogen Replacement Lamp
- 04200-U Halogen Replacement Lamp
- 04400-U Halogen Replacement Lamp For 11470
- 04430-U Halogen Replacement Lamp For LS135 Light
- 04450-U Halogen Replacement Bulb
- 04800-U 2.5V Vacuum Lamp For Larynscope
- 04900-U 3.5V Halogen Replacement Lamp
- 06000-U Halogen Replacement Lamp
- 06200-U Halogen Replacement Lamp
- 06300-U Halogen Replacement Lamp
- 06400-U Halogen Replacement Lamp
- 07200-U Halogen Replacement Lamp For 17200
- 07600-U Halogen Replacement Lamp For 76600
- 07800-U 6V Halogen Replacement Lamp
- 08200-U Halogen Replacement Lamp For 18200
- 08300-U Halogen Replacement Lamp For 18300
- 08800-U 4.6V Halogen Replacement Lamp For 78810 Illumination System
- DA-ARC Lamp Replacement For DLS100
- DA-L400 Lamp Replacement For DLS400 & DCU2000

**ADC® FIBEROPTIC LARYNGOSCOPE BATTERY HANDLES**

Lightsource in handle simplifies blade cleaning, allows for use of brighter, whiter light. Satin chrome plated brass handle with a knurled finish for a positive grip. Fully compatible with all standard hook-on blades that conform to ASTM (F965) and ISO (7376-1). Color coded green to signify compatibility ONLY with other fiberoptic handles. Batteries not included.

- 4065F Laryngoscope Fiberoptic Battery Handle, Medium "C" Size
- 4066F Laryngoscope Fiberoptic Battery Handle, Small "AA" Size
- 4067F Laryngoscope Fiberoptic Battery Handle, Stubby "AA" Size



**ADC® FIBEROPTIC LARYNGOSCOPE BLADES**

Fiberoptic light transmission provides cool obstruction free viewing. 2.5V halogen illumination for true tissue color and consistent long lasting lamp performance. Stainless steel. Satin, non-glare finish. Lightweight design. Fully compatible with all fiberoptic hook on handles that conform to ASTM (F965) and ISO (7376-1). Color coded green to signify compatibility and differentiate from standard hook on models.

- 4070F Size 0 Macintosh Blade, Premie
- 4071F Size 1 Macintosh Blade, Infant
- 4072F Size 2 Macintosh Blade, Child
- 4073F Size 3 Macintosh Blade, Medium Adult
- 4074F Size 4 Macintosh Blade, Large Adult
- 4080F Size 0 Miller Blade, Premie
- 4081F Size 1 Miller Blade, Infant
- 4082F Size 2 Miller Blade, Child
- 4083F Size 3 Miller Blade, Medium Adult
- 4084F Size 4 Miller Blade, Large Adult

**ADC® LARYNGOSCOPE SET**

Choose from either a 4 blade Macintosh or 5 blade Miller set with both AA and C battery handles. Foam lined fitted hard case with carry handle. Batteries included.

- 4079 Standard Macintosh Laryngoscope Set
- 4079F Fiberoptic Macintosh Laryngoscope Set
- 4089 Standard Miller Laryngoscope Set
- 4089F Fiberoptic Miller Laryngoscope Set



**ADC® STANDARD LARYNGOSCOPE BATTERY HANDLES**

Satin chrome plated brass with a knurled finish for a positive grip. Fully compatible with all standard hook-on blades that conform to ASTM (F965) and ISO (7376-1). Batteries not included.

- 4065 Laryngoscope Battery Handle, Medium "C" Size
- 4066 Laryngoscope Battery Handle, Small "AA" Size
- 4067 Laryngoscope Battery Handle, Stubby "AA" Size

**ADC® STANDARD LARYNGOSCOPE BLADES**

Stainless steel construction with satin, non-glare finish. Brass electrical contacts enhance conductivity. Lightweight design. Fully compatible with all standard hook on handles that conform to ASTM (F965) and ISO (7376-1).

- 4070 Size 0 Macintosh Blade, Premie
- 4071 Size 1 Macintosh Blade, Infant
- 4072 Size 2 Macintosh Blade, Child
- 4073 Size 3 Macintosh Blade, Medium Adult
- 4074 Size 4 Macintosh Blade, Large Adult
- 4080 Size 0 Miller Blade, Premie
- 4080-0 Size 00 Miller Blade, Small Premie
- 4081 Size 1 Miller Blade, Infant
- 4082 Size 2 Miller Blade, Child
- 4083 Size 3 Miller Blade, Medium Adult
- 4084 Size 4 Miller Blade, Large Adult
- 4090 Size 0 Wisconsin Blade, Premie
- 4091 Size 1 Wisconsin Blade, Infant
- 40915 Size 1.5 Wisconsin Blade, Small Child
- 4092 Size 2 Wisconsin Blade, Child
- 4094 Size 4 Wisconsin Blade, Large Adult
- 4093 Size 3 Wisconsin Blade, Medium Adult



**WELCH ALLYN KLEENSPEC® DISPOSABLE LARYNGOSCOPES**

Designed for single use. Light source inserts in the handle. After use, remove the light source and discard the blade and handle.

- 56901 Small Laryngoscope, 10/bx
- 56903 Large Laryngoscope, 10/bx

**Accessories**

- 60500 Illuminator For Laryngoscope


**We offer:**

- **Great Service**
- **Great Pricing**
- **Great Products**






68040



68060



68696



69061



69211



69696



69697

## WELCH ALLYN LARYNGOSCOPE BLADES

One-piece stainless steel construction.

|         |  |
|---------|--|
| 05681-U | Fiber Optic - Miller Blade Case  |
| 05691-U | Fiber Optic - MacIntosh Blade Case   |
| 05701-U | Fiber Optic - English MacIntosh Blade Universal Case                                     |
| 05791-U | Fiber Optic - English MacIntosh Blade Standard Case                                      |
| 68040   | Miller Blade, Size 0   |
| 68041   | Miller Blade, Size 1   |
| 68042   | Miller Blade, Size 2   |
| 68043   | Miller Blade, Size 3   |
| 68044   | Miller Blade, Size 4   |
| 68060   | Fiber Optic - Miller Blade, Size 0   |
| 68061   | Fiber Optic - Miller Blade, Size 1   |
| 68062   | Fiber Optic - Miller Blade, Size 2   |
| 68063   | Fiber Optic - Miller Blade, Size 3   |
| 68064   | Fiber Optic - Miller Blade, Size 4   |
| 68470   | Miller Blade, Size 0 (Neonatal)  |
| 68696   | Fiber Optic - Miller Blade, Complete Set w/medium & small handles, 05681 Case            |
| 69041   | MacIntosh Blade, Size 1  |
| 69042   | MacIntosh Blade, Size 2  |
| 69043   | MacIntosh Blade, Size 3  |
| 69044   | MacIntosh Blade, Size 4  |
| 69061   | Fiber Optic - MacIntosh Blade, Size 1  |
| 69062   | Fiber Optic - MacIntosh Blade, Size 2  |
| 69063   | Fiber Optic - MacIntosh Blade, Size 3  |
| 69064   | Fiber Optic - MacIntosh Blade, Size 4  |
| 69211   | Fiber Optic - English MacIntosh Blade, Size 1  |
| 69212   | Fiber Optic - English MacIntosh Blade, Size 2  |
| 69213   | Fiber Optic - English MacIntosh Blade, Size 3  |
| 69214   | Fiber Optic - English MacIntosh Blade, Size 4  |
| 69242   | English MacIntosh Blade, Size 2  |
| 69243   | English MacIntosh Blade, Size 3  |
| 69696   | Fiber Optic - MacIntosh Blade, Complete Set w/medium & small handles, 05691 Case         |
| 69697   | Fiber Optic - English MacIntosh Blade, Complete Set w/medium & small handles, 05791 Case |



60200



60300



60305



60400



60710



60713



60813

## WELCH ALLYN LARYNGOSCOPE HANDLES

Knurled finish for sure grip.

|       |  |
|-------|--|
| 60200 | Standard Laryngoscope Handle D - Size Battery                        |
| 60300 | Standard Laryngoscope Handle C - Size Battery                        |
| 60305 | Standard Stubby Laryngoscope Handle                                  |
| 60400 | Standard Laryngoscope Handle AA - Size Battery                       |
| 60710 | Standard 2.5V Rechargeable Laryngoscope Handle Only                  |
| 60713 | 2.5V Rechargeable Laryngoscope Handle Only For Fiber Optic Blades    |
| 60715 | Standard Rechargeable Battery Handle & Transformer                   |
| 60813 | Laryngoscope Handle For Fiber Optic Blades, Uses C - Size Batteries  |
| 60814 | Laryngoscope Handle For Fiber Optic Blades, Uses AA - Size Batteries |
| 60815 | Stubby Laryngoscope Handle For Fiber Optic Blades                    |

## WELCH ALLYN 3.5V HALOGEN NASAL ILLUMINATOR

When attached to throat illuminator section of No. 20000 otoscope/throat illuminator, provides illumination of the nasal passages through a 9mm speculum.

|       |  |
|-------|--|
| 26530 | 3.5-V Halogen Nasal Illuminator, Complete w/9mm polypropylene speculum |
| 26535 | Nasal Illuminator Section Only   |



**WELCH ALLYN 3.5V HALOGEN TRANSILLUMINATORS**

These instruments use halogen illumination transmitted through fiber optics for brilliant, yet cool illumination. Fits any Welch Allyn 3.5V power source. May be sterilized by steam, gas, or some germicidal solutions.

43300 Curved All Purpose Transilluminator

**Accessories**

41002 Cobalt Blue Filter Only

**WELCH ALLYN 3.5V TONGUE BLADE HOLDER**

New halogen HPX™ lamp that provides 30% more light output, for true tissue color and consistent, long lasting illumination. Fiber optic diffuser for cool distal light. Universal clamp accommodates most tongue blades.

28100 Tongue Blade Holder



**WELCH ALLYN LARYNX & NASOPHARYNX ILLUMINATORS**

Halogen illumination and fiber optics light transmission provide cool, intense light. Wide-angle mirror. Light aperture assists in positioning mirror. Ready for immediate use. Compatible with any Welch Allyn power source.

27010 2.5V or 3.5V Halogen Larynx Mirror Section Only

27050 2.5V or 3.5V Halogen Nasopharynx Illuminator, Complete

**ADC® 2.5V DERMASCOPE ACCESSORIES**

5312-1 Plain Contact Plate

5312-2 Graduated Contact Plate

5312-3 Dermoscope Head

5312-4 Lamp 2.5V Halogen

5312-ZC Black Zipper Case



**WELCH ALLYN EPISCOPE™**

For examining the entire skin surface and making accurate measurements of skin lesions. Bright, 3.5-V halogen illumination and 10x optics system ensure true tissue color and accurate identification of suspicious lesions. Attaches to rechargeable handle. Reticle contact plate allows easy measurement in millimeters. Removable contact plate is autoclavable to eliminate the possibility of cross-contamination. Easy to focus, large knurled focusing ring. Comfortable eyepiece with soft rubber eyewear protector. Fully portable. Soft carrying case accommodates handle and includes pouch for extra front cover. Handles must be ordered separately.

47300 EpiScope Skin Surface Microscope w/reticle contact plate

47351 EpiScope Set (includes EpiScope, Reticle Contact Plate, Rechargeable Handle, Nickel Cadmium Battery & Soft Carrying Case)

**Accessories**

05232-U Soft Carrying Case

47310 Reticle Contact Plate

47320 Plano Contact Plate



65002



65003



65050

**WELCH ALLYN RL-150 RHINOLARYNGOSCOPE**

For upper-respiratory exams. Procedure time is less than 5 minutes and easy to master. High patient acceptance rate. Simple one-step disinfection system. Free lifetime loaner guarantee policy. Short payback period. Reimbursable. Eliminates unnecessary referrals.

65050 RL-150 Rhinolaryngoscope

**Accessories**

65002 Disinfection Tray

65003 Hanger

65006 LX-150 Adaptor For 65051 Rhinolaryngoscope/LX-150 Light Source



AS032L




AS034L

## AMSINO AMSMOOTH™ DISPOSABLE VAGINAL SPECULUM

Made of polystyrene with see through clarity. Provides smooth, one-handed operation for easier insertion, requires no lubricant except water, and allows for secure positioning without slippage. Single use eliminates cross contamination. The self-locking feature speeds examination and saves time.

- AS032L Large Grave Style Vaginal Speculum, 10/bg, 10 bg/cs
- AS032M Medium Grave Style Vaginal Speculum, 10/bg, 10 bg/cs
- AS032S Small Grave Style Vaginal Speculum, 10/bg, 10 bg/cs
- AS034L Disposable Vaginal Specula, Graves Style, Individually Wrapped, Large, w/Light Source Option, 25/bx, 4bx/cs
- AS034M Disposable Vaginal Specula, Graves Style, Individually Wrapped, Medium, w/Light Source Option, 25/bx, 4bx/cs
- AS034S Disposable Vaginal Specula, Graves Style, Individually Wrapped, Small, w/Light Source Option, 25/bx, 4bx/cs



## COOPERSURGICAL SANI-SPEC® SPECULA

Clear plastic material allows 360° visibility and reduces thermal shock. Disposable. Provides protection of the doctor, the nurse and the patient from cross contamination. Plastic material avoids unpleasant temperature shock associated with metal. Smooth pre-lubricated blades with 6 position grooves.

- 82329 Small Specula, 1"W x 4 1/2"L, 10/bg, 10 bg/cs
- 82330 Medium Specula, 1 1/4"W x 4 3/4"L, 10/bg, 10 bg/cs
- 82331 Large Specula, 1 1/3"W x 5 1/3"L, 10/bg, 10 bg/cs
- 82340 Small Specula, Individually Wrapped, 50/cs
- 82341 Medium Specula, Individually Wrapped, 50/cs


## DRG SOFTSPEC VAGINAL SPECULUM

A tampon-like, disposable device that is inserted into the vagina and gently inflated to expand the vaginal walls for guaranteed visualization of the cervix. Allows easy access for an adequate pap smear by inserting the spatula through the donut shape of the SoftSpec. Disposable so you do not have to worry about messy clean up or the spread of bacterial from patient to patient.

- VSI-25-001 SoftSpec Vaginal Speculum Illuminated w/bulb, 25/bx
- VSO-25-001 SoftSpec Vaginal Speculum w/bulb, 25/bx

## COOPER SURGICAL GALENICA GYNE-LITE™ SYSTEM ACCESSORIES

- 9220 Adaptor, 10/bg
- 9230 Sheaths, 100/bx, 4 bx/cs
- 9240 Bulb



## COOPER SURGICAL GALENICA VAGINAL SPECULA

Recommended for routine gynecological exams. Disposable. Allows the examination of the vaginal mucosa at 360° angle without rotation. Wide range of adjustability. Ease of operation and manipulation with one hand. Self-holding in various opening positions. Made from hard, clear and non-toxic plastic resin of high strength. Unit Pack (poly bag) non-sterile. Unit Pack (peeling pouch) sterile.

- 8309 Small Vaginal Specula, Non-Sterile, 10/bg, 10 bg/cs
- 8310 Medium Vaginal Specula, Non-Sterile, 10/bg, 10 bg/cs
- 8311 Large Vaginal Specula, Non-Sterile, 10/bg, 10 bg/cs
- 8615 Small Unit Pack Vaginal Specula, 50/cs
- 8616 Medium Unit Pack Vaginal Specula, 50/cs
- 8617 Large Unit Pack Vaginal Specula, 50/cs

### Accessories

- 9310 Smoke Evacuator For All Specula, 40/cs



**Quality**  
**Value**  
**Availability**





## WELCH ALLYN KLEENSPEC® 580 DISPOSABLE VAGINAL SPECULA WITH BARRIER

Superior halogen illumination for a whiter, brighter light for truer tissue color rendition. Smooth, one-handed angulation and elevation mechanism with positive locking positions. Light handle is completely sealed when inserted into the speculum. Built-in reservoir prevents fluids from contaminating the outside of the speculum handle. Optional disposable sheath covers the exposed sections of the light handle and cord, virtually eliminating the threat of cross-contamination.

58010 Disposable Sheaths, 25/bx, 5 bx/cs



## WELCH ALLYN KLEENSPEC® 580 SERIES DISPOSABLE VAGINAL SPECULA

Halogen light for true tissue color and consistent, long lasting illumination. Disposable for economy, safety and convenience. One-hand operation with angulation and elevation for convenience. Especially narrow blade width available for adolescent and geriatric patients. Individually wrapped specula.

- 58000 Small Specula (Graves Style), 25/bx, 4 bx/cs
- 58000S Small Specula w/sheath, 20/bx, 4 bx/cs
- 58000S-B Small Specula w/sheath, 20/bx
- 58001 Medium Specula (Graves Style), 25/bx, 4 bx/cs
- 58001S Medium Specula w/sheath, 20/bx, 4 bx/cs
- 58001S-B Medium Specula w/sheath, 20/bx
- 58004 Large Specula (Graves Style), 20/bx, 4 bx/cs
- 58004S Large Specula w/sheath, 15/bx, 4 bx/cs
- 58005 Small Specula w/small evacuator tube, 12/bx, 4 bx/cs
- 58006 Medium Specula w/small evacuator tube, 12/bx, 4 bx/cs



## WELCH ALLYN KLEENSPEC® DISPOSABLE VAGINAL SPECULA ILLUMINATION SYSTEM

Permanent light carrier projects intense beam of light through plastic light guide. No external light source necessary.

- 78810 Complete Vaginal Specula Illumination System w/on-off switch & transformer
- Accessories**
- 73305 Transformer w/5 ft cord, 110-130V, 50/60 Hz
  - 73410 Transformer w/5 ft cord, 110-130V, 50/60 Hz, UL Listed, CSA Certified
  - 78800 Permanent 6V Halogen Light Handle w/on-off switch



58600



78600

## WELCH ALLYN KLEENSPEC® VAGINAL SPECULA

Brilliant halogen illumination system allows diagnostic accuracy. Clear, polystyrene construction. Designed for one-hand operation. Single-use prevents cross-contamination. Expendable light carrier for safety. (Illumination system not included.)

- 58600 Small Specula, 25/bx, 4 bx/cs
- 58601 Medium Specula, 25/bx, 4 bx/cs
- 78600 Expendable Illuminator (can be resterilized 50 times before disposal), 6/bx, 6 bx/cs



## SEILER 935 THREE-STEP COLPOSCOPE

- 30935-3 Model 935 Seiler 3-Step Colposcope
- 30935-3-D Model 935 Colposcope w/digital
- 30935-3-V Model 935 Colposcope w/video



## SEILER 955 FIVE-STEP COLPOSCOPE


- 30955-5 Model 955 Colposcope
- 30955-5-D Model 955 Colposcope w/digital
- 30955-5-V Model 955 Colposcope w/video

**Contact your representative  
for information on  
cryosurgical units.**



**SEILER 965Z ZOOM COLPOSCOPE**

- 30965-Z Model 965 Colposcope
- 30965-Z-D Model 965 Colposcope w/digital
- 30965-Z-V Model 965 Colposcope w/video



**WELCH ALLYN IMAGE CAPTURE SYSTEM FOR VIDEO COLPOSCOPE**

A unique software-based documentation system to be used with the Welch Allyn video colposcope. The software can be controlled remotely from the Video Colposcope and enables the user to capture, review, store, print and export images.

- 88801 Image Capture System For Desktop Computer
- 88802 Image Capture System For Laptop Computer



37019



37023



37027

**WELCH ALLYN FIBER OPTIC ANOSCOPES**

Fiber optics provide cool distal ring of light, no reflections, no obstructions. Autoclavable stainless steel. 45° bevel tip (90° available on request). Set includes speculum, obturator, light handle and transformer.

- 37019 19mm x 7cm Speculum Only w/obturator
- 37023 23mm x 7cm Speculum Only w/obturator
- 37027 27mm x 7cm Speculum Only w/obturator



**WELCH ALLYN KLEENSPEC® DISPOSABLE ANOSCOPE**

45° bevel end. Individually wrapped. 19mm x 10cm with obturators.

- 53110 Disposable Anoscope, 25/bx, 4 bx/cs



38108



38614



38900



39614

**WELCH ALLYN STANDARD ANOSCOPES**

Halogen light for true tissue color and consistent, long-lasting illumination. Autoclavable. Comes with obturator.

- 38108 8mm Aperture, Anoscope Only
- 38114 14mm Aperture, Anoscope Only, Standard Speculum, 89mm Length
- 38119 19mm Aperture, Anoscope Only, Standard Speculum, 89mm Length
- 38122 22mm Aperture, Anoscope Only, Standard Speculum, 89mm Length
- 38614 14mm Aperture, Anoscope Only, Operating (Cut-Out) Speculum, 89mm Length
- 38619 19mm Aperture, Anoscope Only, Operating (Cut-Out) Speculum, 89mm Length
- 38622 22mm Aperture, Anoscope Only, Operating (Cut-Out) Speculum, 89mm Length
- 38900 20mm ID, 22mm OD, 68mm Length to Taper, Anoscope Only, Rotating Speculum
- 39614 14mm Aperture, Anoscope Only, Long Speculum, 127mm Length
- 39619 19mm Aperture, Anoscope Only, Long Speculum, 127mm Length
- 39622 22mm Aperture, Anoscope Only, Long Speculum, 127mm Length

**COOPER SURGICAL GALENICA ANOSCOPES**

Easy manipulation for septic anal/rectal examinations. Decreases the risks of cross-infections. 56° bevel at the insertion end. Adaptable to fiber optic lighthouse. Made from hard non-toxic plastic resin of fine medical grade. Disposable and lightweight.

8607 Anoscope, 25/box, 4 bx/cs



32010



32410



32820



35303

**WELCH ALLYN FIBER OPTIC SIGMOIDOSCOPES**

Halogen light. Fiber optics provide cool distal ring of light, no reflections, no obstructions. Autoclavable stainless steel. Hinged window with neoprene seal. Graduated tube with removable tip. Complete set includes: sigmoidoscope speculum, obturator, light handle, insufflation bulb & transformer

- 32010 11mm x 10cm Sigmoidoscope Speculum Only w/obturator
- 32020 11mm x 25cm Sigmoidoscope Speculum Only w/obturator
- 32410 15mm x 15cm Sigmoidoscope Speculum Only w/obturator
- 32420 15mm x 25cm Sigmoidoscope Speculum Only w/obturator
- 32810 19mm x 15cm Sigmoidoscope Speculum Only w/obturator
- 32820 19mm x 25cm Sigmoidoscope Speculum Only w/obturator
- 32830 19mm x 35cm Sigmoidoscope Speculum Only w/obturator
- 33220 23mm x 25cm Sigmoidoscope Speculum Only w/obturator
- 33830 19mm x 30cm Sigmoidoscope Speculum Only w/obturator
- 35303 Set Includes: 32410 Sigmoidoscope, 32820 Sigmoidoscope, 37023 Anoscope, 73210 Light Handle w/cord, 73305 Transformer, 30200 Insufflation Bulb Complete, 07800 Extra Halogen Lamps (2), 30130 Suction Tube, 05391 Case

**Accessories**

05391-U Hard Case


**WELCH ALLYN FX-100 FLEXIBLE FIBER OPTIC SIGMOIDOSCOPE ACCESSORIES**

- 33905 Water Bottle Assembly
- 33950 Biopsy Forceps

**WELCH ALLYN KLENSPEC® FIBER OPTIC DISPOSABLE SIGMOIDOSCOPE**

Speculum attaches to fiber optic light head; no delicate or protruding light carriers to obstruct vision. Brilliant distal illumination is shadow free, without color distortion. For use with any Welch Allyn 3.5V or 2.5V power source and handle adapter or 6.5V handle and transformer.


36103 Complete Illumination System includes Fiber Optic Light Head, Light Handle, Transformer & Insufflation Bulb



30130



30200




34000



34100



36019



38300



53130



73210



73550

**WELCH ALLYN KLENSPEC® FIBER OPTIC DISPOSABLE SIGMOIDOSCOPE ACCESSORIES**

- 30130 30cm Suction Tube
- 30140 40cm Suction Tube
- 30200 Insufflation Bulb, Bladder Net & Tubing
- 34000 Rectal Probe, 7 1/2"
- 34100 Rectal Hook, 7 1/2"
- 36019 Fiber Optic Light Head
- 38300 Incandescent Light Carrier
- 38700 Fiber Optic Light Carrier
- 53130 Disposable Specula (19 mm x 25 cm) w/obturator, 25/box, 4 bx/cs
- 73210 Light Handle w/3 ft cord
- 73500 Halogen 3.5-V Handle Adapter
- 73550 Halogen 2.5-V Handle Adaptor





23200



23210


**WELCH ALLYN AM 232™ MANUAL AUDIOMETER**

Ergonomic panel design and convenient controls allow one-handed operation of this easily portable audiometer. Frequency values and HL levels are easily visible at a glance. The AM 232 also offers selection modes for steady, pulsed or frequency-modulated stimulus for increased reliability of response. +10 dB push-button prevents mistaken presentation of uncomfortable frequency tones. Provides ample storage space inside hard-shell case for standard headphones and optional accessories.

- 23200 AM 232 Audiometer w/AC Only
- 23210 AM 232 Audiometer w/AC Power & Rechargeable Battery




23220



52300

**WELCH ALLYN AM 232™ MANUAL AUDIOMETER ACCESSORIES**

- 23220 Response Hand Switch
- 23221 Patch Cord
- 23222 Audiocups
- 52300 Audiogram Form, 50 sheets/pad, 10 pads/cs
- 72320 Replacement Rechargeable Battery Pack Only
- 72329 Complete Battery Pack Assembly w/rechargeable batteries



23300



92680

**WELCH ALLYN AUDIOSCOPE® 3 SCREENING AUDIOMETER/OTOSCOPE**

Versatile instrument for screening hearing loss combines 20, 25 and 40 dB HL levels in one unit. Tones are presented in random intervals so patients can not second guess their hearing screen. This varied timing means improved results. Built-in 1000 Hz pretone, presented at 20 dB HL above the screening level, allows patients to practice listening before they're actually screened. After pretone, patient is given the easiest tones for enhanced reliability of responses. 110-130V.

- 23300 AudioScope 3 w/AudioSpecs (no transformer)
- 92600 AudioScope 3 Set (AudioScope 3, set of AudioSpecs, charging transformer, recording forms)
- 92680 Complete Set (AudioScope 3, set of AudioSpecs, charging transformer, charging stand, recording forms)



23400



71040



72300

**WELCH ALLYN AUDIOSCOPE® 3 SCREENING AUDIOMETER/OTOSCOPE ACCESSORIES**

- 23303 Small AudioSpec
- 23305 Medium AudioSpec
- 23307 Large AudioSpec
- 23400 Set of AudioSpecs (small, medium, large)
- 55230 Recording Forms, 100/pk, 10 pk/bx
- 71040 Charging Transformer
- 71123 Charging Stand Only
- 72300 3.5V Nickel-Cadmium Rechargeable Replacement Battery
- 94005-0000 Paper, 5.3" x 5.8", 140 shts/pk, 12 pk/cs



**WELCH ALLYN EAR WASH SYSTEM PLUS**

Easy to use by healthcare professionals and patient satisfaction. Chamber regulates the output pressure to prevent damage to the tympanic membrane. Handle indicator will contrast in color when water is within 10 degrees of body temperature. Chamber automatic shut-off will shut off system if water reaches 110 degrees. Multi-functional disposable tip. Device can be moved easily from location to location via a storage tray. One-year warranty.

- 29350 Ear Wash System Plus, English



23640



93690

**WELCH ALLYN MICRO TYMP® 2 PORTABLE TYMPANOMETRIC INSTRUMENT**

Portable, hand-held unit detects otitis media and other middle ear disorders. Measures sound returned from eardrum under positive, negative, and atmospheric pressure to determine status of middle ear. Accurately identifies many middle ear disorders in asymptomatic cases. Printer/charger keeps Micro Tym fully charged and provides hardcopy documentation. Results instantly displayed on LCD before printing. 226 Hz frequency, 85 db HL intensity.

- 23640 Micro Tym II (includes 2 sets of tips)
- 93690 Micro Tym II Set (includes printer/charger, extra paper & tips)



23630



71170

**WELCH ALLYN MICRO TYMP® 2 PORTABLE TYMPANOMETRIC INSTRUMENT ACCESSORIES**

- 05276-U Carrying Case
- 23630 Set of Micro Tympanometry II Tips (4 sizes)
- 23639 Multi-Pack of Micro Tympanometry 2 Probe Tips
- 71170 Printer/Charger (includes paper)
- 72900 2.4V Nickel Cadmium Rechargeable Battery



**WELCH ALLYN TM 262™ AUTO TYMP™**


Offers tympanometry and ipsilateral acoustic reflex testing (with optional manual audiometry). Built-in printer and display. Improved ear seal via lightweight handheld probe. Delivers precision results with each use, and stores up to eight test results.

- 26200 Auto Tympanometry w/Ipsi Reflex
- 26230 Auto Tympanometry w/Audiometer




**WELCH ALLYN TM 262™ AUTO TYMP™ ACCESSORIES**


- 05260-U Carrying Case
- 26100 Eartips (probe), 6 sizes
- 26240 Dust Cover
- 26241 Test Cavity
- 52600 Printer Paper, 5/bx



3050



3051



3052

**TECH-MED EYE CHARTS**

- 3050 Snellen Eye Test Chart, 20 ft, Non-Reflective Matte Finish, 22" x 11"
- 3051 Illiterate Eye Test Chart, 20 ft, Non-Reflective Matte Finish, 22" x 11"
- 3052 Kindergarten Eye Chart 20 ft, Non-Reflective Matte Finish, 22" x 11"
- 3053 Pocket Eye Chart, Use at 14", Provides 20/800 Distance, Laminated Plastic, 6 1/2" x 3 1/2"
- 3060 Illuminated Eye Test Cabinet, Baked Enamel Finish w/handle & power cord, Includes 1 each Illuminated Eye Charts
- 3061 Illuminated Snellen Eye Test Chart, 20 ft
- 3062 Illuminated Snellen Eye Test Chart, 10 ft
- 3063 Illuminated Kindergarten Eye Test Chart, 20 ft
- 3064 Illuminated Illiterate Eye Test Chart, 20 ft
- 3065 Illuminated Illiterate Eye Test Chart, 10 ft

**GRAHAM-FIELD GRAFCO® ILLITERATE PLASTIC EYE CHART**

With green and red color bars. 20' distance.

- 1241 Illiterate Eye Chart, Plastic, 22" x 11"

**GRAHAM-FIELD GRAFCO® KINDERGARTEN PLASTIC EYE CHART**

Shows easy to distinguish objects. 20' distance.

- 1243 Kindergarten Eye Chart, Plastic, 22" x 11"

**GRAHAM-FIELD GRAFCO® SNELLEN TYPE PLASTIC EYE CHART**


Non-reflective, matte finish with green and red color bars. 20' distance.

- 1240 Snellen Eye Chart, Plastic, 22" x 11"
- 1264 Snellen Eye Chart, 10' Distance, Translucent for Use w/Illuminated

**GRAHAM-FIELD GRAFCO® ISHIHARA® TEST CHART BOOKS, FOR COLOR DEFICIENCY**

Accepted by leading authorities as a simple and accurate method for discovering congenital color blindness and red-green blindness, each in two forms, complete and incomplete. Includes four special plates for tests to determine the kind and degree of defect in color vision. Makes use of the peculiarity that in red-green blindness, blue and yellow appear remarkably bright compared with red and green. Color plates encased in specialty designed album-type books for ease of handling. Detailed instructions included.

- 1254 Ishihara Test Chart Book For Color Deficiency, 14 Plate Book
- 1255 Ishihara Test Chart Book For Color Deficiency, 24 Plate Book
- 1256 Ishihara Test Chart Book For Color Deficiency, 38 Plate Book
- 1257 Ishihara Test Chart Book For Color Deficiency, 10 Plate Book
- 1263 Kindergarten Eye Chart, Translucent for Use w/Illuminated



**STEREO OPTICAL OPTEC® VISION TESTER ACCESSORIES**

- 1000-21 Secondary Test Slide For Optec 1000
- 1000-22 Elementary Test Slide For Optec 1000
- 1000-23 Multiple Acuity Test Slide For Optec 1000
- 1000-NPS Near Point Lens For Optec 1000
- 2000-220 Purdue Job Standards Book w/recording forms (5 pads)
- 2000-221 Headrest Tissues (5000 tissues)
- 2000-226 Bulb (6) Color Corrected (120v)
- 2000-230 Plus Lens (+1.75D) For Optec 1000 & 2000
- 2000-231 Plus Lens (+2.25D) For Optec 2000
- 2000-232 Intermediate Lenses for VDT Operator Testing (set of 5)
- 2000-233 Intermediate Lens, 39 1/2", For Optec 2000
- 2000-234 Intermediate Lens, 31 1/2", For Optec 2000
- 2000-235 Intermediate Lens, 26", For Optec 2000
- 2000-236 Intermediate Lens, 22 1/2", For Optec 2000
- 2000-237 Intermediate Lens, 20", For Optec 2000
- 2000-245 Carrying Case For Optec 2000, Hard Sided
- 2000-246 Carrying Case For Optec 1000, Soft Sided
- 2000-247 Carrying Case For Optec 1000, Hard Sided
- 2000-250 Carrying Case For Optec 2000, Soft Sided
- 2000-HO2 H.O.T.V. Slide Set 20/30 (4 slides) w/Binocular Test For Optec 2000
- 2000-HO3 H.O.T.V. Slide Set 20/40 (slides) w/Binocular Test For Optec 2000
- 2000-HO4 H.O.T.V. Slide Set 20/30 (4 slides) Monocular For Optec 2000
- 2000-HO5 H.O.T.V. Slide Set 20/40 (4 slides) Monocular For Optec 2000
- 2500-226 Bulbs (6) For Optec 2500 (24v)



5000



5500

**STEREO OPTICAL OPTEC® VISION TESTERS**

Testers have homogeneous illumination making them the only vision screeners on the market that meet the ANSI standard for illumination in vision. All instruments are micro-processor controlled and continually calibrated. Choose one of the following test packages: 2000MTS - Medical Slide Package with Michigan Preschool; 2000HOTV - Medical Slide Package with H.O.T.V.; 2000ITS - Industrial Slide Package with Job Standards Book; 2000BL - International Slide Package with Job Standards Book; 2000PRS - Ophthalmic Prescreening Slide Package (FAA); 2000RHB - Driver Rehabilitation Slide Package. Accessories include: dust cover, pointer, lens cleaner, headrest tissues, recording forms, and training manual.

- 1000P System w/Peripheral Test & Accessory Kit, Excludes Video, Includes 1000S Test Package - Multiple Acuity Slide, Elementary Test Slide & Nearpoint Lens
- 5000 System w/Manual Control & Accessory Kit, w/o Peripheral Test, Choose a Standard Test Slide Package
- 5000P System w/Manual Control, Peripheral Test & Accessory Kit, Choose a Standard Test Slide Package
- 5500 System w/Remote Control, Peripheral Test & Accessory Kit, Choose a Standard Test Slide Package
- 5500P System w/Peripheral Test, Remote Control & Accessories, Choose a Standard Test Slide Package
- 6500P System w/Microprocessor Controlled Glare & Luminance Levels Peripheral Test, Remote Control, Accessory Kit (Choose a 5000 Standard Slide Pkg or a 6500 Test Pkg)

**STEREO OPTICAL C.A.T.® CONTRAST SENSITIVITY SLIDE PACKAGE**

Package includes: (2000-007) NEAR Point Acuity Right Eye & Left Eye Monocular & Both Eyes Simultaneously (20/100 to 20/20), (2000-003) FAR Point Acuity Right Eye & Left Eye Monocular & Both Eyes Simultaneously (20/100 to 20/20), (2000-037) 100% FAR C.A.T. Contrast Sensitivity Monocular, (2000-038) 50% Contrast Sensitivity Monocular, (2000-039) 25% Contrast Sensitivity Monocular, (2000-040) 12.5% Contrast Sensitivity Monocular, (2000-041) 6.25% Contrast Sensitivity Monocular, (2000-010) FAR Point Color Perception (Pseudo Ishihara), (2000-012) FAR Point Lateral Phoria (1 Diopter Increments), (2000-025) FAR Point Vertical Phoria (1/2 Diopter Increments), (2000-024) FAR Point Stereo Depth Perception (400 to 20 Seconds of Arc), (2000-087) FAR Point Fusion.

2000HS2 C.A.T. Contrast Sensitivity Slide Package

**STEREO OPTICAL DRIVER REHABILITATION SLIDE PACKAGE**

For complete vision screening needs. Ideal for school, athletic, employment physicals, driver licensing exams, D.O.T., F.A.A. and general check-ups. Package includes: (2000-147) Demonstration Highway, (2000-003) FAR Point Acuity Right Eye & Left Eye Monocular & Both Eyes Simultaneously (20/200-20/20), (2000-010) FAR Point Color Perception (Pseudo Ishihara), (2000-007) NEAR Point Acuity Right Eye & Left Eye Monocular & Both Eyes Simultaneously (20/100-20/20), (2000-012) FAR Point Lateral Phoria (1 Diopter Increments), (2000-025) FAR Point Vertical Phoria (1/2 Diopter Increments), (2000-024) FAR Point Stereo Depth Perception (400-20 Seconds of Arc), (2000-087) FAR Point Fusion, (2000-004) FAR Point Acuity Right & Left Eye Monocular Test Also Right & Left Eye Plus Lens Test (20/200-20/20), (2000-065) FAR Point Acuity Numbers Right Eye & Left Eye & Monocular Both Eyes Simultaneously (20/200-20/20), (2000-148) FAR Point Road Sign Recognition & Depth Perception (620-323-145 Seconds of Arc), (2000-149) FAR Point Color Recognition (Red, Yellow, Green, Blue).

2000RHB Driver Rehabilitation Slide Package

**STEREO OPTICAL F.A.C.T.® CONTRAST SENSITIVITY SLIDE PACKAGE**

Package includes: (2000-007) NEAR Point Acuity Right Eye & Left Monocular & Both Eyes Simultaneously (20/100 to 20/20), (2000-003) FAR Point Acuity Right & Left Eye Monocular & Both Eyes Simultaneously (20/200 to 20/20), (2000-171) 1.5 Cycles Per Degree FAR F.A.C.T. Contrast Sensitivity Monocular, (2000-172) 3 Cycles Per Degree Monocular, (2000-173) 6 Cycles Per Degree Monocular, (2000-174) 12 Cycles Per Degree Monocular, (2000-175) 18 Cycles Per Degree Monocular, (2000-010) FAR Point Color Perception (Pseudo Ishihara), (2000-012) FAR Point Lateral Phoria (1 Diopter Increments), (2000-025) FAR Point Vertical Phoria (1/2 Diopter Increments), (2000-024) FAR Point Stereo Depth Perception (400 to 20 Seconds of Arc), (2000-087) FAR Point Fusion.

2000FS3 F.A.C.T. Contrast Sensitivity Slide Package

**STEREO OPTICAL INDUSTRIAL SLIDE PACKAGE**

For complete occupational screening using accurate and easily understood targets to help determine effective vision. Includes Purdue Job Standards Booklet to easily identify pass/fail results for all occupations. Intermediate lenses available for VDT operators which is the fastest growing population in the industry today. Now test at FAR Point, INTERMEDIATE Point and NEAR Point. Package includes: (2000-147) Demonstration Highway, (2000-001) FAR Acuity Both Eyes Binocular (20/200-20/13), (2000-008) FAR Acuity Right Eye Monocular (20/200-20/13), (2000-005) FAR Acuity Left Eye Monocular (20/200-20/13), (2000-024) FAR Stereo Depth Perception (400-20 Seconds of Arc), (2000-010) FAR Color Perception (Pseudo Ishihara), (2000-025) FAR Vertical Phoria (1/2 Diopter Increments), (2000-012) FAR Lateral Phoria (1 Diopter Increments), (2000-002) NEAR Acuity Both Eyes Binocular (20/200-20/13), (2000-009) NEAR Acuity Right Eye Monocular (20/200-20/13), (2000-006) NEAR Acuity Left Eye Monocular (20/200-20/13), (2000-013) NEAR Lateral Phoria (1 Diopter Increments).

2000ITS Industrial Slide Package



# It's OUR Job to Support Yours.



Because people  
depend on you ...  
You need dependable equipment.

Products of Extraordinary Design, Exceptional Quality,  
Superlative Performance with warranties to match.

At ADC®, everything we build, everything we DO is designed  
to help make your job a little easier.

As one of the world's leading manufacturers of value priced  
diagnostic instruments and accessories, ADC® delivers the  
products YOU need to get your job done, at prices that really  
make sense.

After all, what good are great products ... if you can't afford them?



American Diagnostic Corporation





# (say) OSOM



**FIRST TO MARKET** {that's awesome}  
**FIRST IN MARKET** {that's awesome}

**OSOM® Products.** Innovation, quality, accuracy: standards that exemplify all our diagnostic products. Our hCG, Mono, Strep, Trichomonas and Bacterial Vaginosis products are proven choices for doctors and labs world wide. First to market with the dipstick format has made us a leader in Strep A testing. Please call us to find out more about these awesome products from Genzyme Diagnostics.



OSOM® hCG Combo OSOM® Trichomonas OSOM® BVBLUE OSOM® Strep A  
OSOM® Ultra Strep A OSOM® Mono OSOM® hCG Urine OSOM® Card Pregnancy

**genzyme**  
Diagnostics

Alko Enterprises

**STEREO OPTICAL INTERNATIONAL SLIDE PACKAGE**

For complete occupational screening using accurate and easily understood targets to help determine effective vision. Includes Job Standards Booklet to easily identify pass/fail results for all occupations. Intermediate lenses available for VDT operators which is the fastest growing population in the industry today. Now test at FAR Point, INTERMEDIATE Point and NEAR Point. Package includes: (2000-255) FAR Acuity Both Eyes Binocular (20/200-20/17), (2000-256) FAR Acuity Right Eye Monocular (20/200-20/17), (2000-257) FAR Acuity Left Eye Monocular (20/200-20/17), (2000-258) FAR Stereo Depth Perception (362-9.7 Seconds of Arc), (2000-259) FAR Color Perception (Pseudo Ishihara), (2000-253) FAR Vertical Phoria (1/2 Diopter Increments), (2000-254) FAR Lateral Phoria (1 Diopter Increments), (2000-260) NEAR Acuity Both Eyes Binocular (20/200-20/17), (2000-261) NEAR Acuity Right Eye Monocular (20/200-20/17), (2000-262) NEAR Acuity Left Eye Monocular (20/200-20/17), (2000-264) NEAR Lateral Phoria (1 Diopter Increments), (2000-263) NEAR Vertical Phoria (1/2 Diopter Increments).

2000BL International Slide Package

**STEREO OPTICAL MEDICAL SLIDE PACKAGE WITH HOTV**

Complete vision screening needs preschool through adult. Ideal for school, athletic physicals, employment physicals and general check-ups. Package includes: (2000-026) FAR Point Acuity Binocular Test (20/70), (2000-030) FAR Point Acuity Right & Left Eye Monocular Test (20/30), (2000-031) FAR Point Acuity Right & Left Eye Monocular Test (20/30), (2000-032) FAR Point Acuity Right & Left Eye Monocular Test (20/30), (2000-004) FAR Point Acuity Right & Left Eye Monocular Test Also Right & Left Eye Plus Lens Test (20/200-20/20), (2000-020) FAR Point Muscle Balance (Combination Lateral & Vertical Phoria), (2000-019) FAR Point Color Perception (Tumbling "E"), (2000-022) NEAR Point Muscle Balance (Combination Lateral & Vertical Phoria), (2000-003) FAR Point Acuity Right & Left Eye Monocular & Both Eyes Simultaneously (20/200-20/20), (2000-012) FAR Point Lateral Phoria (1 Diopter Increments), (2000-025) FAR Point Vertical Phoria (1/2 Diopter Increments), (2000-007) NEAR Point Acuity Right & Left Eye Monocular & Both Eyes Simultaneously (20/100-20/20)

2000HOTV Medical Slide Package w/HOTV

**GE MEDICAL TUFFSAT™ SENSORS**

OXY-E4-H Ear Clip Sensor  
 OXY-F1-H Finger Clip Sensor (one is included w/TuffSat)  
 OXY-W4-H Wrap Sensor For 6 lbs & above

**KEYSTONE VIEW VISION SCREENER**

The VS-IV Standard is a versatile machine that provides up to 14 tests using both near (16") and far (20") point, as well as a supplemental lens for testing intermediate vision at 26". The VS-IV Universal is an advanced binocular unit that exceeds the requirements of school vision screening programs, and because the tests are completely non-language, it is particularly useful in testing preschool children and children with possible learning disabilities.

818-400 VS-IV Deluxe-Standard Tests  
 818-590 VS-IV Deluxe-Universal Tests

**KEYSTONE VIEW DRIVERS VISION SCREENER**

The DVS-IV is designed for the rapid, yet accurate, evaluation of large groups. Vision testing is both standardized and completely confidential. The DVS-GT contains the same features as the DVS-IV plus tests for contrast sensitivity (simulates snow, fog, rain and night driving conditions) and glare recovery (the ability to adapt to bright lights such as oncoming traffic).

818-700 DVS-IV Drivers Vision Screener  
 818-830 DVS-IV Drivers Vision Screener w/glare

**KEYSTONE VIEW COMPUTER VISION SCREENER AND SOFTWARE**

Working with nearly any PC, the CVS-I both automates the testing process, using a standardized script, and lets you control the test instrument remotely from the computer station. The Visionary Software allows you to easily record data and produce reports.

818-130 CVS-I Screener w/Visionary Software  
 818-150 Visionary Software

**KEYSTONE VIEW TESTING AND TRAINING TARGET SETS FOR USE IN CORRECT-EYE SCOPE**

005-348 Van Orden Transparencies, 16/set  
 006-166 Franzblau Tracing Patterns Forms (Seven Patterns Series)  
 297-170 Red Crayons, 10/pk  
 5348-1 Van Orden Glass Transparency  
 5401 Leavell Exercise Pads Set (A, B, C), 3pads/pk  
 5402 Leavell Exercise Pads Set (AA, BB, CC) for Stereoreader, 3pads/pk  
 725351 Red Rock Games Only, 7 games/bx  
 735-265 Leavell Exercise Pads CC for Stereoreader, 3 pads/pk  
 735-265-1 Leavell Exercise Pads C, 3 pads/pk  
 735-266 Leavell Exercise Pads AA for Stereoreader, 3 pads/pk  
 735-266-1 Leavell Exercise Pads A, 3 pads/pk  
 735-267 Leavell Exercise Pads BB for Stereoreader, 3 pads/pk  
 735-267-1 Leavell Exercise Pads B, 3 pads/pk  
 735-289 Groffman Pattern Cards (2 test cards, 1 demonstration card) Includes Record Form 735-287  
 735-292 Delacato Exercise/Reading Cards, 152/pk  
 735-294 Zweig-Bruno Pads #1 (Series A, B, C), 3 pads/pk  
 735-295 Zweig-Bruno Pads #2 (Series 1, 2, 3), 3pads/pk  
 735-296 Groffman Books #1, 3 books/pk  
 735-297 Groffman Books #2, 3 books/pk  
 735-426 V.O. 2 Star Pattern Pads, 100 sheets/pad, 3 pads/pk  
 735-429 Universal Series Correct-Eye Graphs, 104 drawings/pad, 3 pads/pk  
 735-436 V.O. Base-In Training Charts Pads, 100 sheets/pad, 3 pads/pk  
 735-955 Manas Series Correct Eye Graphs, 110 drawings/pad, 3pads/pk  
 736-690 Peek-A-Boo Red Coloring Books, 10/pk  
 736-831 Brock Technique Transparency Set, 28/set

**KEYSTONE VIEW TESTING AND TRAINING TARGET SETS FOR USE WITH TELEBINOCULAR**

005-100 Visual Skills Test Set  
 005-100-8 Severe Color Blind Test  
 005-100-9 Mild Color Blind Test  
 005-101 Driver's Vision Test Set  
 005-102 IVS (OSV) Rapid Screening Test Set  
 005-103 O.E.P. Test Set (Includes 005-100)  
 005-104 Spache Binocular Reading Test (May also be used in Stereoscope)  
 005-106 Professional Performance Test  
 005-106-11 Gross Stereopsis Test  
 005-106-13 Astigmatic Dial Test  
 005-106-6 Macular Suppression Test  
 005-107 Pilot Test Set  
 005-108 (OVS) (Industrial Screening) Test Set  
 005-109 Snellen Near Point Acuity Test  
 005-110 Snellen Far Point Acuity Test  
 005-111 Near Point Stereopsis Test  
 005-112 D.C. Aciators Test Set  
 005-113 Progress of Fusion Test Set  
 005-115 Industrial Short Tests  
 005-116 Multi-Stereo Test Set  
 005-118 Mark Diplopia Test Set  
 005-119 O'Brian Central Field  
 005-121 AN Series Training Targets  
 005-128 Broward County Test Set  
 005-129 40" Intermediate, Sloan Letter Acuity  
 005-130 Landolt Rings Far Point Test



## KEYSTONE VIEW TESTING AND TRAINING TARGET SETS FOR USE WITH TELEBINOCULAR (Continued)

|           |   |
|-----------|---|
| 005-131   | D.C. Aviators Short Test Set  |
| 005-134   | Intermediate, Landolt Ring Acuity (40")                                   |
| 005-135   | Rapid Screening Tests for Schools   |
| 005-138   | Illinois School Test  |
| 005-138-4 | Illinois Preschool Test   |
| 005-139   | Massachusetts School Test   |
| 005-140   | Michigan School Test  |
| 005-141   | New York School Test  |
| 005-141-1 | Tumbling E Far Point Test   |
| 005-142   | Astigmatism Screening Test  |
| 005-155   | Clinical Fusion Test Set  |
| 005-156   | Walton Acuity/Astigmatism Test  |
| 005-162   | Gesell Test Set   |
| 005-164   | Anisekonia Test Set (Special Order)                                       |
| 005-165   | Near Point Landolt Ring Acuity Test                                       |
| 005-166   | Near Point "Tumbling E" Acuity Test                                       |
| 005-170   | Peek-A-Boo Test Set   |
| 005-171   | Sloan Letter Acuity Test - Near Point                                     |
| 005-178   | Palmer Freed Dissociation Test (May also be used in Stereoscope)          |
| 005-179   | Keystone Binocular Skills Test & Manual (May also be used in Stereoscope) |
| 005-180   | Supplemental NP Test, Sloan Type  |
| 005-184   | Paramacular Suppression Test (May also be used in Stereoscope)            |
| 005-185   | Universal Driver Test   |
| 005-189   | Retinal Rivalry Test (May also be used in Stereoscope)                    |
| 005-192   | Contact Lens Card   |
| 005-193   | Plus Lens Acuity Target   |
| 005-194   | DB-2B Left Eye, DB-3B Right Eye, DB-1B Both Eyes Acuity Test              |
| 005-200   | Illinois College of Optometry Test  |
| 005-202   | Pacific College of Optometry Test   |
| 005-203   | Special Far Point Peek-A-Boo Acuity Test                                  |
| 005-204   | Physicians Visual Rating Tests  |
| 005-209   | Far Point Sloan Letter Acuity Test  |
| 005-210   | Sloan Intermediate Test   |
| 005-213   | Special NP Peek-A-Boo Acuity Target                                       |
| 005-218   | Combined Phoria FP  |
| 005-219   | Combined Phoria NP  |
| 005-220   | Sloan Digit Acuity FP   |
| 005-221   | Sloan Digit Acuity NP   |

## KEYSTONE VIEW TELEBINOCULAR AND ATTACHMENTS

A pioneer in the vision testing field. Offers unmatched testing and training versatility, including diagnosis, visual skills screening and vision training.

|         |   |
|---------|---|
| 1103S   | Visual Survey Telebinocular-Standard  |
| 728-000 | Attachment: Perimeter Attachment  |
| 728-030 | Pilot Plus Visual Survey Telebinocular  |
| 728-044 | Attachment: Plus Lens Attachment, +1.50, Includes 5193 DB-18 Plus Lens Acuity Target  |
| 728-045 | Attachment: Plus Lens Attachment, +1.75D, Includes 5193 DB-18 Plus Lens Acuity Target |
| 728-046 | Attachment: Plus Lens Attachment, +2.00D Includes 5193 DB-18 Plus Lens Acuity Target  |
| 728-047 | Attachment: Plus Lens Attachment, +2.25D Includes 5193 DB-18 Plus Lens Acuity Target  |
| 728-053 | Attachment: Plus Lens Attachment, +2.50D Includes 5193 DB-18 Plus Lens Acuity Target  |
| 728-056 | Ophthalmic Telebinocular  |
| 728-067 | Driver Screening Telebinocular  |
| 728-935 | Industrial Telebinocular  |

## KEYSTONE VIEW EQUIPMENT

|         |  |
|---------|--|
| 004-016 | Colored Circles-Transparent, 10/pk                                   |
| 004-112 | Eccentric Circles-Transparent, 20/pk                                 |
| 297-337 | Red/Green Filters Elastic Band, 1pr                                  |
| 297-449 | Sherman Playing Cards (2 decks), 2 deck/pk                           |
| 725-000 | Correct-Eye Scope  |
| 725-016 | Stereo Reader  |
| 725-019 | Closed Stereoscope - Hand Held                                       |
| 725-025 | Delacato Mask For Stereo Reader                                      |
| 725-026 | Open (725-030) Stereoscope w/pedestal mount                          |
| 725-030 | Open Stereoscope - Hand Held   |
| 725-040 | Closed (725-019) Stereoscope w/pedestal mount                        |
| 725-090 | Spectacles Red/Green Filters, 1pr                                    |
| 725-330 | Franzblau Multi-Sensory Coordinator (includes tracing pattern forms) |
| 725-336 | Van Orden Trainer w/aux lens holder                                  |
| 725-338 | Complete Franzblau Red Rock System                                   |
| 725-366 | Cheiroscopic Attachment for 725-374 & 725-330                        |
| 725-374 | Light Box Transparency Viewer  |
| 725-395 | Colored Circles-Opaque Card, 50/pk                                   |
| 725-396 | Eccentric Circles-Opaque, 25pr/pk                                    |
| 725-924 | Red Rock Face Plate For Use w/725-000 or 725-330                     |
| 728-288 | Sherman TV Training Kit  |

## KEYSTONE VIEW PARTS & ACCESSORIES

|           |   |
|-----------|---|
| 297-390   | VS & DVS Carrying Case  |
| 297-504   | VS II & DSV III Lamps GE400, 10/pk  |
| 297-800   | VS II & DSV III 1A, 250V Fuse, 5/pk   |
| 298-410   | Lamps For Telebinocular, Transparency Viewer, Correct-Eye Scope, Van Orden Trainer, Franzblau Multi-Sensor Coordinator, 25 Watt, 2/pk |
| 298-808   | Telebinocular Large Dust Cover  |
| 298-971   | DVS II Lamps, 10/pk   |
| 725-017   | Pedestal Base for Hand Held Open Stereoscope  |
| 725-087   | Pedestal Base For Hand Held Closed Stereoscope  |
| 725-323   | Plastic Faceplate w/clip for Franzblau  |
| 725-331   | Sensor Pen Assembly For Franzblau Multi-Sensor Coordinator, 1 pair  |
| 728-055   | Auxiliary Lens Holders, 1 pair  |
| 728-241   | Occluders, Flip Type For Telebinocular, 1 pair  |
| 728-310   | Telebinocular Carrying Case   |
| 728-391   | Wooden Pointers, 20/set   |
| 815-249   | DVS II Hand Control   |
| 815-278   | DVS II Transformer  |
| 815-499   | VS II & DVS III Transformer   |
| 815-503   | VS II 14-Pin Hand Control   |
| 815-503-1 | VS II 16-Pin Hand Control   |
| 815-754   | DVS III 14-Pin Hand Control   |
| 815-754-1 | DVS III 16-Pin Hand Control   |
| 818-062   | CVS-I / VS-IV / DVS-IV/GT Target LED  |
| 818-066   | DVS-GT LED Glare Lamp   |
| 818-114   | CVS-I / VS-IV Hand Control  |
| 818-140   | DVS IV - 110V Transformer   |
| 818-145   | DVS IV - 220V Transformer   |
| 818-202   | DVS-GT Hand Control   |
| 818-470   | CVS-I Series 12VDC Universal PS   |
| 818-720   | DVS-IV Hand Control   |



**KEYSTONE VIEW RECORD FORMS**

|           |   |
|-----------|---|
| 5540      | Vinyl Matrix for OVS Form #8 SPECIAL ORDER: Minimum Order 10 Packs, 3 pd/pk                         |
| 727-226   | VS II / VS IV Universal Target Record Form, 3 pd/pk   |
| 735-145   | OVS Rapid Screening Record Form, 3 pd/pk  |
| 735-236   | Orthoscope Refractions Form OTR (5" x 8"), 3 pd/pk  |
| 735-237   | Professional Performance Test Form (8 1/2" x 11"), 3 pd/pk  |
| 735-269   | OVS Test Record Form #8, 3 pd/pk  |
| 735-270   | OVS Test Record Form #9, 3 pd/pk  |
| 735-272   | Doctor's Record Form #3A Visual Skills Profile (8 1/2" x 11"), 3 pd/pk                              |
| 735-272-1 | Doctor's Record Form #3A Visual Skills Profile (5" x 8"), 3 pd/pk                                   |
| 735-273   | Driver Test Set Record Form #10, 3 pd/pk  |
| 735-275   | Visual Survey Short Test Form #4, 3 pd/pk   |
| 735-277   | Delacato Test Summary Sheets, 3 pd/pk   |
| 735-278   | Leavell Test Sheets, 3 pd/pk  |
| 735-287   | Visual Tracing Test Form Groffman, 3 pd/pk  |
| 735-384   | Industrial Short Test Form #12, 3 pd/pk   |
| 735-385   | Professional Performance Test Form (5" x 8"), 3 pd/pk   |
| 735-386   | Doctor's Form #2A VS (8 1/2" x 11"), 3 pd/pk  |
| 735-386-1 | Doctor's Form #2A VS (5" x 8"), 3 pd/pk   |
| 735-388   | Skills Test Form DG-P, 3 pd/pk  |
| 735-389   | Skills Test Form DG-1 / DG-2, 3 pd/pk   |
| 735-390   | Skills Test Form DG-3 / DG-4, 3 pd/pk   |
| 735-394   | Spache Binocular Reading Test Record Form, 3 pd/pk  |
| 735-397   | Vinyl Matrix Overlay for OVS Form #9 SPECIAL ORDER: Minimum Order 10 Packs, 3 pd/pk                 |
| 735-399   | Ready to Read Record Form #6, 3 pd/pk   |
| 735-400   | Multi-Stereo Scoring Form, 3 pd/pk  |
| 735-401   | Mark Diplopia Diagnostic Chart, 3 pd/pk   |
| 735-402   | Visual Survey Parents Notify Form, 3 pd/pk  |
| 735-403   | OVS Referral Form #66, 3 pd/pk  |
| 736-161   | Stereometric Record Form M-DC, 3 pd/pk  |
| 736-386   | School Survey Form #5B (2 Part Form) Use w/005-135, Follow by 005-100, 3 pd/pk                      |
| 736-561   | School Survey Form #5A, 3 pd/pk   |
| 736-698   | Visual Skills Survey Form, Use w/005-135 / 005-100, 3 pd/pk   |
| 736-699   | Peek-A-Boo Series Form (8 1/2" x 11"), 3 pd/pk  |
| 736-744   | Employee Vision Record Form #9-A, 3 pd/pk   |
| 736-776   | Modified Telebinocular Technique Form, 3 pd/pk  |
| 736-804   | Peek-A-Book Series Form (5" x 8"), 3 pd/pk  |
| 736-832   | Palmer-Freed Dissociation Form, 3 pd/pk   |
| 736-878   | Pilot Vision Test Pre-Record Form, 3 pd/pk  |
| 737-017   | Physician's Visual Rating Form, 3 pd/pk   |
| 737-025   | DVS II / DVS III / DVS IV Record Forms, 3 pd/pk   |
| 737-025A  | DVS-GT Standard Record Forms, 3 pd/pk   |
| 737-161   | Walton MVS Form W-T, 3 pd/pk  |
| 737-225-1 | VS II / VS IV Standard Record Form (Replaces 737-225A, 737-131, 737-225 & 737-254) English, 3 pd/pk |
| 737-225B  | VS II / VS IV Standard Target Record Form (2 Part Form) Spanish, 3 pd/pk                            |
| 737-255   | VS II / VS IV Universal Target Record Form (5" x 8"), 3 pd/pk                                       |
| 738-769   | Aniseikonia Record Form R1, 3 pd/pk   |

**KEYSTONE VIEW SERVICES AND SETS**

|         |   |
|---------|---|
| 006-100 | Visual Skills Program 1 (Screen) Includes Instruction Manual, 728-056, 005-100, 735-272 & 735-272-1   |
| 006-101 | Visual Skills Program 1 (Diagnose) Includes Instruction Manual, Rating Cards, 728-056, 005-100 & 735-237  |
| 006-103 | Brock Technique Service, Includes 28 Large Plastic Slides & Lightbox  |
| 006-107 | Combined Programs 1 & 2 Screening & Diagnostic Test Set   |
| 006-116 | Hand-Eye Coordination Test Set, 3 test/set  |
| 006-129 | Keystone 3-D Fusion Games Set, Includes: 20 Red, Green & Black 8 1/2" x 11" Picture Charts, Easy-to-Use Manual, Easel, Pointers, Record Forms & 1 Pair Red/Green Filters  |
| 006-130 | Basic Binocular Test Set, Includes: 12 Special Targets Bound in Flip Chart Form, Non-Glare Plastic-Coated Cards, Red/Green Filters & 1 Pair of Pointers   |
| 006-169 | Economy Color Book Package, Includes: 10 Books, 10 Red Crayons & 1 Pair Red/Clear Filters   |
| 6104    | Complete Van Orden Program, Kit Includes: Keystone View Correct Eye Scope, The Van Orden Trainer, 16 Transparency & Opaque Slides, 1 Set of Split Cards, Red/Green Filters, 2 Pads of Van Orden "Star" Graphs w/VO-1 matrix, 2 Pads of Training Charts & Manual |
| 6106    | School Vision Screening Program: Includes 728-030 w/005-100 & 005-135, Instruction Manual, Familiarization Chart & 736-386  |
| 6193    | Pediatric Vision Screening Service, Comes Equipped with the Visual Survey/School Telebinocular & Peek-A-Boo Test Set, Instruction Manual & 100 Peek-A-Boo Test Forms  |

**WELCH ALLYN SURESIGHT™ VISION SCREENER**

Screen with one-button touch from 14" away. Five-second, automatic test addresses compliance problems and makes testing easier. Minimal cooperation is required, making it ideal for use with young children. The lights and sounds engage the child's attention. Automatically measures and indicates abnormal readings. Addresses eye charts' sensitivity and specificity problems. Automatically screens for common vision problems, including near- and farsightedness, astigmatism, and anisometropia. Hand-held unit operates for three hours of continuous testing. Available battery-operated printer enables you to get printouts anywhere, anytime. Weighs two lbs. and features an ergonomic grip and replaceable, self-adjusting hand strap. Features optics and electronics shock-mounted inside housing to protect unit. The sealed unit helps keep out dirt. Screens babies, children, and adults. You can even test patients while they're wearing glasses or contact lenses.

|       |  |
|-------|--|
| 14000 | SureSight™ Vision Screener Set, Includes 110-130V Charging Transformer & Stand |
|-------|--|

**Accessories**

|         |                      |
|---------|----------------------|
| 05140-U | Carrying Case        |
| 14042   | Thermal Printer      |
| 53600   | Printer Paper, 24/cs |



V50-1ANY



V52-1D

### BURDICK VISION SYSTEM

The software includes: "Search by time" in full disclosure, R-R printed on all strips, Date, day, and time on all pages, Review workstation capability for Vision, Forms Review in Vision, French and Spanish language availability, Availability of a 19" flat panel monitor option for Vision and Vision Premier systems, Availability of Vision 5L recorders with software only configurations, Addition of Rxio to all domestic configurations, distribution of Rxio outside the US is prohibited.

- V50-1ANY Premier System Includes: English Language Software, Recorder NOT included
- V51-1G Premier ACCESSORY: English Language Software with PC Card, Single user license
- V52-1D Vision System with 92515 Holter Recorder, English Language Software
- V53-1DG Vision Software w/PC card interface. Vision 5L recorder and English language Single user license
- V53-1G Vision System ACCESSORY: English Language Software, PC Card Interface

### BCI 3180 PULSE OXIMETER

Provides an exceptional value for everyday oximetry needs. This compact, stand-alone monitor delivers fast, reliable SPO<sub>2</sub>, pulse rate and pulse strength measurements, complete with a plethysmographic waveform and up to 30 hours of on-screen trending. Plus, the bright liquid crystal display, easy-to-use menu and audible and visible alarms enhance clinical efficiency and quality of care for patients ranging from neonate to adult.

- 3180 3180 Pulse Oximeter

### BCI AUTOCORR™ DIGITAL PULSE OXIMETER

Uses refined noise-reducing hardware and a unique patented software algorithm. Instead of relying on a voltage threshold to determine a pulse, it analyzes a digitized signal at a rate of 120 times per second, and compares it to previous pulse data. These pulses are constantly stored and updated, and used as dynamic templates to determine the validity or invalidity for newly-encountered pulsatile signals. Patented Serial Autocorrelation (SAC) technology provides reliable monitoring.

- 3304 Autocorr® Pulse Oximeter
- 3304LS Autocorr® Pulse Oximeter w/loud alarms
- 3404-000 Autocorr® Plus SpO<sub>2</sub> Pulse Oximeter
- 3404-001 Autocorr® Plus SpO<sub>2</sub> Pulse Oximeter/ECG/Respiration

### BCI DIGIT™ FINGER OXIMETER

Combines monitor and sensor into one unit that provides fast, reliable SpO<sub>2</sub>, pulse rate, and pulse strength measurements on patients from pediatric to adult. Ideal for use within emergency settings, in hospital or clinical environments, or for home use. Portable, lightweight pocket-size monitor. One-button keypad. Large LED display. Pulse-strength display bar indicates relative perfusion at the measurement site. Auto power shutdown after 8 seconds. Low battery indicator flashes when 30 minutes of battery time remains.

- 3420G Finger Oximeter, Gray
- 3420Y Finger Oximeter, Yellow



### BCI FINGERPRINT® HAND-HELD PULSE OXIMETER

Portable pulse oximetry monitoring system with built-in printer includes: 3401 Monitor (with 3 rolls of paper); 1884 Operation Manual; 3418 Protective Rubber Boot; 3044 Oximetry Finger Sensor; (4) Alkaline AA Cell Batteries.

- 3401-000 FingerPrint® Hand-Held Pulse Oximeter
- 3403-000 FingerPrint® Sleep Hand-Held Pulse Oximeter/Summary Printer



### BCI HAND-HELD PULSE OXIMETER

Portable pulse oximetry monitoring system includes: 3301 Monitor; 1818 Operation/Service Manual; (3) Alkaline C Cell Batteries.

- 3301 Hand-Held Pulse Oximeter, with Oximetry Finger Sensor
- 3301-EP Hand-Held Pulse Oximeter, with 3078 Ear Probe



### BCI HAND-HELD PULSE OXIMETER WITH ALARMS

Portable continuous or spot check pulse oximetry monitoring system with alarms and internal rechargeable battery includes: 3303 Monitor; 3311 Oximetry Cable, 5 ft.; 3044 Oximetry Finger Sensor; 1850 Clinician's Operation Manual; 1854 Home-Use instruction Book; 8210 AC Charger, 105-125V 60Hz; 3353 Protective Rubber Boot with Carrying Strap and Tilt Stand.

- 3303 Hand-Held Pulse Oximeter w/alarms

### BCI MINICORR DIGITAL HAND-HELD PULSE OXIMETER

Designed to provide fast, reliable measurements on challenging patients. The device rejects spontaneous motion and other artifacts. Informs the operator of excess artifacts and improper sensor positioning and processes low perfusion signals. Designed for spot-checking or bedside monitoring. This portable unit is ideal for use in home care and clinical environments, as well as during sleep screening. Easy to use, 5-key operation. Infrared link to printer. Battery or charger. Boot is available. System includes 3402 Monitor; 3311 Oximetry Cable, 5 ft.; 3044 Oximetry Finger Sensor; 1896 Operation Manual.

- 3402 MiniCorr Digital Hand-Held Pulse Oximeter

**BCI PULSE OXIMETER ACCESSORIES**

1100 Airway Adapter, Straight w/o filter (8200, 8400), 10/pkg  
 1101 Adult O2 Nasal Cannula + Sample Line, (8200, 8400) 10/pkg  
 1102 Pediatric O2 Nasal Cannula + Sample Line (8200, 8400), 10/pkg  
 1105 Airway Adapter, Straight w/filter (8400), 10/pkg  
 1121 Adult Oral/Nasal CO2 Cannula + Sample Line (8200, 8400), 10/pkg  
 1123 Adult Oral/Nasal CO2 Sample Line (8200, 8400), 10/pkg  
 1124 Pediatric Oral/Nasal CO2 Sample Line (8200, 8400), 10/pkg  
 1129 Adult Nasal CO2 Sample Line (8200, 8400), 10/pkg  
 1130 Pediatric Nasal CO2 Sample Line (8200, 8400), 10/pkg  
 1131 Infant Nasal CO2 Sample Line (8200, 8400), 10/pkg  
 1140 15 ft CO2 Sample Line Extension (8200, 8400), 10/pkg  
 1151 Airway Adapter, Pediatric w/o filter (8200, 8400), 10/pkg  
 1152 Airway Adapter, Pediatric w/filter (8200, 8400), 10/pkg  
 1178 Moisture Trap (9004), 5/pkg  
 1179 External CO2 Absorber (9004), 2/pkg  
 1187 Moisture Trap Adapter (1178)  
 1300 Adult Disposable Oximetry Finger Sensor, 10/bx  
 1301 Pediatric Disposable Oximetry Finger Sensor, 10/bx  
 1302 Neonate Disposable Oximetry Finger Sensor, 10/bx  
 1303 Infant Disposable Oximetry Finger Sensor, 10/bx  
 1310 D.O.T. Reusable Oximetry Sensor, 3/bx  
 1311 D.O.T. Sensor Replacement Tape, 24/pkg  
 1606 Oximeter/ECG Simulator (3303, 8400)  
 1606HH Oximeter Simulator & 5 ft Cable (3301, 3301, 3401, 3402)  
 1611 115V AC Power Supply (3402)  
 1614 AC Charger, 115V, 60Hz (6004, 9004)  
 1617 105-125V AC Power Supply (3404, 3408)  
 1621 Oximetry Data Management Program  
 1818 3301 Operation/Service Manual  
 1818R 3301 Operation/Service Manual  
 1835 8200 Operation Manual  
 1836 8200 Service Manuals  
 1850 3303 Clinician Operation Manual  
 1851 3303 Service Manual  
 1854 3303 Home-Use Instruction Book  
 1857 Autocorr® Oximeter Clinician Operation Manual (3304)  
 1858 Autocorr® Oximeter Service Manual (3304)  
 1861 Autocorr® Oximeter Home-Use Instruction Book (3304)  
 1872E Capnocheck® Plus/Capnocheck Sleep Operation Manual (9004)  
 1873 Capnocheck® Plus/Capnocheck Sleep Service Manual (9004)  
 1876E Mini-Torr Plus® NIBP Monitor Operation Manual (6004)  
 1877 Mini-Torr Plus® NIBP Monitor Service Manual (6004)  
 1884 3401 FingerPrint® Operation Manual  
 1884SS FingerPrint® Sleep Oximeter Operation Manual (3403)  
 1884SSHU FingerPrint® Sleep Oximeter Home-Use Instruction Note (3403)  
 1885 3401 FingerPrint® Service Manual  
 1885SS FingerPrint® Sleep Oximeter Service Manual (3403)  
 1886 Advisor® Vital Signs Monitor Operation Manual (9200)  
 1887 Advisor® Vital Signs Monitor Service Manual (9200)  
 1888 Autocorr® Plus Oximeter/ECG Monitor Clinician Operation Manual (3404)  
 1889 Autocorr® Plus Oximeter/ECG Monitor Service Manual (3404)  
 1892 Graphics Printer Operation Manual (3408)  
 1893 Graphics Printer Service Manual (3408)  
 1894 Capnocheck® II Monitor Operation Manual (8400)  
 1895 Capnocheck® II Monitor Service Manual (8400)  
 1896 MiniCorr Pulse Oximeter Operation Manual (3402)  
 1897 MiniCorr Pulse Oximeter Service Manual (3402)  
 1898 MiniCorr Pulse Oximeter Home Use Instruction Book (3402)  
 1900 Digit® Operation Manual (3420)  
 1904 3180 Oximeter Clinician Operation Manual  
 1905 3180 Oximeter Service Manual  
 3005 Power Cord (3180, 9100, 9200)  
 3025 Infant Oximetry Wrap Sensor  
 3026 Neonate Oximetry Wrap Sensor  
 3043 Universal "Y" Oximetry Sensor  
 3044 Standard Oximetry Finger Sensor  
 3049 Adhesive Strip Tape, 40/pkg  
 3078 Ear Oximetry Probe Sensor  
 3106 ECG Shielded Safety Type Cable, 3-Lead (3404, 9100, 9200)

3107 ECG Shielded Leads, 40", Safety Type Wires, 3-Lead (3404, 9100, 9200), 3/pkg  
 3134 Neonate Attachment Tape, 50/pkg  
 3135 Infant Attachment Tape, 50/pkg  
 3136 Neonate Attachment Tape, 100/pkg  
 3137 Infant Attachment Tape, 100/pkg  
 3138 Universal "Y" Posey Wrap Cover, 10/pkg  
 3178 Pediatric Finger Sensor  
 3310 Seiko Printer Model #DPU-201GSU (3301, 3303, 8200)  
 3311 Oximetry Cable, 5 ft (3301, 3303, 3401, 3402, 8400)  
 3311L Oximetry Extension Cable, 15 ft (3403, 9004-051)  
 3311N Nellcor Oximetry Extension Cable 4 ft  
 3313 Printer Interface Cable, 1 ft (3301)  
 3315 3301 Carrying Case w/belt clip/shoulder strap, No Logo  
 3318 3301 Protective Rubber Boot w/carrying strap  
 3322 3303 Carrying Case w/belt clip/shoulder strap  
 3339 PC Adapter Cable (3301, 3303, 3401, 3402, 8200)  
 3350 Printer Interface Cable  
 3351 Analog Output Adapter (3304)  
 3352NC Digital Output Adapter (Nurse Call - Closed Circuit) (3304, 3404)  
 3352NO Digital Output Adapter (Nurse Call - Open Circuit) (3304, 3404)  
 3353 3303 Protective Rubber Boot w/carrying strap/tilt stand  
 3354 Pole Mount Bracket (3303, 8200)  
 3356 Clarity® Series IV Pole Bracket (3304, 3404, 6004, 9004)  
 3361 Printer Interface Cable (3304)  
 3362 PC Interface (3304/3404)/Graphics Printer Interface (6004)  
 3365 Printer Interface Cable (9004)  
 3366 PC Interface (9004)/Graphics Printer Interface (9004)  
 3369 Clarity® Series Carrying Case w/shoulder strap (3304, 3404, 6004)  
 3371 Graphics Printer Interface (3404)  
 3389 Interface Cable (1621)  
 3406 ECG Shielded Safety Type Cable, 5-Lead (3404, 9200)  
 3407 ECG Shielded Leads, 40" Safety Type Wires, 5-Lead (3404, 9200), 5/pkg  
 3408 Graphics Printer  
 3409 Printer Paper (3402, 8400), 4/bx  
 3410A DC Power Cable (3404, 3408)  
 3414 Analog Output Cable (3404)  
 3415 3401 FingerPrint Carrying Case  
 3418 3401 FingerPrint Protective Rubber Boot w/carrying strap  
 3419 MiniCorr Oximeter Protective Rubber Boot w/carrying strap (3402)  
 3426 Digit® Lanyard (3420)  
 3427 Digit® Carrying Case (3420)  
 3444 Comfort-Clip® Oximetry Finger Sensor (3303, 8400)  
 5093 ETCO2 Calibration Gas (8200, 8400)  
 6007 Mini-Torr Plus® NIBP Monitor Printer Upgrade Kit (6004)  
 6008 Mini-Torr Plus® NIBP Monitor Oximetry Upgrade Kit (6004)  
 6012 3401 FingerPrint® Printer Paper, 4/pkg  
 6014 Electronic Thermometry Temperature Probe Covers (6004), 1000/ctn  
 6015L Electronic Thermometry Oral Temperature Probe w/9 ft cable (6004)  
 6015R Electronic Thermometry Rectal Temperature Probe w/9 ft cable (6004)  
 6016 Electronic Thermometry Calibration Key (6004)  
 6110 NIBP Cuff Supply Line (Luer-Lock), 10 ft (may need Transition Supply Line)  
 6110S NIBP Cuff Supply Line (Luer-Lock), 5 ft (may need Transition Supply Line)  
 6111 Printer Paper (9100, 9200), 4/bx  
 6120A Adult (25.3-34.3cm) Non-Disposable NIBP Cuff  
 6120C Child (15.8-21.3cm) Non-Disposable NIBP Cuff  
 6120I Infant (9.8-13.3cm) Non-Disposable NIBP Cuff  
 6120LA Large Adult (32.1-43.4cm) Non-Disposable NIBP Cuff  
 6120N Newborn (7.7-10.5cm) Non-Disposable NIBP Cuff  
 6120SA Small Adult (20.0-27.0cm) Non-Disposable NIBP Cuff  
 6120SC Small Child (12.4-16.8cm) Non-Disposable NIBP Cuff  
 6121 10 ft NIBP Cuff Transition Supply Line  
 6121S 5 ft NIBP Cuff Transition Supply Line  
 6122 10 ft NIBP Cuff Supply Line (Quick Connect)  
 6122S 5 ft NIBP Cuff Supply Line (Quick Connect)  
 6123 5 ft NIBP Cuff Neonate Transition Supply Line  
 6124 5 ft NIBP Cuff Neonate Supply Line (Quick Connect)  
 6140 Neonate (3-6cm) Single Use NIBP Cuff, 3/pkg  
 6141 Neonate (4-8cm) Single Use NIBP Cuff, 3/pkg  
 6142 Neonate (6-11cm) Single Use NIBP Cuff, 3/pkg  
 6143 Neonate (7-13cm) Single Use NIBP Cuff, 3/pkg



**BCI PULSE OXIMETER ACCESSORIES (Continued)**

- 6144 Neonate (8-15cm) Single Use NIBP Cuff, 3/pkg
- 6145 Neonate Variety Pack Single Use NIBP Cuff (2 of each size), 10/pkg
- 7108 Clarity® Series/9200 Roll Stand w/mounting bracket & basket (3304, 3404, 6004, 9004)
- 8044 8 ft CO2 Sample Line (8200, 8400), 10/pkg
- 8061 Calibration Gas Flow Regulator w/gauge (8200, 8400)
- 8208 Gas Manifold Filters (8200, 8400), 10/pkg
- 8210 AC Charger, 105-125V, 60Hz (3303, 8200)
- 8211 4 ft CO2 Sample Line (8200, 8400), 10/pkg
- 8213 Printer Interface Cable, 1 ft (8200)
- 8214 Patient Attachment Kit: 1100, 8211, 8208 (8200, 8400), 10/pkg
- 8217 Capnocheck Calibration Kit (8200, 8400)
- 8223 Calibration Adapter (9004)
- 8403 Battery Eliminator 105-15VAC, 60Hz (8400)
- 8404 Battery Charger (8400)
- 8408 Battery 7.4V Lithium-Ion Rechargeable (8400)
- 8409 Pole Mount Bracket (8400)
- 8411 Infrared Printer Model #HP82240B (3402, 8400)
- 8412 Capnocheck® II Monitor Protective Rubber Boot (8400)
- 8414 Capnocheck® II Capnograph/Oximeter Carrying Case (8400, 8401)
- 9007 Capnocheck® Plus/Capnocheck Sleep Oximetry Upgrade Kit (9004)
- 9014 3 Channel Analog Output Cable (polysomnography) (9004)
- 9015 4 Channel/Distance Analog Output Cable (polysomnography) (9004)
- 9047 Capnocheck® Plus Capnograph Filter Upgrade Kit (9004)
- 9048 External 3.5 Micron Filter (9004), 2/pkg
- 9189 FIO2 Adapter T (9004), 10/pkg
- 9190 FIO2 Sensor (9004)
- 9191 FIO2 Cable (9004)
- 9193T Capnocheck® Plus Capnograph FIO2 Upgrade Kit (9004)



**BURDICK OXY 100 PULSE OXIMETER**

Fast, accurate SpO2 levels whenever and wherever needed. Displays SpO2 level, pulse rate and pulse strength on an easy to read LED screen. Uses 3 "C" cell batteries and automatically switches off when not in use. Compact, portable and durable. Up to 99 individual tests can be stored for printout using the optional printer.

- 92521-1A Burdick Oxy 100 w/finger sensor
- 92522-1A Burdick Oxy 200 w/finger sensor



**GE MEDICAL TUFFSAT® PULSE OXIMETER**

Ideal oximetry spot-check solution for physicians' offices and clinics. The hand-held oximeter fits in your pocket or is conveniently worn using the included carrying case with belt clip. Includes a Plr® (relative perfusion index measurement) for easy quantification of signal quality. Available in four colors; yellow, blue, pink and white.

- 6051-0000-160 TuffSat® Yellow w/sensor
- 6051-0000-185 TuffSat® Blue w/sensor
- 6051-0000-186 TuffSat® Pink w/sensor
- 6051-0000-187 TuffSat® White w/sensor

**GE MEDICAL MASIMO PULSE OXIMETRY SENSORS & CABLES**

Accurate under patient motion or movement. Accurate for low perfusion patients. Latex free.

- 2009743-001 Masimo PCO8 Cable
- 2009745-001 Masimo Finger Sensor, Adult, Reusable
- 2002799-001 Masimo Finger Sensor, Pediatric, Reusable
- 2002800-001 Masimo Finger Sensor, Adult, Reusable

**GE MEDICAL NELLCOR DISPOSABLE PULSE OXIMETRY SENSORS (PROBES)**

Single patient use reduces potential for cross-contamination. Sensor placement indexes printed on wrap for easy alignment. Latex free.

- 407705-001 Nellcor OxiSensor Neonatal Disposable Probes, Patient Range
- 407705-002 Nellcor OxiSensor Infant Disposable Probes, Digital Probe, Patient Range 3-20 KG, 24/cs
- 407705-003 Nellcor OxiSensor Pediatric Disposable Probes, Digital Probe, Patient Range 10-50 KG, 24/cs
- 407705-004 Nellcor OxiSensor Adult Disposable Probes, Digital Probe, Patient Range > 30 KG, 24/cs

**GE MEDICAL NELLCOR REUSABLE PULSE OXIMETRY SENSORS & CABLES**

Hinged shell accommodates a wide range of digits. Soft vinyl padding holds probe in place. Latex free.

- 2008773-001 Nellcor Interface Cable, OxiSmart DOC10 Cable
- 407705-006 Nellcor DuraSensor Reusable Finger Probe
- 2006644-001 12 ft Nellcor Interface Cable OxiSmart
- 4033CAX Intermediate Cables For Nellcor Sensors
- DS100A Nellcore Reusable SpO2 Fingercip Sensor

**GE MEDICAL NONIN PULSE OXIMETER SENSORS & CABLES**

Latex free.

- 412926-001 8" Apex S Nonin Pulse Oximetry Interface Cable
- 412926-002 12" Apex S Nonin Pulse Oximetry Interface Cable
- 418497-001 8" ApexPro Nonin Pulse Oximetry Interface Cable
- 420970-001 Reusable Adult Articulated Finger Clip, 36", Patient Range > 30 Kg
- 420970-004 Reusable Ear Finger Clip, 18", Patient Range > 40 Kg

**GE MEDICAL OHMEDA PULSE OXIMETRY SENSORS & CABLES**

Latex free.

- 9587-002 Ohmeda FingerClip, Reusable Finger Probe, 12 ft Integrated Cable
- 9587-011 Ohmeda EarProbe, Adult Reusable Ear Probe, 8 ft Integrated Cable
- 9587-022 Ohmeda Flex II Probe, Reusable Wrap Probe, 12 ft Integrated Cable
- 9587-208 Ohmeda Adult/Pediatric Oxy-Tip Disposable Probes, Use w/9587-209 Interconnect Cable, 25/cs
- 9587-209 Ohmeda Interconnect Cable For Oxy-Tip Disposable Probes, 12 ft
- 9587-210 Ohmeda Neonatal Oxy-Tip Disposable Probes, Use w/9587-209 Interconnect Cable, 10/cs

**GE MEDICAL SENSORMEDICS PULSE OXIMETRY PROBES & CABLES**

Latex free.

|        |                                  |
|--------|----------------------------------|
| 008731 | Sensormedics Interconnect Cable  |
| 009084 | OMNI-Sat Finger Sensor           |
| 009120 | OMNI-Sat Ear Clip Sensor         |
| 009130 | Neonatal Sensor, Semi-Disposable |
| 009136 | Adult Sensor, Semi-Disposable    |

**MALLINCKRODT N-45 HANDHELD PULSE OXIMETER**

Compatible with the complete line of OxiMax® pulse oximetry sensors, so clinicians have the convenience of Nellcor's most advanced sensor technology with a lightweight, battery-operated, handheld device. Designed with an artifact-reducing algorithm to deliver highly accurate SpO<sub>2</sub> and pulse rate values, even during challenging monitoring conditions. The N-45 oximeter can be used for spot-checks or continuous monitoring. It offers the security of audible and visible alarms, and features 14 or 30 hour battery life, depending on battery type. The 24-hour downloadable trend memory enhances patient management and enables overnight oximetry studies in the hospital or home. The N-45 oximeter is compatible with ProFox software.

|     |   |
|-----|---|
| N45 | N-45 Handheld Pulse Oximeter, Includes: 9 V Lithium Battery, Cable Clip & Operator's Manual |
|-----|---|

**Accessories**

|               |                              |
|---------------|------------------------------|
| N45ACADPT     | N-45 110-230 V A/C Adapter   |
| N45BOOTCORAL  | N-45 Protective Boot, Coral  |
| N45BOOTPURPLE | N-45 Protective Boot, Purple |
| N45DTCBL      | N-45 Data Transfer Cable     |



**MALLINCKRODT N20 SERIES HANDHELD PULSE OXIMETERS**

Rugged, durable construction with protective boot. 32 to 37 hours long battery life. Easy-to-see day/night display. Three-year warranty. Comes with four C alkaline batteries, a finger clip sensor, a sensor extension cable, protective boot and moisture-resistant jacket. All models have a constant-pitch beep tone that sounds with each detected pulse. Displays oxygen saturation, pulse rate, vertical bar graph pulse amplitude indicator/battery meter, pulse search indicator, low battery indicator.

|        |   |
|--------|---|
| N20PA  | N-20PA Handheld Pulse Oximeter, Includes: 4-C Size Batteries, 1-Finger-Clip Sensor (DS-1000A), DEC-4 Sensor Extension Cable, Integral Thermal Printer w/12 rolls paper (901624), Protective Boot (024605), Moisture-Resistant Jacket & Manual |
| N20PAV | Pulse Oximeter w/continuous operation, Printer & Alarms (Vet Only)  |
| N20PE  | Pulse Oximeter w/printer, 1/cs  |

**Accessories**

|        |  |
|--------|--|
| 024603 | Protective boot for N-20, N-20E & N-20V            |
| 024605 | Protective Boot, N-20P Series                      |
| 024607 | Soft Carrying Case for N-20 Series Pulse Oximeters |
| 029043 | Dual Soft Carrying Case for Pulse Oximeters        |
| 901624 | Printer Paper, 12 rolls/box                        |




**MALLINCKRODT NPB40 HANDHELD PULSE OXIMETER**

Both compact and affordable. It is simple to operate and provides spot-check monitoring of oxygen saturation and pulse rate. The durable monitor is small enough to fit in your pocket. Designed for monitoring neonatal through adult patients in almost any setting, from hospital to home. Features an easy-to-read display and printing capabilities via optional Hewlett-Packard printer. Operates on AA batteries that provide approximately 19 hours of battery life.

|        |  |
|--------|--|
| NPB40  | NPB40 Handheld Pulse Oximeter, 1/bx      |
| NPB40V | NPB40 Handheld Pulse Oximeter (Vet Only) |

**Accessories**

|        |  |
|--------|--|
| 035183 | Soft Carrying Case, 1/bx               |
| 035328 | Protective Boot, 1/bx                  |
| 902114 | Printer, HP 82240B HP Infrared Printer |
| 902115 | Printer Paper, HP 82175, 6/bx          |




**MALLINCKRODT OXIMAX®/OXISENSOR® II ADHESIVE SENSORS**

|         |   |
|---------|---|
| ASP-3N  | Oxisensor II Neonatal Sensor, Assorted, 1/bg  |
| ASP3    | Oxisensor II Adult Sensor, Assorted, 1/bg     |
| ASP3P   | Oxisensor II Pediatric Sensor, Assorted, 1/bg |
| D20     | Oxisensor II Sensor, Pediatric, 24/cs         |
| D25     | Oxisensor II Sensor, Adult, 24/cs             |
| D25L    | Oxisensor II Adult Sensor w/36" cable, 24/cs  |
| I20     | Oxisensor II Sensor, Infant, 24/cs            |
| MAXA    | Oximax Sensor, Adult, 24/cs                   |
| MAXFAST | Oximax Sensor, Adult Forehead, 24/cs          |
| MAXI    | Oximax Sensor, Infant, 24/cs                  |
| MAXN    | Oximax Sensor, Neonate/Adult, 24/cs           |
| MAXP    | Oximax Sensor, Pediatric, 24/cs               |
| N25     | Oxisensor II Sensor, Neonate/Adult, 24/cs     |
| R15     | Oxisensor II Sensor, Adult Nasal, 24/cs       |

**Accessories**

|       |  |
|-------|--|
| DEC4  | Oximax 4 ft Extension Cable, 1/bx        |
| DEC8  | OxiMax 8 ft Sensor Extension Cable, 1/bx |
| DOC10 | Oximax 10 ft Sensor/Patient Cable        |



D-YS



DS100A-1

## MALLINCKRODT REUSABLE SENSORS

For spot checks or short-term monitoring. Non-sterile and easy-to-use design. Built-in shielding protects signal from electronic noise. High quality LEDs maximize tracking capabilities. Sensor head design optimizes signal and shields detector from ambient light.

|          |  |
|----------|--|
| A        | OxiCliq Sensor, Adult, 24/cs                                     |
| D-YS     | Dura-Y Multisite Sensor, 1/bx                                    |
| D-YSE    | Dura-Y Ear clip For Use with the Dura-Y Sensor                   |
| D-YSPD   | Pedi Check Pediatric Spot-Check Clip for the Dura-y Sensor, 1/bg |
| DS100A-1 | DuraSensor Adult Finger Clip Sensor, 1/bx                        |
| I        | OxiCliq Sensor, Infant, 24/cs                                    |
| N        | OxiCliq Sensor, Neonatal/Adult, 24/cs                            |
| OXI-A/N  | Oxiband Adult/Neonatal Sensor, 1/bx                              |
| OXI-P/I  | Oxiband Pediatric/Infant Sensor, 1/bx                            |
| P        | OxiCliq Sensor, Pediatric, 24/cs                                 |
| PDSLV    | Pedi Check Replacement Sleeve                                    |
| RS10     | Reflectance Sensor, 6/cs   |
| V-SAT    | Veterinary Sensor  |
| VSC-L    | Veterinary Sensor Clip, Large                                    |
| VSC-S    | Veterinary Sensor Clip, Small                                    |

### Accessories

|          |  |
|----------|--|
| ADH-A/N  | Adhesive Wrap for Reusable Sensors, Adult/Neonatal, 100/bx   |
| ADH-P/I  | Adhesive Wrap for Reusable Sensors, Pediatric Infant, 100/bx |
| FOAM A/N | Foam Wrap for Reusable Sensors, Adult/Neonatal, 12/pk        |
| FOAM P/I | Foam Wrap for Reusable Sensors, Pediatric/Infant, 12/pk      |
| OC-3     | OxiCliq Sensor Cable, 3 ft, 1/bx                             |
| POSEY    | Velcro Wrap for Reusable Sensors, 12/pk                      |



## PEPPER MEDICAL PULSE OXIMETRY SENSORS

Cost effective alternative pulse oximetry sensors that will plug ("Plug and Play") directly into your customer's existing monitors. Compatible with 8 different manufactures including AvalonAire, NPB, BCI, and Ohmeda. Available in 4 sizes: adult pediatric, infant, and neonate. 2 tape styles available: Microfoam, Transpore, and Medplast. FDA approved.

|      |   |
|------|---|
| 2211 | Disposable Pulse Oximetry Sensor, 24/bx |
|------|---|



## BCI CAPNOCHECK® II CAPNOGRAPH

Hand-held capnograph is a full-featured monitor providing ETCO<sub>2</sub>, Inspired CO<sub>2</sub>, Respiration Rate SpO<sub>2</sub>, and Heart Rate measurements. Ideal for use in and out of the hospital environment. A back-light feature on the display enables easy viewing in difficult lighting. Battery operated with an optional battery bypass offers flexibility for stationary and transport applications. Includes: 8400 Monitor, 1894 Operation Manual, 3044 Oximetry Finger Sensor, 8412 Protective Rubber Boot, 8214 Patient Attachment Kit includes the following: 8208 Capnocheck Filters (10), 8211 Capnocheck Sample Lines, 4 ft (10), 1100 Airway Adapter, straight without filter (10).

|        |  |
|--------|--|
| 8400   | Capnograph/Oximeter w/battery & charger  |
| 8400E  | Capnograph/Oximeter w/battery eliminator |
| 8400NP | Capnograph/Oximeter w/o Power Supply     |
| 8401   | Capnograph w/battery & charger           |
| 8401E  | Capnograph w/battery eliminator          |
| 8401NP | Capnograph w/o power supply              |

## BCI CAPNOCHECK® PLUS CAPNOGRAPH

Offers crisp, clear waveforms with large numerics for maximum visibility. High resolution vacuum fluorescent display. Simplified menu structure. Monitoring of both intubated and non-intubated patients. Flexibility to incorporate SPO<sub>2</sub> and FIO<sub>2</sub> technology. Patented moisture removal system. Four-channel analog output for polysomnography interfacing.

|          |   |
|----------|---|
| 9004-000 | Capnocheck® Plus Capnograph                           |
| 9004-001 | Capnocheck® Plus Capnograph/Oximeter                  |
| 9004-002 | Capnocheck® Plus Capnograph/FIO <sub>2</sub>          |
| 9004-003 | Capnocheck® Plus Capnograph/Oximeter/FIO <sub>2</sub> |

## BCI CAPNOCHECK® SLEEP CAPNOGRAPH/OXIMETER

Specially designed for sleep screenings and lab-based sleep studies. Its 4-channel analog output capability lets you select ETCO<sub>2</sub> numeric value, ETCO<sub>2</sub> waveform, inspired CO<sub>2</sub>, respiratory rate, heart rate, as well as SPO<sub>2</sub> numeric waveform displays. Patented Serial Autocorrelation (SAC) pulse oximetry technology provides reliable SPO<sub>2</sub> measurements.

|          |                                       |
|----------|---------------------------------------|
| 9004-050 | Capnocheck® Sleep Capnograph          |
| 9004-051 | Capnocheck® Sleep Capnograph/Oximeter |





**MALLINCKRODT HANDHELD CAPNOGRAPH/PULSE OXIMETER**

Delivers two highly advanced monitoring technologies in a convenient, portable monitor: Microstream® CO<sub>2</sub> technology and Nellcor's latest OxiMax pulse oximetry technology. Microstream technology and broad selection of breath sampling accessories ensure effective capnography monitoring for intubated and nonintubated patients alike. A low sample flow rate of 50 mL/min allows monitoring on the widest range of patients, from neonates to adults. OxiMax pulse oximetry technology provides exceptional performance, even during challenging monitoring conditions. The NPB-75 monitor is compatible with Nellcor's growing line of OxiMax sensors, including the Max-Fast™ Forehead Sensor and SoftCare™ Nonadhesive Sensors.

- NPB70 NPB70 Handheld Capnograph Includes: AC Adapter, Filterline Accessory Kit, Rechargeable NiMH Battery Pack, Operator's Manual & Quick Guide
- NPB75MAX Capnograph/Pulse Oximeter, Includes: AC Adapter, Accessory Kit, MAXPAC Oximetry Sensor Asst. Pack, DEC-4 Sensor Extension Cable, Rechargeable NiMH Battery Pack, Manual & Quick Guide

**Accessories**

- 047165 Carrying Case
- 047166 Extra Battery Pack Pouch
- 047184 Pole Mount Clamp
- 047203 External Battery Charger, 110-250 VAC
- 048127 Communication Adapter Kit: Printer/Computer Interface



**BURDICK EK-10 SINGLE-CHANNEL ELECTROCARDIOGRAPH**

Combines the cost advantages of a single-channel electrocardiograph with the increased productivity of automatic operation. Features a thermal array printer, liquid crystal display, 'soft-touch' membrane pad and 50 ECG acquisition/print on battery (optional). Compact and portable.

- 101070 EK-10 Single-Channel Electrocardiograph
- Accessories**
- 007175 1V Output Cable
  - 086027 Service Manual For EK-10
  - 086063 Operators Manual For EK-10

**BURDICK SINGLE CHANNEL ECG MOUNTS**

- 007036 ECG Mounts, Elite Folder, Pressure-Sensitive, Full Information, 100/bx, 5 bx/cs
- 007037 ECG Mounts, Card, Pressure-Sensitive, EK10 or Abridged Elite, 100/bx, 5 bx/cs
- 007857 ECG Mounts, Card, Pressure-Sensitive, 12-Lead, 100/bx, 5 bx/cs
- 007909 ECG Mounts, Card, Pressure-Sensitive, 3 Sec., 12-Lead, 100/bx, 5 bx/cs
- 007911 ECG Mounts, Folder, Pressure-Sensitive, 3 Sec., 12-Lead, 100/bx, 5 bx/cs



**BURDICK ELITE II INTERPRETIVE ECG RECORDER**

Simultaneously acquires 10 seconds of ECG data for all 12 leads. Prints interpretive reports in 4-channel format. Menu displays guide user through all functions. Stores up to 30 ECGs in memory. 11" x 13" x 3 1/4". 120/240 VAC, 50/60 Hz. Includes 10 lead accessory group and 12.5-V nickel cadmium battery. Compact and portable.

- 101158 Elite II, Single Channel Interpretative Electrocardiograph
- Accessories**
- 007697 Elite II Factory Installed Fax Option
  - 086136 Elite II Operators Manual
  - 086175 Elite II Service Manual
  - 086245 ECG Interpretation Criteria, GRI



010-1639-00



A30-1DIL01

**BURDICK ATRIA™ 3000 RESTING ECG SYSTEM**

Simple to use, compact and reliable. The lightweight and compact device operates either manually or automatically, printing 3-, 4-, 6- or 12-channel ECG tracings. The combination of a standard QWERTY keyboard and easy-to-navigate menu makes it easy to use. Features advanced functionality, including a stable baseline filter for greater accuracy, pacemaker algorithm, user-definable QTc calculations and Burdick's exclusive AccuPrint™ functionality, which restricts ECG printing if lead wires are not functioning properly. Variety of upgrade paths.

- A30-1DIL01 Atria™ 3000 ECG w/interpretation, extended storage 40 patients & no communication
- A30-1DKL01 Atria™ 3000 ECG w/measurements, 40 patient storage & no communication

**Accessories**

- 010-1639-00 Cart For The Atria 3000 ECG System
- 040-1457-00 Atria 3000 Kit Upgrade Interpretation
- 040-1458-00 Atria 3000 Kit Upgrade Extra Storage
- 040-1474-00 Atria 3000 Kit Upgrade Measurements



**BURDICK ECLIPSE™ 850 ELECTROCARDIOGRAPHS**

Multi-channel electrocardiograph with preview screen. Thermal printer produces crisp accurate 3, 4, 6 or 12 channel ECG reports and continuous rhythm reports on 8 1/2" paper. Acquires and prints on battery and/or AC power. Optional internal modem with fax and SCP-ECG communications. Includes: 10-Lead patient cable with 90 day warranty, power cord, one pad assurance 50 paper, cardiosens/ultra sensors & adapter clips, manual, handbook and video. Three year warranty on parts and labor.

- 92300-1AK09 Interpretive Eclipse 850i w/storage for 40 ECG
- 92300-1AL09 Non-Interpretive Eclipse 850 w/storage for 40 ECG
- 92300-1BKP09 Interpretive Eclipse 850i w/SCP-ECG & fax communications, internal modem & storage for 1 ECG
- 92300-1BLP09 Non-Interpretive Eclipse 850 w/SCP-ECG & fax communications, internal modem & storage for 1 ECG

**BURDICK EQUIPMENT CARRYING CASE**

- 007039 Soft Shelled Carrying Case for ELITE, EK-10 & Presto Spirometer #100231
- 007160 Soft Shell Carry Case For Use w/Extol LE & LEI
- 007732 Carrying Case For Eclipse 4
- 007771 Carry Case For Use w/Eclipse 8



**BURDICK UNIVERSAL ECG**

Portable and easy to use. Patient data is instantly and accurately obtained, viewing up to 12 channels. Unlimited data analysis including the ability to reanalyze data if system parameters are modified. Compatible with your favorite off-the-shelf PC, laptop, tablet PC or Compaq iPAQ.

- 92309-01 Universal ECG Non-Interpretive, AHA, English
- 92309-02 Universal ECG Interpretive, AHA, English



**BURDICK ASSURANCE 50™ ECG RECORDING PAPER**

Guaranteed to hold a trace for 50 years. Single and multi-channel formats. Special thermal paper resists effects of alcohol, sunlight, plastics and other common contaminants that degrades image quality on standard thermal papers.

- 716-0233-00 ECG Recording Paper, Red Grid Use w/ EK10, Elite, Elite II, Medic 5, SUB FOR 007957 & 007975, 100rl/cs
- 716-0234-00 ECG Recording Paper, Z-Fold Pad Red Grid, Perforated No Header use w/ Eclipse LE, LEII, E550, E560, E600, ExTOL 200 sheets/pd, 10pd/cs
- 716-0235-00 ECG Recording Paper, Z-Fold Pad Red Grid, Perforated w/blank header use w/ Eclipse8, 800, 300 sheets/pd, 10pd/cs
- 716-0237-00 EKG Thermal Paper, 200 sheet/pad, 10 pad/cs
- 716-0239-00 ECG Recording Paper, Z-Fold Pad Red Grid, Perforated w/blank header use w/Eclipse 850, 300 sheets/pd, 10pd/cs

**BURDICK HEARTLINE THERMAL PAPER**

- 007868 Chemical/Thermal Heavyweight Paper, Perf, Multi-Channel w/blank header, Matrix, Red Grid/Black Trace, 200 shts/pd, 10 pd/cs
- 007940 Single Channel Thermal Paper, Grid, 45mm, Red Grid/Black Trace, 100'/roll, 100 rl/cs
- 007941 Thermal Paper, Orange Grid, For Use w/EK4, EK5/5A, Ek8, E100, E200, 100rl/cs
- 007944 Thermal Paper, Orange Grid, Perforated Header, For Use w/EK670 S/N 6000 & Above, EK7 S/N 7000 & Above, E310, E320, 20rl/cs
- 007958 Single Channel Thermal Paper, Dot Matrix, 45mm, Red Grid/Black Trace, 150 ft/rl, 80 rl/cs (Sub for 007957)
- 007966 Thermal Paper, Z-Fold Pad, Red Grid, Perforated w/header, Use w/E350, 200 sheets/pd, 10pads/cs
- 007971 Thermal Paper, Red Grid, Perforated w/blank header, Use w/Eclipse 4/400, 200 sheets/pd, 40pd/cs
- 007979 Chemical/Thermal Paper, Perf, Multi-Channel, No Header, 8.5" x 5.5", Red Grid/Black Trace, 300 shts/pd, 10 pd/cs
- 007983 Chemical/Thermal Paper, Perf, Multi-Channel w/blank header, Matrix, Red Grid/Black Trace, 200 shts/pd, 10 pd/cs
- 007988 Chemical/Thermal Paper, Reinforced Perf, Multi-Channel, Blank Header, 8.5" x 5.5", Red Grid/Black Trace, 300/pd, 10 pd/cs



92400-01



B24-1AMP01

**BURDICK QUEST® EXERCISE STRESS SYSTEM**

Features easy, intuitive touch screen lets you focus on the patient. Performs both resting 12-lead ECGs and stress tests on one system. Accurate ECG data acquired at the patient eliminating non-cardiac electrical noise. Flexible programming means you choose the protocol. Generous storage, enhanced online data editing and several data management options. Full-range reporting offers over 15 different report configurations. Optional treadmill not included, See TM55 or TM65. Optional T-Wave Alternans helps predict Sudden Cardiac Death (SCD).

- 92400-01 TM55 Treadmill Option For Quest Exercise Stress System
- 214-1022-00 Stress Belt (Black) w/pouch and quick release buckle, 48"
- 214-1023-00 Belt Extender w/quick release buckle no pouch, 24"
- B24-1AMP01 Quest Exercise stress system, English language, with AHA accessories kit, TM55 treadmill and network communications option (120V), ea
- B24-1AP Quest System w/o treadmill. Includes standard US Thermal Paper and AHA Patient Cable, English Language (120V)
- B24-1AP01 Quest Exercise Stress System, English Language w/AHA accessories kit and TM55 Treadmill
- B24-1APT Quest System w/T-Wave Option w/o treadmill. Includes Standard US Thermal Paper and AHA Patient Cable, English Language and 120V
- B24-1APT01 Quest Exercise Stress System, English Language w/AHA accessories kit, TM55 Treadmill and T-Wave Option

**Accessories**

- 004115 Emergency Stop Switch For T600
- 007772 Patient Module
- 004219 Patient Data Acquisition Cable

**BURDICK ELECTROCARDIOGRAPH ACCESSORIES**

- 007001 Patient Cable, Resting, 10 Lead, Non-replaceable (E500)
- 007026 Cable, E600/ELITE/E500tw To Modem
- 007028 Interconnect Cable, E600/E500tw/ELITE To E-Scribe
- 007053 Junction Box For Modem On Stress Cart
- 007059 Interconnect Cable, Cart To Cart
- 007153 Plastic Cardiogram Mount Inserter
- 007157 Adult Limb Plate Sensor, 4 electrodes/pk
- 007158 Limb Strap For Sensor, Single Hole, 18", Adult
- 007159 Limb Strap For Sensor, 2 Holes, 20", Adult
- 007161 Infant Limb Plate Sensor, 4/pk
- 007162 Limb Strap For Sensor, Infant, Single Hole
- 007163 Welsh Suction Cup Sensor, Self-Retaining, Adult
- 007211 Patient Cable, Resting, 5 Lead, Non-replaceable (EK-5A)
- 007264 Alligator (Veterinary) Adapter Clip
- 007364 Welsh Suction Cup Sensor, Spring-Clip, Infant
- 007374 Stress Cable Belt, Black Buckle, 5 ft
- 007463 Lead Adapter For Snap End
- 007469 Lead Adapter For Grabber End
- 007517 Patient Cable, Stress, 10 Lead, Non-replaceable (ELITE, EK-10,E350, E350i, E550, E560)
- 007519 Patient Cable, Stress, Non-replaceable, EK-8, E100, E200, E310, E320, E600
- 007526 Stress Belt, 4 ft, White, Velcro
- 007550 Adult Limb Plate Sensor, 4mm, Set of 4
- 007557 How To Operate The Burdick EK-10
- 007558 A Guide To Good ECG Technique
- 007597 External Modem w/cable
- 007655 How To Operate The Burdick ELITE
- 007656 How To Operate The Burdick E350i
- 007679 How To Operate The Burdick E350
- 007690 Welsh Suction Cup Sensor 4mm, Set of 6
- 007704 Patient Cable, Resting, 10 Lead, Replaceable, Eclipse 4, 400, 8, 800, 850, Plus
- 007705 Patient Cable, Resting, 10 Lead, Replaceable, ELITE, E50, E350i, E550, E560, Eclipse LE/Lei
- 007706 Patient Cable, Stress, 10 Lead, Replaceable, E350, E350i, E560, Eclipse LE/Lei
- 007710 Patient Cable, Three Lead Wires (Replacement Kit) For 007704, 007704, Includes: one 20", one 30", one 40" Lead Wire & 10 Lead Identifiers
- 007711 Patient Cable, Three Lead Wires (Replacement Kit) For 007706, 007725, Includes: One 20", One 30", One 20" Lead Wire & 10 Lead Identifiers
- 007725 Patient Cable, Stress, 10 Lead, Replaceable, Eclipse 4, 400, 8, 800, 850, Plus
- 007726 Patient Cable, Three Lead Wires (Replacement Kit) For 007727, 007709, Includes: One 20", One 30", One 20" Lead Wire & 10 Lead Identifiers
- 007728 Patient Cable, Three Lead Wires (Replacement Kit) For 007707, 007708, Includes: One 20", One 30", One 40" Lead Wire & 10 Lead Identifiers
- 007751 EKS-95 Electronic Patient
- 007797 Lead Adapter, 0.158" (4mm) diameter to snap end
- 007861 Retaining Strap w/weight for 107155 Set
- 008003 How To Operate The Burdick Eclipse
- 010-1608-00 ECG Replacement Cables Interpretive
- 010-1609-00 ECG Replacement Cables Non-Interpretive
- 012-0700-00 Patient Cable, Resting, 10 Lead, Non-replaceable (E350, E350i, E550, E560)
- 012-0701-00 Patient Cable, Resting, Leads, Non-replaceable, (EK-8, E100, E200, E310, E320, E600)
- 012-0702-00 Patient Cable, Resting, 5 Lead, Non-Replaceable, EK10, Elite
- 012-0844-00 Nonreplaceable Lead ECG Cable, 10 Lead, AHA, Combs Located Between Lead Wires to Decrease Tangling, Use With Atria, Eclipse Premier & Eclipse 850
- 013-0150-00 30" Replacement Stress Leads, Green and Brown
- 047224 3-Lead Patient Cable
- 047236 5-Lead Patient Cable
- 047246 DC Input Cable (1V to 1mV)
- 047254 Patient Cable, Twist-Lock Connector For PA-101 Preamp
- 086061 Electrocardiography: A Better Way Workbook
- 097152 8 ft Power Cord- 3 Cond For EK-5A
- 097159 Power Cord For EK-7, EK670, DC-190, 10 ft
- 107067 Welsh Sensor Kit, Infant Plates/Straps & Bulbs
- 107068 Welsh Sensor Kit, Adult Plates/Straps, Infant Bulbs
- 107546 4mm Sensor Kit, Adult Plates/Straps, Infant Bulbs
- 147234 3-Lead Wire Kit For 047224 & 047254 Cables Only
- 147235 5-Lead Wire Kit For 047236 Cable Only, Snap
- 700-0277-00 USB to Serial Cable Adapter For Use w/Eclipse, CardioVault, Pyramid & AED Products
- 860062 Retaining Strap Only For 107155 Set
- 860403 Bulb Only For Welsh Infant Sensor (007364)
- 860404 Bulb Only For Welsh Adult Sensor (007163)
- 862278 Battery Pack For ELITE, EK-10



**GE MEDICAL MAC 1200 STANDARD RESTING ECG**

A 12-lead electrocardiograph with a full-size paper, alphanumeric keypad and waveform display. Optional Marquette® 12SL™ ECG analysis program, data storage, and bi-directional data transmission capability. AC line or battery operation from NiCd battery. All models listed include: Printed measurements and starter kit (Silver Mactrode® electrodes, grabber adapters, paper, patient cable, power cord).

- MAC1200-AAX-CXX-XXX      MAC 1200 'Plus' Interpretive Electrocardiograph w/12SL ECG Analysis Program & Waveform Display

**Accessories**

- 901168-001      MAC 1200 Office Cart

**GE MEDICAL DEFIBRILLATOR SUPPLIES**

- 009/00310/24V      Battery Pack, NICAD For 360 Defib
- 009/09956/USA      Charger Lead-Acid Battery USA
- 099/01270/000      Cable Vehicle Charger For Responder
- 210/34223/000      Battery For Data Event Card
- 414034-002      Data Event Card, 80 Minutes
- 850/40156/005      Adult Defibrillator Pads Responder, 5 pr/cs
- 850/40156/025      Adult Defibrillator Pads Responder, 25 pr/cs
- 850/40518/005      Pediatric Defibrillator Pads Responder, 5 pr/cs
- 850/40518/025      Pediatric Defibrillator Pads Responder, 25 pr/cs
- 900634-001      Custom Setup Card

**GE MEDICAL ECG CABLES & LEADWIRES**

- 107230      Leadwire Set, Snap, 3 LD, Gray, 100cm AHA, 3/set
- 107231      Leadwire Set, Grab, 3 LD, Gray, 100cm AHA, 3/set
- 107235      Leadwire Set, Grab, 6 LD, Gray, 100/130cm AHA, 6/set
- 107289      Leadwire Set, Grab, 6 LD, Gray, 100 AHA, 6/set
- 107290      Leadwire Set, Snap, 6 LD, Gray, 100 AHA, 6/set
- 107314      Leadwire Set, Grab, 3 LD, Gray, Pediatric, AHA, 3/set
- 107320      ECG Cable, 6 LD, AHA
- 107326      ECG Cable, 3 LD, AHA
- 107328      Leadwire Set, Snap, 3 LD, AHA, 3/set
- 107331      ECG Cable, 3 LD w/ESU Filter For OR, AHA
- 107333      ECG Cable, 6 LD w/ESU Filter For OR, AHA
- 107384      ECG Cable, 3 LD, Safety, AHA



**GE MEDICAL ECG & DEFIB RECORDING PAPERS**

Consolidate recording paper purchases. Reduce time and expense to generate multiple purchase orders.

- 9801-001 Thermal Paper, Spacelabs Equivalent 307438-001, All White, 1.97" x 100 ft, For Spacelabs ECG Systems PC Bedside, 90600 Series, 90731, 90323, 10/pk, 10 pk/cs
- 9801-002 Thermal Paper, Quinton Equivalent 15022-001, All White, Black Trace, 8.5" x 11", Queue Hole, Z-Fold, For Quinton Stress Systems Q3000, Q4000, 216 sheets/pk, 20 pk/cs
- 9801-003 Thermal Paper, Quinton Equivalent 16211-002, All White, Queue Mark, 8.5" x 11", Z-Fold, For Quinton Stress Systems Q5000 & QCATH, 200 sheets/pk, 20 pk/cs
- 9801-004 Thermal Paper, Hewlett-Packard Equivalent M1707A, Red Grid, 8.5" x 11", Z-Fold, For HP M1700A, M1701A Pagewriter XL Systems, 200 sheets/pk, 8 pk/cs
- 9801-010 Thermal Paper, Hewlett-Packard Equivalent 40453A/B, Single Channel, Black Grid, 1.97" x 100 ft, For HP Defibrillators, 10/pk, 10 pk/cs
- 9801-016 Thermal Paper, Burdick Equivalent 007961, Red Grid, Black Trace, 8.5" x 11", Z-Fold, For Burdick E500TW, E350, E-600 ECG Systems, 200 sheets/pk, 10 pk/cs
- 9801-028 Thermal Paper, General Scanning Equivalent AR-200, All White, Black Grid, 8.5" x 11", Z-Fold, Queue Marks, For General Scanning Systems AR-200, 70206, 300 sheets/pk, 10 pk/cs
- 9801-057 Thermal Paper, Hewlett-Packard Equivalent M2481A, Red Grid, Black Trace, 8.5" x 5.71", 8 pk/cs
- 9801-059 Thermal Paper, Spacelabs Equivalent 006-0196-001, Black Trace, 1.97" x 97 ft, For Spacelabs Systems Using 006-0196-001 Paper, 10/bx, 10 bx/cs

**GE MEDICAL LASER PRINTER PAPERS**


Quality paper stock minimizes wear on laser printers. Save time printing traces on preprinted grid paper. Red grids improve trace legibility.

- 9402-053 Laser Printer Paper, Red Grid w/blank interpretation area, 8.5" x 11", 500 sheets/pk, 5 pk/cs
- 9402-054 Laser Printer Paper, Full Page Red Grid, 8.5" x 11", 500 sheets/pk, 5 pk/cs

**GE MEDICAL LEG PLATE ADAPTER**

An alternative to the disposable attachment pads (2464AAO). Smooth surface against patient's skin. Easy to clean or sterilize. Latex free.

- 1594AAO 2" x 2" Leg Plate Strap Adapter



**GE MEDICAL ULTRA-ARCHIVIST® THERMAL RECORDING PAPERS**

100 year warranty for trace legibility. Eliminates need for special storage precautions. Heavyweight paper to minimize potential for paper jams. Clear, precise traces.

- 2007854-001 Ultra-Archivist Thermal Paper, Red Grid w/blank interpretation area, 8.5" x 11", Z-Fold, 250 sheets/pk, 8 pk/cs
- 2007854-002 Ultra-Archivist Thermal Paper, Red Grid A-4 Modified, 8.26" x 11", Z-Fold, 250 sheets/pk, 8 pk/cs
- 2007854-003 Ultra-Archivist Thermal Paper, Full Page Red Grid, 8.5" x 11", Z-Fold, 250 sheets/pk, 8 pk/cs

**GE MEDICAL UNIVERSAL THERMAL RECORDING PAPERS**

Bright white paper offers excellent trace legibility. Best choice for MUSE customers. Electronic record reduces need for more expensive, long-term trace retention.

- 9902-020 Universal Thermal Paper, Full Page Red Grid, 8.5" x 11", Z-Fold, Queue Hole, 300 sheets/pk, 8 pk/cs
- 9902-021 Universal Thermal Paper, Full Page Red Grid, A-4 Modified, 8.26" x 11", Z-Fold, Queue Hole, 300 sheets/pk, 8 pk/cs
- 9902-024 Universal Thermal Paper, Red Grid w/blank interpretation area, 8.5" x 11", Z-Fold, Queue Hole, 300 sheets/pk, 8 pk/cs
- 9902-024GY Universal Thermal Paper, Gray Grid w/blank interpretation area, 8.5" x 11", Z-Fold, Queue Hole, 300 sheets/pk, 8 pk/cs
- 9902-030 Universal Thermal Paper, All White, Black Queue Mark, 8.5" x 11", Z-Fold, 300 sheets/pk, 8 pk/cs
- 9902-031 Universal Thermal Paper, Green Grid w/blank interpretation area, 8.5" x 11", Z-Fold, Queue Hole, 300 sheets/pk, 8 pk/cs
- 9902-324 Universal Thermal Paper, Red Grid w/blank interpretation area, 3-Hole Punch, 8.5" x 11", Z-Fold, Queue Hole, 300 sheets/pk, 8 pk/cs
- 9902-920 Universal Thermal Paper, Full Page Red Grid, 8.5" x 11", Z-Fold, Queue Hole, 300 sheets/pk, 10 pk/cs
- 9902-924 Universal Thermal Paper, Red Grid w/blank interpretation area, 8.5" x 11", Z-Fold, Queue Hole, 300 sheets/pk, 10 pk/cs
- 9902-931 Universal Thermal Paper, Green Grid w/blank interpretation area, 8.5" x 11", Z-Fold, Queue Hole, 300 sheets/pk, 10 pk/cs

**GE MEDICAL PREMIUM THERMAL RECORDING PAPERS**

Exclusive chemical composition resists trace fading for 5-7 years following recommended storage procedures. Clear, precise trace. Designed for reading ease.

- 22616606 Premium Thermal Paper, Full Page Brown Grid, Queue Mark, 8.5" x 11", Z-Fold, 150 sheets/pk, 10 pk/cs
- 9402-060 Premium Thermal Paper, Full Page Red Grid, Queue Hole, 8.5" x 11", Z-Fold, 150 sheets/pk, 16 pk/cs
- 9402-061 Premium Thermal Paper, Red Grid w/blank interpretation area, Queue Hole, 8.5" x 11", Z-Fold, 150 sheets/pk, 16 pk/cs
- 22613001 Premium Thermal Paper, 55mm x 40M
- 22616802 Premium Thermal Paper, Brown Dot Grid, 3.54" x 115 ft, 10/cs
- 385-1000-015 Premium Thermal Paper, Brown Grid, 8.5" x 11", Z-Fold, 400 sheets/pk, 10 pk/cs
- 9402-018 Premium Thermal Paper, Red Grid, 4" DDW 4.25" x 150 ft, 32/cs
- 9402-020 Premium Thermal Paper, Full Page Red Grid, 8.5" x 11", Z-Fold, Queue Hole, 300 sheets/pk, 8 pk/cs
- 9402-021 Premium Thermal Paper, Full Page Red Grid, A-4 Modified, 8.26" x 11", Z-Fold, Queue Hole, 300 sheets/pk, 8 pk/cs
- 9402-023 Premium Thermal Paper, Red Grid, 4.25" x 75 ft, 2/pk, 9 pk/cs
- 9402-024 Premium Thermal Paper, Red Grid w/blank interpretation area, 8.5" x 11", Z-Fold, Queue Hole, 300 sheets/pk, 8 pk/cs
- 9402-028 Premium Thermal Paper, All White, 2" x 100 ft, 48/cs
- 9402-029 Premium Thermal Paper, All White, Queue Hole, 8.5" x 11", Z-Fold, 300 sheets/pk, 8 pk/cs
- 9402-030 Premium Thermal Paper, All White, Black Queue Mark, 8.5" x 11", Z-Fold, 300 sheets/pk, 8 pk/cs
- 9402-031 Premium Thermal Paper, Green Grid w/blank interpretation area, 8.5" x 11", Z-Fold, Queue Hole, 300 sheets/pk, 8 pk/cs
- 9402-046 Premium Thermal Paper, 2" x 100 ft, 48/cs
- 9402-050 Premium Thermal Paper, Red Grid w/blank interpretation area, 8.5" x 11", Z-Fold, Queue Hole, 500 sheets/pk, 5 pk/cs
- 9402-324 Premium Thermal Paper, Red Grid w/blank interpretation area, 3-Hole Punch, 8.5" x 11", Z-Fold, Queue Hole, 300 sheets/pk, 8 pk/cs
- 9472-103 Premium Thermal Paper, All White, 2" x 100 ft, 4/bx

**GE MEDICAL CARDIOLOGY LEADWIRE ADAPTERS**

Latex free. 90 day warranty.

|          |  |
|----------|--|
| 9490-210 | Set Adapter GE Wide-Mouth Universal Alligator/Banana Socket/Snap, 10/set |
| 38401804 | Set Leadwire Adapters Grabber/Banana, 10/set                             |
| 9459-901 | Set Adapters Grabber/Banana, 6" Wire, 10/set                             |
| 9490-105 | Set Adapter Grey/Black Alligator/Banana Socket, 4mm, 10/set              |
| 9490-106 | Set Adapter Grey/Blue Alligator/Banana Pin, 3mm, 10/set                  |
| 9490-111 | Set Adapter Clip Black Alligator/Banana Socket, 4mm, 10/set              |
| 9490-112 | Set Adapter Clip Black Alligator/Banana Pin, 3mm, 10/set                 |
| 9490-114 | Set Adapter Clip Black Alligator/Banana Socket, 4mm, 14/set              |
| 9490-115 | Set Adapter Grey/Black Alligator/Banana Socket, 4mm, 14/set              |

**GE MEDICAL CARDIOSMART, CARDIOSYS, MAC CABLES & LEADWIRES**

Easy to clean, wipe with damp cloth. Built for durability. Latex free. One year warranty. Coaxial cable reduces potential noise.

|          |  |
|----------|--|
| 22341809 | 10 ft ECG Cable For MAC 1100/1200/1200ST, AHA, Uses CAM 14 Leadwires           |
| 22341802 | 7 ft ECG Cable for CardioSmart, MAC 1000, CardioSys, AHA Uses CAM 14 Leadwires |
| 38401742 | Leadwire Set, Banana/Cardiosmart, 0.156" Banana, AHA, 10/set                   |



**WELCH ALLYN CP 10 & CP 20 ECG UNITS**

Includes the following accessory items standard: 10-lead resting AHA patient cable; 500 disposable tab electrodes; 10 resting electrode adaptors; physician's guide to the analysis program (interpretive systems) on CD; in-servie video on CD; 200 sheets of paper; operator manual on CD; power cord.

|           |   |
|-----------|---|
| 11102-200 | CP 10 ECG Non-Interpretive 12 Lead Multichannel ECG w/ECG measurements, 8.25" x 11" Paper, Alphanumeric Keyboard & AHA Patient Cable  |
| 11102-500 | CP 10 Interpretive 12-Lead Multichannel ECG, 8.25" x 11" Paper, Alphanumeric Keyboard & AHA Patient Cable   |
| 11103-200 | CP 20 Non-Interpretive 12-Lead Multichannel ECG w/ECG measurements, ECG Preview Display, 8.25" x 11" Paper, Alphanumeric Keyboard & AHA Patient Cable                       |
| 11103-210 | CP 20 Non-Interpretive 12-Lead Multichannel ECG w/ECG measurements, Communication Module, ECG Preview Display, 8.25" x 11" Paper, Alphanumeric Keyboard & AHA Patient Cable |
| 11103-500 | CP 20 Interpretive 12-Lead Multichannel ECG, ECG Preview Display, 8.25" x 11" Paper, Alphanumeric Keyboard & AHA Patient Cable  |
| 11103-510 | CP 20 Interpretive 12-Lead Multichannel ECG w/communication module, ECG Preview Display, 8.25" x 11" Paper, Alphanumeric Keyboard & AHA Patient Cable                       |

**Accessories**

|            |   |
|------------|---|
| 80185-0000 | Data Transfer Cable   |
| 96152-0000 | CP 10 Interpretation Upgrade (must provide system serial #) |

**WELCH ALLYN ECG ACCESSORIES**

|             |   |
|-------------|---|
| 08250-0000  | CP 10/CP 20 Cart w/cable arm                                |
| 08265-0000  | ECG Utility Cart  |
| 08285-0000  | ECG Carrying Case   |
| 25254-0000  | Cable Separators, 10/set                                    |
| 45008-0000  | Resting Tab Electrodes, 1000/bx                             |
| 58548-0000  | Pinch Clip Adapters, 10/set                                 |
| 65159-0000  | Cable Arm   |
| 65160-0000  | CP 10/CP 20 Cart w/o cable arm                              |
| 80141-0000  | Hospital Grade Power Cable                                  |
| 9293-032-60 | 10-Lead Patient Cable, Resting ECG, AHA CP 10/CP 20         |
| 9293-032-61 | 10-Lead Patient Cable, Resting ECG, IEC CP 10/CP 20         |
| 94016-0000  | 8.25" x 11" Paper For CP 10 & CP 20, 200 sheets/pk, 5 pk/cs |
| 58581-0000  | ACCESSORY: Universal Adapter Clip, 10/pk                    |



**MIDMARK IQMARK™ DIGITAL ECG (SERIAL VERSION)**

Includes: IQmark™ Digital ECG Acquisition Module (serial version), IQmark™ Digital ECG Software Package, IQmark™ Diagnostic Workstation - Patient Management Software, 10-lead patient cable, 10-pack universal clips, 100-pack ECG electrodes, 4-pack AA batteries, mouse pad, operation manuals and quick reference guide.

|            |                                      |
|------------|--------------------------------------|
| 4-000-0000 | IQmark™ Digital ECG (Serial Version) |
|------------|--------------------------------------|

**MIDMARK IQMARK™ ACCESSORIES**

|            |  |
|------------|--|
| 1-370-0004 | USB External Compact Flash Card Reader |
| 3-370-0005 | PCMCIA Compact Flash Adapter           |
| 1-136-0005 | USB to Serial Port Adapter             |

**MIDMARK IQMARK™ DIGITAL RECORDER KIT**

Kit includes: IQmark™ Digital Recorder, lead wire set, recorder pouch and belt, compact flash card and operation manual.

|            |   |
|------------|---|
| 4-000-0170 | IQmark™ Digital Recorder Kit For EZ & Advanced Holter |
|------------|---|

**MIDMARK IQMARK™ DIGITAL RECORDER KIT - 3 CHANNEL**

Kit includes: IQmark™ Digital Recorder, lead wire set, recorder pouch and belt, 16MB PCMCIA flash card and operation manual.

|            |   |
|------------|---|
| 4-000-0105 | IQmark™ Digital Recorder Kit - 3 Channel For 7700/8800 Holter |
|------------|---|

**MIDMARK IQMARK™ EZ HOLTER KIT WITH DIGITAL RECORDER**

Kit includes: IQmark™ Digital Recorder Kit (lead wire set, recorder pouch and belt, compact flash card, operation manual), IQmark™ EZ Holter Software Package, IQmark™ Diagnostic Workstation - Patient Management Software, compact flash card reader, PCMCIA compact flash adapter, sentinel data protection key, (4) holter prep kits, mouse pad, operation manual and quick reference guide..

|            |  |
|------------|--|
| 4-000-0150 | IQmark™ EZ Holter Kit w/digital recorder |
|------------|--|



**MIDMARK IQMARK™ ADVANCED HOLTER KIT WITH RECORDER & PACEMAKER OPTION**

Includes: leadwire set, recorder pouch and belt, compact flash card, operation manual, IQmark™ Advanced Holter with Pacemaker Analysis Software Package, IQmark™ Diagnostic Workstation - Patient Management Software, Compact Flash Card Reader, PCMCIA Compact Flash Adapter, Sentinel Data Protection Key for Pacemaker (parallel port key shipped standard, USB key available upon request), 4 Holter Prep Kits, mouse pad, and quick reference guide.

4-000-0158 IQmark™ Advanced Holter Kit w/recorder & Pacemaker Option



**BURDICK BASELINE STRESS PREP KIT**

Each kit contains: 10 Stress Electrodes, Foam, 1 Skin Scrub, 1 Instruction Sheet, 3 Alcohol Prep Pads, 1 Bacitracin Ointment, 1 Moist Wipe, 1 Surgical Prep Razor.

015-0629-00 Blue Max Stress Prep Kit, 30 kit/cs  
097008 Baseline Stress Prep Kit, 30/cs



**BURDICK TM55 AND TM65 TREADMILLS - DIGITALLY CONTROLLED**

Thousands of clinicians worldwide rely on the TM55 and TM65 to deliver the accuracy and reliability required for exercise stress tests. Fewer moving parts reduces wear and tear and maximizes system reliability. Controlled "soft start" for patient safety. Digitally controlled AC-drive system provides high torque for accurate speed and grade, regardless of weight. Quiet Operation through entire speed range eliminates office disruption. Fully integrated into the Burdick Quest exercise stress system.

TM55HXXRBE TM55, 230V US & Canada  
TM55LXXRBE TM55, 120V



**BURDICK TM55 & TM 65 TREADMILLS**

Proven design maximizes system reliability. Controlled soft start improves patient safety. Accurate speed and grade provides better test results. Quiet operation calms patients and eliminates office disruption. Fully integrated into the Burdick Quest® exercise stress system

000488-010 TM65 Treadmill, 120V  
000488-011 TM65 Treadmill, 230V US/Canada

**Accessories**

032758-001 Short Handrail  
032759-001 Long Handrail  
037288-001 Emergency Stop

**BURDICK TWAVE ALTERNANS OPTION FOR QUEST**

015-0631-00 PM-3 Alternans Lead Wire Set Used For Microvolt T-Wave Alternans Test on Quest Exercise Stress Systems AHA  
015-0635-00 Micro-V Alternans Sensors, 7/pk, 20 pk/bx  
015-0637-00 PM-3 Patient Module & System Cable For Microvolt T-Waves Alternans Test Option  
040-1219-00 Microvolt T-Wave Alternans Test S/W & H/W Upgrade Kit For Quest Systems AHA  
92400-T Microvolt T-Wave Alternans Option AHA

**BURDICK HEARTWORKS™ DIAGNOSTIC WORKSTATION**

A suite of diagnostic devices conveniently integrated on one PC: Vision Holter, Vision 5L Recorder, ABP Monitor and software, Universal ECG, SpiroCard and Office Medic software.

D20-1ADH HeartWorks™ Diagnostic Workstation with Laptop Computer and Printer  
D20-1D HeartWorks™ Diagnostic Workstation w/out Laptop and Printer

**BURDICK HEARTSTRIDE EXERCISE STRESS SYSTEM**

Designed to meet the unique needs of primary care physicians looking to add cardiac stress testing capabilities to their practices. With its simple, graphical user controls, plus optional network connectivity capability (including ability to interface with EMRs), integrating cardiac stress testing into the office has never been simpler. Designed to deliver value, the PC-based HeartStride allows smaller offices to leverage their PC hardware investment while increasing patient care. HeartStride is built on 50 years of diagnostic experience and reliably delivers clinical accuracy and diagnostic confidence. Plus it's backed by the Cardiac Science Corporation field service network. Features: Superior ECG trace quality; Easy, intuitive operation; Simple to learn and use; Comprehensive Burdick on-site training; In-office Practice Support Program; Fast, accurate reporting.

H20-1AEPT01 Network Configuration: English Language w/AHA Accessories, Includes: CPU, Printer, LCD Flat Panel Monitor, Laser Printer, Network Connectivity, TM55 Treadmill & Custom Cart.  
H20-1AKPT01 Base Configuration: English Language w/AHA Accessories, Includes: CPU, Printer, LCD Flat Panel Monitor, Laser Printer, TM55 Treadmill & Custom Cart.  
H21-1AKT01 SW Only: English Language w/AHA Accessories, Includes: TM55 Treadmill & Custom Cart.






**MIDMARK IQMARK™ DIGITAL ECG (USB VERSION)**

Includes: IQmark™ Digital ECG Acquisition Module (USB Version), IQmark™ Digital ECG Software Package, IQmark™ Diagnostic Workstation - Patient Management Software, 10-lead patient cable, 10-pack universal clips, 100-pack ECG electrodes, mouse pad, operation manuals and quick reference guide.

4-000-0056 IQmark™ Digital ECG (USB Version)

**MIDMARK IQMARK™ EZ ADVANCED HOLTER SUPPLIES**

- 2-100-0060 Lead-Lok Electrodes, 50/bx
- 2-100-0080 Prep Kit - IQmark™ Holter, 5-Lead, 2AAs
- 2-100-0090 Prep Kit - IQmark™ Holter, 7-Lead, 2AA's
- 3-370-0001 Lead Wire, 5-Lead, IQmark™ Digital Recorder
- 3-370-0002 Lead Wire, 7-Lead, IQmark™ Digital Recorder
- 3-370-0003 Pouch w/belt for IQmark™ Digital Recorder
- 1-370-0002 Compact Flash Card
- 2-100-0081 Prep Kit - IQmark™ Holter, 5-Lead 2 AAs, 30/cs
- 2-100-0091 Prep Kit - IQmark™ Holter, 7-Lead, 2AA's, 30/cs



**MIDMARK IQMARK™ PACEMAKER OPTION FOR EXISTING ADVANCE HOLTER**

IQmark™ Advanced Holter Pacemaker Analysis Software, Sentinel Data Protection Key for Pacemaker (parallel port key shipped standard, USB key available upon request), operation manual, and quick reference guide.

4-000-0140 IQmark™ Pacemaker Option For Existing Advanced Holter

**MIDMARK IQMARK™ 7700/8800 HOLTER SUPPLIES**

- 3-053-0001 Lead Wire, 7-Lead, IQmark™ Digital Recorder
- 3-079-0001 Pouch w/belt for IQmark™ Digital Recorder
- 3-513-0005 Junction Box For IQmark™ Digital Recorder
- 1-079-0125 16MB Flash Card

**MIDMARK IQMARK™ MANUALS**

- 3-100-1006 IQmark™ Digital Recorder Manual
- 3-100-1007 IQmark™ Holter 7700 Manual
- 3-100-1008 IQmark™ Holter 8800 Manual
- 3-100-1021 IQmark™ Digital ECG Manual
- 3-100-1025 IQmark™ Digital Spirometer Manual
- 3-100-1038 IQmark™ Diagnostic Workstation Manual
- 3-100-1044 IQmark™ EZ Holter Manual
- 3-100-1045 IQmark™ Digital Recorder Manual
- 3-100-1049 IQmark™ Digital ECG PDA & SpiroPDA Manual
- 3-100-1050 IQmark™ Advanced Holter Manual

**MIDMARK IQMARK™ ECG SUPPLIES**

- 2-100-0205 Disposable ECG EZ-Trodes, 500/bx
- 3-047-0001 Universal ECG Clips (3mm, 4mm & Snap)
- 3-100-0200 IQmark™ Patient Cable
- 3-999-0011 24-pack, AA Batteries
- 2-100-0206 Disposable ECG EZ-Trodes, 500/bx, 10 bx/cs

**KENDALL-LTP MEDI-TRACE® CHART PAPERS**

Kendall-LTP manufactures one of the largest offerings of recording chart papers in the industry. With more than 40 years of manufacturing experience, our customers can expect consistent trace quality as well as a product that is designed to meet or exceed the instrument specifications.

- 01091735 Fetal Monitoring Paper, Compatible to Philips 9270-0678-030, w/label sets, 30 pk/cs
- 20608594 Recording Chart Paper, Compatible to Spacelabs/Burdick 007868, 10 pads/cs
- 30015682 Recording Chart Paper, Compatible to Spacelabs 307438, 50 r/cs
- 30540499 Recording Chart Paper, Compatible to Spacelabs 307568, 20 pad/cs
- 30561354 Recording Chart Paper, Compatible to Quinton 16211-002, 20 pad/cs
- 30597226 Recording Chart Paper, Compatible to GE/Corometrics 4305AAO, 40 pad/cs
- 30659745 Recording Chart Paper, Compatible to Nihon Kohden FQW50-3-100, 100 pad/cs
- 30718970 Recording Chart Paper, Compatible to Physio-Control 804700, 50 r/cs
- 30725389 Recording Chart Paper, Compatible to Philips 40477A, 80 r/cs
- 30729225 Recording Chart Paper, Compatible to Quinton 18416-001, 20 pad/cs
- 30748696 Recording Chart Paper, Compatible to Philips M1910A, 40 pad/cs
- 30763711 Recording Chart Paper, Compatible to Spacelabs 006-0196-00, 100 pad/cs
- 30768678 Recording Chart Paper, Compatible to Philips M1707A, 8 pad/cs
- 30781747 Recording Chart Paper, Compatible to GE/Marquette 9402-046, 50 r/cs
- 31002176 Recording Chart Paper, Compatible to GE/Marquette 9402-020, 8 pad/cs
- 31002184 Recording Chart Paper, Compatible to GE/Marquette 9402-024, 8 pad/cs
- 31091427 Recording Chart Paper, Compatible to Physio-Control 805319-03, 18 r/cs

**JONES-ZYLON HOLTER MONITOR KIT**

Contains: patient diary, audio cassette tape, 5 electrodes, surgical prep razor, 2 AA alkaline batteries and 6" x 9" minigrip bag.

568803 Holter Monitor Kit, 50/cs



**BURDICK BASELINE HOLTER CASSETTE PREP KIT**

Each kit contains: 1 Hook-Up Guide, 1 Surgical Prep Razor, 1 Skin Prep Scrub, 3 Alcohol Prep Pads, 12 Stress Loop Tape, 1 60 Minute Audio Cassette, 1 9V Battery, 1 Patient Diary.

- 043084 2 Channel Cassette Prep Kit, w/5 Holter Electrodes, 30/cs
- 043085 3 Channel Cassette Prep Kit, w/7 Holter Electrodes, 30/cs



**BURDICK BASELINE HOLTER DISC PREP KIT**

Each kit contains: 1 Hook-Up Guide, 1 Surgical Prep Razor, 1 Skin Prep Scrub, 3 Alcohol Prep Pads, 12 Stress Loop Tape, 1 Patient Diary, 1 Patient ID Card, 2 9V Batteries.

- 043101 2 Channel 6532 Prep Kit, w/5 Holter Electrodes, 30/cs
- 043104 3 Channel 6532 Prep Kit, w/7 Holter Electrodes, 30/cs



043253



043273

**BURDICK HOLTER MONITORING SYSTEM ACCESSORIES**

- 012-0141-00 Universal Test Box For 90205 & 90208
- 012-0439-00 5-Lead Bonded Cable/Lead Set For 90205 & 90208
- 012-0440-00 7-Lead Bonded Cable/Lead Set For 90205 & 90208
- 016-0475-00 Nylon Carrying Pouch For 90205 & 90208
- 016-5066-00 Nylon Adjustable Shoulder Strap For 90205 & 90208
- 043092 Test Box For 6532, 6632, 6300
- 043113 Interface Kit For 6300 Cassette Recorder
- 043138 2 Channel Safety Cable (cable only) For 6632
- 043139 2 Channel Lead Set, 5 Safety Wires For 6632
- 043140 3 Channel Lead Set, 7 Safety Wires For 6632
- 043229 Signal Averaging/Late Potential Software
- 043253 3 Channel Lead Set, 7 Safety Wires, 6732 Recorder
- 043255 2 Channel Lead Set, 5 Safety Wires, 6732 Recorder
- 043273 Compact Flash Card, 30Mb For Use w/6732 Compact Digital Recorder
- 086418 Altair PC Holter System Operating Manual
- 086419 Altair Simplicity Holter Operating Manual
- 086421 6500 DISC Recorder Operators Manual
- 086422 Altair 6300 Cassette Recorder Operators Manual
- 143207 Interface Kit For 6732 & 6632 PC Card Recorder
- 175-1401-00 3 Channel Safety Cable (cable only) For 6632
- 685-0200-00 Compact Flash Card, 64 MB For 6732
- 92514-01 Digital Holter Recorder, 3-Channel, 128 MB

**BURDICK PC CARD RECORDER PREP KIT**

Each kit contains: 1 Surgical Prep Razor, 3 Alcohol Prep Pads, 7 Holter Electrodes, 1 Patient ID Card, 1 Hook-Up Guide, 1 Skin Prep Scrub, 12 Stress Loop Tape, 1 Patient Diary, 4 "AA" Batteries

- 043202 PC Card Recorder Prep Kit, 30/cs



**BURDICK THREE-CHANNEL COMPACT DIGITAL RECORDER PREP KIT**

Each kit contains: 1 Surgical Prep Razor, 3 Alcohol Prep Pads, 7 Holter Electrodes, 1 Patient ID Card, 1 Hook-Up Guide, 1 Skin Prep Scrub, 12 Stress Loop Tape, 1 Patient Diary, 2 "AA" Batteries

- 043272 3-Channel Compact Digital Recorder Prep Kit, 30/cs

**BURDICK TWO-CHANNEL COMPACT DIGITAL RECORDER PREP KIT**

Each kit contains: 1 Surgical Prep Razor, 3 Alcohol Prep Pads, 5 Holter Electrodes, 1 Patient ID Card, 1 Hook-Up Guide, 1 Skin Prep Scrub, 12 Stress Loop Tape, 1 Patient Diary, 2 "AA" Batteries

- 043250 2-Channel Compact Digital Recorder Prep Kit, 30/cs

**GE MEDICAL HOLTER 12-LEAD AMBULATORY CABLES**

Ribbon cable design with strategic electrode snap placement. No need for separate leadwires; one piece design. Cable can be cleaned or disinfected. Lightweight wire can be used on a large cross section of patients. Latex free. 90 day warranty.

- 421825-001 Cable SEER MC Recorder, 12-Ld Ambulatory, AHA, 36", Includes Electrode Snap Connections
- 421825-002 Cable SEER MC Recorder, 12-Ld Ambulatory, IEC, 36", Includes Electrode Snap Connections

**GE MEDICAL HOLTER 5 & 7 LEAD SNAP-LOCK CABLES & LEADWIRES**

Locks leadwires in place to ensure a good connection. Release button allows the leadwires to be removed when depressed. Snap end of leadwire provides secure connection to snap style electrodes. Latex free. 90 day warranty.

- 409106-001 Cable 8500 Holter, Snap-lock, 5-Ld, DIN, 42"
- 409106-002 Cable 8500 Holter, Snap-lock, 7-Ld, DIN, 42"
- 409106-003 Cable SEER Holter, Snap-lock, 5-Ld, DIN, 42"
- 409106-004 Leadwire Set, 8500/SEER/SEER MC Holter, 5-Ld, 18", 5/set
- 409106-005 Leadwire Set, 8500/SEER MC Holter, 7-Ld, 18", 7/set
- 409106-006 Leadwire Set, 8500/SEER/SEER MC Holter, 5-Ld, 24", 5/set
- 409106-201 Cable 8500 Holter, Snap-lock, 5-Ld, DIN, 15"
- 419598-005 Cable SEER MC, Snap-lock, 5-Ld, DIN, 42"
- 419598-007 Cable SEER MC, Snap-lock, 7-Ld, DIN, 42"

**GE MEDICAL HOLTER DOWNLOAD CABLES**

Download patients data quickly and easily from SEER to Holter scanner. Durable flex cable. Latex free. 90 day warranty.

- 80784-001 20" Cable Download SEER Recorder

**GE MEDICAL HOLTER MULTI-LINK 12-LEAD CABLES & LEADWIRES**

Shielded leadwires with multi-link connectors. Provides 12-lead monitoring capability. Latex free. One year warranty.

- 411200-001 Multi-Link Leadwire Set, 5-Ld, Snap, AHA, 29", RA, LA, RL, LL, V, 5/set
- 411200-003 Multi-Link Leadwire Set, 5-Ld, Snap, IEC, 29", R, L, N, F, C, 5/set
- 419598-010 Cable SEER MC, Multi-Link, 12-Ld, AHA, 36"
- 419598-011 Cable SEER MC, Multi-Link, 12-Ld, IEC, 36"

**GE MEDICAL HOLTER SUPPLIES**

- 3801-902 9V Battery Energizer Alkaline, 12/bx
- 3801-904 1.5V AA Battery Energizer Alkaline, 144/bx
- 3801-905 9V Battery Energizer Alkaline, 72/bx
- 400068-002 Holter Kits 5-Lead Kit Contents: 5 Vinyl Tape Electrodes, 24-Hour Cassette, 9V Battery, 2 Prep Pads & Patient Diary, 24 kits/cs
- 400068-003 Holter Kits 7-Lead Kit Contents: 7 Vinyl Tape Electrodes, 24-Hour Cassette, 9V Battery, 2 Prep Pads & Patient Diary, 24 kits/cs
- 400068-004 Holter Kits 7-Lead Kit Contents: 7 Vinyl Tape Electrodes, 48-Hour Cassette, 9V Battery, 2 Prep Pads & Patient Diary, 24 kits/cs
- 400068-009 Holter Kits 5-Lead Kit Contents: 5 Vinyl Tape Electrodes, 4 AA Batteries, 2 Prep Pads & Patient Diary, 24 kits/cs
- 400068-016 Holter Kits 5-Lead Kit Contents: 5 Foam Oval Electrodes, 4 AA Batteries, 2 Prep Pads & Patient Diary, 24 kits/cs
- 400068-021 Holter Kits 5-Lead Kit Contents: 5 Vinyl Tape Electrodes, 48-Hour Cassette, 9V Battery, 2 Prep Pads, Patient Diary & Prep Razor, 24 kits/cs
- 400068-024 Holter Kits 7-Lead Kit Contents: 7 Vinyl Tape Electrodes, 4 AA Batteries, 2 Prep Pads & Patient Diary, 24 kits/cs
- 400068-025 Holter Kits 12-Lead Kit Contents: 10 Vinyl Tape Electrodes, 4 AA Batteries, 2 Prep Pads & Patient Diary, 24 kits/cs
- 57995-002 8500 Holter Recorder Belt For Pouch 57995-005
- 57995-003 8500 Holter Recorder Shoulder Pouch w/strap
- 57995-005 8500 Holter Recorder Pouch w/o straps & belt
- 57995-901 8500 Holter Recorder Pouch w/strap & belt
- 90026-902 Holter Patient Diary, 50/pk
- 9435-001 24-Hour Capacity Cassette Tapes, 10/bx
- 9435-002 48-Hour Capacity Cassette Tapes, 10/bx
- 9435-003 24-Hour Capacity Cassette Tapes, 20/bx
- 9435-004 48-Hour Capacity Cassette Tapes, 20/bx

**GE MEDICAL HOLTER TEST CABLES**

Verify ECG signals on 8500 recorder. Positive fit connector to the recorder; banana socket for ECG cart leadwires. Accepts 3mm or 4mm leadwires. Latex free. 90 day warranty.

- 900248-001 8500 Holter Test Cable, 3 Channel, Fits Recorders With Mini-9 Pin Connectors
- 9401-901 8500 Holter Test Cable, 2 Channel, Fits Recorders With 3 Pin Round Connectors

**GE MEDICAL HOLTER VALUE LINE CABLES & LEADWIRES**

Easy insertion of leadwire into cable. Latex free. 90 day warranty.

- 2014788-001 Set Leadwire Holter, 5-Ld, DIN, 20", Value Line, 5/set
- 2014789-001 Set Leadwire Holter, 7-Ld, DIN, 20", Value Line, 7/set
- 2014790-001 Cable 8500 Holter Recorder, 5-Ld, DIN, 42", Value Line
- 2014791-001 Cable 8500 Holter Recorder, 7-Ld, DIN, 42", Value Line
- 2014792-001 Cable SEER Recorder, 5-Ld, DIN, 42", Value Line
- 2014793-001 Cable 8500 Holter Recorder, 5-Ld, DIN, 15", Value Line
- 2014794-001 Cable 8500 Holter Recorder, 7-Ld, DIN, 15", Value Line
- 2014795-001 Cable SEER Recorder, 5-Ld, DIN, 15", Value Line

**WELCH ALLYN HOLTER SYSTEM 2000 ACCESSORIES**

- 08239-0000 Template Matching, Software Option
- 08240-0000 Pouch For MT-100
- 08241-0000 Holter Hook-Up Kit, Includes 1.5V AA Batteries, Patient Diary, Razor, Abrading Pad, Alcohol Prep Pads, 5 Disp. Holter Electrodes
- 80188-0000 Patient Cable, Short, MT-100
- 80189-0000 Fiber Optic Transfer Cable, MT-200

**KENDALL-LTP MEDI-TRACE® TELEMETRY/HOLTER POUCHES AND KITS**

Kendall-LTP offers a large variety of Holter and Telemetry pouches, kits, and accessories. All of our pouches are designed for patient comfort, transmitter protection and improved monitor performance. Our kits are offered in many variations to accommodate almost any requirement.

- 30951514 Telemetry Kit, MPK-MOD5, Contains 535 Electrode, Prep Pad, (2) 1.5 Alkaline Batteries and Pouch, 50/cs
- 30951555 Telemetry Pouch, TMY-BLUE-50, Small, 50/cs
- 30951746 Holter Kit, HMPK-HOLTR1, Contains 535 & 533 Electrodes, Tape, Prep Pad, Abrasive Pad, 9V Battery, Razor and Diary, 50/cs
- 30951811 Telemetry Pouch, TMY-BLUE-2-100, Large, 100/cs
- 31027942 Telemetry Kit, MPK-MOD3, Contains 533 Electrode, Prep Pad, 8.4V Battery & Pouch, 50/cs
- 31027983 Telemetry Kit, MPK-MOD2, Contains 535 Electrode, Prep Pad, 8.4V Zinc-Air Battery & Pouch, 50/cs
- 31098448 Telemetry Kit, MPK-MOD4, Contains 535 Electrode, Prep Pad, 9V Battery & Pouch, 50/cs
- MK05140 Holter Kit, SK7828, Contains 535 Electrode, 9V Battery, Tape & Diary, 20/cs
- MK06200 Holter Kit, Digial, Contains 535 & 533 Electrodes, (2) AA Batteries, Razor, Prep Pad, Abrasive Pad & Diary, 20/cs
- MK07740 Telemetry Pouch, TP11 Large w/window, 50/cs
- MK07860 Telemetry Pouch, TP3 Non-woven, 100/cs
- MK07900 Telemetry Pouch, TP5 Non-woven w/2 straps, 100/cs
- MK07980 Telemetry Pouch, TP7S Large w/2 straps, 50/cs
- MK08000 Telemetry Pouch, TP8 w/slit & loop, 50/cs
- MK08080 Telemetry Pouch, TPV Vinyl w/clip, 50/cs



**BURDICK BI-PHASIC MEDIC 6 PORTABLE DEFIBRILLATOR**

Compact, streamlined design that provides easy to use, comprehensive defibrillation therapies in a compact 13 1/2 lb package. A radiant 6.4", non-glare screen makes functions easy to use in a variety of light settings. The easy to follow function keys, make system feature selections quick and easy to navigate. Includes a powerful rechargeable battery, which delivers 40+ consecutive discharges at 360j or 2 hours of continuous monitoring operation. Sync and a broad range of joule (energy) selections are available for emergency resuscitation, internal electrode procedures or other cardiac therapeutic modalities. Offers ECG capabilities via 3-lead (included) patient cable. Comprehensive pre/post treatment summary documentation which is provided on 50mm paper and stored into internal memory.

- 92450-1APV Medic 6 Defibrillator w/ext pacing & advisory
- 92450-1AV Medic 6 Defibrillator w/advisory



## BURDICK DEFIBRILLATOR ACCESSORIES

|             |  |
|-------------|--|
| 015-0605-00 | Patient Cable, 3-Lead, Medic 6                                     |
| 015-0608-00 | Battery Pack For Medic 6   |
| 015-0617-00 | Standard Thermal Paper, 150 ft x 45mm For Medic 6, 3 r/bx          |
| 015-0610-00 | AC Power Cord For Medic 6  |
| 047262      | Power Cord For DC-200, Medic 4 & Medic 5                           |
| 086144      | Operators Manual For DC-200  |
| 086183      | Service Manual For DC-200  |
| 015-0612-00 | Defibrillation/Pacer/Monitoring Pads, Adult For Medic 6, 10/bx     |
| 015-0621-00 | Defibrillation/Pacer/Monitoring Pads, Pediatric For Medic 6, 10/bx |
| 067006      | External Paddle Set For Medic 5, Adult                             |
| 067015      | Adult Internal Paddles w/handles For DC-190                        |
| 067016      | Pediatric Internal Paddles w/handles For DC-190                    |
| 067019      | Adult Internal Paddles w/handles For Medic 5                       |
| 067107      | 3-Lead Patient Cable For Medic 4 & Medic 5                         |
| 067112      | 5-Lead Patient Cable For Medic 4 & Medic 5                         |
| 850711      | Battery Pack For Medic 5   |



## WELCH ALLYN PIC30 DEFIBRILLATOR

|            |  |
|------------|--|
| 973039C00E | PIC30 Manual Defibrillator                           |
| 973039C0AE | PIC30 Manual Defibrillator w/shock advisory          |
| 973039CPOE | PIC30 Manual Defibrillator w/pacing                  |
| 973039CPAE | PIC30 Manual Defibrillator w/pacing & shock advisory |



## WELCH ALLYN PIC40 DEFIBRILLATOR

|         |                                     |
|---------|-------------------------------------|
| 971016  | EtCO2 Capnography                   |
| 973092E | PIC40 Defibrillator                 |
| 973093E | PIC40 Defibrillator w/NIBP          |
| 973094E | PIC40 Defibrillator w/NIBP & pacing |
| 973095E | PIC40 Defibrillator w/pacing        |



## WELCH ALLYN PIC50 DEFIBRILLATOR

|         |                                     |
|---------|-------------------------------------|
| 971081E | PIC50 Defibrillator, 3 Lead, Basic  |
| 971082E | PIC50 Defibrillator, 3 Lead, Deluxe |
| 971083E | PIC50 Defibrillator, 12 Lead, Basic |

## WELCH ALLYN PIC™ RESCUER AED/MONITOR/DEFIBRILLATOR ACCESSORIES

|              |   |
|--------------|---|
| 001537-U     | PIC Pediatric Electrode                 |
| 001647-U     | SuperPak Rechargeable Battery           |
| 001726-U     | 10-Pack Disposable Electrodes, 10 pk/bx |
| 001739-U     | Chart Plain White Paper, 40mm, 3 r/bx   |
| 001790-U     | Cable: 3-Lead ECG (AHA) 6 ft            |
| 001795-U     | PIC 5-Lead Patient Cable                |
| 001934-U     | CO2 Sensor Kit                          |
| 001950-U     | Adult Airway Adapter, 5/bx              |
| 560639       | PIC Tinted Lens Package                 |
| 560640       | PIC Clear Lens Package                  |
| 601215       | PIC Hanger, Black                       |
| 90021-4      | PIC Soft Carrying Case                  |
| 900323       | PIC30 Carrying Case                     |
| 971029-E     | PIC Paddle Tray                         |
| 971091       | Reliable Integral External Pacing       |
| 971107       | PIC Hands Free Adapter, Grey            |
| 971108-E     | PIC Deluxe External Paddle Set, Grey    |
| 980132       | PIC 12V Power Adapter                   |
| 999895       | PIC EMS Service Plus Program (per year) |
| 971104-PIC-E | PIC Quick Charger                       |



## 3M™ DEFIB PADS

Ready-to-use, conductive gel pads for use in defibrillation. Self-contained to eliminate inadvertent gel bridging which can cause skin trauma to patient or operator. Bright orange solid gel pads are easy to see, require no cleanup, and leave no slippery residue that can interfere with CPR.

|       |  |
|-------|--|
| 2345N | 4 1/2" x 4 1/2" Defib Pads, 2/pk, 10 pk/bx, 10 bx/cs |
| 2346N | 4 1/2" x 6" Defib Pads, 2/pr, 10 pr/bx, 10 bx/cs     |



## ZOLL AED PRO

Rugged construction and exceptional durability support the needs of Basic Life Support providers and tiered systems. Superior design and rubber over-molds are built to withstand harsh environments. ECG Monitoring with 3-lead cable and high-quality display. Manual override option satisfies Advanced Life Support professionals. Unrivaled battery compatibility and options that meet your specific needs. Internal memory data storage and IrDA or USB data transfer.

|                    |  |
|--------------------|--|
| 90110820499992010N | AED Pro Includes: (2) Sets STAT Padz II, (1) Charger, (2) Rechargeable Lead Acid Batteries |
|--------------------|--|

**ZOLL D 900 DEFIBRILLATOR**

The Zoll D 900 offers advanced defibrillation technology in an easy-to-use unit ideal for “first response” defibrillation in hospital areas where complete ACLS-response capabilities are not required. It provides ECG monitoring with a large, easy-to-read monitor. The Zoll D 900 can be used either with Zoll Multi-Function Electrodes for “hands-off” defibrillation or with standard defibrillation paddles. A significant advantage is that pacing functions may be added to upgrade the unit at any time. The Zoll D 900 operates on both AC line and battery power and is completely compatible with other models in the Zoll cardiac resuscitation system.

8009-0006 1 Set Defib Paddles, 1 ECG Cable, Recorder Tape, 2 Guides & 1 Tube Defib Gel, Refurbished

**ZOLL PD 1200 DEFIBRILLATOR  
PACEMAKER/DEFIBRILLATOR**

The Zoll PD 1200, combining defibrillation, pacing and monitoring, offers complete ACLS-response capabilities in a single, easy-to-use unit that is ideal for critical care specialty care areas. It can be configured with the Zoll “hands-off” Multi-Function Electrodes, with standard defibrillation paddles and pacing electrodes, and with internal defibrillation paddles and pacing electrodes. The Zoll PD 1200 operates on both AC line and battery power and is completely compatible with other models in the Zoll cardiac resuscitation system.

8002-0011 1 Set Standard Defib Paddles, 1 ECG Cable, 1 Pacer Output Cable, 2 Sets Electrodes (Adult, Pedi), Recorder Paper, 2 Guides, 1 Tube Gel, 1 Inservice Video, 2 Yr. Warranty, Refurbished

**ZOLL MEDICAL AED PLUS DEFIBRILLATOR**

Combines superior biphasic defibrillation with fast and simplified electrode placement. The AED Plus uses a brilliant new graphical/voice user interface to help the rescuer every step of the way. The electrode pad helps the rescuer perform CPR. A sensor on the CPR-D pads tracks the rescuer's compression depth and rate as soon as CPR begins. The electrode comes complete with rescue essentials and tamper-evident seals and it lasts four years. Features include audio and visual coaching and Rectilinear Biphasic™ Waveform. Requires only one electrode pad. Runs on inexpensive, consumer lithium batteries.

21000010102011010 AED Plus Defibrillator

**ZOLL MEDICAL AED DEFIBRILLATOR ACCESSORIES**

1000-0113-01 Replacement Trainer Remote Control  
 8000-0800-01 Simulator/Tester For Functional Testing  
 8000-0804-01 Universal Simulator Cable  
 8000-0811-01 Flush Mount Wall Cabinet  
 8000-0813-01 Zoll Data Review  
 8000-0817-01 Surface Mount Wall Cabinet  
 8000-0821 Replacement Trainer Phone Cord  
 8000-0822 Replacement Trainer AC Adapter  
 8008-0104-01 AED Trainer w/8 pre-programmed AED scenarios & remote control  
 8900-2055 Adult Multi Function Electrode, 12/cs  
 8900-2065 Pediatric Multi Function Electrode, 6/cs  
 8900-0800-01 CPR-D Padz w/compression  
 9310-0737 8 1/2" x 11" Wall Sign  
 9310-0738 V-Shape Wall Sign  
 9650-0301-01 Administration Guide  
 9650-0304-01 Replacement Trainer Operator's Guide  
 8000-0301 80mm Fanfold Paper Use Only w/M-Series Defibrillator, 20/cs  
 8000-0802-01 Replacement Soft Case  
 8000-0807-01 Type 123 Lithium Batteries, 10/pk  
 8000-0808-01 Replacement PASS Cover  
 8000-0809-01 Wall Mounting Bracket  
 8000-0814-01 Recessed Mount Wall Cabinet  
 8000-0815-01 USB IRDA Adapter  
 8000-0816-01 RS-232 IRDA Adapter  
 8900-0801-01 Stat-Padz II HVP Multi-Function Electrodes  
 8900-0802-01 Stat-Padz II, 1 pr  
 8900-0803-01 Replacement Adhesive Gels for CPR-D Padz - Training Electrode replacements, 5pair/bx  
 8900-0804-01 CPR-D Padz Training Electrodes, reusable "z-design" electrode with CPR Hand Placement Indicator and 1 pair of Disposable Adhesive Gels  
 8900-0805-01 Stat Padz II (reusable for Training)  
 8900-0810-01 ACCESSORY: Pedi-padz II, Pediatric Multi-Function Electrodes, 1 pair  
 8900-4003 HVP Multi Function Electrodes, stat-padz, 24 month shelf life, 12pr/cs  
 8900-4004 Stat-padz HVP Multi-Function Electrodes, Shelf Life 24 Months, 1 pair  
 8900-4005 Solid Gel Multi Function Electrodes, pro-padz, 24 month shelf life, 12pr/cs  
 8900-4012 Multi-Function Electrode, Shelf-Life 9 mo, Sterile, 6 pr/cs  
 9650-0300-01 Operator's Guide/Wall Poster  
 9656-0156 AED Plus Brochure Literature, 25/pk  
 9656-0158 AED Plus Automated External Defibrillator w/CPR feedback literature, 50/pk  
 9656-0161 AED Plus Die-Cut Literature, 25/pk  
 9656-0162 AED Plus Accessories & Training Products Literature, 25/pk

**BURDICK FIRSTSAVE AED DEFIBRILLATOR**

A consumer-grade AED. It includes many of the innovations found on the CardioVive AED at a value price. Audible alarms indicate system failure. Patented RHYTHMx® analysis software. Patented STAR® biphasic defibrillation energy waveform. Patented technology performs automated daily, weekly, and monthly self-tests on the battery, electrode presence and internal electronics. Pediatric capability. Simple, easy-to-use, one-button operation. Instructive voice prompt guides.

92531-201 FirstSave AED Defibrillator

**BURDICK FIRSTSAVE AED DEFIBRILLATOR ACCESSORIES**

010-1564-00 Adult Electrodes, 1 set/pk  
 010-1572-00 Wire Rack  
 010-1596-00 Pediatric Electrodes, 1 set/pk  
 010-1649-00 Replacement Battery  
 010-1651-00 Soft Carrying Case



## BURDICK CARDIOVIVE AED DEFIBRILLATOR

Automatic daily self-test ensures safe, reliable operation. Pre-connected, non-polarized, disposable pads eliminates confusion by saving critical time. Easy to use: open cover, open electrode package, place electrodes on bare chest & push button as directed. Patented biphasic waveform measures impedance and adjusts to deliver appropriate defibrillation level. Advanced technology differentiates between shockable life-threatening cardiac rhythms and non-shockable arrhythmia events.

92532-201 CardioVive Defibrillator



## WELCH ALLYN AED-10 AUTOMATED EXTERNAL DEFIBRILLATOR

Comes with an illuminated defib pad placement illustration on the device front panel, backlit display to provide the user with test prompts - battery status - number of shocks and energy settings, concise voice prompts and CPR timer. Includes the following standard accessory package: (1) non-rechargeable lithium battery (90+ shocks at 360 joules or 5 hours of continuous use), user's manual, quick reference card, defibrillation pads and five-year limited warranty program.

970300E AED-10 Device & Standard Accessories Only  
 970302E AED-10 w/carrying case & standard accessories  
 970308E AED-10 w/lead out pads  
 970311E AED-10 RAPID RESPONSE KIT incl: AED-10 w/new expanded carrying case, 2 AED-10 Lithium Batteries, 2 sets of Defib Pads, 1 Pediatric Energy Reducer, 1 1st Responder Prep Kit, Quick Ref Card, In-Svc Video, and Operator Manual

### Accessories

900432 AED 10 Carrying Case, 1 Year Limited Warranty

## WELCH ALLYN AED AUTOMATED EXTERNAL DEFIBRILLATOR ACCESSORIES


00172-0 Electrode: Bi-Trode Limb Electrode Set (2 wrist clamp electrodes)  
 00183-0E AED-20 Non-Rechargeable (Lithium) Battery  
 00185-2 AED-10 Non-Rechargeable (Lithium) Battery  
 00185-3 AED Defibrillation Pads, 2/pouch  
 001851-U AED Trainer, Adult Pads  
 001855-U JumpStart Adult Defibrillation Pads, 10 pr/bx  
 001857-E Lead Out Pads For AED10 & 20  
 001962-U SmartLink Software Program  
 00211-9 AED-20 Parallel Printer Communications Kit  
 00213-0 Cable: AED-20 ECG (AHA)  
 00213-7 AED Quick Response Prep Kit  
 00214-3 AED IR Communication Kit  
 002155 AED Wire Mount  
 002156-U AED Compact Wall Cabinet  
 90021-6E Rugged AED-20 Carrying Case w/quick reference instruction placard, storage space for 1 battery & storage pockets for other supplies  
 90042-1 Rugged AED-10 Soft Shell Carrying Case, Includes Quick Reference Instruction Placard & Storage Area For Supplies  
 970205 AED-20 Upgrade ECG to Manual  
 980136 Cardiology Datacard, 4MB "Linear"  
 980143 AED-20 Trainer  
 981125E PowerStick Battery Charging Station, 1 Bay  
 990178 AED-20 In Service Video (VHS)



## WELCH ALLYN AED-20 AUTOMATED EXTERNAL DEFIBRILLATOR

Includes a diagnostic quality backlit display which provides the user with text prompts, battery status, heart rate, event timer, number of shocks, volume control, display brightness control and energy settings. Also included are concise voice prompts and a CPR timer. Can be upgraded in the future for expanded functions. Comes complete with the standard accessory package: (1) battery (either a non-rechargeable long life lithium for infrequent use or rechargeable nickel metal hydride for frequent use/monitoring with charging station), operator's manual, quick reference card, defibrillation pads and five-year limited warranty program.

972211E AED 20 Basic w/non-rechargeable battery  
 972212E AED 20 Basic w/rechargeable battery  
 972213E AED 20 Basic Plus w/non-rechargeable battery  
 972214E AED 20 Basic Plus w/rechargeable battery  
 972215E AED 20 Deluxe w/non-rechargeable battery  
 972216E AED 20 Deluxe w/rechargeable battery



## 3M™ RED DOT™ CLEAR PLASTIC MONITORING ELECTRODES

A tear-drop style monitoring electrode with clear plastic backing that is ideal for situations when immediate skin assessment is required without the need to remove the electrode. The gel's low chloride content minimizes skin irritation for increased patient comfort.

2235 5.1cm Dia. Monitoring Electrodes w/o abrader, 50/bg, 20 bg/cs  
 2235-3 5.1cm Dia. Monitoring Electrodes w/abrader, 3/bg, 20 bg/bx, 10 bx/cs  
 2235-5 5.1cm Dia. Monitoring Electrodes w/abrader, 5/bg, 20 bg/bx, 10 bx/cs  
 9642 5.1cm Dia. Monitoring Electrodes w/o abrader, 50/bg, 20 bg/cs



## 3M™ RED DOT™ DIAPHORETIC FOAM MONITORING ELECTRODES

A rectangular-shaped electrode that offers a high-tack adhesive and foam backing designed for diaphoretic clinical situations. Patented solid gel penetrates the skin quickly allowing high quality. Low chloride content minimizes skin irritation.

2230 5.1cm x 5.5cm Monitoring Electrodes w/o abrader, 50/bg, 20bg/cs  
 2270 5.1cm x 5.5cm Monitoring Electrodes w/abrader, 25/bg, 40bg/cs  
 2270-3 5.1cm x 5.5cm Monitoring Electrodes w/abrader, 3/bg, 20 bg/bx, 10 bx/cs  
 2270-5 5.1cm x 5.5cm Monitoring Electrodes w/abrader, 5/bg, 20 bg/bx, 10 bx/cs







### 3M™ RED DOT™ DIAPHORETIC SOFT CLOTH MONITORING ELECTRODES

A rectangular-shaped electrode that offers a high-tack adhesive and soft cloth backing designed for diaphoretic clinical situations. Patented solid gel penetrates the skin quickly allowing high quality. Low chloride content minimizes skin irritation.


|        |   |
|--------|---|
| 2231   | 5.1cm x 5.5cm Monitoring Electrodes w/o abrader, 50/bg, 20 bg/cs        |
| 2271   | 5.1cm x 5.5cm Monitoring Electrodes w/abrader, 25/bg, 40 bg/cs          |
| 2271-3 | 5.1cm x 5.5cm Monitoring Electrodes w/abrader, 3/bg, 20 bg/bx, 10 bx/cs |
| 2271-5 | 5.1cm x 5.5cm Monitoring Electrodes w/abrader, 5/bg, 20 bg/bx, 10 bx/cs |



2269T



2280



2284

### 3M™ RED DOT™ ECG MONITORING ELECTRODES WITH PRE-ATTACHED LEAD WIRE

|       |   |
|-------|---|
| 2269T | Neonatal, 2cm x 4cm, Pre-Wired Radiolucent Electrode w/clear tape, 3/bg, 10 bg/bx, 10 bx/cs |
| 2280  | Neonatal, 22mm x 22mm, Pre-Wired Non-Radiolucent Electrode w/soft cloth, 3/bg, 100 bg/cs    |
| 2282  | Neonatal, 22mm x 22mm, Pre-Wired Radiolucent Electrode w/soft cloth, 3/bg, 100 bg/cs        |
| 2283  | Infant, 33mm x 22mm, Pre-Wired Radiolucent Electrode, 3/bg, 100 bg/cs                       |
| 2284  | Neonatal, 17mm x 22mm, Pre-Wired Radiolucent Limb Band Electrode, 3/bg, 100 bg/cs           |



### 3M™ RED DOT™ FOAM MONITORING ELECTRODES

Tear-drop style monitoring electrodes with foam backing. Patented solid gel penetrates the skin quickly allowing high quality, uninterrupted traces. Low chloride content minimizes skin irritation.

|        |  |
|--------|--|
| 2237   | 4.4cm Dia. Monitoring Electrodes w/o abrader, 50/bg, 20 bg/cs          |
| 2237-3 | 4.4cm Dia. Monitoring Electrodes w/o abrader, 3/bg, 20 bg/bx, 10 bx/cs |
| 2237-5 | 4.4cm Dia. Monitoring Electrodes w/o abrader, 5/bg, 20 bg/bx, 10 bx/cs |
| 2259   | 5.1cm Dia. Monitoring Electrodes w/abrader 25/bg, 40 bg/cs             |
| 2259-3 | 5.1cm Dia. Monitoring Electrodes w/abrader, 3/bg, 20 bg/bx, 10 bx/cs   |
| 9640   | 5.1cm Dia. Monitoring Electrodes w/o abrader, 50/bg, 20 bg/cs          |



### 3M™ RED DOT™ INFANT SOFT CLOTH MONITORING ELECTRODES

2258-3 Infant, 3.2cm Dia., Soft Cloth Electrode w/o abrader, 3/bg, 20 bg/bx, 10 bx/cs



### 3M™ RED DOT™ MONITORING ELECTRODES WITH FOAM TAPE & STICKY GEL

Multi-purpose monitoring electrodes with foam backing and sticky gel are ideally suited for the OR, ER and diagnostic stress testing applications. Conveniently packaged in individual strips of ten electrodes for easy use and less waste. Remains fresh for 30 days out of the bag.

|        |  |
|--------|--|
| 2560   | Monitoring Electrode w/o abrader, 10/strip, 5 strips/bg, 20 bg/cs              |
| 2570   | Monitoring Electrode w/abrader, Radiolucent Stud, 50/bg, 20 bg/cs              |
| 2570-3 | 4cm x 3.5cm Monitoring Electrodes w/abrader, Radiolucent Stud, 3/bg, 200 bg/cs |
| 2570-5 | 4cm x 3.5cm Monitoring Electrodes w/abrader, Radiolucent Stud, 5/bg, 200 bg/cs |



### 3M™ RED DOT™ MONITORING ELECTRODES WITH MICROPORE™ TAPE BACKING

A round-style monitoring electrode with 3M™ Micropore™ Surgical Tape backing. Patented solid gel penetrates the skin quickly allowing high quality. Low chloride content minimizes skin irritation. Radiolucent to minimize the need to replace or remove electrodes in most radiological imaging applications.

|        |   |
|--------|---|
| 2239   | 6cm Dia. Monitoring Electrodes w/o abrader, 50/bg, 20 bg/cs                     |
| 2249   | 6cm Dia. Monitoring Electrodes w/abrader, 25/bg, 40 bg/cs                       |
| 2249-3 | 6cm Dia. Monitoring Electrodes w/abrader, 3/bg, 20 bg/bx, 10 bx/cs              |
| 2248   | Pediatric, 4.4cm Dia., Micropore Paper Tape Electrode, 25/bg, 40 bg/cs          |
| 2248-3 | Pediatric, 4.4cm Dia., Micropore Paper Tape Electrode, 3/bg, 20 bg/bx, 10 bx/cs |

### 3M™ RED DOT™ REPOSITIONABLE MONITORING ELECTRODES

|        |  |
|--------|--|
| 2660-3 | Repositionable Monitoring Electrodes, 1.56" x 1.25", 3/bg, 200 bg/cs |
| 2660-5 | Repositionable Monitoring Electrodes, 1.56" x 1.25", 5/bg, 200 bg/cs |
| 2670-3 | Repositionable Monitoring Electrodes, 1.56" x 1.25", 3/bg, 200 bg/cs |
| 2670-5 | Repositionable Monitoring Electrodes, 1.56" x 1.25", 5/bg, 200 bg/cs |



**3M™ RED DOT™ RESTING MONITORING ELECTRODES**

Resting electrodes for 12-lead EKG applications.

- 2330 2.2cm x 2.2cm Monitoring Electrodes , Tab Style, 10/strip, 10 strips/bg, 10 bg/bx, 4 bx/cs
- 2352 1 3/4" x 7/8" Monitoring Electrodes , Snap Style, 10/bg, 100 bg/cs
- 9650 2.2cm x 2.2cm Monitoring Electrodes , Tab Style, 100/bg, 10 bg/bx, 4 bx/cs
- 2360 Monitoring Electrode w/conductive adhesive, Tab Style, 10/card, 10 cards/bg, 40 bg/cs



**3M™ RED DOT™ SOFT CLOTH MONITORING ELECTRODES**

A round-style monitoring electrode with soft cloth backing for general ECG monitoring. Patented solid gel penetrates the skin quickly allowing high quality. Low chloride content minimizes skin irritation.

- 2238 6cm Dia. Monitoring Electrodes w/o abrader, 50/bg, 20 bg/cs
- 2238-3 6cm Dia. Monitoring Electrodes w/o abrader, 3/bg, 20 bg/bx, 10 bx/cs
- 2238-5 6cm Dia. Monitoring Electrodes w/o abrader, 5/bg, 20 bg/bx, 10 bx/cs
- 2255 6cm Dia. Monitoring Electrodes w/abrader, 25/bg, 40 bg/cs
- 2255-3 6cm Dia. Monitoring Electrodes w/abrader, 3/bg, 20 bg/bx, 10 bx/cs
- 9641 6cm Dia. Monitoring Electrodes w/o abrader, 50/bg, 20 bg/cs



**3M™ RED DOT™ TRACE PREP**

A gentle abrasive material designed to remove the non-conductive skin layer to achieve a better trace quality.


- 2236 3/4" x 196" Trace Prep w/dispenser, 36/cs



**3M™ RED DOT™ WET GEL FOAM MONITORING ELECTRODES**

Tear-drop style monitoring electrodes with wet gel. Low chloride content minimizes skin irritation.

- 2246 7.2cm Dia. Monitoring Electrodes w/abrader & Micropore Tape, 25/bg, 40 bg/cs
- 2256 5.1cm Dia. Monitoring Electrodes w/abrader, 25/bg, 40 bg/cs



**BURDICK BLUEMAX ELECTRODES, HOLTER/STRESS**

Offset snap.

- 097001 40mm Dia. Electrodes, Vinyl, 25/pch, 40 pch/cs
- 097004 50mm Dia. Electrodes, Foam, 10/pch, 100 pch/cs
- 097010 68mm Dia. Electrodes, Micropore, 25/pch, 40 pch/cs
- 097014 30mm Dia. Electrodes, Foam, 10/pch, 200 pch/cs
- 015-0630-00 Blue Max SUPATAB Resting ECG Electrode, Wet Gel, 100/pk, 24 pk/cs
- 097002 50mm Dia. Electrodes, Foam, 25/pch, 40 pch/cs



**BURDICK CARDIOSENS™ DISPOSABLE ELECTRODE SENSORS**

Three styles of disposable sensors: CardioSens™/Ultra II has aggressive adhesive. Fast signal stabilization with rapid defibrillator recovery and superior conductivity. CardioSens™/K has an exceptionally high conductive performance with natural Karaya adhesive. CardioSens™/Advantage has a narrow design, special formulation adhesive and unique, high conductivity tin backing.

- 047029 CardioSens/Ultra II, 500/bx, 10 bx/cs
- 047866 CardioSens/K, 500/bx, 10 bx/cs
- 047881 CardioSens/Advantage (replaces 047057), 500/bx, 10 bx/cs

**Accessories**

- 704-0102-00 Clear Choice Lead Adapter Pin to Universal Clip, 10/pk



### BURDICK HEARTLINE® MONITORING ELECTRODES

Designed to deliver high quality, easy-to-read tracings and long-lasting adhesion for consistent, dependable performance. Choose from give versatile electrodes, each offering a variety of performance possibilities.

|        |   |
|--------|---|
| 007778 | HeartLine Comfort Tape Wet Gel Electrodes, 48mm, 30/pouch, 600/cs |
| 007781 | HeartLine DuraFoam Wet Gel Electrodes, 48mm, 30/pouch, 600/cs     |
| 007782 | HeartLine Radiolucent Universal Electrode, 30/pouch, 600/cs       |
| 007784 | HeartLine SensaCloth Electrodes, 30mm, 30/pouch, 600/cs           |

### GE MEDICAL DISPOSABLE RESTING ECG ELECTRODES

Latex free. Rapid, stable baseline. Reduces potential for cross-contamination. Secure adhesion. Silver/silver chloride construction.

|             |   |
|-------------|---|
| 9623-003P   | Silver Mactrode® Plus Disposable Resting ECG Electrodes, Tab, 10/card, 10 card/pk, 10 pk/cs                           |
| 2014778-001 | SilverTRACE® Plug Disposable Resting ECG Electrodes, Plug B4MO Foam Rectangle, Banana Socket, 30/pch, 10 pch/bx       |
| 402207-210  | BabyMAC® Pediatric Disposable Resting ECG Electrodes, Tab, 10/card, 10 card/pk, 5 pk/cs                               |
| 900703-310  | SilverTRACE® Disposable Resting ECG Electrodes, Adult Foam Rectangle, Wet Gel, Snap, 5/strip, 2 strip/pch, 160 pch/cs |
| 9623-003    | Silver Mactrode® Disposable Resting ECG Electrodes, Tab, 10/card, 10 card/pk, 10 pk/cs                                |
| 9623-103P   | Silver Mactrode® Plus Disposable Resting ECG Electrodes, Tab, 14/card, 10 card/pk, 10 pk/cs                           |

### GE MEDICAL EXERCISE/STRESS, HOLTHER, TELEMETRY SNAP ELECTRODES

Latex free foam and vinyl tape electrodes Two-piece snap connectors. Silver/silver chloride electrodes with wet or solid get.

|             |   |
|-------------|---|
| 2014768-001 | SilverTRACE® Multi P21MO, Disposable Adult ECG Electrodes, Foam Oval w/tab, Solid Gel, 30/pch, 10 pch/bx      |
| 2014775-001 | SilverTRACE® First P28MO, Disposable Adult ECG Electrodes, Wet Gel, 30/pch, 10 pch/bx                         |
| 2014777-001 | SilverTRACE® Soft P55MO, Disposable Adult ECG Electrodes, Foam Round, Wet Gel, 30/pch, 10 pch/bx              |
| 2014780-001 | SilverTRACE® Soft Stress PS50MO, Disposable Adult ECG Electrodes, Foam Round, Wet Gel, 30/pch, 10 pch/bx      |
| 2014783-001 | SilverTRACE® Window P50TR, Disposable Adult ECG Electrodes, Vinyl Tape Round, Solid Gel, 30/pch, 10 pch/bx    |
| 2014784-001 | SilverTRACE® Window P50TR, Disposable Adult ECG Electrodes, Vinyl Tape Round, Solid Gel, 5/strip, 60 strip/bx |
| 2014785-001 | SilverTRACE® Window P50TR, Disposable Adult ECG Electrodes, Vinyl Tape Round, Solid Gel, 7/pch, 50 pch/bx     |
| 900703-205  | SilverTRACE® Disposable Adult ECG Electrodes, Foam Oval, Wet Gel, 5/strip, 30 strip/bx, 4 bx/cs               |
| 900703-230  | SilverTRACE® Disposable Adult ECG Electrodes, Foam Oval, Wet Gel, 30/pch, 10 pch/cs                           |
| 9653-305    | SilverTRACE® Disposable Adult ECG Electrodes, Clear Vinyl Tape, Solid Gel, 5/pch, 60 pch/cs                   |
| 9653-507    | SilverTRACE® Disposable Adult ECG Electrodes, Clear Vinyl Tape, Solid Gel, 7/pch, 50 pch/cs                   |

### GE MEDICAL MONITORING/TELEMETRY ELECTRODES

Latex free. Conductive gel.

|             |   |
|-------------|---|
| 2014779-001 | SilverTRACE® Micro P55MI Disposable Adult ECG Electrodes, Paper Round, Wet Gel, 30/pch, 10 pch/bx           |
| 2014781-001 | SilverTRACE® Soft Stress PS50MO Disposable Adult ECG Electrodes, Foam Round, Wet Gel, 5/strip, 120 strip/bx |
| 2014782-001 | SilverTRACE® Tissair P40CL Disposable Adult ECG Electrodes, Cloth Round, Solid Gel, 30/pch, 10 pch/bx       |
| 900690-030  | SilverTRACE® Disposable Adult ECG Electrodes, Foam Round w/tab, Solid Gel, 30/pch, 10 pch/cs                |
| 900690-403  | SilverTRACE® Disposable Pediatric ECG Electrodes, Cloth Rectangle, Solid Gel, 3/bg, 15 bg/pch, 10 pch/cs    |
| 900690-404  | SilverTRACE® Disposable Neonatal ECG Electrodes, Cloth w/preattached leadwires, Solid Gel, 3/bg, 100 bg/cs  |
| 9431-004    | SilverTRACE® Disposable Adult ECG Electrodes, Foam Round, Wet Gel, 30/pch, 10 pch/cs                        |

### GE MEDICAL QWIK CONNECT® PLUS DISPOSABLE ATTACHMENT PADS

Reduces potential for cross-contamination. Nickel-plated stud insures good conductivity. Eliminates the need for uncomfortable leg plate adapter straps. Latex free. Secure adhesive.

|         |  |
|---------|--|
| 2464AAO | Spiral Leg Plate Attachment Pads, 1" x 3", 50/cs |
|---------|--|

### GE MEDICAL QWIK CONNECT® PLUS FETAL SPIRAL ELECTRODES

Locking mechanism secures the electrode during attachment. Electrode wires stabilized inside the guide tube until unlocked by the clinician. Latex free. Sterile.

|         |                               |
|---------|-------------------------------|
| 7000AAO | Fetal Spiral Electrode, 50/cs |
|---------|-------------------------------|

### GE MEDICAL RADIOTRANSLUCENT ELECTRODES

Latex free polyethylene foam base. Carbon snap is x-ray translucent and non-magnetic. Conductive wet gel provides clear ECG traces.

|              |  |
|--------------|--|
| 2014786-001  | SilverTRACE® First Carbon P28MO, Disposable Radiotranslucent ECG Electrodes, Foam Rectangle, Wet Gel, 30/pch, 10 pch/bx    |
| 2014787-001  | SilverTRACE® First Carbon P28MO, Disposable Radiotranslucent ECG Electrodes, Foam Rectangle, Wet Gel, 5/strip, 60 strip/bx |
| 320-1534-000 | SilverTRACE® Disposable Adult ECG Electrodes, Cloth w/preattached leadwires, Solid Gel, 3/pch, 10 pch/pk                   |
| 9408-401     | SilverTRACE® Disposable Adult Radiotranslucent ECG Electrodes, Foam Oval, Wet Gel, 5/strip, 60 strip/cs                    |

### GE MEDICAL REUSABLE RESTING ECG ELECTRODES & ACCESSORIES

German silver electrodes minimize baseline drift and reduce contact resistance. GE suction electrode perfect fit for 0.156" (4mm) banana plug leadwires. Limb clamps that are latex free, easy to apply, and accepts both banana plug and snap leadwires.

|          |   |
|----------|---|
| 5608-003 | GE Gray Suction Bulb                            |
| 5623-202 | Metal Limb Plate Electrode                      |
| 5880-004 | Metal Snap Electrodes, 5/pk                     |
| 5881-002 | Rubber Boots For Metal Snap Electrodes, 5/pk    |
| 9033-006 | Rubber Limb Strap, Latex Free                   |
| 9033-014 | GE Suction Electrode                            |
| 9033-015 | GE Suction Electrode & Gray Bulb                |
| 9033-401 | Limb Clamps, 4/set                              |
| 9235-001 | Adhesive Tapes For Metal Snap Electrodes 100/bx |






**LEONHARD LANG SKINTACT® DIAGNOSTIC RESTING ECG ELECTRODES**

Designed for resting 12 lead diagnostic ECG's. Easitabs are disposable eliminating the spread of cross infection associated with older style suction cups. The C-LINE technology ensures that every electrode will be fresh providing exceptional consistency in quality of recordings at very affordable prices.

- FSWA00 Diagnostic ECG Electrodes, Foam, Ultra-Tack Adhesive, Wet Gel, 500/bx
- RT-24 Diagnostic ECG Electrodes, All-Purpose, 1000/bx
- RT-41 Diagnostic ECG Electrodes, Ultra-Tack Adhesive, Oversize, 1000/bx, 4 bx/cs

**Accessories**


- 58100 EasiClip (Universal Alligator Clip), 10/bx
- 58727 EasiPrep (FingerTip Abrasion Strips), 500/bx



**LEONHARD LANG SKINTACT® STRESS TEST & HOLTER ECG ELECTRODES**

Designed to provide the highest quality tracings for stress test procedures. The offset design and Aqua-Set gel significantly reduce any motion artifacts ensuring a very clear stable trace even under the most strenuous conditions. The C-LINE technology ensures that every electrode will be fresh providing exceptional consistency in quality of recordings at very affordable prices.


- FS-40 40mm Stress Test ECG Electrode, Aqua-Wet Gel, Foam Backing, Diaphoretic, Hypoallergenic, 1200/bx
- FS-50 50mm Stress Test ECG Electrode, Aqua-Wet Gel, Snap Connector, Foam Backing, Enhanced Adhesive, 1200/bx
- FS-50C 50mm Stress Test ECG Electrode, Radio Translucent, Foam w/carbon stud, Diaphoretic, Hypoallergenic, 1200/bx
- FS-RC-4 35 x 45mm Stress Test ECG Electrode, Aqua-Wet Gel, Foam Backing, Diaphoretic, Hypoallergenic, 600/bx
- FS-TBI 50mm Stress Test ECG Electrode, Solid Gel, Lift Tab, Foam Backing, Diaphoretic, Hypoallergenic, 1500/bx
- FSVB01 36 x 53 mm Stress Test ECG Electrode, Aqua-Wet Gel, Offset Snap Connector, Foam Backing, Enhanced Adhesive, 1200/bx
- FSVB011 35 x 53 mm Stress Test ECG Electrode, Solid Gel, Offset Snap Connector, Foam Backing, Enhanced Adhesive, 1500/bx
- FS-RC-3 35 x 45mm Stress Test ECG Electrode, Aqua-Wet Gel, Foam Backing, Diaphoretic, Hypoallergenic, 3/strip, 1200/bx
- FS-RC-5 35 x 45mm Stress Test ECG Electrode, Aqua-Wet Gel, Foam Backing, Diaphoretic, Hypoallergenic, 5/strip, 1200/bx



**NIKOMED DIAGNOSTIC TAB ELECTRODES**

The Nikotab System with Nikomed diagnostic ECG is easy. Unique design allows for quick and easy prep. The highest quality adhesive and innovative shape prevent twisting and deliver the greatest percentage of first time recordings. Wide-body design holds Nikotab absolutely flat, securing it in place to eliminate the twist. This product is latex free.


- 0315 Economy Size Electrodes, 100/pk, 500/bx, 5000/cs
- 0515 Nikotabs E Synthetic Electrodes, Foil, 100/pk, 500/bx, 5000/cs
- 0715 Nikotabs II Synthetic Electrodes, Foil, Larger Size, 100/pk, 500/bx, 5000/cs
- 0815 Pediatric Tab Electrodes, 100/pk, 500/bx, 5000/cs



**NIKOMED TRACE1™ SOLID GEL MONITORING ELECTRODES**

This product is latex free.

- 1015 Solid Gel Foam Monitoring Electrode, Compact Design, Super Tack Adhesive, 32 x 38mm, 6 per strip, 60/pk, 300/bx
- 2001 Specialty Nonwoven Monitoring Electrodes, 7/8" x 7/8", 4/pk, 400/bx
- 2002 Specialty Nonwoven Monitoring Electrodes, 1 7/8" x 7/8", 4/pk, 400/bx
- 2010 Solid Gel Stress/Holter Monitoring Foam Electrodes, Eyeball Shape, 50mm, Helps Eliminate Artifacts, 30/pk, 300/bx
- 2012 Solid Gel Monitoring Electrodes, Clear Tape, Radio Translucent X-Ray Opaque, 38 x 38mm, 30/pk, 300/bx
- 2013 Solid Gel Monitoring Electrodes, Non-Woven Cloth, Helps Eliminate Artifacts, 55mm, 30/pk, 300/bx
- 2014 Solid Gel Monitoring Electrodes, Clear Tape, Multi-Purpose, 38 x 38mm, 30/pk, 300/bx
- 2015 Solid Gel Foam Monitoring Electrodes, 35 x 45mm, Helps Eliminate Artifacts, 30/pk, 300/bx
- 2016 ICU/CCU Monitoring Cloth, 55mm, 30/pk, 300/bx
- 2054 Solid Gel, Clear Tape Monitoring Electrode, Kidney Shape Design, Super Tack Adhesive, 34 x 50mm, 30/pk, 300/bx
- 5520 Solid Gel Foam Monitoring Electrode, Multi-Snap Design, Super Tack Adhesive, 32 x 45mm, 30/pk, 300/bx
- 5521 Solid Gel Foam Monitoring Electrode, Snap-Tab Design, Super Tack Adhesive, 32 x 45mm, 30/pk, 300/bx
- 5525 Solid Gel Foam Monitoring Electrode, Dynamic Off Center Design, Super Tack Adhesive, 32 x 45mm, 30/pk, 300/bx



**NIKOMED TRACE1™ WET GEL MONITORING ELECTRODES**

This product is latex free.

- 4500 Wet Gel, Foam Stress/Holter Monitoring Electrode, Super Tack Adhesive, 50mm, 60/pk, 300/bx
- 4530 Wet Gel, Foam Direct Connect Snap Monitoring Electrode, Patient Friendly Adhesive, 32 x 45mm, 5 per strip, 60/pk, 300bx
- 4533-5 Wet Gel, Foam Monitoring Electrode, Super Tack Adhesive, 32 x 45mm, 5 per strip, 60/pk, 300/bx
- 4535 Wet Gel, Foam Direct Connect Snap Monitoring Electrode, Patient Friendly Adhesive, 32 x 45mm, 3 per strip, 60/pk, 300bx
- 4540 Wet Gel, Foam Direct Connect Snap Monitoring Electrode, Patient Friendly Adhesive, 32 x 45mm, 5 per strip, 60/pk, 300bx
- 4555 Wet Gel, Foam Monitoring Electrode, Super Tack Adhesive, 45mm, 5 per strip, 10/pk, 100/bx
- 4560 Wet Gel, Foam Stress Monitoring Electrode, Super Tack Adhesive, 38 x 60mm, 60/pk, 300/bx
- 5500 Stress/Holter Monitoring Electrodes, 50mm, Super-Tack Disphoretic Adhesive, Eyeball Shape, 30/pk, 300/bx
- 5511 General Purpose Monitoring Electrodes, Clear Tape, Diaphoretic Adhesive, 50mm, 30/pk, 300/bx
- 5533-3 General Purpose Monitoring Electrodes, 32 x 44mm, Super-Tack Diaphoretic Adhesive, 3 Per Strip, 30/pk, 300/bx
- 5533-5 General Purpose Monitoring Electrodes, 32 x 44mm, Super-Tack Diaphoretic Adhesive, 5 Per Strip, 30/pk, 300/bx
- 5535 Diagnostic Resting ECG Monitoring Electrodes, Pediatric Monitoring, 32 x 44mm, Easy-Tack Patient-Friendly Adhesive, 30/pk, 300/bx
- 5539-3 Thallium Stress Monitoring Electrodes, Cardiac Cath., 28 x 44mm, Super-Tack Diaphoretic Adhesive, Radio-Translucent X-Ray Opaque Snap/Sensor System, 3 Per Strip, 30/pk, 300/bx
- 5539-5 Thallium Stress Monitoring Electrodes, Cardiac Cath., 28 x 45mm, Super-Tack Diaphoretic Adhesive, Radio-Translucent X-Ray Opaque Snap/Sensor System, 5 Per Strip, 30/pk, 300/bx
- 5550 Wet Gel, Foam Monitoring Electrode, Kidney Shape Design, Super Tack Adhesive, 34 x 50mm, 30/pk, 300bx
- 5559 Wet Gel, Foam Monitoring Electrode, Kidney Shape Design, Radio-Translucent, Super Tack Adhesive, 34 x 50mm, 30/pk, 300bx
- 5560 Stress Monitoring Electrodes, 38 x 59mm, Super-Tack Disphoretic Adhesive, 30/pk, 300/bx



**CHESTER GEL WARMER**

Chester Labs Gel Warmers are the perfect companion for any of the Transmission Media Gel products. Available in 2 and 3 bottle warming units, these durable warmers are made from stainless steel construction with a white powder-coat exterior. The adjustable temperature control maintains contents just above body temperature. Equipped with an illuminated on/off rocker switch and a six foot power cord, these warmers operate for just pennies a day and can be used on table tops or mounted to walls.

- 1205-2 2 Bottle Gel Warmer w/temperature adjustment, White Steel Case
- 1205-3 3 Bottle Gel Warmer w/temperature adjustment, White Steel Case



**CHESTER LIQUASONIC® ULTRASOUND GEL**

LiquaSonic Ultrasound Transmission Gel has a thicker, no-slip formula that remains viscous even when applied to salty skin. Our salt-free, formaldehyde free formula eliminates the air between the transducer and skin, offering exceptional coupling efficiency for more effective treatment/diagnosis. LiquaSonic Gel has been formulated to be acoustically correct for the broad range of frequencies used in medical ultrasound applications. Our formula is water soluble and non-greasy preventing damage to expensive transducers and clothing.

- 001205 8 oz bottle w/mustard cap, Ultrasound Gel, 12/cs
- 001222 5 Liter Collapsible Cubitainer w/empty bottle, Ultrasound Gel, 4/cs
- 001297 128 oz Multi-Purpose, Medium Viscosity Ultrasound Gel, 4/cs
- 001298 8 oz Bottle w/mustard cap, Multi-Purpose, Medium Viscosity Ultrasound Gel, 12/cs
- 001305 8 oz Bottle w/mustard cap, Clear Ultrasound Gel, 12/cs
- 001333 128 oz Clear Ultrasound Gel, 4/cs



**CHESTER LIQUASONIC® ULTRASOUND LOTION**

Similar in use to the ultrasound gel, the lotion is often preferred by medical professionals. The aqueous formulation is pleasing to the patient and is acoustically correct for the large spectrum of frequencies used in medical ultrasound applications. The lotion is hypoallergenic, bacteriostatic, salt free, formaldehyde free and will not harm clothing or equipment.

- 001267 128 oz Ultrasonic Lotion, 4/cs
- 001268 8 oz Bottle w/mustard cap, Ultrasonic Lotion, 12/cs



## CHESTER LIQUASONIC® ULTRASOUND SCANNING GEL

This gel has been designed for diagnostic ultrasound scanning procedures requiring a less viscous formulation. Unlike other scanning materials, Ultrasound Scanning Gel possesses characteristics of a gel and fluid. Like a gel, it is applied with ease and like a fluid it has high-quality scanning abilities. The scanning gel is an acceptable alternative to mineral oil, is acoustically correct, completely aqueous and hypoallergenic.

- 001277 128 oz Scan Gel, 4/cs
- 001278 8 oz Bottle w/mustard cap, Scan Gel, 12/cs

## GE MEDICAL SKIN PREP PRODUCTS & ELECTRODE ACCESSORIES

- 3704-901 Prep Razors, 100/pk
- 418195-100 GE Electrode Clips For Stress Loops, 100/pk
- 4828-004 Alcohol Pumice Prep Pads, 100/bx
- 4828-005 Sandpaper Prep Strips, 100/bx
- 4829-002 Pre-Cut Tape Strips For Stress Loops, 6/card, 100 card/pk
- 9386-001 Dry Skin Prep Pads, 200/pk

## PARKER LABS REDUX® ELECTROLYTE CRÈME

For routine ECG. Cosmetic quality, pleasing to patients. Non-irritating. Highly conductive electrode cream.

- 66-04 Electrolyte Creme, 114g (4 oz) Tube, 10/bx, 5 bx/cs

## PARKER LABS REDUX® ELECTROLYTE GEL

Non-abrasive, water-based gel for routine ECG and defibrillation. Highly conductive. Water soluble, non-staining, non-gritty.

- 65-04 Electrolyte Gel, 114g (4 oz) Tube, 10/bx, 5 bx/cs

## PARKER LABS REDUX® ELECTROLYTE PASTE

Water-based, abrasive gel for short term monitoring, ECG and defibrillation. Highly conductive and water soluble.

- 67-05 Electrolyte Paste, 142g (5 oz) Tube, 10/bx, 5 bx/cs



## PARKER LABS SIGNACREME® ELECTRODE CREAM

Highly conductive cream electrolyte for use in electromedical procedures. Bacteriostatic, non-irritating and hypoallergenic. Cosmetic quality; pleasing to patients.

- 17-05 Electrode Cream, 5 oz Dispenser, 12/bx, 6 dz/cs
- 17-20 Electrode Cream, 2 Liter (70 oz) w/2 dispensers & 1 Dispenser Pump, 4 bx/cs



## PARKER LABS SIGNAGEL® ELECTRODE GEL

Highly conductive, multi-purpose electrolyte is an ideal saline electrode gel. Recommended for all electromedical procedures. Bacteriostatic, non-irritating, non-gritty. Water-soluble and non-staining.

- 15-25 Electrode Gel, 250g Tube, 6 dz/cs
- 15-60 Electrode Gel, 60g Tube, 12 dz/cs



## PARKER LABS SIGNASPRAY® ELECTRODE SOLUTION & SKIN PREP

Highly conductive, multi-purpose spray electrolyte and non-gritty skin prep which provides superior conduction between skin and electrodes. For ECG and muscle stimulation. Assures clean, artifact-free ECG tracings. Bacteriostatic, hypoallergenic, non-staining. No residue build-up.

- 18-04 Signaspray Pac® Electrode Solution w/4 gallon containers, 2 Spray Dispenser Bottles, 1 Pump, 1/cs
- 18-25 Electrode Solution, 250-ml Spray Dispenser Bottle, 12/bx, 6 dz/cs
- 18-28 Electrode Solution, Gallon w/spray dispenser bottle, 4/cs




## TRIAD® ULTRASOUND GEL

Conducts a broad range of frequencies, does not allow air to intrude between the transducer and skin and requires only minimal pressure to make contact. It is water soluble, non-irritating and hypoallergenic. It will not stain skin or clothing.

- 10-9208 8 oz Ultrasound Gel, 12/cs
- 10-9228 Gallon Ultrasound Gel w/dispenser, 4/cs





**ADC® EKG CALIPER**

A precision crafted instrument of hardened steel for accurate measurement of tracing intervals. Point protector has built-in 40mm metric scale. Measures 4 3/4"


- 394 EKG Ruler
- 395 EKG Caliper



2121



EZ-123



NIK-10

**NIKOMED ECG ACCESSORIES**

- 2121 Niko-Prep, 10/card, 500/cs
- 2121-E Niko-Prep, Individual, 500/cs
- EZ-123 Diagnostic Enhancement Kit (1 pk NIK-10 [choice of 0315, 0515, 0715], 1 pk NIKOtabs, 1 card 2121)
- NIK-10 Niko-Clip, Green, 10/pk
- NIK-10W Niko Clip Adapter, Fits 3-4mm, Modified, White, 10/pk



**WESTERN SURGILAST® STRESS VEST**

Pre-cut stretchable net material holds electrodes in place during electrocardiograms and stress tests. Does not restrict patient's circulation or breathing. No tapes or adhesives required.

- GL750 Small/Medium Stress Vest, 20/bx
- GL752 Large/X-Large Stress Vest, 20/bx



**GE MEDICAL COROMETRICS® 170 SERIES FETAL MONITORS**

Represents a new generation of fetal monitors designed to meet prenatal monitoring needs for the physician's office. The monitor features large LED displays allowing easy, at-a-glance confirmation of fetal heart rate and uterine activity. Offers 9-Crystal Pulsed Doppler Ultrasound technology providing a wide-beam profile for high-quality fetal heart rate detection. Includes standard Z-fold paper and an interface for the optional Model 146 Fetal Acoustic Stimulator. Available in single or twin models.

- 0171WAT-B Model 171 External Fetal Monitor w/Nautilus Ultrasound & Toco Transducers (loop style)
- 0172WAT-B Model 172 Twins External Fetal Monitor w/Nautilus Ultrasound & Toco Transducers (loop style)

**GE MEDICAL FETAL CHART PAPERS**

Designed for GE Corometrics monitors. Serialized pages for case documentation. Thermal array provides accurate, vivid tracings.

- 4305CAO Fetal Chart Paper, Red Grid, 30-240 BPM, Serialized Pages, Z-Fold, 160 sheets/pk, 40 pk/cs
- 4305AAO Fetal Chart Paper, Red Grid, 30-240 BPM, Serialized Pages, Z-Fold, 160/pk, 40 pk/cs
- 4305BAO Fetal Chart Paper, Red Grid, 50-210 BPM, Serialized Pages, Z-Fold, 160/pk, 40 pk/cs
- 4305DAO Fetal Chart Paper, Red Grid, 50-210 BPM, Serialized Pages, Z-Fold, 160/pk, 40 pk/cs
- 4305FAO Fetal Chart Paper, Red Grid, 50-210 BPM, Serialized Pages, Additional 1 Minute Time Marks, Z-Fold, 160/pk, 40 pk/cs
- 4306CAO Fetal Chart Paper, Green Grid, 50-210 BPM, Serialized Pages, Z-Fold, 150/pk, 40 pk/cs
- 4482AAO Fetal Chart Paper, 30-240 BPM, Serialized Pages, Pre-Loaded on Cartridges, Z-Fold, 17 cart/cs
- 4483AAO Fetal Chart Paper, Light Red Grid, 30-240 BPM, Serialized Pages, 350 sheets/rl, 6 rl/cs
- 4483BAO Fetal Chart Paper, Light Red Grid, 50-210 BPM, Serialized Pages, 350 sheets/rl, 6 rl/cs
- 4528AAO NST Patient Record Forms, For Recording From GE Corometrics Fetal Monitor Model 140, Compatible w/4482AAO recording paper, 300/cs
- 4528BAO Antepartum NST Patient Record Forms, For Recording From GE Corometrics Fetal Monitors Models 145 & 155, Compatible w/4483AAO & 4483BAO recording paper, 150/cs
- 4914AAO Chart Guard Label Packer For Model 115, 116, 118, 170, 150, 151, 120 Series Fetal Monitor Recordings, 200/cs
- 5231AAO Fetal Chart Paper, Red Grid, Serialized Pages, Z-Fold, 20 pk/cs
- 5235AAO Fetal Chart Paper, Red Grid, Serialized Pages, Z-Fold, 20 pk/cs
- 5236AAO Fetal Chart Paper, Red Grid, Serialized Pages, Z-Fold, 10 pk/cs
- 9801-012 Hewlett-Packard Fetal Chart Paper Equivalent M1910A, 2 Fetal Channels, Red Grid, Black Trace, 5.97" x 3.81", Z-Fold, 40 pk/cs
- 9801-058 Hewlett-Packard Fetal Chart Paper Equivalent M1910A, 2 Fetal Channels, Red Grid, Black Trace, 5.96" x 46 ft, Z-Fold, 40 pk/cs

**GE MEDICAL FETAL BLOOD SAMPLING KITS**

Latex free. Designed for the GE Corometrics Model 220 pH monitor. Sterile. Non-returnable and non-refundable.

- 4520AAO Fetal Blood Sampling Kit Contains: Endoscope, 6 Swabs, Lancet w/handle, 7 Heparinized Capillary Tubes, Wax Sealant, Capillary Tube Holder, 3ml Silicone Gel



**MALLINCKRODT RENAISSANCE® II SPIROMETRY SYSTEM**

Offers a range of sophisticated, yet user-friendly features for diagnosing respiratory function in a variety of clinical settings. Clinicians will appreciate the intuitive design and technological capabilities, including a true, real-time graphic display that enable tests to be viewed as patients perform them. Measures a broad selection of physiological parameters and provides color printout reports. It also has a 1,000-patient memory. Five-year warranty exchange program.

- P-000700-00 PB700 Portable Spirometer
- P-000710-00 Base Station
- P-000715-00 Basic Starter Kit, Includes: PB700 Spirometer, PB710 Base Station, AC Adapter, (10) FSII Flow Sensors & Syringe Adapter
- P-000720-00 Basic Starter Kit w/printer, Includes: PB700 Spirometer, PB710 Base Station, AC Adapter, (10) FSII Flow Sensors & Syringe Adapter
- P-000722-00 Kit w/printer, Includes: Basic Starter Kit & Data Flow Data Management Software, Plus Compatible Printer
- P-000725-00 Kit w/syringe, Includes: Basic Starter Kit Plus VS300 Calibration Syringe
- P-000730-00 Complete Kit, Includes: Basic Starter Kit Plus VS300 Calibration Syringe, Compatible Printer
- P-000732-00 Complete Kit, Includes Basic Starter Kit & Data Flow Data Management Software, Plus VS300 Calibration Syringe, Compatible Printer



**MALLINCKRODT FLOW SENSORS**

A low cost, disposable, single-patient-use flow sensor that minimizes the possibility of cross-contamination and reduces staff clean-up time.

- P-000200-00 FS200 Flow Sensor, 50/cs
- P-000201-00 FS200 Flow Sensor, 250/cs
- P-000250-00 BD250 Flow Sensor, 25/cs
- P-000251-00 BD250 Flow Sensor, 100/cs
- P-000263-00 FSII Flow Sensors, 50/cs
- P-000265-00 FSII Flow Sensors, 250/bx



**MALLINCKRODT SIMPLICITY® SPIROMETRY SYSTEM**

Offers a simple, affordable solution for performing basic respiratory function screening. The lightweight, handheld unit is portable and easy to use, yet performs all basic spirometry measurements. It also provides automated comparison to predicted values, pre- and post-medication comparison, interpretation of test results and comprehensive test quality control. Equipped with a 999-patient memory, and a blow and go option that allows clinicians to conduct quick, on-the-spot screening without entering patient data.

- P-000500-00 NPB-500 Spirometer
- P-000510-00 Spirometer Base Station
- P-000515-00 Starter Kit, Includes: NPB-500 Spirometer, Base, Flow Sensors, AC Adapter, Syringe Adapter & (2) Nose Clips
- P-000525-00 Starter Kit w/syringe, Includes: NPB-500 Spirometer, Base, (10) flow Sensors, Calibration Syringe, AC Adapter, Syringe Adapter, (2) Nose Clips
- P-000530-00 Complete Kit, Includes: NPB-500 Spirometer, Base, (10) Flow Sensors, Printer, Calibration Syringe, AC Adapter, Syringe Adapter, (2) Nose Clips

**MALLINCKRODT SPIROMETRY ACCESSORIES**

- 067100 Data Flow Data Management Software
- P-000150-01 Parallel Printer, 1/bx
- P-006548-00 Nose Clips, 25/pk
- P-010701-00 Syringe Adapter
- P-025513-SP Pressure Tube for REN II
- P-062033-00 NPB510/PB-710 Serial Cable
- P-062522-00 Renaissance Battery
- P-062596-00 Parallel Printer Cable, 1/bg
- P-495208-00 AC Adapter, 1/bx
- P-495219-00 NICAD Battery Pack
- P-000300-00 3 Liter Calibration Syringe

**WELCH ALLYN PNEUMOCHECK SPIROMETER ACCESSORIES**

Easy to use. PnuemoCheck will automatically calculate and display the best effort and five most needed values on its LCD screen. For your records, the printer/charger will print the data in either flow/volume or volume/time formats in as little as 20 seconds. The printer/charger recharges the PnuemoCheck.

- 56100 Thermal Paper, 5 rl/bx, 12 bx/cs
- 56120 Pediatric Mouthpieces, 100/bg, 5 bg/cs
- 56121 Adult Mouthpieces, 100/bg, 5 bg/cs
- 56130 Disposable Nose Clips, 100/bg, 4 bg/cs
- 61010 Hose for 3-L Syringe.
- 61020 3-L Calibration Syringe Only
- 61030 3-L Calibration Syringe & Hose
- 72500 7.2-V Rechargeable Batteries (2) - Replacement

**WELCH ALLYN SP-2 SPIROMETER ACCESSORIES**

- 26006-0001 Disposable Flow Tubes For SP-2 Only, 100/bx

**WELCH ALLYN SPIROMETER FLOW TUBES**

- 26001-0000 Disposable Flow Tubes, SP-150/SP-160, 10/bx
- 26001-0001 Disposable Flow Tubes, SP-150/SP-160, 100/bx

**Can't find what  
you are looking for?  
Give us a call for fast,  
personal service.**

**NDD EASYONE™ DIAGNOSTIC SPIROMETRY SYSTEM**

Ultrasonic flow measurement assures fast, accurate, reliable results. No moving parts, No codes to enter. Single patient use spirette has hygienic mouthpiece that reduces risk of cross infection. Alpha numeric keypad/graphic display makes navigation and data entry fast and easy. Graphic curve display allows to review curves without waiting for a printout. Meets ATS 1994 diagnostic standards. Meets NIOSH/OSHA and SSD guidelines. Customizable configuration. Automated coaching messages and QC. Automated interpretation. Two AA batteries. Stores up to 700 test sessions. Optional 8.5 x 11 or A4 printout. PC data transfer. Toll free technical support. Three-year exchange warranty. If your spirometer fails, we will replace if overnight without charge. Extended warranties are activated upon receipt of end-user information and must be received within six months of invoice date.

|        |  |
|--------|--|
| 1000-1 | EasyOne™ Diagnostic Spirometry System II Includes: #2001-2 Diagnostic Spirometry System I, #3104-4 Four Year Extended Warranty, #2040-2 EasyWare Software, #2050-1 Case of 50 Spirettes™ |
| 2000-1 | EasyOne™ Frontline Spirometer Includes: Calibration Adapter & 10 Spirettes™  |
| 2001-1 | EasyOne™ Diagnostic Spirometer Includes: Calibration Adapter & 10 Spirettes™   |
| 2001-2 | EasyOne™ Diagnostic Spirometry System I Includes: #2001-1 Diagnostic Spirometer, #2010-4 Cradle, #2020-1Color InkJet Printer   |
| 2000-2 | EasyOne™ Frontline Spirometry System Includes: #2000-1 Frontline Spirometer, #2010-4 Cradle, #2020-1Color InkJet Printer   |
| 2010-4 | EasyOne™ Cradle Includes: USB Cable  |

**NDD EASYONE™ DIAGNOSTIC SPIROMETRY SYSTEM ACCESSORIES**

|          |   |
|----------|---|
| 2020-1   | Color InkJet Printer  |
| 2030-1   | Calibration Adapter for Spirometer Serial #'s Below 46000   |
| 2030-2   | 3-Liter Calibration Syringe   |
| 2040-1   | EasyWare Software for use with Parallel/Serial Cradle #2010-1   |
| 2050-1   | Spirettes™, 50/cs   |
| 2050-5   | Spirettes™, 200/cs  |
| 2060-2   | Parallel Cable  |
| 3101-1   | 1 Year Warranty Extension EasyOne™ Diagnostic Spirometer (cradle included at no charge) *Extended Warranties are activated upon receipt of end-user information and must be received within 6 months of invoice date. |
| 3101-2   | 2 Year Warranty Extension EasyOne™ Diagnostic Spirometer (cradle included at no charge) *Extended Warranties are activated upon receipt of end-user information and must be received within 6 months of invoice date. |
| 3101-3   | 3 Year Warranty Extension EasyOne™ Diagnostic Spirometer (cradle included at no charge) *Extended Warranties are activated upon receipt of end-user information and must be received within 6 months of invoice date. |
| 3101-4   | 4 Year Warranty Extension EasyOne™ Diagnostic Spirometer (cradle included at no charge) *Extended Warranties are activated upon receipt of end-user information and must be received within 6 months of invoice date. |
| 3105-1   | 1 Year Warranty Extension EasyOne™ Frontline Spirometer (cradle included at no charge) *Extended Warranties are activated upon receipt of end-user information and must be received within 6 months of invoice date.  |
| 3105-2   | 2 Year Warranty Extension EasyOne™ Frontline Spirometer (cradle included at no charge) *Extended Warranties are activated upon receipt of end-user information and must be received within 6 months of invoice date.  |
| 3105-3   | 3 Year Warranty Extension EasyOne™ Frontline Spirometer (cradle included at no charge) *Extended Warranties are activated upon receipt of end-user information and must be received within 6 months of invoice date.  |
| 3105-4   | 4 Year Warranty Extension EasyOne™ Frontline Spirometer (cradle included at no charge) *Extended Warranties are activated upon receipt of end-user information and must be received within 6 months of invoice date.  |
| 2000-2NP | EasyOne™ Frontline Spirometry System Includes: #2000-1 Frontline Spirometer, #2010-4 Cradle, Does not include Printer   |
| 2001-2NP | EasyOne™ Diagnostic Spirometry System I Includes: #2001-1 Diagnostic Spirometer, #2010-4 Cradle, Does not include Printer   |
| 2001-2S  | EasyOne™ Screen Spirometry System Includes: #2001-1 Diagnostic Spirometer, 2010-6S, EasyOne screen  |
| 2010-6S  | EasyOne™ Screen, Includes: Screen Connector & #2040-2 EasyWare  |
| 2030-3   | Calibration Adapter for Spirometer Serial #'s 46000 and above   |
| 2031-1   | EasyOne™ Carry Case   |
| 2040-2   | EasyWare Software for use with USB Cradle p/n 2010-4  |

**BURDICK SENAIRE HANDHELD SPIROMETER**

Full functioning ATS compliant spirometer that weighs less than 13 ounces. User friendly touch screen provides quick navigation and easy interpretation. Base station allows easy battery recharging and data transfer. Print reports by connecting the spirometer or base station to a computer.

92524-01 Sensaire Handheld Spirometer, English Version, 120V

**BURDICK SPIROCARD SPIROMETER**

Features patented, pre-calibrated, disposable patient filters. Performs all the tests you would expect from a conventional device. Interface to Office Medic software for real-time data access, reporting and archival. Compatible with your favorite off-the-shelf PC, laptop, tablet PC or Compaq iPAQ.

92526-01 SpiroCard PC Card Spirometer, English

**BURDICK SPIROMETRY ACCESSORIES**

|             |  |
|-------------|--|
| 010-1600-00 | Sensaire Rechargeable Battery Pack   |
| 010-1601-00 | Sensaire User's Manual   |
| 010-1604-00 | Sensaire Mouthpieces, 100/bx   |
| 010-1605-00 | SpiroCard/SpiroCard Mouthpieces, 100/bx  |
| 010-1606-00 | SpiroCard/SpiroCard Pressure Tubes, 25/bg  |
| 010-1607-00 | OxiCard Adult Fingercip Sensors  |
| 010-1615-00 | Sensaire Pressure Tube For Calibration Use Only, 10/pkg  |
| 010-1616-00 | Sensaire Stylus, 10/pkg  |
| 010-1628-00 | OxiCard Pediatric Fingercip Sensors  |
| 043231      | Video "How To Operate Flash"   |
| 043232      | Video "How To Operate SpiroSense"  |
| 043245      | Video: How To Operate Spirotouch   |
| 086613      | Spirometry Book: A Better Way Mouthpiece Adapter For Calibration Syringe SpiroTouch, Sensaire, SpiroCard & SpiroCard |
| 100022      | 3.0L Calibration Syringe   |
| 100049      | Presto™ Thermal Paper, 6 r/cs  |
| 100052      | RJ-11 Cable For Use w/Presto   |
| 10005201    | DB-9 Cable For Use w/SpiroSense  |
| 100160      | Presto Operator Manual   |
| 100161      | Plastic Mouthpieces, 100/pkg   |
| 100162      | Noseclips, 5/pkg   |
| 100163      | 12" Calibration Hose For Use w/SpiroSense  |
| 10021718    | 36" Hose & Adapter For Use w/Presto & SpiroSense   |
| 100219      | Flow Sensor Bracket For Use w/Presto & SpiroSense  |
| 100220      | Flow Sensor Insert   |
| 100243      | Cardboard Mouthpieces (Presto & SpiroSense), 100/pkg   |
| 100245      | SpiroSense Operator Manual   |
| 100327      | Flash Memory Card, 256K Memory (40 Tests)  |
| 100338      | Flash Operator Manual  |
| 100340      | Flash Battery Recharger  |
| 100351      | Flash Carry Case   |
| 100362      | Flash Cardboard Mouthpieces, 100/pkg   |
| 100363      | Flash Flow Sensor Hose Replacement   |
| 100381      | Flash Flow Sensor Replacement  |
| 100414      | SpiroSense Battery Booster   |
| 100611      | Flash Compatible Printer w/cable   |
| 100612      | Mouthpiece, 25/pkg   |
| 100613      | Flow Sensor Insert For Spirotouch  |
| 103-0208-00 | Mouthpiece Adapter for Calibration Syringe SpiroTouch, Sensaire, SpiroCard & SpiroCard                               |





**BURDICK SPIROX CARD SPIROMETER/PULSE OXIMETER**

Combines a spirometer and a pulse oximeter on a single PCMCIA card providing two diagnostic tests in one portable package. Performs all measurements performed on conventional devices. Features patented, disposable and biodegradable patient filters. Devices interface to Office Medic software for real-time data access, reporting and archival. Compatible with your favorite off-the-shelf PC, laptop, tablet PC or Compaq iPAQ.

92525-01 SpirOx Card PC Card Spirometer/Pulse Oximeter, English



**MIDMARK IQMARK™ DIGITAL SPIROMETER (SERIAL VERSION)**

Includes: IQmark™ Digital Spirometer (serial version), IQmark™ Digital Spirometer Software Package with trending, IQmark™ Diagnostic Workstation - Patient Management Software, 10-pack disposable pneumotach mouthpieces, 10-pack disposable nose clips, 4-pack AAA batteries, mouse pad, operation manuals and quick reference guide.

4-000-0001 IQmark™ Digital Spirometer (Serial Version)



**MIDMARK IQMARK™ DIGITAL SPIROMETER (SERIAL VERSION) W/3-LITER CALIBRATION SYRINGE**

Includes: IQmark™ Digital Spirometer (Serial Version), IQmark™ Digital Spirometer Software Package with Trending, IQmark™ Diagnostic Workstation - Patient Management Software, IQmark™ 3-Liter Calibration Syringe, 10-pack disposable pneumotach mouthpieces, 10-pack disposable nose clips, 4-pack AAA batteries, mouse pad, operation manuals and quick reference guide.

4-000-0002 IQmark™ Digital Spirometer (Serial Version) w/IQmark™ 3-liter calibration syringe



**MIDMARK IQMARK™ DIGITAL SPIROMETER (USB VERSION)**

Includes: IQmark™ Digital Spirometer (USB Version), IQmark™ Digital Spirometer Software Package with trending, IQmark™ Diagnostic Workstation - Patient Management Software, 10-pack disposable pneumotach mouthpieces, 10-pack disposable nose clips, mouse pad, operation manuals and quick reference guide.

4-000-0004 IQmark™ Digital Spirometer (USB Version)



**MIDMARK IQMARK™ DIGITAL SPIROMETER (USB VERSION) W/3-LITER CALIBRATION SYRINGE**

Includes: IQmark™ Digital Spirometer (USB Version), IQmark™ Digital Spirometer Software Package with Trending, IQmark™ Diagnostic Workstation - Patient Management Software, IQmark™ 3-Liter Calibration Syringe, 10-pack disposable pneumotach mouthpieces, 10-pack disposable nose clips, mouse pad, operation manuals and quick reference guide.

4-000-0003 IQmark™ Digital Spirometer (USB Version)

**MIDMARK IQMARK™ DIGITAL SPIROPDA (PDA VERSION) W/3-LITER CALIBRATION SYRINGE**

Includes: IQmark™ Digital Spirometer (PDA Version), IQmark™ SpiroPDA Digital Spirometer Software Package, IQmark™ Digital Spirometer Software Package with Trending, IQmark™ PDA/IQmark™ Sync Software Package, IQmark™ Diagnostic Workstation - Patient Management Software, IQmark™ 3-Liter Calibration Syringe, 10-pack disposable pneumotach mouthpieces, 10-pack disposable nose clips, PDA Cable (PDA brand and model number required), USB to serial port adapter, 4-pack AAA batteries, mouse pad, operation manuals and quick reference guide (Does NOT include a Pocket-PC device, i.e. an iPAQ).

4-000-0007 IQmark™ Digital SpiroPDA (PDA Version) w/IQmark™ 3-Liter Calibration Syringe

**MIDMARK IQMARK™ SPIROMETRY SUPPLIES**

- 1-100-0003 IQmark™ 3-Liter Calibration Syringe
- 2-100-1200 Disposable Pneumotach Mouthpieces, 25/bx
- 2-244-0001 Disposable Nose Clips, 25/bx
- 3-100-8500 Syringe Adapter
- 3-999-0005 16-pack AAA Batteries
- 2-100-1201 Disposable Pneumotach Mouthpieces, 25/bx, 8 bx/cs

**MIDMARK IQMARK™ MANUALS**

3-100-1025 IQmark™ Digital Spirometer Manual



**PERMAXIM REDISCREEN™ STREP A TEST THROAT SWAB**

97% sensitivity. Equivalent results to a traditional swab culture. Add sample, wait five minutes and read results. Kit includes reagents, controls, tubes, swabs, and test strips. Remove strip from foil pouch. Add reagents to tube. Swab throat and add to reagent solution. Place test strip into solution. Read results in five minutes.

2033 Strep A Test Throat Swab Dipstick, CLIA Waived, 25/bx



**QUIDEL QUICKVUE® DIPSTICK STREP A TEST**

A rapid qualitative test for the detection of Strep A antigen directly from a throat swab. Results in 5 minutes. Room temperature storage CLIA waived.

20108 Dipstick Strep A Test, 50 test/kit



**QUIDEL QUICKVUE® IN-LINE® ONE-STEP GROUP A STREP KIT**

Allows for the rapid detection of group A streptococcal antigen directly from patient throat swab specimens. The unique in-line extraction and rapid procedure allow testing to take place at the point of care. Two-color results give you clear readability, and the proven accuracy of QuickVue provides dependable results. The kit stores at room temperature and includes positive and negative external controls. CLIA waived.

0343 Strep A Test, 25 test kit, 12 kits/cs  
 0345 Strep A Control Swabs, 1/kit, 10 kits/cs  
 0347 Strep A Swab Pack, 75 swabs/pk



**PERMAXIM REDISCREEN™ URINE HCG PREGNANCY CASSETTE DEVICE**

99% clinical accuracy. 25 mIU/ml sensitivity. Some positive results can be seen in as soon as one minute. One step, no reagents. Remove cassette from foil pouch and lay on flat surface. Add three drops of urine to the sample well. Read results in three minutes for urine.

2036 Urine hCG Pregnancy Cassette Device, CLIA Waived, 25/bx  
 2036-100 Urine hCG Pregnancy Cassette Device, 25/bx, 4bx/cs



**PERMAXIM REDISCREEN™ URINE HCG PREGNANCY TEST STRIPS**

99% clinical accuracy. 25 mIU/ml sensitivity. Some positive results can be seen as soon as one minute. One step, no reagents. Remove strip device from canister. Immerse the strip in urine sample. Read results in three minutes.

2032 Urine hCG Pregnancy Test Strips, CLIA Waived, 100/bx



**PERMAXIM REDISCREEN™ URINE/SERUM HCG PREGNANCY CASSETTE DEVICE**

99% clinical accuracy. 25 mIU/ml sensitivity. Some positive results can be seen in as soon as one minute. One step, no reagents. Remove cassette from foil pouch and lay on flat surface. Add three drops of urine or serum to the sample well. Read results in three minutes for urine and five minutes for serum.

2012 Urine/Serum hCG Pregnancy Cassette Device, CLIA Waived, 25/bx



**QUIDEL RAPIDVUE® HCG KIT**

Easy to interpret dipstick. Results in 3 minutes. Each test is individually foil wrapped. Room Temperature storage. CLIA waived.

0257 HCG Kit, 25 test/kit, 10 kits/cs



0317



20183

**QUIDE L QUICKVUE® INFLUENZA A/B TESTS**

Designed specifically for exam room use providing test results in ten minutes or less. Test accurately detects Influenza type A and B viral antigens in three easy steps from a nasal swab or nasal aspirate/nasal wash specimen. The A+B test identifies type A or type B in positive samples. The easy-to-use, easy-to-read kit can be stored at room temperature. No instrumentation is required. The 10-test kit is conveniently packaged for single-use testing and disposal. The 25-test kit offers bulk packaging that is ideal for batch testing. Both the 10-test and 25-test kits include everything to perform the test including positive and negative controls.

- 0317 Influenza A/B Test, Identifies Positive or Negative, 25 test/kit, 10 kits/cs
- 0324 Influenza A/B Test, Identifies Positive or Negative, 10 test/kit, 10 kits/cs
- 20171 Transport Swab Set, 25 swabs/set, 40 set/pk
- 20183 QuickVue® Influenza A+B Test, Dipstick Format, Identifies Type A, Type B, or Both, Two-Color Endpoint, 25 test/kit, 10 kits/cs



2080



2081

**PERMAXIM REDISCREEN™ ULTRA GUAIAIC OCCULT BLOOD TEST/ACCESSORIES**

- 2080 Patient Pack, comes with 100 Triple Slide Tests, 6 x 10mL Developers, Applicators, Foil Lined Envelopes, CLIA Waived, 50/pk, 2 pk/bx
- 2081 Lab Pack, comes with 100 Single Slide Tests, 2 x 10mL Developers, Applicators, CLIA Waived, 100 test/bx
- 2084 Office Pack, comes with 34 Triple Slide Tests, 2 x 10mL Developers, Applicators, 34/bx

**Accessories**

- 2082 8 x 10mL Developers and Applicators, CLIA Waived, 1/bx
- 2083 2 x 100 Test Tape, 2 x 10mL Developers, Applicators, CLIA Waived, 200/bx



**QUIDE L QUICKVUE® IFOB TEST**

- 20194 Kit-Cassette Format, 20 Test Kit, Includes: 1 Pkg Insert, 1 Procedure Card, 20 Specimen Collection & return Kits Each Containing 1 Collection Tube w/2mL FOB Buffer, 1 Collection Paper w/adhesive, 1 Pouch, 1 Absorbent Sleeve, 1 Return Mail Pouch & 1 Patient Instructions
- 20196 Specimen and Collection Kit, 10 Specimen & Return Kits Each Containing 1 Collection Tube w/2mL FOB Buffer, 1 Collection Paper w/adhesive, 1 Pouch, 1 Absorbent Sleeve, 1 Return mail Pouch & 1 Patient Instructions
- 20197 External Controls, Positive & Negative



555001




555050

**HYOGUARD ASSURE™ 3 BLOOD GLUCOSE MONITORING SYSTEM**

Makes testing three times faster. Uses biosensor technology. Calibration is easy - simply insert the code chip and the meter is ready to test. Testing is easy - insert a test strip and the meter turns on. Apply sample and obtain plasma-referenced result in 10 seconds. Uses a specially designed test strip. The test strip with ergonomic. Guide-me-Curve™ facilitates blood application. Capillary action wicks the sample into the test chamber. Features: smaller sample - uses only 3ul blood sample, automatic on/off with insertion and removal of test strip, and wide test range 30-550 mg/dL.

- 555001 Assure 3 Meter Only
- 555006 Assure 3 Control Solution Level 1&2
- 555008 Assure 3 Check Strip
- 555050 Assure 3 Test Strips w/code chip, 50/btl
- 555100 Assure 3 Test Strips w/code chip, 100/btl



**SURGILANCE™ SAFETY LANCETS**

Needle is concealed before and after use. It is automatically rendered inoperative after a single use. No arming is required and the lancets are color-coded for easy recognition. The lancet is triggered by a gentle push on the test site area, giving optimum blood flow at shallow depths. High-speed penetration minimizes pain. The hidden needle and quiet operation helps increase patient comfort.

- SLB200 Blade, 1.8mm Depth, 75 to 100uL Blood Sample, Green, 100/bx
- SLB250 Blade, 2.3mm Depth, 150 to 200uL Blood Sample, Blue, 100/bx
- SLN200 Needle, 1.8mm Depth, 10 to 20uL Blood Sample, Gray, 100/bx
- SLN240 Needle, 2.2mm Depth, 20 to 40uL Blood Sample, Orange, 100/bx
- SLN300 Needle, 2.8mm Depth, 40 to 60uL Blood Sample, Pink, 100/bx

See Laboratory Section for a complete listing of physician's office lab diagnostic products.